




JAY Your Way is a modification program designed to individualize and customize JAY products. From clinical improvements to colored covers, JAY Your Way can be used to create cushions and backs manufactured to users' unique individual needs. This catalog includes information on customizing JAY bases, fluid inserts, covers and backrests, and using our simple online ordering process.


Cushion Base Modifications

Size Modifications



Reducing Size

Mod Detail: Removing foam from the length or width of the cushion base and sewing the cover to fit



Increasing Size

Mod Detail: Adding dummy foam to the length or width of the cushion base and sewing the cover to fit


Clinical Application: To match the size of the client and/or wheelchair



Increasing/Reducing Height

Mod Detail: Adding/Removing foam to/from cushion base to increase/decrease height

Clinical Application: To accommodate wheelchair fitting, environment needs, and/or client needs, e.g. transfers, ground clearance, foot propulsion



Well Sizing

Mod Detail: Increasing or decreasing the standard, factory set size of the Pelvic Loading Area (PLA)

Clinical Application: To allow for anatomically optimized loading areas when considering proper immersion, pressure reduction and lateral stability


Base and Positioning Modifications



Notch for Rails

Mod Detail: Removing foam along the bottom, lateral edges of the foam base

Clinical Application: To lower the cushion between seat rails to achieve the lowest possible seat-to-floor height



Notch for Posts

Mod Detail: Notching the back corners of the foam base


Clinical Application: To slide cushion between back posts when loading surface is posterior to back posts, to lock in cushion in the case of frequent transfers and to build in some depth growth opportunities



Leg Length Cut (Right/Left)

Mod Detail: Decreasing the length of the foam base and fluid pad (if applicable) on the left or right side and sewing the cover to fit

Clinical Application: To accommodate varying femoral lengths



Contracture Cut

Mod Detail: Undercutting the front of the foam base at an angle


Clinical Application: To clear the calves in cases of tight hamstrings, large calves, or close front thigh to lower leg angles, to accommodate chair configurations, and to accommodate foot propulsion



Scrotal Relief

Mod Detail: Scooping out the front and center of the Pelvic Loading Area (PLA)


Clinical Application: To accommodate organ tissue, reduce discomfort, and increase sitting tolerance



Custom Drop Seat

Mod Detail: Fitting a drop seat to the cushion base


Clinical Application: To minimize seat-to-floor height for foot propelling and transfers



Angle Cut Posterior/Anterior Slope

Mod Detail: Removing foam at an angle from underneath the cushion base


Clinical Application: **Anterior Slope:** To accommodate hip flexion limitations and assist with foot propelling and transfers. **Posterior Slope:** To close seat to back angle (e.g. for postural stability)



Wedge Glue Posterior/Anterior Slope

Mod Detail: Adding a foam wedge underneath the cushion base

Clinical Application: **Anterior Slope:** To accommodate hip flexion limitations and assist with foot propelling and transfers. **Posterior Slope:** To close seat to back angle (e.g. for postural stability)



Medial Support Cut

Mod Detail: Removing the pre-contoured medial support


Clinical Application: To improve positioning when femoral alignment is not neutral (e.g. in windswept postures) or when transfers are inhibited by the medial thigh support



Femoral Relief

Mod Detail: Shaving away the front of PLA to transition from PLA to femoral support area


Clinical Application: To reduce pressure at front of PLA area to accommodate tissue and increase sitting tolerance for users with sensation



Leg Trough Cut

Mod Detail: Scooping away 1" of base material from the front of the PLA forward to the front edge of the femoral support area

Clinical Application: To improve medial and lateral positioning and support of the femurs and to accommodate hip limitations by opening hip angles




Trochanteric Relief

Mod Detail: Removing material at the lateral borders of the PLA, beneath the trochanters

Clinical Application: To reduce pressure at or surrounding the trochanters and evenly distribute load across the whole sitting surface

Positioning Components

Positioning pieces, such as adductors, abductors, and hip guides, can be added to several foam cushion bases to improve sitting posture and positioning. The below positioning components are shown on the JAY J3 Cushion base.



Lateral Thigh



Medial Thigh



Lateral Pelvic

Fluid Insert Modifications



Underfill

Mod Detail: Reducing factory set fluid volume in the JAY Flow fluid pad

Clinical Application: To increase immersion of the pelvis into the fluid to achieve appropriate skin protection and/or lateral stability when bottoming out is not a concern



Overfill

Mod Detail: Increasing factory set fluid volume in the JAY Flow fluid pad


Clinical Application: To redistribute load and reduce pressure on the bony prominences when a client experiences "bottoming out" or discomfort on a standard JAY Flow fluid pad



Accommodating Obliquities

Mod Detail: Applying an overfill to either the right or left side of the JAY Flow fluid pad

Clinical Application: To optimize immersion (underfill), prevent bottoming out (overfill), or substitute for bone or tissue asymmetry (overfill) on one side of the client in the case of a flexible obliquity



Sacral or Scrotal Seal

Mod Detail: Adding a heat seal in the sacral or scrotal area of the JAY Flow fluid pad

Clinical Application: To control the migration of fluid in order to reduce pressure in areas that are at-risk for skin breakdown, pain, or discomfort

Cover Modifications

Each JAY cushion and back comes with a carefully designed cover manufactured from top-of-the-line materials that enhance their clinical effectiveness and comfort. JAY Your Way can be used to add some color or better fit users' unique lifestyles.



Sealed Inner Cover

Mod Detail: The water-resistant, stretchable Sealed Inner Cover reduces surface tension and protects from dirt, moisture and odor build-up.

Clinical Application: To reduce surface tension, to protect foam base from dirt and moisture build-up, and to prevent odor



Custom Outer Covers

Specialty outer covers are available on select JAY cushions and backs. Options include:

- Incontinent Shower Cap Cover
- Incontinent Zipper Cover
- JAY Stretch Zipper Cover
- X-static® Cover with Breathable Spacer Fabric
- Dartex™ Cover with Aquaguard Zipper
- Reverse Dartex™ Cover with Aquaguard Zipper



Front Zipper Pocket

The front zippered pocket makes it convenient to access and store small items, including keys, credit cards, or cell phones.



Colored Outer Cover with Zipper

Select JAY cushions are available with colored outer covers or colored semi-stretch sidebands for added flair! Available color options are red, blue pink, purple and toxic green.

Backrest Modifications



Military Service Patches

To show past service in the military or support for the armed forces, Military Service Patches can be added to the JAY J3 and GO backs. Military Service Patches are also available on Quickie wheelchairs' tension adjustable upholstery!



JAY J3 Back Cover


A rugged desert camouflage cover can be chosen at no charge to replace the JAY J3's standard cover.



Spinal Fluid Pad

Mod Detail: Adding a JAY Flow insert over spinal area of foam back

Clinical Application: To reduce pressure along the spine, which is often at high risk of skin breakdown due to its bony structure



Visco Foam Insert

Mod Detail: Adding a Visco memory foam insert with a cover sewn to match

Clinical Application: To add a comfort layer over the standard supportive foam for increased sitting tolerance and pressure reduction



Painted Back Shell

Select JAY backrests are available with painted back shells to match or complement the wheelchair! Choose from any of the standard or matte frame colors available on Quickie and Zippie manual wheelchairs.

Arm Trough Modifications



Adding Pommel to Hand Pad

Mod Detail: Add ergonomic pommel mount to hand portion of arm trough

Clinical Application: To improve comfort and pressure distribution for improved sitting tolerance

JAY Your Way Cushions Modification Matrix

Cushions																							
Modifications Offered	Additional Details	MOD Codes	Upcharge Y or N	Catalog Reference Page	JAY Basic Cushion	JAY BasicPRO Cushion	JAY GO Cushion	JAY Zip Cushion	JAY Ion Cushion	JAY Union Cushion	JAY LITE Cushion	JAY Soft Combi P Cushion	JAY Care Cushion	JAY Duo Cushion	JAY Easy Cushion	JAY X2 Cushion	JAY J2/J2 Recline Cushion	JAY J2DC Cushion	JAY J2 Plus Cushion	JAY GS Cushion	JAY Fusion Cushion	JAY J3 Cushion	
Size Modifications																							
Increase Length (added to front of cushion)	1"	IL1	N	2								•					•	•		•	•	•	
	2"	IL2	N	2								•								•	•	•	
	3"	IL3	N	2								•								•	•	•	
	4"	IL4	N	2								•									•	•	•
Reduce Length (removed from front of cushion)	1"	RL1	N	2	•	•						•								•	•	•	
	2"	RL2	N	2	•	•						•								•	•	•	
	3"	RL3	N	2								•								•	•	•	
	4"	RL4	N	2	•							•	•							•	•	•	
Increase Width (added to sides of cushion)	1" - 1/2" added to each side	IW1	N	2								•								•	•	•	
	2" - 1" added to each side	IW2	N	2								•								•	•	•	
	3" - 1 1/2" added to each side	IW3	N	2								•								•	•	•	
	4" - 2" added to each side	IW4	N	2								•								•	•	•	
Reduced Width (removed from sides of cushion)	1" - 1/2" removed from each side	RW1	N	2	•	•						•	•							•	•	•	
	2" - 1" removed from each side	RW2	N	2	•							•	•							•	•	•	
	3" - 1 1/2" removed from each side	RW3	N	2								•	•							•	•	•	
	4" - 2" removed from each side	RW4	N	2								•	•							•	•	•	
Increase Height (added to bottom of cushion)	1"	IH1	N	2						•	•						•	•		•	•	•	
	2"	IH2	N	2						•	•						•	•		•	•	•	
	3"	IH3	N	2							•						•	•		•	•	•	
Reduce Height	1/2"		N	2																			
	5/8"	RP	N	2				•															
Increase PLA / Well Sizing	1" - 1/2" increase to each side	IWS1	N	2													•	•	•	•	•	•	
	2" - 1" increase to each side	IWS2	N	2													•	•	•	•	•	•	
	3" - 1 1/2" increase to each side	IWS3	N	2													•	•	•	•	•	•	
	4" - 2" increase to each side	IWS4	N	2													•	•	•	•	•	•	
Reduce PLA / Well Sizing	1" - 1/2" decrease to each side	RWS1	N	2													•	•	•	•	•	•	
	2" - 1" decrease to each side	RWS2	N	2													•	•	•	•	•	•	
	3" - 1 1/2" decrease to each side	RWS3	N	2													•	•	•	•	•	•	
	4" - 2" decrease to each side	RWS4	N	2													•	•	•	•	•	•	
Move PLA / Well	1" Forward	CPLAL1	N																		•	•	
	2" Forward	CPLAL2	N																		•	•	
	3" Forward	CPLAL3	N																		•	•	
Base Modifications																							
Notch for Rails	1.25" from side, 1.5" from bottom	NR	N	2	•	•	•	•	•	•	•	•				•	•	•	•	•	•	•	
Notch for Posts	1.25" from side, 1.5" from back	NP	N	2	•	•	•	•	•	•	•	•				•	•	•	•	•	•	•	
Leg Length Cut Right (removed from front, right of cushion)	1"	LLCR1	N	2	•							•					•	•	•	•	•	•	
	2"	LLCR2	N	2	•	•						•					•	•	•	•	•	•	
	3"	LLCR3	N	2	•							•					•	•	•	•	•	•	
	4"	LLCR4	N	2								•					•	•	•	•	•	•	
Leg Length Cut Left (removed from front, left of cushion)	1"	LLCL1	N	2	•	•						•					•	•	•	•	•	•	
	2"	LLCL2	N	2	•	•						•					•	•	•	•	•	•	
	3"	LLCL3	N	2	•							•					•	•	•	•	•	•	
	4"	LLCL4	N	2								•					•	•	•	•	•	•	
Contracture Cut	(1" x 45°)	CCT	N	2				•	•	•	•	•	•	•	•		•	•	•	•	•	•	
Scrotal Relief Cut		SR	N	2	•	•						•	•	•	•		•	•	•	•	•	•	
Custom Drop Seat		CDS	Y	3					•	•	•	•	•	•	•		•	•	•	•	•	•	
Angle Cut Posterior Slope (front to back slope)	1/2"	ACP12	N	3	•												•	•	•	•	•	•	
	1"	ACP1	N	3								•								•	•	•	
	1 1/2"	ACP112	N	3								•								•	•	•	
Angle Cut Anterior Slope (back to front slope)	1/2"	ACA12	N	3								•					•	•	•	•	•	•	
	1"	ACA1	N	3	•							•					•	•	•	•	•	•	
	1 1/2"	ACA112	N	3								•					•	•	•	•	•	•	
Wedge Glue	Posterior Slope	WGP	N	3				•	•	•	•					•	•	•	•	•	•	•	
Medial Support Cut	Anterior Slope	WGA	N	3				•	•	•	•					•	•	•	•	•	•	•	
Femoral Relief Cut		MSC	N	3		•						•	•	•	•		•	•	•	•	•	•	
		FR	N	3								•	•	•	•		•	•	•	•	•	•	
Leg Trough Cut 1" Deep	Right	LTC1R	N	3	•	•						•					•	•	•	•	•	•	
	Left	LTC1L	N	3	•	•						•					•	•	•	•	•	•	
Trochanteric Relief Cut		TR	N	3																			
Lateral Thigh Positioning Pieces	1" Pair	LT1P	Y	3							•						•	•	•	•	•	•	
	1" Right only	LT1R	Y	3							•						•	•	•	•	•	•	
	1" Left only	LT1L	Y	3							•						•	•	•	•	•	•	
	2" Pair	LT2P	Y	3							•						•	•	•	•	•	•	
	2" Right only	LT2R	Y	3							•						•	•	•	•	•	•	
	2" Left only	LT2L	Y	3							•						•	•	•	•	•	•	
Medial Thigh Positioning Piece	1"	MT1	Y	3							•						•	•	•	•	•	•	
	2"	MT2	Y	3							•						•	•	•	•	•	•	
Lateral Pelvic Positioning Pieces	1" Pair	LP1P	Y	3							•						•	•	•	•	•	•	
	1" Right only	LP1R	Y	3							•						•	•	•	•	•	•	
	1" Left only	LP1L	Y	3							•						•	•	•	•	•	•	
	2" Pair	LP2P	Y	3							•						•	•	•	•	•	•	
	2" Right only	LP2R	Y	3							•						•	•	•	•	•	•	
	2" Left only	LP2L	Y	3							•						•	•	•	•	•	•	
Visco Foam Overlay		VFO																				•	
DO NOT Install Positioning Pieces (pieces ship separately)		DNIP	N															•	•				



Bullet indicates that the modification is available on the product



Gray box indicates that the modification is not available on the product.

Please consult the JAY Price List for upcharges.

		Cushions																						
Modifications Offered	Additional Details	MOD Codes	Upcharge Y or N	Catalog Reference Page	JAY Basic Cushion	JAY BasicPRO Cushion	JAY GO Cushion	JAY Zip Cushion	JAY Ion Cushion	JAY Union Cushion	JAY LITE Cushion	JAY Soft Combi P Cushion	JAY Care Cushion	JAY Duo Cushion	JAY Easy Cushion	JAY X2 Cushion	JAY J2/J2 Recline Cushion	JAY J2DC Cushion	JAY J2 Plus Cushion	JAY GS Cushion	JAY Fusion Cushion	JAY J3 Cushion		
Fluid Insert Modifications																								
Underfills	Up to 20%	UF	N	4																				
Overfills	Up to 40%	OF	N	4																				
Obliquity Underfills	Left or Right Side, up to 20%	OUF	N	4																				
Obliquity Overfills	Left or Right Side, up to 40%	OOF	N	4																				
Sacral Seal		SAS	N	4																				
Scrotal Seal		SCS	N	4																				
Cover Modifications																								
Custom Outer Cover	Sealed Inner Cover (Ion Style)	ZIC	Y	4																				
	Incontinent Showercap Outer Cover (J2 Style)	JISC	N	4																				
	Incontinent Zipper Outer Cover (Fusion Style)	JIOC	Y	4																				
	JAY Stretch Zipper Outer Cover (Fusion Style)	JSOC	Y	4																				
	X-static® Outer Cover with Breathable Spacer Fabric (Ion Style)	JBOC	Y	4																				
	Dartex® Outer Cover with Aquaguard Zipper (GO Style)	JDOC	Y	4																				
	Reverse Dartex® OuterCover with Aquaguard Zipper (GO Style)	JRDOC	Y	4																				
Embroidery (Arial Font up to 32 Characters - Single color)		EMBROIDERY	Y	4																				
Embroidery (Custom Artwork from Digital File)		EMBDIG	Y	4																				
Colored Outer Cover with Zipper	Blue (semi stretch)	COCB	Y	4																				
	Red (nylon)	COCR	Y	4																				
	Purple (nylon)	COCPU	Y	4																				
	Pink (semi stretch)	COCPI	Y	4																				
	Blue (semi stretch)	COCBB	Y	4																				
Colored Outer Cover Sidebands	Red (nylon)	COCBR	Y	4																				
	Purple (nylon)	COCBPU	Y	4																				
	Pink (semi stretch)	COCBPU	Y	4																				
	Toxic Green (nylon)	COCBTG	Y	4																				
	Add Front Pocket		FZP	Y	4																			

JAY Your Way Backs Modification Matrix

Backs												
Modifications Offered	Additional Details	MOD Codes	Upcharge Y or N	Catalog Reference Page	JAY Focus/Tail Back	JAY Basic Back	JAY J2 Back	JAY J2DC Back	JAY J2 Plus Back	JAY GO Back	JAY J3 Back	JAY Zip Back
VA Service Patch	Purple Heart	VAPH	N	5								
	Marines	VAM	N	5								
	Army	VAA	N	5								
	Air Force	VAAF	N	5								
	Navy	VAN	N	5								
	Coast Guard	VACG	N	5								
Spinal Fluid Pad		SFPXT	Y	5								
Visco Foam Insert		VF	Y	5								
Scapular Relief Cut		SR	Y	5								
Upside Down Back Shell w/ Head rest		TBSH	N									
Upside Down Back Cover		USDC	N									
Custom Outer Cover	Sealed Inner Cover	ZIC	Y	4								
	Incontinent Zipper Outer Cover	JIOC	Y	4								
	X-static® Cover with Breathable Spacer Fabric	JBOC	Y	4								
	Dartex® Cover with Aquaguard Zipper (Go style)	JDOC	Y	4								
	Reverse Dartex® Cover with Aquaguard Zipper	JRDOC	Y	4								
Embroidery (Arial Font up to 32 Characters - single color)		EMBROIDERY	Y	4								
Embroidery (Custom Artwork from digital file)		EMBDIG	Y	4								
Colored Outer Cover with Zipper (semi stretch)	Blue (semi stretch)	COCB	Y	4								
	Red (nylon)	COCR	Y	4								
	Purple (nylon)	COCPU	Y	4								
	Pink (semi stretch)	COCPI	Y	4								
	Toxic Green (nylon)	COCTG	Y	4								
	Camouflage (Ballistic)	BBCC	N									
Painted Back Shell	All Quickie & Zippie Colors available	CBS	N	5								

JAY Your Way Arm Trough and Foot Box Modification Matrix

Arm Troughs & Foot Boxes					
Modifications Offered	Additional Details	Upcharge Y or N	Catalog Reference Page	Arm Troughs	Foot Boxes
Add Pommel to Hand Pad		N	5	•	•
Gel Overlay		Y		•	•
Visco Overlay		Y		•	•
Add Additional Strap		Y		•	•
Reduce Height of Side Wall	1/2"	N		•	•
	1"	N		•	•
Remove Side Walls		N		•	•
Colored Side Bands	Blue	Y		•	•
	Red	Y		•	•
	Purple	Y		•	•
	Pink	Y		•	•
	Toxic Green	Y		•	•
	Camouflage	Y		•	•
Contact Surface	Fleece	Y		•	•
	Naugahyde	Y		•	•
	Sheep Skin	Y		•	•
	Microclimatic	Y		•	•
Reduce Width	1"	N		•	•
Increase Leg Trough Size	1/2"	N		•	•
	1"	N		•	•
Reduce Leg Trough Size	1/2"	N		•	•
	1"	N		•	•

JAY Your Way Online Ordering Process for Dealers

Select JAY cushions and backs are now easier than ever to order with JAY Your Way modifications using our online ordering process! Follow these simple steps to ensure you receive a unique JAY product that is manufactured to your exact specifications. For an updated list of cushions available to order online with modifications, please visit the JAY Your Way product page on www.sunrisemedical.com.

- Step 1:** Consult the JAY Your Way Catalog and JAY Your Way Price List to select the modifications you need.
- Step 2:** Visit www.sunrisemedical.com
- Step 3:** Log on to your My Sunrise account.*
- Step 4:** Begin an order or quote.
- Step 5:** Enter an "M", followed by the item number of the cushion or back you would like to modify.
- Step 6:** Select from drop-down menus the modifications you require.



JAY Your Way
Drop Down Menus

*To set up a My Sunrise account, visit www.sunrisemedical.com, click "Sign up for an Online Account," and follow the prompts.