



GT / GTI

## Quickie® GT & GTi

**Supplier:** This manual must be given to the rider of this wheelchair.

**Rider:** Before using this wheelchair read this entire manual and save for future reference.

Owner's Manual

## Quickie® GT y GTi Titanium

**Distribuidor:** Este manual debe ser entregado al pasajero de esta silla de ruedas.

**Pasajero:** Antes de usar esta silla de ruedas, lea este manual en su totalidad y guárdelo para futura referencia.

Manual de Instrucciones

Cada una de las sillas se envía con un manual de instrucciones en inglés. El manual en español o francés está disponible en formato PDF en nuestra página en Internet: [www.SunriseMedical.com](http://www.SunriseMedical.com). Ingrese a la página del producto específico para descargar el manual, o comuníquese con el proveedor autorizado de Sunrise Medical.

## Quickie® GT et GTi Titanium

**Au Fournisseur :** Ce manuel doit être remis à l'utilisateur de ce fauteuil roulant.

**À l'Usager :** Veuillez lire tout le manuel et le conserver pour référence ultérieure.

Chaque fauteuil est livré avec un manuel d'instructions en anglais. Les versions en espagnol et en français sont à votre disposition en format PDF sur le site : [www.SunriseMedical.com](http://www.SunriseMedical.com). Veuillez vous rendre à la page de votre produit pour télécharger le manuel dans la langue souhaitée, ou contactez un fournisseur agréé Sunrise Medical.

Manuel d'instructions



**SUNRISE MEDICAL LISTENS**

Thank you for choosing a Quickie wheelchair. We want to hear your questions or comments about this manual, the safety and reliability of your chair, and the service you receive from your Sunrise supplier. Please feel free to write or call us at the address and telephone number below:

**Sunrise Medical  
Customer Service Department  
2842 Business Park Ave  
Fresno, CA 93727  
(800) 333-4000**

Be sure to return your warranty card, and let us know if you change your address. This will allow us to keep you up to date with information about safety, new products and options to increase your use and enjoyment of this wheelchair. If you lose your warranty card, call or write and we will gladly send you a new one.

**FOR ANSWERS TO YOUR QUESTIONS**

Your authorized supplier knows your wheelchair best, and can answer most of your questions about chair safety, use and maintenance. For future reference, fill in the following:


Supplier: \_\_\_\_\_

Address: \_\_\_\_\_  
\_\_\_\_\_

Telephone: \_\_\_\_\_

Serial #: \_\_\_\_\_ Date/Purchased: \_\_\_\_\_

|  |    |  |           |
|--|----|--|-----------|
| I. INTRODUCTION .....                            | 2  | I. Push Handles Option.....                      | 13        |
| II. TABLE OF CONTENTS .....                      | 3  | J. Quick-Release Axles .....                     | 13        |
| III. YOUR CHAIR AND ITS PARTS .....              | 4  | K. Rear Wheels .....                             | 13        |
| IV. NOTICE - READ BEFORE USE .....               | 5  | L. Rear Wheel Locks .....                        | 13        |
| A. Choose the Right Chair & Safety Options ..... | 5  | M. Modified Seat Systems .....                   | 13        |
| B. Review This Manual Often .....                | 5  | N. Seat and Back Upholstery .....                | 13        |
| C. Warnings .....                                | 5  | O. Rear Suspension .....                         | 13        |
| V. GENERAL WARNINGS.....                         | 5  | <b>VIII. USE AND MAINTENANCE .....</b>           | <b>14</b> |
| A. Weight Limits .....                           | 5  | A. Introduction .....                            | 14        |
| B. Intended Use .....                            | 5  | B. Critical Maintenance Tips.....                | 14        |
| C. Attendants and Caregivers.....                | 5  | C. Maintenance Chart .....                       | 14        |
| D. Accessories .....                             | 5  | D. Cleaning Tips.....                            | 14        |
| E. Know Your Chair .....                         | 5  | E. Troubleshooting Chart.....                    | 14        |
| F. Reduce The Risk of an Accident .....          | 6  | F. To Mount & Remove Rear Wheels .....           | 15        |
| G. Safety Checklist .....                        | 6  | G. Wheel Locks.....                              | 15        |
| H. Changes & Adjustments .....                   | 6  | H. Anti-tips Standard & Swing-away .....         | 15        |
| I. Environmental Conditions .....                | 6  | I. Padded Swing-away Armrests .....              | 16        |
| J. Terrain .....                                 | 6  | J. Height Adjustable Armrests.....               | 16        |
| K. Street Use .....                              | 6  | K. Back Angle Performance .....                  | 16        |
| L. Motor Vehicle Safety .....                    | 6  | L. Cushion Installation.....                     | 16        |
| M. When You Need Help .....                      | 7  | M. Storage Tips.....                             | 16        |
| VI. SAFETY WARNINGS: FALLS & TIP-OVERS .....     | 7  | N. Check-Out .....                               | 16        |
| A. Center of Balance .....                       | 7  | <b>IX. DEALER SERVICE &amp; ADJUSTMENT .....</b> | <b>17</b> |
| B. Dressing or Changing Clothes .....            | 7  | A. Dealer Service Introduction.....              | 17        |
| C. Obstacles .....                               | 7  | B. Critical Maintenance Tips.....                | 17        |
| D. Front Caster Lift.....                        | 7  | C. Cleaning.....                                 | 17        |
| E. Reaching or Leaning.....                      | 8  | D. Rear Axle .....                               | 17        |
| F. Moving Backward.....                          | 8  | E. GT With XTR Suspension Rear Axle .....        | 19        |
| G. Escalators.....                               | 8  | F. Rear wheel Axle Nut Adjustment.....           | 19        |
| H. Ramps, Slopes & Sidehills .....               | 9  | G. Rear Wheel Spacing.....                       | 19        |
| I. Transfer .....                                | 9  | H. Wheel Lock Adjustment .....                   | 20        |
| J. Curbs & Single Steps.....                     | 10 | I. Casters .....                                 | 20        |
| K. Climbing a Curb or Single Step.....           | 10 | J. Standard Anti-tip Tubes.....                  | 21        |
| L. Descending a Curb or Single Step.....         | 10 | K. Anti tip Swing-away .....                     | 21        |
| M. Stairs .....                                  | 11 | L. GT with XTR Suspension Anti-tip Tubes .....   | 21        |
| N. Climbing Stairs .....                         | 11 | M. Padded Swing-away Armrests .....              | 21        |
| O. Descending Stairs .....                       | 11 | N. Height Adjustable Armrests.....               | 22        |
| VII. WARNINGS: COMPONENTS & OPTIONS .....        | 12 | O. Armrest receiver Attachment Clamp .....       | 22        |
| A. Anti-Tip Tubes .....                          | 12 | P. Footrest Height Adjustment.....               | 22        |
| B. Armrests .....                                | 12 | Q. Back Adjustments .....                        | 22        |
| C. Caster Pin Locks.....                         | 12 | R. Seat Sling Adjustment.....                    | 23        |
| D. Cushion & Sling Seats .....                   | 12 | S. Rear Seat Height Adjustment.....              | 24        |
| E. Fasteners .....                               | 12 | T. Rear Axle Plate Suspension .....              | 24        |
| F. Footrests.....                                | 12 | U. XTR Suspension.....                           | 24        |
| G. Pneumatic Tires .....                         | 12 | <b>X. SUNRISE LIMITED WARRANTY .....</b>         | <b>25</b> |
| H. Positioning Belts Option .....                | 13 |  |           |



1. Backrest
2. Tire
3. Seat cushion
4. Seat sling
5. Footplate
6. Caster Leg
7. Caster Fork
8. Caster Wheel
9. Caster Housing
10. Hand rim
11. Wheel rim
12. Rear Wheel hub
13. Quick-release axle
14. Axle Plate
15. Camber Tube
16. Spokes
17. Folding Backrest Bracket

Note: All features may not be available with some chair setups or in conjunction with another chair feature. Please consult your authorized supplier for more information.

### GT/GTi

**Weight**

Chair weight will vary per component options

**Adjustments**

Caster angle, CG, rear seat height, back angle, footrest length, rear wheel spacing

**Color**

Natural Titanium Finish, Quickie colors

**Frame Dimensions**

Frame width: Std 12" - 20"  
(XTR Suspension: 14" - 20")  
Seat depth: 12" - 20"  
(XTR Suspension: 14" - 20")

**Backrest**

Non-folding  
Folding Lock Down Angle Adjustable  
Height-adjustable & angle: 11"-14", 12"-16", 16"-20"  
Fixed Height: 10"-20"  
Opt- adjustable upholstery, vented adjustable upholstery, folding push handles

**Footrest**

Std- Tubular  
Opt- Tubular w/Cover,  
Angle Adjustable Platform  
Aluminum & Titanium

**Casters**

Std- 5" Low-profile Polyurethane  
Opt- 4" Low-profile Poly,  
4" Micro Light Up  
4"x1" Aluminum Poly  
4"x1.5" Aluminum soft Roll  
3" Micro,  
3" Micro Light Up  
4", 5" & 6" x 1-1/2" wide poly  
6" Poly  
Frog Legs

**Armrests**

Opt- Padded Swing-Away  
Opt- Single-Post Height-Adjustable

**Axle Bracket**

Std- Fore/Aft CG Adjustable Bracket  
Opt Camber- 0°, 3°, 6°, 9°  
XTR Camber- 0°, 4°, 8°, amputee

**Axles**

Std- stainless steel  
Opt- quad-release axle nuts

**Anti-tip tubes**

Opt- Rear anti-tip tubes, swing-away anti-tips

**Standard rear wheels**

Std- Lite Spoke  
Opt- Mag

**Quickie Performance wheels**

Opt- Spinergy  
Opt- Ultra lightweight  
Opt- Mountain

**Tire size**

Std- 24"  
Opt- 20", 22", 25", 26"

**Tire types**

Std - Pneumatic  
Opt - Airless insert, Poly, Kevlar (24" only),  
High-pressure clinchers (24",26" only),  
Mountain bike (24" only), Full Poly, Kik,  
Knobby V-TRAK, Schwalbe Marathon  
Plus, Urban Assault

**Hand rims**

Std- aluminum  
Opt- plastic-coated  
vertical/oblique projections  
Natural Fit

**Wheel locks**

High-push, high-pull, ergo scissor  
Opt- compact

**A. CHOOSE THE RIGHT CHAIR & SAFETY OPTIONS**

Sunrise provides a choice of many wheelchair styles to meet your needs. This product is intended for single person use only. Final selection of the type of wheelchair, options and adjustments rests solely with you and your health care professional. Choosing the best chair for your unique mobility needs, will depend on such things as:

1. Your disability, strength, balance and coordination.
2. The types of hazards you must overcome in daily use (where you live and work, and other places you are likely to use your chair).
3. The need for options that will improve your positioning, safety and comfort (such as anti-tip tubes, positioning belts, or special seating systems).

**B. REVIEW THIS MANUAL OFTEN**

Before using this chair you, and each person who may assist you, should read this entire manual and make sure to follow all instructions. Review the warnings often, until they are second nature to you.

**C. WARNINGS **

The word “WARNING” refers to a hazard or unsafe practice that may cause severe injury or death to you or to other persons. The “Warnings” are in three main sections, as follows:

**I. V — GENERAL WARNINGS**

Here you will find a safety checklist and a summary of risks you need to be aware of before you ride this chair.

**2. VI— SAFETY WARNINGS: FALLS & TIP-OVERS**

Here you will learn about practices for the safe use of your chair, and how to avoid a fall or tip-over while you perform daily activities in your chair.

**3. VII — WARNINGS — COMPONENTS & OPTIONS**

Here you will learn about the components of your chair and options you can select for safety. Consult your authorized supplier and your health care advisor to help you choose the best set-up and options for safe use.

**NOTE—** Where they apply, you will also find “Warnings” in other sections of this manual. Heed all warnings in these sections. If you fail to do so a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

## V. GENERAL WARNINGS

**A. WEIGHT LIMITS**** WARNING**

NEVER exceed the weight limit specified by Sunrise Medical. The weight capacity provided by your manufacturer is for the combined weight of a rider and items carried using on-board storage. If you do exceed the weight limit, damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

Weight Capacities

Quickie GT - 265lbs/(120Kg).  
Quickie GTi - 265lbs/(120kg).

**B. INTENDED USE**

The Quickie GT Series of wheelchair's intended use is to provide mobility to persons limited to a sitting position.

** WARNING**

DO NOT use this device for purposes other than what is intended by the manufacturer

1. The wheelchair is not designed for weight training and is unsafe for use as a seat while weight training. Weight training from the wheelchair substantially changes the stability of the chair and may cause tipping.
2. DO NOT stand on the frame of the wheelchair.
3. NEVER allow someone to stand on your chair or use it as a step ladder.
4. This chair is designed for a single rider only.
5. Unauthorized modifications and use of parts or accessories not supplied or approved by Sunrise Medical may change the chair structure. This will void the warranty and may cause a safety hazard. If the warning is ignored, damage to your chair, and the potential severe injury of the person using the chair for unintended purposes can occur.

**C. ATTENDANTS AND CAREGIVERS**** WARNING**

Before you assist a rider, be sure to read all warnings contained in this manual, and follow all instructions that apply. Be aware that after consulting a healthcare advisor, you will need to learn safe and proven body mechanics to use and create assistive methods best suited to your abilities.

**D. ACCESSORIES**** WARNING**

Unauthorized modifications or use of parts, or accessories not supplied or approved by Sunrise Medical may change the chair structure. This will void the warranty and may cause a safety hazard.

Some problems that may occur, but are not limited to:

1. Incorrect Wheels and/or tires that put the rider at risk of a fall or tip-over.
2. Adding a component to the frame, changing the structural integrity of the chair.
3. Any modification or disassembly can potentially create an unsafe situation where rider and/or attendant are put at risk.

**E. KNOW YOUR CHAIR**** WARNING**

Every wheelchair is different. Take the time to learn the feel of this chair before you begin riding. Start slowly, with easy, smooth strokes. If you are used to a different chair, you may use too much force and tip over. If you use too much force, damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**F. REDUCE THE RISK OF AN ACCIDENT****⚠ WARNING**

1. BEFORE you begin riding, you should be trained in the safe use of this chair by your health care advisor.
2. Practice bending, reaching and transfers until you know the limit of your ability. Have someone help you until you know what can cause a fall or tip-over and how to avoid doing so.
3. Be aware that you must develop your own methods for safe use best suited to your level of function and ability.
4. NEVER try a new maneuver on your own. Ask the advice of your health care provider to lower the risk of a fall or tip over.
5. Get to know the areas where you plan to use your chair. Look for hazards and learn how to avoid them.

**G. SAFETY CHECKLIST****⚠ WARNING**

Before Each Use Of Your Chair:

1. Make sure the chair rolls easily and that all parts work smoothly. Check for noise, vibration, or a change in ease of use. (This may indicate low tire pressure, loose fasteners, or damage to your chair).
2. Immediately repair any problem. Your authorized supplier can help you find and correct the problem.
3. Check to see that both quick-release rear axles are locked. When locked, the axle button will "pop out" fully. If not locked, the wheel may come off and cause you to fall.
4. If your chair has anti-tip tubes, lock them in place.

**H. CHANGES & ADJUSTMENTS****⚠ WARNING**

1. See your healthcare advisor and have them adjust seating components any time a change or adjustment needs to be made.
2. Unauthorized modifications or use of parts not supplied or approved by Sunrise may change the chair structure. This will void the warranty and may cause a safety hazard.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**I. ENVIRONMENTAL CONDITIONS****⚠ WARNING**

1. Use extra care if you must ride your chair on a wet or slick surface. If you are in doubt, ask for help.
2. Contact with water or excess moisture may cause your chair to rust or corrode. Avoid all extreme weather situations if possible.
  - a. Do not use your chair in a shower, pool or other body of water. The chair tubing and parts are not water-tight and may rust or corrode from the inside.
  - b. Avoid excess moisture (for example, do not leave your chair in a damp bathroom while taking a shower).
  - c. Dry your chair as soon as you can if it gets wet, or if you use water to clean it.



If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**J. TERRAIN****⚠ WARNING**

1. Based on ANSI/RESNA testing, Sunrise Medical recommends the use of a caster wheel with a minimum diameter of 5", if the wheelchair will be overcoming obstacles up to 1/2" on a regular basis.
2. Your chair is designed for use on firm, even surfaces such as concrete, asphalt, indoor flooring, and carpets.
3. Do not operate your chair in sand, loose soil, or over rocky terrain.
4. If you use your chair on terrain that is rougher than described above there is a danger that screws and bolts will loosen prematurely, and that damage to wheels or axles could put the rider at risk of a fall, tip-over, or loss of control.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**K. STREET USE****⚠ WARNING**

1. This product is not intended for street use.
2. Avoid streets whenever possible.
3. Obey and follow all legal pedestrian pathways, and laws that apply to pedestrians.
4. Be alert to the danger of motor vehicles in parking lots, or if you must cross a road.

If you fail to heed this warning, damage to your chair, and a possible collision may occur and cause severe injury to the rider or others.

**L. MOTOR VEHICLE SAFETY**

If possible and feasible, the rider should transfer to the Original Equipment Manufacturer vehicle seat and use the OEM vehicle restraint.

**⚠ WARNING**

If your chair is **NOT** equipped with the Transit Option:

1. NEVER let anyone sit in this chair while in a moving vehicle.
  - a. ALWAYS move the rider to an approved vehicle seat.
  - b. ALWAYS secure the rider with proper motor vehicle restraints.
2. In an accident or sudden stop the rider may be thrown from the chair. Wheelchair seat belts will not prevent this, and further injury may result from the belts or straps.
3. NEVER transport this chair in the front seat of a vehicle. It may shift and interfere with the driver.
4. ALWAYS secure this chair so that it cannot roll or shift.
5. Do not use any chair that has been involved in a motor vehicle accident. The frame and/or components may have been changed due to the accident. Such items could be, but are not limited to: bent, loosened, and/or broken components that were subjected to an impact.

Failure to heed these warnings puts you and others that may assist you at a high risk of injury.

**M. WHEN YOU NEED HELP****⚠ WARNING**

**For The Rider:** Make sure that each person who helps you reads and follows all warnings and instructions that apply.

**For Attendants:**

1. Work with the rider's doctor, nurse or therapist to learn safe methods best suited to your abilities and those of the rider.
2. Tell the rider what you plan to do, and explain what you expect the rider to do. This will put the rider at ease and reduce the risk of an accident.
3. Make sure the chair has push handles. They provide secure points for you to hold the rear of the chair to prevent a fall or tip-over. Check to make sure push handle grips will not rotate or slip off.

4. To prevent injury to your back, use good posture and proper body mechanics. When you lift or support the rider or tilt the chair, bend your knees slightly and keep your back as upright and straight as you can.
5. Remind the rider to lean back when you tilt the chair backward.
6. When you descend a curb or single step, slowly lower the chair in one easy movement. Do not let the chair drop the last few inches to the ground. This may damage the chair or injure the rider.
7. To avoid tripping, unlock and rotate anti-tip tubes up, out of the way.
8. Whenever you aren't attending the wheelchair, ALWAYS use the wheel-lock to secure the rear wheels, and lock the anti-tip tubes in place.

If you fail to ask for help when in doubt, you run a high risk of a fall, tip-over or loss of control that may occur and cause severe injury to the rider or others.

**VI. SAFETY WARNINGS: FALLS & TIP-OVERS****.A. CENTER OF BALANCE****⚠ WARNING**

The point where this chair will tip forward, back or to the side depends on its center of balance and stability. How your chair is set up, the options you select and the changes you make may affect the risk of a fall or tip-over.

1. The Most Important Adjustment Is:  
The position of the rear wheels. The more you move the rear wheels forward, the more likely your chair will tip over backward.
2. The Center Of Balance Is Also Affected By:
  - a. A change in the set-up of your chair, including:
    - The distance between the rear wheels.
    - The amount of rear wheel camber.
    - The seat height and seat angle.
    - Backrest angle.
  - b. A change in your body position, posture or weight distribution.
  - c. Riding your chair on a ramp or slope.
  - d. A back pack or other options and the amount of added weight.
3. To Reduce The Risk Of An Accident:
  - a. Consult your doctor, nurse or therapist to find out what axle and caster position is best for you.
  - b. Consult your authorized supplier BEFORE you modify or adjust this chair. Be aware that you may need to make other changes to correct the center of balance.
  - c. Have someone help you until you know the balance points of your chair and how to avoid a tip-over.
  - d. Use anti-tip tubes.

If you fail to heed these warnings, you are at a high risk of a fall, tip-over or loss of control that could cause severe injury to the rider or others.

**B. DRESSING OR CHANGING CLOTHES****⚠ WARNING**

Your weight may shift if you dress or change clothes while seated in this chair.

To reduce the risk of a fall or tip-over:

1. Rotate the front casters until they are as far forward as possible. This makes the chair more stable.
2. Lock anti-tip tubes in place. (If your chair does not have anti-tip tubes, back it up against a wall and lock both rear wheels).

If you fail to heed this warning, damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**C. OBSTACLES****⚠ WARNING**

Riding over curbs or obstacles can cause tipping and serious bodily harm. If you have any doubt that you can safely cross any curb or obstacle, ALWAYS ASK FOR HELP. Be aware of your riding skills and personal limitations. Develop new skills only with the help of a companion.

To avoid the risk of traversing obstacles:

1. Keep a lookout for danger – scan the area well ahead of your chair as you ride.
2. Make sure the floor areas where you live and work are level and free of obstructions.
3. Remove or cover threshold strips between rooms.
4. Install a ramp at entry or exit doors. Make sure there is not a drop off at the bottom of the ramp.
5. To Help Correct Your Center Of Balance:
  - a. Lean your upper body FORWARD slightly as you go UP over an obstacle.
  - b. Press your upper body BACKWARD as you go DOWN from a higher to a lower level.
6. If your chair has anti-tip tubes, do not go over an obstacle without help.
7. Keep both of your hands on the handrims as you go over an obstacle.
8. Never push or pull on an object (such as furniture or a doorjamb) to propel your chair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**D. FRONT CASTER LIFT****⚠ WARNING**

Front caster lift can occur when the front wheels, ordinarily in contact with the ground, are either intentionally, or unintentionally caused to lift from the ground while the rear wheels remain in contact. Pitch control (partial or full) should NEVER be attempted without consulting your healthcare advisor and making sure that anti-tips are installed.

If you fail to heed this warning, you are at a high risk of damage to your chair, a fall, tip-over or loss of control that may occur and cause severe injury to the rider or others.

**E. REACHING OR LEANING**

**⚠ WARNING**

If you reach or lean it will affect the center of balance of your chair. This may cause you to fall or tip over. When in doubt, ask for help or use a device to extend your reach.

1. NEVER reach or lean if you must shift your weight sideways or rise up off the seat.
2. NEVER reach or lean if you must move forward in your seat to do so. Always keep your buttocks in contact with the backrest.
3. NEVER reach with both hands (you may not be able to catch yourself to prevent a fall if the chair tips).
4. NEVER reach or lean to the rear unless your chair has anti-tip tubes locked in place.
5. DO NOT reach or lean over the top of the seat back. This may damage one or both backrest tubes and cause you to fall.
6. If You Must Reach Or Lean:
  - a. Do not lock the rear wheels. This creates a tip point and makes a fall or tip-over more likely.
  - b. Do not put pressure on the footrests.
  - c. Move your chair as close as you can to the object you wish to reach.
  - d. Do not try to pick up an object from the floor by reaching down between your knees. You are less likely to tip if you reach to the side of your chair.
  - e. Rotate the front casters until they are as far forward as possible. This makes the chair more stable.
  - f. Firmly grasp a rear wheel or an armrest with one hand. This will help to prevent a fall if the chair tips.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**F. MOVING BACKWARD**

**⚠ WARNING**

Use extra care when you move your chair backward. Your chair is most stable when you propel yourself forward. You may lose control or tip over if one of the rear wheels hits an object and stops rolling.

1. Propel your chair slowly and smoothly.
2. If your chair has anti-tip tubes, make sure to lock them in place.
3. Stop often and check to be sure your path is clear.

If you fail to heed this warning, you are at a high risk of damage to your chair, a fall, tip-over or loss of control that may occur and cause severe injury to the rider or others.

**G. ESCALATORS**

**⚠ WARNING**

NEVER use this chair on an escalator, even with an attendant. If you do, a fall or tip-over is likely.

**H. RAMPS, SLOPES & SIDEHILLS** 3 4 5 6


**⚠ WARNING**

Whenever possible, avoid riding on a slope, which includes a ramp or side-hill. This will change the center of balance of your chair. Your chair is less stable and more difficult to maneuver when it is at an angle. When moving up a hill, anti-tip tubes may not prevent a fall or tip-over.


1. Do not use your chair on a slope steeper than 10%. (A 10% slope means: one foot in elevation for every ten feet of slope length.)
2. Always go as straight up and as straight down as you can. (Do not "cut the corner" on a slope or ramp.)
3. Do not turn or change direction on a slope.

4. When you have to use a ramp, always stay in the CENTER of the ramp. Make sure ramp is wide enough that you have no risk of going off the edge.
5. Do not stop on a steep slope. If you stop, you may lose control of your chair.
6. NEVER use rear wheel locks to try to slow or stop your chair. This is likely to cause your chair to veer out of control.
7. ALWAYS be aware of:
  - a. Wet or slippery surfaces.
  - b. A change in grade on a slope (or a lip, bump or depression). These may cause a fall or tip-over.
  - c. A drop-off at the bottom of a slope or ramp. A drop-off as small as 3/4 inch can stop a front caster and cause the chair to tip forward.
8. To Reduce The Risk Of A Fall Or Tip-Over:
  - a. Lean or press your body UPHILL. This will help adjust for the change in the center of balance caused by the slope or sidehill.
  - b. Keep pressure on the handrims to control your speed on a down slope. If you go too fast you may lose control.
  - c. Ask for help any time you are in doubt.
9. Ramps At Home & Work— For your safety, ramps at home and work must meet all legal requirements for your area:
  - a. Avoid A Drop Off. You may need a section at the top or bottom to smooth out the transition.
  - b. Always stay in the center of the ramp and control your speed.


If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and can possibly cause severe injury to the rider or others.




- IF you need to go down a hill:
- a. The downslope should be less than 10°
  - b. GO SLOW
  - c. ALWAYS go straight down
  - d. ALWAYS control your speed
  - e. DO NOT Turn on the hill.



- IF you need to go up a hill:
- a. The Upslope should be less than 10°
  - b. ALWAYS go straight up
  - c. NEVER stop on the upslope
  - d. NEVER turn on a hill.
  - e. ALWAYS lean forward.



- IF you need to traverse a sidehill:
- a. The sidehill must be less than 6° or you risk a tip-over
  - b. DON'T STOP
  - c. ALWAYS go straight down
  - d. DO NOT Turn on the hill.



- IF you need to traverse a ramp:
- a. Go slowly
  - b. Watch out for Drop-offs.
  - c. ALWAYS go straight down.
  - d. Stay centered on the ramp.
  - e. ALWAYS lean forward.

I. TRANSFERS 7 8

**⚠ WARNING**


It is dangerous to transfer on your own. It requires good balance and agility. Be aware that there is a point during every transfer when the wheelchair seat is not below you. To avoid a fall:

1. Work with your health care advisor to learn safe methods for transfers.
  - a. Learn how to position your body and how to support yourself during a transfer.
  - b. Have someone help you until you learn safe transfer methods.
2. Lock the rear wheels before you transfer.
3. Be aware that the chair can still slide and/or tip. The wheel lock keeps the rear wheels from rolling while you are performing the transfer.
4. Make sure that the pneumatic tires are properly inflated. Low tire pressure may allow the rear wheel locks to slip. (see table in Section G "Pneumatic Tires")
5. Move your chair as close as you can to the seat you are transferring to. If possible, use a transfer board.
6. Rotate the front casters until they are as far forward as possible.
7. If you can, remove the footrests, or swing them out of the way.
  - a. Make sure your feet do not catch in the space between the footrests.
  - b. Avoid putting weight on the footrests as this may cause the chair to tip.
8. Make sure armrests are removed, or out of the way and do not interfere with the transfer.
9. Transfer as far back onto the seat surface as you can. This will reduce the risk that the chair will tip or move away from you.

Optimum Transfer position



Before transferring, remove the armrests, and footrests, and make sure the wheel lock is set.




**J. CURBS & SINGLE STEPS** 9**⚠ WARNING**

Before riding over curbs, and negotiating even a single step, ALWAYS ask for assistance first. Curbs and steps can cause tipping and serious bodily harm. When in doubt as to your ability to avoid, or traverse any obstacle, always ask for help. Be aware of your riding skills and personal limitations. Develop new skills only with the help of a companion.

**For Attendant:** Each person who assists the rider with curbs and steps should read and follow all instructions and warnings pertaining to attendants, and caregivers.

1. Do not try to climb a high curb or step (more than 4 inches high) UNLESS you have help. Doing so may cause your chair to exceed its balance point and tip over.
2. Go straight up and straight down a curb or step. If you climb or descend at an angle, a fall or tip-over is likely.
3. Be aware that the impact of dropping down from a curb or step can damage your chair or loosen fasteners.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and can possibly cause severe injury to the rider or others.

**K. CLIMBING A CURB OR SINGLE STEP** 9**⚠ WARNING**

**For Attendant:** follow these steps to help the rider climb a curb or single step going FORWARD:

1. Stay behind the chair.
2. Face the curb and tilt the chair up on the rear wheels so that the front casters clear the curb or step.
3. Move forward, placing the front casters on the upper level as soon as you are sure they are past the edge.
4. Continue forward until the rear wheels contact the face of the curb or step. Lift and roll the rear wheels to the upper level.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**L. DESCENDING A CURB OR SINGLE STEP** 9**⚠ WARNING**

**For Attendant:** Follow these steps to help a rider descend a curb or single step going BACKWARD:

1. Stay at the rear of the chair.
2. Several feet before your reach the edge of the curb or step, turn the chair around and pull it backward.
3. While looking over your shoulder, carefully step back until you are off the curb or stair and standing on the lower level.
4. Pull the chair toward you until the rear wheels reach the edge of the curb or step. Then allow the rear wheels to slowly roll down onto the lower level.
5. When the rear wheels are safely on the lower level, tilt the chair back to its balance point. This will lift the front casters off the curb or step.
6. Keep the chair in balance and take small steps backward. Turn the chair around and gently lower front casters to the ground.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**M. STAIRS** 10**⚠ WARNING**

1. NEVER use this chair on stairs UNLESS you have someone to help you. Doing so is likely to cause a fall or tip-over.
2. Negotiating stairs with a wheelchair always requires at least 2 attendants for safety.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**N. CLIMBING STAIRS** 10**⚠ WARNING**


1. Use at least two attendants to move a chair and rider up stairs.
2. Move the chair and rider BACKWARD up the stairs.
3. The person at the rear is in control. He or she tilts the chair back to its balance point.
4. A second attendant at the front firmly grasps a non-detachable part of the front frame and lifts the chair up and over one stair at a time.
5. The attendants move to the next stair up. Repeat for each stair, until you reach the landing.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**O. DESCENDING STAIRS** 10**⚠ WARNING**

1. Use at least two attendants to move a chair and rider down stairs.
2. Move the chair and rider FORWARD down the stairs.
3. The person at the rear is in control. He or she tilts the chair to the balance point of the rear wheels and rolls it to the edge of the top step.
4. A second attendant stands on the third step from the top and grasps the chair frame. He or she lowers the chair one step at a time by letting the rear wheels roll over the stair edge.
5. The attendants move to the next stair down. Repeat for each stair, until you reach the landing.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.



**A. ANTI-TIP TUBES**



Anti-tip tubes can help keep your chair from tipping over backward in normal conditions.

1. Sunrise recommends the use of anti-tip tubes:
2. When locked in place (in the “down” position) anti-tip tubes should be BETWEEN 1 1/2 to 2 inches off the ground.
  - a. If set **higher** than 2 inches, they may not prevent a tip-over.
  - b. If set **lower** than 2 inches, they may “hang up” on obstacles and cause a fall or tip over.
3. If you have to climb or descend a curb, or overcome an obstacle it may be necessary to have an attendant make sure the anti-tip tubes are rotated up, and out of the way, so that the chair and rider do not get stuck and/or become unstable.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**B. ARMRESTS**



Armrests detach and will not bear the weight of this chair.

1. NEVER lift this chair by its armrests. They may come loose or break.
2. Lift this chair only by non-detachable parts of the main frame.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**C. CASTER PIN LOCKS**



Always use both pin locks when transferring in and out of the chair.

Caster pin locks only keep the front caster wheels from changing direction, they do not lock-out the front wheels from rotation. To prevent the wheelchair from rolling, you must engage the rear wheel locks.

If you propel your wheelchair with the caster pin locks engaged, the chair will be prevented from turning.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**D. CUSHIONS & SLING SEATS**



1. Quickie sling seats and standard foam cushions are not designed for the relief of pressure.
2. If you suffer from pressure sores or if you are at risk that they will occur, you may need a special seat system or a device to control your posture. Consult your doctor, nurse or therapist to find out if you need such a device for your well-being.
3. Seat slings are not intended to be used as a direct seating surface. A cushion or other seating surface should be placed on the sling before use.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**E. OFFSET SEATING (OPTIONAL)**



1. A fore/aft change in the backrest position will affect the center of balance of your chair.
2. Only use this option when installing a forward offset mounted seating solution, not using the standard back sling.
3. Each time you adjust or modify the backrest plane location, the change may make it easier to tip backwards. Use anti-tips

**F. FASTENERS**



Many of the screws, bolts and nuts on this chair are special high-strength fasteners. Use of improper fasteners may cause your chair to fail.

1. ONLY use fasteners provided by an authorized supplier (or ones of the same type and strength, as indicated by the markings on the heads).
2. Over- or under-tightened fasteners may fail or cause damage to chair parts.
3. If bolts or screws become loose, tighten them as soon as you can.

**F. FOOTRESTS**



1. At the lowest point, footrests should be AT LEAST 2 inches off the ground. If set too LOW, they may “hang up” on obstacles you can expect to find in normal use. This may cause the chair to stop suddenly and tip forward.
2. To Avoid A Trip Or Fall When You Transfer:
  - a. Make sure your feet do not “hang up” or get caught in the space between the footrests.
  - b. Avoid putting weight on the footrests, as the chair may tip forward.
  - c. Footrests should be swung out of the way or removed whenever entering or exiting the wheelchair.
3. NEVER lift this chair by the footrests. Footrests detach and will not bear the weight of this chair. Lift this chair only by non-detachable parts of the main frame.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**G. PNEUMATIC TIRES**



Proper inflation extends the life of your tires and makes your chair easier to use.

1. Do not use this chair if any of the tires are under- or over-inflated. Check weekly for proper inflation level, as listed on the tire sidewall.
2. Low pressure in a rear tire may cause the wheel lock on that side to slip and allow the wheel to turn when you do not expect it.
3. Low pressure in any of the tires may cause the chair to veer to one side and result in a loss of control.
4. Over-inflated tires may burst.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**Recommended Tire Pressures**

| Tires                        | T | Recommended Pressure or Pressure Range |        |
|------------------------------|---|--|--------|
|                              |   | kPa                                    | PSI    |
| Pneumatic - 20"              | x | 448                                    | 65     |
| Pneumatic - 22"              | x | 448                                    | 65     |
| Pneumatic - 24"              | x | 517                                    | 75     |
| Pneumatic - 26"              | x | 448                                    | 65     |
| High Pressure Clincher (22") | x | 690                                    | 100    |
| High Pressure Clincher (24") | x | 758                                    | 110    |
| High Pressure Clincher (25") | x | 758                                    | 110    |
| High Pressure Clincher (26") | x | 758                                    | 110    |
| Knobby V-Track (24")         | x | 448                                    | 65     |
| Schwalbe Marathon (24")      | x | 586 - 1000                             | 85-145 |
| Schwalbe Marathon (25")      | x | 586 - 1000                             | 85-145 |
| Schwalbe Marathon (26")      | x | 586 - 1000                             | 85-145 |

**H. POSITIONING BELTS (OPTIONAL)****⚠ WARNING**

Use positioning belts **ONLY** to help support the rider's posture. Improper use of these belts may cause severe injury to or death to the rider.

1. Make sure the rider does not slide down in the wheelchair seat. If this occurs, the rider may suffer chest compression or suffocate due to pressure from the belts.
2. The belts must be snug, but must not be so tight that they interfere with breathing. You should be able to slide your open hand, flat, between the belt and the rider.
3. A pelvic wedge or a similar device can help keep the rider from sliding down in the seat. Consult with the rider's doctor, nurse or therapist to find out if the rider needs such a device.
4. Use positioning belts only with a rider who can cooperate. Make sure the rider can easily remove the belts in an emergency.
5. NEVER Use Positioning Belts:
  - a. As a patient restraint. A restraint requires a doctor's order.
  - b. On a rider who is comatose or agitated.
  - c. As a motor vehicle restraint. In an accident or sudden stop the rider may be thrown from the chair. Wheelchair seat belts will not prevent this, and further injury may result from the belts or straps.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**I. PUSH HANDLES (OPTIONAL)****⚠ WARNING**

When you have an attendant, make sure that this chair has push handles.

1. Push handles provide secure points for an attendant to hold the rear of this chair, to prevent a fall or tip-over.
2. Check to make sure push handle grips will not rotate or slip off.

**J. QUICK-RELEASE AXLES****⚠ WARNING**

1. Do not use this chair **UNLESS** you are sure that both quick-release rear axles are locked.
2. An axle is not locked until the quick-release button pops out fully.
3. An unlocked axle may come off during use, resulting in a fall, tip-over or loss of control and cause severe injury to the rider or others.
  - Quick-Release Axles should be periodically cleaned and inspected for function and signs of wear or bending. Replace as necessary.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**K. REAR WHEELS****⚠ WARNING**

A change in set-up of the rear wheels will affect the center of balance of your chair.

1. The farther you move the rear axles **FORWARD**, the more likely it is that your chair will tip over backward.
2. Consult your healthcare advisor to find the best rear axle set-up for your chair. Do not change the set-up **UNLESS** you consult your healthcare advisor first.
3. Adjust the rear wheel locks after you make any change to the rear axles.
  - a. If you fail to do so, the wheel locks may not work.
  - b. Make sure lock arms embed in the tires at least 1/8 inch when locked.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**L. REAR WHEEL LOCKS****⚠ WARNING**

Rear wheel locks are **NOT** designed to slow or stop a moving wheelchair. Use them only to keep the rear wheels from rolling when your chair is at a complete stop.

1. **NEVER** use rear wheel locks to try to slow or stop your chair when it is moving. Doing so may cause a fall or tip-over
2. To keep the rear wheels from rolling, always set both rear wheel locks when you transfer to or from your chair.
3. Low pressure in a rear tire may cause the wheel lock on that side to slip and may allow the wheel to turn when you do not expect it.
4. Make sure lock arms embed in tires at least 1/8 inch when locked. If you fail to do so, the locks may not work.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**M. MODIFIED SEAT SYSTEMS****⚠ WARNING**

Use of a seat system not approved by Sunrise may alter the center of balance of this chair. This may cause the chair to tip over.

1. Do not change the seat system of your chair **UNLESS** you consult your authorized supplier first.
2. Use of a seating system not approved by Sunrise may affect the folding mechanism of this chair.
3. Use of a seating system not provided by Sunrise is prohibited for transit use.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**N. SEAT AND BACK UPHOLSTERY****⚠ WARNING**

1. Replace worn or torn fabric of seat and seat back as soon as you can. If you fail to do so, the seat or seat back may fail.
2. Sling fabric will weaken with age and use. Look for fraying, thin spots, or stretching of fabrics especially at edges and seams.
3. "Dropping down" into your chair will weaken fabric and result in the need to inspect and replace the seat more often.
4. Be aware that laundering or excess moisture will reduce the flame retardation qualities of the fabric.
5. If Tension Adjustable upholstery is present, loosen straps or remove upholstery prior to folding or unfolding the chair.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

3. Adjust the rear wheel locks after you make any change to the rear axles.
  - a. If you fail to do so, the locks may not work.
  - b. Make sure lock arms embed in tires at least 1/8 inch when locked.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**O. REAR SUSPENSION (OPTIONAL)****⚠ WARNING**

When replacing the suspension bushings, always use genuine Quickie parts. Use of non-Quickie replacement parts voids the warranty and could cause structural failure.

If you fail to heed these warnings damage to your chair, a fall, tip-over or loss of control may occur and cause severe injury to the rider or others.

**⚠ WARNING**

The owner of this chair is responsible for making sure that it has been setup and adjusted by a trained service professional under the advice of a healthcare advisor. The chair may require periodic maintenance or certain in-use adjustments that may be performed by the owner or caregiver. Always use parts and/or accessories that have been recommended or approved by Sunrise Medical when servicing this chair.

**A. INTRODUCTION**

1. Proper maintenance will improve performance and extend the useful life of your chair.
2. Clean your chair regularly. This will help you find loose or worn parts and make your chair easier to use. You will need a mild detergent solution and plenty of cleaning rags.
3. If discovered, repair or replace loose, worn, bent or damaged parts **before using the chair.**
4. To protect your investment, have all major maintenance and repair work done by your authorized dealer.
5. **Inspect and maintain this chair strictly per the maintenance chart.**
6. **If you detect a problem, make sure to order parts, or have service, and repair work done at your authorized dealer before use.**
7. **At least once per year, have a complete inspection, safety check, and service of your chair made by an authorized dealer.**

**B. CRITICAL MAINTENANCE TIPS**

**1. Tire Air Pressure:**

Check air pressure in pneumatic tires at least **ONCE PER WEEK**. The wheel locks will not grip properly if you fail to maintain the air pressure shown on tire sidewall.

**C. MAINTENANCE CHART**

| Weekly | 3 Months | 6 Months | Annually | Check...                              |
|--------|----------|----------|----------|---------------------------------------|
| ✓      |          |          |          | Tire inflation level                  |
| ✓      |          |          |          | Wheel locks                           |
|        | ✓        |          |          | All fasteners for wear and tightness  |
|        |          | ✓        |          | Armrests                              |
|        |          | ✓        |          | Axle and axle sleeves                 |
|        | ✓        |          |          | Quick-release axles                   |
|        | ✓        |          |          | Wheels, tires and spokes              |
|        | ✓        |          |          | Casters                               |
|        | ✓        |          |          | Anti-tip tubes                        |
|        |          | ✓        |          | Frame and crossbrace                  |
|        | ✓        |          |          | Upholstery                            |
|        |          |          | ✓        | <b>Service by Authorized Supplier</b> |

**2. Axles & Axle Sleeves:**

When checking axles and axle sleeves every six months, make sure they are clean and tight. Loose sleeves will damage the axle plate and will affect performance.

**3. Camber Tubes and Camber Plugs**

Inspect the camber tubes and plugs every 6 months ensuring that all fasteners are tight and no parts have moved from their preset position.

**D. CLEANING**

**1. Paint Finish**

- a. Clean the painted surfaces with mild soap or detergent at least once a month.
- b. Protect the paint with a coat of non-abrasive auto wax every three months.

**2. Axles and Moving Parts**

- a. Clean around axles and moving parts WEEKLY with a slightly damp (not wet) cloth.
- b. Wipe off or blow away any fluff, dust or dirt on axles or moving parts.

**3. Upholstery**

- a. Hand-wash (machine washing may damage fabric).
- b. Drip-dry only. DO NOT machine dry as heat will damage fabric.

**4. Basic Maintenance Materials available at local store**

- a. Tire Pump
- b. Clean rags or cotton cloth.
- c. Teflon-based Lubricant
- d. Mild Detergent.

**NOTE**– You do not need to grease or oil the chair.

**E. TROUBLESHOOTING CHART**

| Symptom | Left turn in chair | Right Turn in chair | Looseness in chair | Sluggish turning | Squeaks and rattles | Caster flutter | Solution   |
|---------|--------------------|---------------------|--------------------|------------------|---------------------|----------------|--|
|         |                    |                     |                    |                  |                     |                | <p>Make sure tire pressure is correct and equal in both rear tires and front caster tires, if pneumatic.</p> |
|         |                    |                     | ✓                  | ✓                | ✓                   | ✓              | Make sure all nuts and bolts are tight.  |
|         |                    |                     | ✓                  |                  | ✓                   |                | Make sure all spokes and nipples are tight on radial spoke wheels.   |
|         |                    |                     |                    |                  | ✓                   |                | Use Tri-Flow Lubricant (Teflon®-based) between frame connections and parts.                                  |
| ✓       | ✓                  |                     | ✓                  |                  |                     | ✓              | Check for proper caster fork adjustment. See instructions for caster fork adjustment.                        |

**F. TO MOUNT & REMOVE REAR WHEELS** 11 12

The rear wheel axle (B) attaches the rear wheel (D) to the camber plug (C). When the detent balls (E) engage into the camber plug the rear wheel effectively becomes locked onto the chair. Pushing the quick-release button (A) disengages the detent balls and allows mounting, or removal of the axle and wheel assembly.

**⚠ WARNING**

Do not use this chair UNLESS you are sure both quick-release axles are locked. An unlocked axle may come off during use and cause a fall.



**NOTE**— Setup for this step should include putting the wheelchair on a workbench or table to make the process easier.

**1. To install the Wheel**

- a. Depress quick-release button (A) fully
- b. Insert axle (B) through hub or rear wheel.
- c. Keep button (A) depressed as you slide axle (B) into camber plug (C).
- d. Release button to lock axle in the camber plug. Adjust the axle if it does not lock.
- e. Repeat these steps for the other side of your wheelchair.

**2. To Remove the Wheel**

- a. Depress quick-release button (A) fully.
- b. Remove wheel by sliding axle (B) out of camber plug (C)
- c. Repeat steps on the other side of your wheelchair.

**G. WHEEL LOCKS****⚠ WARNING**

Wheel locks are installed at Sunrise and should be adjusted by your qualified service person. Inspect wheel locks weekly per the maintenance chart. Do not use your chair UNLESS you are sure both wheel-locks can fully engage. A wheel-lock that is not correctly adjusted may allow your chair to roll, or turn unexpectedly. Wheel-locks must be adjusted after making sure the tires have the correct air pressure. When fully engaged, the arm should be imbedded into the tire at least 1/8" to be effective. If you find the wheel locks have slipped or are not working correctly contact your service provider for proper adjustment.

**H. ANTI-TIPS STANDARD AND SWING AWAY** 13 14 15

Sunrise Medical recommends anti-tip tubes for all wheelchairs.

**1. Inserting Anti-tip Tubes Into Receiver** 13

- a. Press the rear anti-tip release button (D) on the anti-tip tube so that both release pins are drawn inside.
- b. Insert the anti-tip tube (E) into the anti-tip receiver (F).
- c. Turn the anti-tip tube down until release pin is through the receiver mounting hole.
- d. Insert second anti-tip tube the same way.

**2. Adjusting Anti-Tip Tube Wheel** 13

The anti-tip wheels may have to be raised or lowered to achieve the proper clearance of 1 1/2" to 2".




Press the anti-tip wheel release button (G) so that both release pins are drawn inside. Raise or lower to one of the three predrilled holes (H). Release button. Adjust the second anti-tip tube wheel the same way. Both wheels should be the same height.

**3. Turning Anti-Tip Tubes Up** 14

Turn anti-tip tubes up when being pushed by attendant, overcoming obstacles or climbing curbs. Press the rear anti-tip tube release button (D). Hold lever in and turn anti-tip tube up. Release the lever and repeat with second anti-tip tube. Remember to return anti-tip tubes to down position after completing maneuver.


**4. Positioning Swing-Away Anti-Tip Tubes** 15

The Swing-Away Anti-Tip can be positioned in the forward position when the wheelchair is being pushed by the attendant, overcoming obstacles or climbing curbs. Press release tab (I) in a downward motion and swing the anti-tip (J) to either the forward or rearward position. The anti-tip will lock back into place by moving upward. Remember to return the anti-tip to the rearward position after completing the maneuver.



**I. PADDED SWING-AWAY ARMRESTS (Optional)** 16

Swing-away armrests can be detached, or can swing away to allow lateral transfers. The Swing-away armrests are installed by sliding armrest (A) into receiver (B) on back frame, and aligning the detent button with the hole in the receiver making sure that the button pops out and armrest is locked.



**J. HEIGHT-ADJUSTABLE ARMRESTS (Optional)** 17

**1. Installation**


- a. Slide the outer armpost into the receiver mounted to the wheelchair frame.
- b. The armrest will lock into place if pressed down firmly.

**2. Height Adjustment**

- a. Rotate height release lever (2) to second stop.
- b. Slide armrest pad up or down to desired height.
- c. Return lever to locked position against armpost.
- d. Push arm pad (4) down until upper armpost locks firmly into place.

**3. Removing Armrest**

- a. Push on the release lever (3) and remove the armrest.



1. Outer armpost
2. Height Release Lever
3. Release lever
4. Armrest pad
5. Transfer bar
6. Side panel
7. Receiver
8. Clamp

**K. BACK ANGLE PERFORMANCE** 18 19

**Folding and Unfolding the Backrest**


1. Release the backrest by pulling the release cord (C) with a simple, short force.
2. While maintaining tension on the release cord, rotate the backrest to the desired angle or to the folded down position.
3. Lock the backrest in place by releasing the cord and allowing the pins (D) to engage in the nearest hole (E) of the backrest bracket (F).

**⚠ WARNING**

Be sure the pins (D) are fully engaged through the backrest bracket (F) on both sides of the chair. **Test the backrest to be sure it is locked before use.**

**L. CUSHION INSTALLATION**

Place cushion on seat sling with VELCRO® side down. The beveled edge of the cushion should be in front. Press firmly into place.




**M. STORAGE TIPS**

1. Store your chair in a clean, dry area. If you fail to do so, parts may rust or corrode.
2. Before using your chair, make sure it is in proper working order.
3. Inspect and service all items on the "Maintenance Chart".
4. If stored for more than three months, have your chair inspected by an authorized supplier before use.

**N. CHECK-OUT**

After the wheelchair is assembled and adjusted, it should roll smoothly and easily. All accessories should also perform smoothly. If you have any problems, follow these procedures:

1. Review the set-up and check-out sections and operating guide to make sure chair was properly prepared.
2. Review the troubleshooting chart at the beginning of this section.
3. If your problem persists, contact your authorized supplier. If you still have a problem after contacting your authorized supplier, contact Sunrise customer service. See the introduction page for details on how to contact your authorized supplier or Sunrise customer service.



## ⚠ WARNING

The owner of this chair is responsible for making sure that it has been setup and adjusted by a trained service professional under the advice of a healthcare advisor. Service and/or adjustments should only be done with the advice of a healthcare professional. Always use parts and/or accessories that have been recommended and approved by Sunrise Medical when servicing this chair.

### A. DEALER SERVICE INTRODUCTION 20

1. At least once per year, this chair should have a complete inspection, safety check, and regular service made by an authorized dealer. The tool list at right is a good starting point for proper maintenance.
2. Torque settings—A torque setting is the optimum tightening which should be made on a particular fastener. It is important to use proper torque settings where specified. When not specified, torque settings should be 60 in-lbs
3. If you have discovered a worn, bent, or damaged part, repair or replace them with recommended parts **before returning this chair to service.**
4. All major maintenance and repair work should be done by the authorized dealer.

### B. CRITICAL MAINTENANCE TIPS

#### 1. Torque settings:

A torque setting is the optimum tightening which should be made on a particular fastener. It is important to use proper torque settings where specified. When not specified, torque settings should be 60 in-lbs

#### 2. Dealer Service and Adjustment Reference materials:

Go to [www.sunrisemedical.com](http://www.sunrisemedical.com) for parts manuals, instruction sheets, and instructional videos that will aid in the repair of the Quickie GT Family of wheelchairs.

### C. CLEANING

#### 1. Paint Finish

- a. Clean the painted surfaces with mild soap or detergent.
- b. Protect the paint with a coat of non-abrasive auto wax.

#### 2. Axles and Moving Parts

- a. Clean around axles and moving parts with a slightly damp (not wet) cloth.
- b. Wipe off or blow away any fluff, dust or dirt on axles or moving parts.
- c. **DO NOT USE 3 in 1 oil, or WD-40 ® for lubrication.** Only use Teflon based Lubricant when working on this wheelchair.

### D. REAR AXLE

The most important adjustment on this Quickie wheelchair is the position of the rear axle. The center of gravity and wheel camber are determined by the axle adjustment.

#### 1. Center of Gravity Adjustment 21

## ⚠ WARNING


The more you move your rear wheels forward, the more likely your chair will tip over backwards. Always make adjustments in small increments, and check the stability of your chair with a spotter to prevent a tip-over. We recommend that you use anti-tip tubes. Refer to additional Warnings in Section VI "Falls and Tip Over".

**NOTE—** Changes to the center of gravity may affect the rear seat height (Section H), toe-in/toe-out of the rear wheels (Section J-4) and the squareness of the casters (Section M). If you change your center of gravity position, re-adjust all of these settings if necessary.

**NOTE—** Adjusting your chair's center of gravity will require re-adjusting the location of the wheel locks (if provided). See Section O for instructions on adjusting the wheel locks.


You can adjust the center of gravity position by moving the two axle plate clamps (A) forward or rearward on the seat tube (B). Moving the axle plate clamps forward shortens the wheelbase and lightens the front end, making your chair more maneuverable. Moving the axle plates rearward makes the chair more stable and less likely to tip over rearward. To adjust the center of gravity location, remove both rear wheels. Loosen the 4 screws (C) (2 per side) that secure the axle plates (D) to the seat tubes. Slide both axle clamps forward or rearward along the seat tube to the desired position. Make sure both sides are adjusted symmetrically before tightening screws. Tighten screws to 144 in-lbs. (16.3 Nm).

Once the axle plate clamps are secured, attach the rear wheels, occupy the chair and maneuver it with a spotter to get a feel for the new adjustment.



#### TOOLS YOU WILL NEED


1. Imperial Hex key set
2. Metric Hex key set
3. Imperial Open-end Wrench set
4. Metric Open-end Wrench set
5. Torque wrench
6. Imperial Hex bit Socket set
7. Metric Hex bit Socket set
8. Phillips and Flat screwdrivers



**2. Wheel Camber** 22


Wheel camber, shown as angular relationship (A), provides greater side-to-side stability due to the increased width and angle of the wheelbase. It also allows for quicker turning and greater access to the top of the hand-rims.

Wheel camber is determined by pairs of interchangeable camber plugs which are available from your authorized supplier in 0°, 3°, 6° and 9° angles.



**3. Changing the Camber Plug** 23

Before changing the camber plug, note the camber tube bolt positions. To change the camber of your chair, loosen bolt (D) on both camber tube clamps (E). Then loosen the two bolts (F). Repeat on other side. Now the camber tube (C) and camber tube clamps (E) can be removed from the axle plate (G). Once removed from the axle plate, the camber tube clamps can be removed from the camber tube. You can now remove the camber plug (B). Insert the desired camber plug into each end of the camber tube. Then replace the camber tube clamps on the ends of the camber tube. Be sure to space camber plugs equal distance from the ends of the camber tube. Tighten bolt (D). Slide camber tube assembly onto axle plate. Realign axle plate holes to their original position. Replace and tighten the 2 bolts (F).



**4. Setting Toe-in toe-out to zero** 24 25 26 27



**NOTE**— A wheelchair equipped with 0° camber plugs cannot have a toe-in toe-out condition. This adjustment is only required when using 3°, 6° and 9° camber plugs.

Toe refers to how well the rear wheels of the chair are aligned relative to the ground. It affects how well the chair will roll. Drag or rolling resistance is optimally minimized when the wheel toe is set to zero.

**5. Setting the Toe to Zero** 24 25

Loosen the 2 cap screws (D) (1 per side) that secure the camber tube clamp. Observe the ball in the level (I), and rotate the camber tube (H) until the ball is centered in the level. The toe is now set at zero.



Before tightening the screws (D), make certain that the camber tube is centered left-to-right relative to the wheelchair frame. There should be an equal gap on both sides or none at all. Torque fasteners (D) to 144 in-lbs. (16.3 Nm).



**6. Setting the Toe to Zero— Alternative Method** 24 25 26 27

Place the entire wheelchair on a flat horizontal table or ground surface. Loosen the 2 cap screws (D) (1 per side) that secure the camber tube clamp. Locate the flat surfaces on the front and rear of the camber plugs (J). Place an object that is known to have an accurate 90° corner (such as a carpenter's square, drafting triangle, etc.) down on the flat horizontal surface and up against the flat of the camber plug. Rotate the camber tube and plug assembly until the flat surface of the camber plug is parallel to the measuring tool.

Before tightening the screws (D), make certain that the camber tube is centered left-to-right relative to the wheelchair frame. There should be an equal gap on both sides of the wheelchair or none at all. Torque the fasteners to 144 in-lbs. (16.3 Nm).



**E. GT WITH XTR SUSPENSION REAR AXLE** 28 29

The most important adjustment on your Quickie wheelchair is the position of the rear axle. The center of gravity, wheel camber and wheelbase are all determined by the axle adjustment.

**NOTE:** Use a torque setting of 100 in.-lbs. when adjusting rear axle.

**1. Center of Gravity**

Moving the center of gravity to the front will lighten the force required to turn the chair. The further back the axle and clamp, the more stable the chair becomes.

The center of gravity is adjusted by moving the rear axle and clamp forward or backward on the frame. Moving the axle and clamp forward will increase your turning speed and lighten the front end.

**IMPORTANT NOTE:** Changing the center of gravity increases the risk of a fall or tip over. Be very careful to assure the chair will not tip over backwards when tilted.

**2. To Move Center of Gravity** 28

Loosen bolts (A). Slide camber tube assembly forward or rearward along frame tube (B). Retighten bolts. Be sure to measure the new position for evenness on both sides of the chair.

Once the axle plate clamps are secured, attach the rear wheels, occupy the chair and maneuver it with a spotter to get a feel for the new adjustment.

**3. Wheelbase Width Adjustment** 29

Adjusting the wheelbase width allows the rider the option to move the wheels closer or further away from the hips. It also compensates for camber adjustment and gives the proper wheel spacing to maximize pushing efficiency.

Adjust the wheelbase width by loosening the nut (C) and repositioning the axle sleeve (D) in its corresponding hole. Make sure the axle sleeve (C) spacing is equivalent on both sides of the chair. Retighten nuts.

**F. REAR WHEEL AXLE NUT ADJUSTMENT** 30 31

To adjust the axle you will need a 3/4" wrench to turn the adjustment nut (E). You will also need a 1/2" wrench to securely hold the ball detent end of the axle (F) to prevent it from turning.

If the wheel and axle (G) will not lock into the camber plug then the axle requires adjustment. Turn the nut (E) counter-clockwise approximately 1/4 revolution and try to lock the axle into the camber plug. If it doesn't lock, continue making small nut adjustments until it securely locks.

If the wheel is locked on the chair but there is excessive wheel play (the wheel hub can be pushed back and forth on the axle) then adjust the nut clockwise until there is no perceptible gap between the wheel and camber tube and the axle is securely locked onto the chair.







**CAUTION-** Quick-release button (H) must be flush with the edge of the wheel hub for detent balls (F) to be engaged.

**G. REAR WHEEL SPACING** 32 33

Rear wheel spacing is measured as the gap between the top of the rear wheels and the backcanes, and is shown as dimension (M) Factory settings will be as narrow as possible depending on the options chosen.

**NOTE-** When setting the rear wheel spacing only make adjustments to one side of the chair at a time. Loosening both sides will undo the toe setting.

To adjust the rear wheel spacing, the camber inserts (I) telescope in and out of the camber tube (K) with hard stops at each end of travel. On the left side of the chair, loosen screw (L) (the one closest to the camber tube). Slide the camber insert in or out to establish the required wheel spacing. Tighten screw to 144 in.-lbs. (16.3 Nm). Repeat on the right side of the chair, matching the wheel spacing set on the left side.



**H. WHEEL LOCK ADJUSTMENT**

Quickie GT and Quickie GTi wheelchairs are shipped with one of five types of wheel locks. Wheel locks are installed at the factory unless you have requested otherwise.

Use a torque setting of 100 in.-lbs. when setting up wheel locks.







**1. High-Mount Push-to-Lock or Pull-to-Lock 34, and the One-Hand/Unilateral 35**

- a. Using a 5mm Allen wrench, turn one of the screws (A) in the clamp counterclockwise one-quarter turn.
- b. Repeat the same process with the second of the two screws.
- c. Alternately loosen the screws (B) (two turns each) until both screws are removed.
- d. Slide clamp toward the rear wheel until the wheel lock is embedded into the tire to prevent wheel movement, when in the locked position.
- e. Tighten screws (B) to 144 in.-lbs. (16.3 Nm).

**2. Ergo Scissor Wheel Locks 36 or Compact Wheel Locks**

- a. Using a 5mm Allen wrench, turn one of the screws (C) in the clamp counterclockwise one-quarter turn.
- b. Repeat the same process with the second of the two screws.
- c. Alternately loosen the screws (two turns each) until both screws are removed.
- d. Slide clamp toward the rear wheel until the wheel lock is embedded into the tire to prevent wheel

**NOTE**– Clamp and wheel lock may need to be rotated to clear frame tubing.



**I. CASTERS**

**NOTE**– Setting the wheelchair on a very flat table or a workbench makes this setting more accurate.

**1. Adjusting Caster Angle 37 38**

**Checking for Caster Squareness:**

Your Quickie wheelchair leaves the factory with the casters adjusted square. Any time you change your wheelchair rear axle height, rear tire size, center of gravity position or camber tube, you should check that the casters are still square. A caster that is not square will result in the caster lifting off the floor as it rotates around the stem bolt and may also lift the front of the wheelchair up and down. Check for squareness by placing your wheelchair on a flat, horizontal surface. Rotate both forks so that the casters are towards the rear of the wheelchair. Place a square or pocket level against the flat surface on the rear of the fork and note whether this surface is vertical (square). If this surface is more than 1 degree away from vertical, your casters require re-squaring.


**2. Re-squaring Your Casters 38 39**

Gently pry cap (D) off of both caster tubes. Locate the set screw (E) inside the right caster tube and loosen using an 8mm Allen wrench. Do not remove the set screw. Simply loosen it enough to allow the teeth on the upper and lower splines (F) to clear one another. With the stem bolt loosened, the caster fork should pivot freely. Place the wheelchair on the floor so that the two rear wheels and left caster contact the floor. Place some weight in the seat towards the rear to help stabilize the chair. Place a square or pocket level against the flat surface of the fork and pivot the fork until it is vertical (square). With the fork vertical, engage the teeth of the two splines to the closest mating position. Make certain that the teeth are properly engaged and not crossed. Shine a bright light behind the mating spline teeth and check for proper teeth engagement. Tighten the set screw (G) to secure the spline mating position. The set screw should be tightened to 192 in.-lbs. (21.7 Nm) to ensure that it does not loosen during use.

Record the teeth engagement position of the splines on the right fork, and adjust the left fork in a similar fashion to the same position. Reassemble both casters, and re-check that they are square. Replace caps (D).

**3. 3-D Perpendicular Alignment** 40

The allen bolt adjustment (A) provides another dimension of adjustability in addition to the spline system. The bolt may be loosened to relieve tension in the stem assembly (B). The caster angle may then be fine-tuned and the bolt retightened.



**J. STANDARD ANTI-TIP TUBES**


Sunrise Medical recommends anti-tip tubes for all wheelchairs. Use torque setting of 106 in.-lbs. (12 Nm) when installing anti-tip tubes.

**K. ANTI-TIP SWING AWAY** 41

**I. Swing-Away Anti-Tip**

The anti-tip wheel may have to be raised or lowered to achieve proper clearance of 1 1/2" to 2".

To adjust, using an Allen key and wrench, remove bolt (C) from current position. Adjust to proper height, reinsert bolt, and tighten back into place.





**L. GT WITH XTR SUSPENSION ANTI-TIP TUBES (Option)** 42 43

Sunrise Medical recommends anti-tip tubes for all wheelchairs.

**NOTE-** Use torque setting of 100 in./lbs. when setting-up anti-tip tubes.

**I. Adjusting Anti-tip Receiver**

- a. Remove anti-tip tubes (D) from anti-tip receivers (E).
- b. Remove button head bolt (F) from anti-tip tube receiver.
- c. Reposition bracket to desired height.
- d. Reinsert button-head bolt and tighten.
- e. Repeat with second anti-tip tube receiver. Both brackets should be at exactly the same height.




**M. PADDED SWING-AWAY ARMRESTS (Option)** 44

Swing-away armrests can be detached or can swing away to allow lateral transfers. They are height adjustable (2") by moving bolts (G) up or down in predrilled holes on armrest bracket.

Swing-away, removable armrests are installed by sliding armrest into receiver (H) on back frame.

**Swing-away, angle adjustable**

To change the angle, remove bolts (I), reposition to needed angle. Reinsert bolt and tighten.



**N. HEIGHT-ADJUSTABLE ARMRESTS (Option)** 45

**1. Installation**


- a. Slide the outer armrest into the receiver mounted to the wheelchair frame.
- b. The armrest will automatically lock into place.

**2. Height Adjustment**

- a. Rotate height release lever to second stop.
- b. Slide armrest pad up or down to desired height.
- c. Return lever to locked position against armrest.
- d. Push down on arm pad until upper armrest locks firmly into place.

**3. Removing Armrest**

- a. Press on release lever to unlock and remove the armrest.



**O. ARMREST RECEIVER ATTACHMENT CLAMP** 46


**Adjusting Armrest Receiver Fit**

To tighten or loosen the fit of the outer armrest in the receiver:

1. Loosen the four receiver adjustment bolts (A) on the sides of the receiver.
2. With the armrest in the receiver, squeeze the receiver to achieve the desired fit.
3. Tighten the four bolts (A).


**Position Adjustment**

1. Loosen the three clamp bolts (C) until clamp (B) is loose.
2. Slide armrest receiver to desired position along the chair frame.
3. Tighten all 3 bolts to 212 in.-lbs. (24Nm).



**P. FOOTREST HEIGHT ADJUSTMENT** 47

To adjust the height of your footrest, loosen the screw (D) on both footrest clamps. Slide the footrest up or down to the new desired height and tighten the two set screws. Use a torque setting of 144 in.-lbs (16.3 Nm). It is recommended that you maintain a minimum of 2.0 inches (5 cm) between the lowest point on the footrest and the floor. This will provide adequate clearance for uneven surfaces and prevent damage to your footrest.




**Q. BACK ADJUSTMENTS**

**CAUTION** – When changing the back angle of the wheelchair both the left and right angle adjustable back bracket must be in the same hole position before the wheelchair is to be used.

**Adjusting Backrest Options** 48

The backrest bracket has two positions, 1 or 2, that can be achieved by removing bolt (E) and pivoting the bracket forward or rearward. With the bracket in the desired position, reinsert and tighten the bolt. Then follow the instructions for “Folding and Unfolding the Backrest”.


**1= Forward Mount. 2= Rearward Mount**



**Backrest Angle Preset Plugs** 49

The preset plugs allow you to unemploy the backrest bracket holes you do not wish to use. This allows the backrest to automatically lock into your desired angle position after unfolding.

1. Determine which backrest bracket holes you do not want to use.
2. Insert plugs (A) into the predetermined holes from the chair side of the bracket (B).
3. To ensure the backrest locks in the folded position, do not insert a plug into the lowest hole (C).

**Depth Adjustable Back Plugs** 49

The DAB Frame insert plugs (D) allow you to compensate for a thicker or thinner after-market backrest without having to change the Center of Gravity of the person in the wheelchair.

1. The backrest may be adjusted by first removing the JAY style quick-release lanyard cord (not shown) from the backrest release pins.
2. Using a 5mm and 10mm Allen Wrench, remove the bolt holding the backrest bracket to the chair.
3. Remove the quick-release pins, frame bushing, and springs.
4. If the desired back depth is 1" greater, move the plunger bushing back one hole. If the desired back depth adjustment is 2", the plunger bushing must be moved to the most rearward hole.
5. To reassemble the back, reinsert the springs into the bushings. Then reinsert the quick-release pins into the bushings with the lanyard attachment holes first.
6. Reposition the backrest in place in the folded position. Insert the quick-release pin into its appropriate placement hole.
7. Insert the large beveled washer in between the backrest bracket and the frame. Insert the backrest plate pivot bolt into the chair and secure using the washer and nylock nut.
8. Using the 5mm and 10mm Allen Wrench, tighten the nylock nut so that the end of the bolt reaches the end of the nut. Repeat on the opposite side of the chair.

**NOTE-** Do not over-tighten the backrest plate pivot bolt as it may hinder the backrest from folding properly.


9. Reattach the JAY®-style quick-release lanyard to the quick-release pins.

**NOTE-** This adjustment is intended to compensate for differences in after-market backrest thickness and/or contour. The back-canes can be adjusted fore and aft relative to the seat so that the backrest support surface does not change the location of the user, regardless of the backrest mounting location. This ensures that the center of gravity of the user relative to the axle location, is not altered due to after-market backrest thickness or contour. There is no recommendation for Depth Adjustable Backrest use for transit.

**Non-folding back instructions** 50


1. Remove backrest bracket bolt (A).
2. Rotate the backrest to the desired angle.
3. Align desired backrest bracket hole with seat frame opening.
4. Reinsert bolt and tighten.

**CAUTION** – When changing the back angle of the wheelchair both the left and right angle adjustable back bracket must be in the same hole position before the wheelchair is to be used.

**R. SEAT SLING ADJUSTMENT** 51

Remove the bolts (B) and rear bolt nut (C) to allow seat upholstery (D) to be removed from frame. Retain plastic saddles (not shown) for placement between seat sling and frame. Repeat steps for other side.

To reattach seat sling, place seat sling (C) over chair frame. Line up seat sling and plastic saddles. Replace the screws (A) and rear nut (B) and tighten accordingly. If there is any difficulty in reattaching the screws, try using a probe to help align the holes. Ensure that the plastic saddles are under the seat rail in all fastener positions.




**S. REAR SEAT HEIGHT ADJUSTMENT** 52

To adjust rear seat height, loosen and remove the 6 cap screws (E) (3 per side) and nut plate (F) that secure the camber tube clamps (G) to the axle plates (H). Move both camber tube clamps (I) to the desired height and reinstall the 4 cap screws. Prior to tightening, follow the instructions for setting the toe-in/out to zero.

Tighten screws to 144 in-lbs. (16.3 Nm).

**NOTE**— Adjustment to rear seat height may also require caster angle adjustment.




**T. REAR AXLE PLATE SUSPENSION** 53

The performance of your suspension system is determined by the compression elastomers (A). There are several different elastomers available from Quickie which are selected based on the user's weight.

To change out elastomers, remove the 2 screws (C) (1 per side). Rotate swingarms (B) downward and the compression elastomer (A) can be freely removed from its top and bottom engagement sockets.

Reverse the order to install new bushings, making sure the compression elastomers are securely engaged into the top and bottom sockets.



**U. XTR SUSPENSION (optional)** 54

- Tuned to perform in conjunction with Quickie GT.
  - Adjustable spring pre-load
1. If the shock ever loses oil or if it makes sounds of excessive topping out, stop use immediately and have the shock inspected by a dealer. Continued use with the shock in either of these conditions could result in loss of control with possible serious and/or fatal injuries.
  2. Always use genuine Quickie parts. Use of after-market replacement parts voids the warranty and could cause structural failure to the shock. Structural failure could result in loss of control with possible serious and /or fatal injuries

**Tuning Your Shock**

By matching the proper spring rate and damping adjustment, your shocks can be tuned to your weight.

When tuning suspension, always make one change at a time and write it down. This takes patience, but allows you to understand how each change affects your ride. Keeping a record lets you know what changes you have tried that work and suggests what changes you might try.

**Washing Your Deluxe Shock**

Clean your shock with mild soap and a toothbrush. Keep the shaft and body eye bushings cleaned and greased.

**IMPORTANT: NEVER USE A HIGH-POWERED WASHER FOR CLEANING THE SHOCK!**


**Adjusting Your Shock Tension**

Locate the spring preloading adjuster ring (4). Rotate adjuster ring clockwise to decrease tension. Rotate adjuster ring counter-clockwise to increase tension.



**WARNING**

Test your new chair settings with a spotter before proceeding with use.



1. Body eye
2. Reservoir
3. Body threads
4. Spring adjuster ring
5. Coil spring
6. Seal head
7. Shock shaft
8. Bottom out bumper
9. Shaft eye

**A. FOR LIFETIME**

Sunrise warrants the frame and cross brace if applicable of this wheelchair against defects in materials and workmanship for the life of the original purchaser.

The expected life of this frame is five years.

**B. FOR ONE (1) YEAR**

We warrant all Sunrise-made parts and components of this wheelchair against defects in materials and workmanship for one year from the date of first consumer purchase.

**C. LIMITATIONS**

1. We do not warrant:
  - a. Tires and tubes, upholstery, pads and push-handle grips.
  - b. Damage from neglect, accident, misuse, or from improper installation or repair.
  - c. Products modified without Sunrise Medicals express written consent.
  - d. Damage from exceeding the weight limit.
2. This warranty is VOID if the original chair serial number tag is removed or altered.
3. This warranty applies in the USA only. Check with your supplier to find out if international warranties apply.
4. This warranty is not transferable and only applies to the first consumer purchase of this wheelchair through an authorized Sunrise Medical dealer.

**D. WHAT WE WILL DO**

Our sole liability is to repair or replace covered parts. This is the exclusive remedy for consequential damages.

**E. WHAT YOU MUST DO**

1. Obtain from us, while this warranty is in effect, prior approval for return or repair of covered parts.
2. Return the wheelchair or part(s), freight pre-paid, to Sunrise Mobility Products Division at: 2842 Business Park Avenue, Fresno, CA 93727-1328.
3. Pay the cost of labor to repair, remove or install parts.

**F. NOTICE TO CONSUMER**

1. If allowed by law, this warranty is in place of any other warranty (written or oral, express or implied, including a warranty of merchantability, or fitness for a particular purpose).
2. This warranty gives you certain legal rights. You may also have other rights that vary from state to state.

Record your serial number here for future reference:

---





Sunrise Medical • 2842 Business Park Ave. • Fresno, CA 93727 • USA  
In Canada (800) 263-3390



Customer Service: 800.333.4000  
[www.SunriseMedical.com](http://www.SunriseMedical.com)

©2014 Sunrise Medical (US) LLC  
6.14 MK-100067 Rev.B