JAY[®]

J3 CUSHION

JAY[®] J3 Cushion

IMPORTANT CONSUMER INFORMATION

NOTICE: This manual contains important instructions that must be passed on to the user of this product. Please do not remove this manual before delivery to the end user.

SUPPLIER: This manual must be given to the user of this product.

USER: Before using this product, read this entire manual and save for future reference.

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NOTE: Check all parts from shipping damage. In case of damage DO NOT use. Contact Carrier/Sunrise for further instructions.

Sunrise Medical recommends that a clinician such as a doctor or therapist experienced in seating and positioning be consulted to determine if a J3 Cushion is appropriate. Cushions should only be installed by an authorized Sunrise Medical supplier.

GENERAL WARNINGS

WARNING

DO NOT install this equipment without first reading and understanding this manual. If you are unable to understand the Warnings, Cautions and Instructions, contact a qualified clinician or supplier - otherwise injury or damage may occur.

WARNING

Hook and Loop velcro should not come into prolonged contact with moisture as this may degrade the adhesive and lead to a failure of the velcro.

WARNING

Installing a cushion on a wheelchair may affect the center of gravity of the wheelchair and may cause the wheelchair to tip backwards, potentially resulting in injury. Always assess for the potential need for antitippers or amputee axle adapter brackets to be added to the wheelchair to help increase stability.

WARNING

Prior to prolonged sitting, any cushion should be tried for a few hours at a time while a clinician inspects your skin to ensure that red pressure spots are not developing. You should regularly check for skin redness. The clinical indicator for tissue breakdown is skin redness. If your skin develops redness, discontinue the use of the cushion immediately and see your doctor or therapist.

The J3 cushion is designed to help reduce pressure. However, no cushion can completely eliminate sitting pressure or prevent pressure sores. The J3 cushion is not a substitute for good skin care including: proper diet, cleanliness, and regular pressure reliefs.

WARNING

Do not leave the fluid cushion outside overnight at temperatures below 40° F (5° C). Allow the cushion to warm to room temperature before using. Sitting on extremely cold or hot surfaces can cause skin damage.

CAUTION- AVOID SHARP OBJECTS OR EXPOSURE TO EXCESSIVE HEAT OR OPEN FLAME.

WARNING

PROPER Fluid Level: The amount of fluid contained in the insert is determined by the clinician during a seating evaluation. This should not be adjusted without a clinician directive as it can alter the pressure relieving properties and lead to injury or damage.

WARNING

OBSTRUCTIONS: DO NOT place any obstructions between the user and the cushion because it will reduce product effectiveness.

WARNING

CUSHION AND COVER ORIENTATION: Product must be used with the fluid insert facing up. If the cover is not used correctly it may reduce or eliminate the cushion's benefits and could increase risk to the skin and soft tissue.

WARNING

ROHO® DRY FLOATATION® PAD: refer to applicable warnings in the ROHO® DRY FLOATATION® Pad Supplement.

INTRODUCTION

The Jay J3 Cushion is designed for clinicians to use with clients who need superior pressure distribution and stability. The J3 Cushion offers a complete solution to meet a large range of skin protection and positioning needs.

Maximum user weight: 500 lbs.

CUSHION ADJUSTABILITY

Fluid and Air Adjustability

The J3 Cushion offers volume adjustability to accommodate for users needs through four models, the “Factory Filled” or the “Field Variable” model.

The determination of Air Single or Air Dual Models is made by your clinician during evaluation. The inserts should only be changed at the direction of your clinician.

The fluid field variable insert has the ability to adjust to the user's need over time. Your supplier will be involved in the adjustment once it is determined to be needed for better clinical efficacy.

Adjusting through the “Field Variable Volume” J3 Cushion

The Field Variable version of the J3 Cushion allows the appropriate healthcare provider to apply a specified amount of fluid at the time of prescription and to modify the amount of fluid if needs change through the use of discreet Fluid Supplement Pads. This method of adjustment is useful where there is no concern of unqualified fluid adjustments occurring after prescription.

Adjusting through the ROHO® DRY FLOATATION® Pad Insert

The air insert version of the JAY J3 Cushion can be adjusted at any time through the use of the hand pump (provided) and release valve on the ROHO® DRY FLOATATION® insert.

WARNING

The Field Variable Fluid Pad comes standard with 2 Fluid Supplement Pads and has a Blue Urethane Cover to differentiate it versus the standard grey urethane on the Factory Filled Fluid Pad. This color difference identifies that the Blue Fluid Pad comes underfilled in comparison to the Factory Filled Pad, and has 4 velcro attachment points to attach the standard Fluid Supplement pads to handle obliquities, and overfills. To determine how to use the Fluid Supplement Pads, please reference the “Fitting a Pressure Distributing Cushion” section.

FITTING A PRESSURE DISTRIBUTING CUSHION

Using Fluid Supplements pads

Fluid Supplement Pads are to be only used with the Field Variable Fluid Pad. To determine if you have a Field Variable Fluid Pad, unzip the cushion cover, and look at the pad in the Pelvic Loading Area. If the pad has blue outer cover, it is a Field Variable Fluid pad and designed to work with Fluid Supplement Pads.

An authorized clinician or therapist should use Fluid Supplement Pads by attaching or removing to the Fluid Pad Velcro Attachment points located between the folds on the upper Fluid Pad surface to increase or decrease the amount of fluid needed.

Fitting the cushion

Select a cushion to match your hip width when in the seated position. The cushion length should extend to within 1" - 2" (2.5cm - 5cm) from the back of the knees. Place the cushion to the back of the wheelchair with the material identification tags and pressure relieving pad toward the rear. When fit correctly, your hips should be to the back of the chair. Your ischials (seat bones) should be centered on the insert in the pelvic loading area of the cushion. Proper footrest adjustment may enhance sitting comfort and help lower peak sitting pressures. Adjust the footrests so your legs rest firmly but comfortably on the cushion.

Use the cushion indexing notches located in the rear of the foam base to ensure proper alignment of the cushion to the backrest. Misalignment of the insert and pelvic loading area could lead to injury and skin breakdown.

Avoid using the cushion on overly stretched seat upholstery. Stretched upholstery does not provide adequate support for the cushion and improper positioning may result. Replace stretched upholstery or purchase the Jay Adjustable Solid Seat.

Reference the ROHO® DRY FLOATATION® Pad Owner's Manual Supplement for proper fitting and set-up instructions.

EASY MAINTENANCE & CLEANING

Monthly cleaning and regular maintenance may help extend the life of your cushion. During cleaning, component inspection is recommended.

Check the cover for tears and excessive wear and replace if ripped, torn, or otherwise not fully functional. Inspect the fluid pad for punctures or any other abnormalities.

Check the foam base to ensure foam consistency. While checking the fluid pad if you ever feel the fluid is firmer in one area, simply knead the fluid back to its original consistency and/or call your local authorized supplier.

To clean cover

1. Remove the cover from the foam base and turn inside out.
2. Machine wash in warm water and drip or tumble dry on low heat.

No ironing required.

CAUTION- DO NOT DRY CLEAN COVERS OR USE INDUSTRIAL WASHERS AND DRYERS TO CLEAN COVER. DO NOT STEAM AUTOCLAVE. DO NOT BLEACH.

To refit cover after washing, reinstall the cover on the foam base. Ensure the back of the cover matches up to the back of the foam base.

NOTE: If utilized in an institutional setting, write the resident's name with a permanent marker on the content label of the covers. This will help to ensure the return of the covers after washing.

Use only Jay-designed covers on Jay cushions. The cover is an important part of allowing the immersion into the fluid insert which redistributes the pressure. Changing the cover will alter this design and its efficacy.

To clean foam base, PLA reducer ring and positioning components

Remove the cover, fluid pad, and PLA reducer ring (if applicable) from the foam base. Wipe lightly with a damp cloth. Do not use soap. Do not submerge in water. Wipe off with a clean cloth. Let it completely air-dry before reassembly.

NOTE: Avoid submersion of the foam base in water.

Reinstall the cover over foam base. Ensure the back of the cover matches up to the back of the foam base.

To clean fluid pad and fluid inserts (if applicable)

Remove the cover from the foam base. Remove fluid pad from its foam base and wipe with warm water and soap. Rinse with a clean, damp cloth. Then wipe dry with a clean cloth.

CAUTION- NEVER SUBMERGE THE FLUID PAD IN ANY LIQUID.

CAUTION- AVOID HARSH CLEANING OR ROUGH HANDLING AS THIS MAY LEAD TO A DEGRADATION OF THE FLUID INSERT.

If at any time, the cover, Pelvic Loading Ring, Foam base, fluid insert or any other component appears to be wearing or you notice a degradation in the texture, contact your supplier/Sunrise for evaluation and further instructions.

To clean the ROHO® DRY FLOATATION® Pad (if applicable),

Reference the ROHO® DRY FLOATATION® Pad Owner's Manual Supplement for proper cleaning instructions.

Laundry Care Symbols	
	Machine wash in warm water (60° C)
	Only non-chlorine bleach, when needed
	Drip dry
	Do not iron
	Do not dry clean


REASSEMBLING THE JAY J3 CUSHION AFTER CLEANING

Follow cleaning and maintenance instructions as described in "Easy maintenance and cleaning" then reassemble as follows.

⚠ WARNING

Insert the JAY Fluid or ROHO® DRY FLOATATION® Pad on the foam base making sure that all Velcro® strip locations are aligned. A misalignment or creation of an edge could cause an unwanted pressure point which could lead to injury.

Insert the foam base into the outer cover. Ensure that the zipper is at the rear of the foam base (the part number tag should be to the rear of the foam base) where the pelvic loading area is located.

WHAT IS BOTTOMING OUT?

Bottoming out may occur on a J3 Cushion if you displace the fluid or air underneath your pelvic bones or buttocks, which would leave you sitting on the foam base. It sometimes occurs on very thin individuals, people using recliner wheelchairs, people who have lost weight or have other changes in body size or shape, people who slouch when sitting, or a person using a P.L.A size inappropriately sized for their pelvis.

WARNING

When bottoming out occurs, increased pressure is placed onto the ischials and coccyx increasing the risk for skin breakdown. Immediately discontinue use of the pad. See your healthcare professional.

HOW TO CHECK FOR BOTTOMING OUT

JAY Flow Fluid Pad

To check for bottoming out, sit on the cushion, with the cover on and the zipper opened, for seven to ten minutes. Transfer up and off the cushion (or have someone help you transfer), trying not to disturb the fluid underneath you. Fold the back portion of the cover forward over the front, exposing the fluid pad. Push down in the depressions on the pad where your ischials (seat bones) and coccyx (tailbone) were. You should have to push through at least 1/2" (1.3 cm) of fluid before you feel the firm cushion base below.

1/2" (1.3 cm) is approximately this deep.

If the cushion is properly positioned, the footrests are properly adjusted, and there is not at least the minimum 1/2" (1.3 cm) of fluid, the cushion is bottoming out and should not be used. If you are bottoming out, discontinue use of the cushion and see your clinician. Usually bottoming out is easily solved by ordering additional fluid in the form of an Overfill on the Factory Filled model. On the Field Variable version, you can use the Fluid Supplement pads to add additional fluid. Call your local authorized supplier to see if this is appropriate for you.

ROHO® DRY FLOATATION® Pad

Reference the ROHO® DRY FLOATATION® Owner's Manual Supplement for instructions on preventing and checking for bottoming out.

WHEN TO CHECK

With JAY Flow Fluid, check every month when you clean it.

With ROHO® DRY FLOATATION®, check daily.

JAY J3 CUSHION WARRANTY

Each Jay cushion is carefully inspected and tested to provide peak performance. Every Jay cushion is guaranteed to be free from defects in materials and workmanship for a period of 24 months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within 24 months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge. This warranty does not apply to punctures, tears or burns, nor to the cushion's removable cover.

The removable cover is guaranteed to be free from defects in materials and workmanship for a period of six (6) months from the date of purchase, provided normal use. Should a defect in materials or workmanship occur within six (6) months from the original date of purchase, Sunrise Medical will, at its option, repair or replace it without charge.

Claims and repairs should be processed through the nearest authorized supplier. Except for express warranties made herein, all other warranties, including implied warranties of merchantability and warranties of fitness for a particular purpose are excluded.

There is no implied warranty beyond what is contained herein. Remedies for breach of express warranties herein are limited to repair or replacement of the goods. In no event shall damages for breach of any warranty include any consequential damages or exceed the cost of non-conforming goods sold.

ADDENDUM TO WARRANTY USE AND CARE GUIDE

JAY Flow Fluid Pad

Consult your healthcare professional for proper fitting and use, and read the following statements carefully:

- The fluid in this cushion may lose volume over time.
- Volume loss in this cushion may result in bottoming out the cushion.
- Monthly, during routine cleaning and inspection, check this cushion for bottoming out (see instructions under "checking for bottoming out" in the warranty use and care guide) and any fluid inconsistencies. If the fluid is firmer in one area simply knead it back to its original consistency.
- If bottoming out occurs discontinue use of this cushion and contact your healthcare professional. If further assistance is necessary, contact Sunrise Medical Customer Service at 1-800-333-4000.

ROHO® DRY FLOATATION® Pad

Reference the ROHO® DRY FLOATATION® Owner's Manual Supplement for additional use and care guidelines.

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©2009 Sunrise Medical Inc. 6.09
ROHO® and DRY FLOATATION® are registered trademarks of ROHO, Inc.
Sunrise Medical AB
Sunrise Medical B.V.