

PLEASE SEE THE DISCLAIMER INFORMATION PROVIDED AT THE BOTTOM OF THIS DOCUMENT

Guide to Memory Seat Programming

Step 1: Move the seating to desired position.

- Change each actuator and element of the seating system to achieve your desired position.



Step 2: Using the joystick, display the desired memory position you would like to change.

- EX: *"M3 Pressure Relief" or "relax"*



Step 3: Plug a single switch (ex: egg switch) into inhibit #6 on the CxSM.



Step 4: Press and hold the single switch, then pull the joystick back and hold for 3 seconds.



Step 5: Once you hear the joystick beep the position is set.

- Release the joystick then release the button.



Reset: To reset a memory position back to factory settings repeat steps 4-5 but hold for 6 seconds. (second tone will sound)

Note: Only Sunrise Medical power wheelchairs utilizing the CxSM control module have access to the memory seating option. If you have questions as to if a specific power wheelchair has this option, contact your local Sunrise Medical representative or Customer Service at 800-333-4000 with the serial number.

DISCLAIMER: FOR PROFESSIONAL USE ONLY. THIS DOCUMENT (AND THE WEBSITE REFERENCED HEREIN) DO NOT PROVIDE MEDICAL ADVICE. Sunrise Medical (US) LLC does not provide clinician services. The information contained on this document (and the website referenced herein), including, but not limited to, the text, graphics, images, and descriptions, are for informational purposes only and should be utilized as a general resource for clinicians and suppliers to then use clinical reasoning skills to determine optimal seating and mobility solutions for individual patients. No material on this document (or on the website) is intended to be used as (or a substitute for) professional medical advice, diagnosis or treatment. Clinicians should adhere to their professional medical training and their institution's or practice's recommended guidelines. Reliance on this document (and the information contained herein) is solely at your own risk.