

Elegance ACCESSORIES



Elegance Accessories



❑ Anti-tips

Inserted into the rear frame of the chair, anti-tips will give you confidence and extra stability when climbing steep hills or curbs.

PN: 129412



❑ Step Tube

This easily attachable insert can be "stepped" on by your caregiver to help boost the chair up over a curb.

PN: 064931



❑ Back Pack

This ample sized carry-all attaches quickly to the back of your Elegance and features a Velcro® closure to protect your valuables.

PN: 128177



❑ Elevating Legrests

For clients requiring the elevation of their legs for various conditions. Available individually.

PN: Pair-129415*, Left-129413, Right-129414

*Footplates with extensions are sold separately.



❑ Crutch Holder

An attachment that allows crutches to be carried on the back of your Elegance chair.

PN: 067708



❑ Extra Cushion

For times when your cushion becomes soiled or damaged, having an extra cushion prevents any delays in adventure!