



Quickie TS

EIQT, Effective January 2018

ACCOUNT INFORMATION

Account No: _____	Drop Ship/Ship To: _____
Date: _____	Name: _____
P.O. #: _____	Address: _____
Buyer: _____	City: _____ State: _____
Marked For: _____	Zip Code: _____ Phone Number: _____

The HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payor's coverage policy.

QUICKIE TS

The Quickie TS is NA w/ seat or back upholstery nor does it include seating
To order seating, see "Standard Conventional Seating System" on page 4.

<input type="checkbox"/> 043WM35 E1161 Quickie TS	\$4,000.00
---	------------

FRAMES

Base Frame		
<input type="checkbox"/> 043FS4	Growing Base	NC
<input type="checkbox"/> 043FS5	Flared Front Growing Base	NC
Required for Contracture Footrest System		

Seat Frame		
<input type="checkbox"/> 043FT49	Growing Seat Frame	STD
<input type="checkbox"/> 043FT46	Contracture Growing Seat Frame	\$200.00
Manufactured w/o footrest hanger tubes Use with Contracture Footrest System.		

Heavy Duty Option

Use 24" Pneumatic rear tires and 8" x 2" pneumatic front tires only
Max Wt. Capacity 350 lbs

<input type="checkbox"/> 043WT2 Heavy-Duty Option	\$145.00
---	----------

Seat Frame Position (See Page 4 for Chart)

<input type="checkbox"/> 043SP1	Std. Seat Frame Position	NC
<input type="checkbox"/> 043SP2	-1" Seat Frame Position	NC
<input type="checkbox"/> 043SP3	-2" Seat Frame Position	NC

Tilt Option

<input type="checkbox"/> 043TT1	Gas Cylinders	NC
<input type="checkbox"/> 043TT2	Dual Mech Locks	NC
NA w/ reclining back		

FRAME COLOR

<input type="checkbox"/> 043C140	Aztec Gold	NC
<input type="checkbox"/> 043C2	Black (NA as base w/Zebra)	NC
<input type="checkbox"/> 043C141	Black Cherry	NC
<input type="checkbox"/> 043C103	Black Opal	NC
<input type="checkbox"/> 043C38	Blue Opal	NC
<input type="checkbox"/> 043C36	Candy Blue	NC
<input type="checkbox"/> 043C14	Candy Purple	NC
<input type="checkbox"/> 043C35	Candy Red	NC
<input type="checkbox"/> 043C142	Electric Blue	NC
<input type="checkbox"/> 043C143	Evergreen	NC
<input type="checkbox"/> 043C126	Glow	NC
<input type="checkbox"/> 043C043	Green Apple	NC
<input type="checkbox"/> 043C147	Hot Sparkle Pink	NC
<input type="checkbox"/> 043C148	Mauve Pink	NC
<input type="checkbox"/> 043C103	Rootbeer	NC
<input type="checkbox"/> 043C114	Sparkle Silver	NC
<input type="checkbox"/> 043C145	Sunrise Orange	NC
<input type="checkbox"/> 043C113	Titanium Color Paint	NC
<input type="checkbox"/> 043C8	Yellow	NC

MATTE FRAME COLOR

<input type="checkbox"/> 043C116	Matte Black	NC
<input type="checkbox"/> 043C151	Matte Black Cherry	NC
<input type="checkbox"/> 043C150	Matte Electric Blue	NC
<input type="checkbox"/> 043C149	Matte Evergreen	NC
<input type="checkbox"/> 043C152	Matte Purple	NC

CUSTOM FRAME COLOR

<input type="checkbox"/> 043C112	Desert Camouflage	\$295.00
<input type="checkbox"/> 043C43	Mossy Oak Camouflage	\$295.00
<input type="checkbox"/> 043C42	Stars & Stripes	\$295.00
<input type="checkbox"/> 043C108	Zebra (Choose above base color)	\$295.00
Std with White unless above base paint is chosen NA with Black.		

Frame Width

<input type="checkbox"/> 043W5	14"	NC
<input type="checkbox"/> 043W6	15"	NC
<input type="checkbox"/> 043W7	16"	NC
<input type="checkbox"/> 043W8	17"	NC
<input type="checkbox"/> 043W9	18"	NC
<input type="checkbox"/> 043W10	19"	NC
<input type="checkbox"/> 043W11	E2201 20"	\$500.00

Frame Depth & Back Post Location (See Diagram on Page 4)

Dual post armrests NA w/ 14" short, 16" regular, or 19" long options

043FL5 Short		
<input type="checkbox"/> 043BL1	14"	NC
<input type="checkbox"/> 043BL2	15"	NC
<input type="checkbox"/> 043BL3	16"	NC
<input type="checkbox"/> 043BL4	17"	NC
043FL11 Regular		
<input type="checkbox"/> 043BL3	16"	NC
<input type="checkbox"/> 043BL4	17"	NC
<input type="checkbox"/> 043BL5	18"	NC
<input type="checkbox"/> 043BL6	19"	NC
043FL21 Long		
<input type="checkbox"/> 043BL6	19"	NC
<input type="checkbox"/> 043BL7	E2203 20"	\$550.00
<input type="checkbox"/> 043BL8	E2203 21"	\$550.00
<input type="checkbox"/> 043BL9	E2204 22"	\$900.00

SEATING**Attach Order Form for Seating For:**

- Jay Fit Seating
- Standard Conventional Seating (See Page 4)
- Custom Conventional Seating

043PN1 E2231 Aluminum Solid Seat Base \$75.00

Back Type

043BT8 Angle Adjustable Backrest NC
 043BT5 E1225 Reclining Backrest \$700.00
NA w/ low, SA adj stroller handles, or dual mech locks

Stroller Handle Back Posts (For Adaptive Seating)

043BP9 Low (25-3/4" Ht.) NC
 043BP10 Medium (28-3/4" Ht.) NC
 043BP11 Tall (31-3/4" Ht.) NC
 043BP12 Removable Stroller Handle Ext. \$85.00
NA w/ swing-away adjustable stroller handles
 043BP15 Swing Away Adjustable Stroller \$100.00
NA w/ removable stroller handle ext.

ARMREST**Armrest**

Dual post armrests NA w/ 14" short, 16" regular, or 19" long frame depth/backpost options

043AA26 Dual Post NC
 043AA39 E0973 Dual Post Height Adjustable \$230.00
 043AA43 Dual Post Flip Back \$95.00
 043AA48 E0973 Dual Post Flip Back Height Adjust \$295.00
 043AA8 E0973 Height Adjustable Adult \$285.00
 043AA9 E0973 Height Adjustable Kids \$285.00
 043AA1 E0973 Adjustable Locking Flip-Up \$255.00
NA w/ Armpads

Armpad

043AM1 Desk Length (10") NC
 043AM2 Full Length (14") NC

FOOTREST

Heel loop or leg strap not included, See Footplate Options Section

1. Hanger Mounted Footrest (Footrests Mount to Hanger Tubing)

043H7 70° Swing-Away Front Mount NC
 043H8 80° Swing-Away Front Mount NC
 043H6 90° Swing-Away Front Mount NC
 043H9 70° Heavy-Duty Front Mount Lift-Off NC
 043H10 80° Heavy-Duty Front Mount Lift-Off NC
 043H23 E0990 90° Elevating \$280.00
NA w/ widths over 16"
Available w/ adult and kids angle adj footplates only

Extensions

043E23 2" \$55.00
 043E16 High Mount \$55.00
Only available w/ adult and kids angle adjustable footplates
NA w/ ELRs

Footplate

043F4 K0040 Kids Angle Adjustable \$185.00
Only available on 14"-16" width frames
 043F3 K0040 Adult Angle Adjustable \$185.00

2. Extension Tube Mounted Footrests (Footrests Mount to Ext. Tubes)

043H1 60° Swing-Away NC
 043H2 70° Swing-Away STD
 043H14 E0990 Elevating \$295.00
Only available w/ aluminum and angle adj footplates
 043H24 K0053 Articulating Adult \$350.00
NA w/ high-push wheel locks and short footrest ext tubes
or 14" width. Only available w/ comp and adult angle adj footplate:

Extension Tube

043E1 Short (11.5" - 14") NC
Only available w/ aluminum and angle adj footplates
 043E2 Medium (14" - 17") NC
 043E5 Long (17" - 19") NC

Footplate

043F1 E0951 Composite Std
Only available with 60° or 70° hangers and ALRs
 043F5 Aluminum (ELRs only) NC
 043F3 K0040 Adult Angle Adjust. \$185.00
 043F48 K0040 Platform Flip-Up \$190.00

Footrest Options

043HL1 E0951 Heel Loop \$30.00
 043HL5 K0038 Leg Strap \$35.00
 043IG2 Impact Guard Neoprene \$85.00
Available w/ 60° or 70° hangers with ext tubes only
 043HL5 K0038 Leg Strap \$35.00
NA w/ ELR, ALR, 90° ELR or contracture footrest.
 043H32 Omit Hanger -\$15.00
 043F7 Omit Hanger & Extension Tube -\$15.00
 043HL4 Omit Heel Loop NC
 043HL6 Omit Single Leg Strap NC

3. Contracture Footrest System (CFS)

Multi-Adjustable Padded Hangers & Angle Adjustable Aluminum Footrests

Best when used with Flared Front End

Length measurement taken from top of seat rail to top of footrest plate when CFS is set at 90°. Footrest plate width determined by seat frame width. CFS may interfere with front caster rotation when Quickie TS is in non-tilt position.

043H26 K0108 Contracture Footrest System \$325.00

043E25 Short

043EP1 10.5" Extension Height NC
 043EP2 11.0" Extension Height NC
 043EP3 11.5" Extension Height NC
 043EP4 12.0" Extension Height NC
 043EP5 12.5" Extension Height NC
 043EP6 13.0" Extension Height NC

043E26 Regular

043EP5 12.5" Extension Height NC
 043EP6 13.0" Extension Height NC
 043EP7 13.5" Extension Height NC
 043EP8 14.0" Extension Height NC
 043EP9 14.5" Extension Height NC
 043EP10 15.0" Extension Height NC

043E27 Long

043EP9 14.5" Extension Height NC
 043EP10 15.0" Extension Height NC
 043EP11 15.5" Extension Height NC
 043EP12 16.0" Extension Height NC
 043EP13 16.5" Extension Height NC
 043EP14 17.0" Extension Height NC

CASTERS

<input type="checkbox"/>	043CW5	5" Low-Profile Polyurethane	Std
<input type="checkbox"/>	043CW7	E2219 5" x 1.5" Semi-Pneumatic	\$85.00
<i>Caster wheel interferes with caster fork in the upper hole position.</i>			
<input type="checkbox"/>	043CW33	E2214 6" Pneumatic	\$75.00
E2215			
<input type="checkbox"/>	043CW34	6" Polyurethane	NC
<input type="checkbox"/>	043CW35	E2219 6" x 1.5" Semi-Pneumatic	\$85.00
<i>Caster wheel interferes with caster fork in the upper hole position.</i>			
<input type="checkbox"/>	043CW21	E2214 8" Pneumatic	\$75.00
E2215 NA w/ 90° hangers and 90° ELR			
<input type="checkbox"/>	043CW32	8" Polyurethane	NC
NA w/ 90° hangers and 90° ELR			
<input type="checkbox"/>	043CW39	E2219 8" x 1.5" Semi-Pneumatic	\$85.00
NA w/ 90° hangers and 90° ELR			
<i>Caster wheel interferes with caster fork in the upper hole position.</i>			
<input type="checkbox"/>	043CW36	E2214 8" x 2" Pneumatic	\$145.00
E2215 NA w/ 90° hangers and 90° ELR			
<input type="checkbox"/>	043CW37	E2214 8" x 2" Pneumatic w/Airless Ins.	\$225.00
K0108 NA w/ 90° hangers and 90° ELR			

Caster Options

<input type="checkbox"/>	043CS6	3/4" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	043CS10	1-1/2" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	043CW12	Omit caster wheel and tire	NC
<input type="checkbox"/>	043CF30	Omit Caster Wheels and Forks	NC
<input type="checkbox"/>	043CS4	Quick-Release Caster Stems	\$160.00
NA w/ 3/4" or 1-1/2" fork stem bolts or 8" x 2" front casters			

AXLE AND WHEEL**Axle Plate**

Self-Propelling available w/ 24" and 26" rear wheels only. Tilt is limited to 33°. Rear seat to floor height 18" with 24" rear wheels, 19" with 26" rear wheels.

<input type="checkbox"/>	043AL1	Standard Axle Plate	NC
<input type="checkbox"/>	043AL21	Self-Propelling 0° Axle	\$150.00
<input type="checkbox"/>	043AL22	Self-Propelling 4° Axle	\$150.00
<input type="checkbox"/>	043AL17	E0958 One Arm Drive Axle Plate	\$850.00
Only available w/ Mag Wheels			
<input type="checkbox"/>	043OA2	One Arm Drive Standard Handrim	NC
<input type="checkbox"/>	043HM1	Right Hand Mount	NC
<input type="checkbox"/>	043HM2	Left Hand Mount	NC
<input type="checkbox"/>	043OA1	One Arm Drive Small Handrim	NC
<input type="checkbox"/>	043HM1	Right Hand Mount	NC
<input type="checkbox"/>	043HM2	Left Hand Mount	NC

Axle

<input type="checkbox"/>	043AX1	Quick Release Axles	NC
<input type="checkbox"/>	043RR25	Omit Rear Wheel and Axle	-\$55.00

Rear Wheel

12", 16", and 18" wheels only available w/ Pneumatic and Pneumatic w/ Airless Inserts.

20", 22" only available w/ Pneumatic, Pneumatic w/ Airless Inserts and Full Poly

12" and 16" Mags NA w/ handrims

<input type="checkbox"/>	043RW48	12" Mag Wheel	\$40.00
<input type="checkbox"/>	043RW52	16" Mag Wheel	\$40.00
<input type="checkbox"/>	043RW3	20" Mag Wheel	NC
<input type="checkbox"/>	043RW8	22" Mag Wheel	NC
<input type="checkbox"/>	043RW1	24" Mag Wheel	NC
<input type="checkbox"/>	043RW12	18" Spoke Wheel	\$100.00
<input type="checkbox"/>	043RW4	20" Spoke Wheel	\$200.00
<input type="checkbox"/>	043RW9	22" Spoke Wheel	\$200.00
<input type="checkbox"/>	043RW51	24" Spoke Wheel	\$100.00

Rear Wheel Tire

<input type="checkbox"/>	043RT3	Full Polyurethane	NC
NA w/ 12", 16", and 18"			
<input type="checkbox"/>	043RT1	Pneumatic	NC
<input type="checkbox"/>	043RT2	E2213 Pneumatic Airless Insert	\$80.00
E2211			
<input type="checkbox"/>	043RT6	E2211 High Pressure Clincher	\$95.00
E2212 Available with 24" spoke only			

Handrim

NA for 12" or 16" wheels

<input type="checkbox"/>	043HR1	Aluminum Anodized	NC
<input type="checkbox"/>	043HR4	Plastic Coated	\$95.00
<input type="checkbox"/>	043HR7	Omit Handrim	NC
Does not apply to 12" & 16" wheels			

Projections

NA on plastic coated

<input type="checkbox"/>	043HR8	8 Oblique	\$325.00
--------------------------	--------	-----------	----------

Wheel Lock

<input type="checkbox"/>	043WL1	High Mount Push Wheel Lock	NC
<input type="checkbox"/>	043WL2	High Mount Pull Wheel Lock	NC
<input type="checkbox"/>	043WL11	Foot Wheel Lock (Only 12" whl)	\$50.00
<input type="checkbox"/>	043WL10	Omit Wheel Lock	NC

Grade Aid

<input type="checkbox"/>	043GA1	E0974 Grade Aid	\$150.00
NA w/ poly or high-pressure clincher			

Extension Handle

<input type="checkbox"/>	043EH1	E0961 6" Pair	\$75.00
--------------------------	--------	---------------	---------

Anti-Tip

<input type="checkbox"/>	043AO1	E0971 Front Anti-Tip	\$150.00
<input type="checkbox"/>	043AT5	E0971 Rear Anti-Tip	\$140.00

ACCESSORIES

Backpack

- 043A1 Black Utility Backpack \$65.00
- 043A2 Black/Charcoal Kids \$60.00

Seat Pouch

- 043A8 Black \$50.00

Extra Leg Strap

- 043A54 K0038 Black \$30.00

Miscellaneous

- 043A82 Caddy \$50.00
Only available w/ footrests with ext tubes

Positioning Belt

- 043A141 E0978 1.5" Auto Buckle w/Clamps \$70.00
- 043A12 E0978 2" Aircraft Buckle \$110.00
- 043A13 E0978 2" Aircraft Buckle Padded \$120.00
- 043A127 E0978 2" Dual Pull Side Squeeze Padded \$85.00

Spoke Guard

NA w/ 12", 16", and 18" wheels

- 043A24 K0065 Clear \$115.00
- 043A23 K0065 Black \$115.00

Wheelchair Tray Table *(Must include armrest)*

- 043A176 E0950 Medium (14") \$255.00
- 043A177 E0950 Large (15"-17") \$255.00
- 043A178 E0950 X-Large (18"-20") \$255.00

Side Guard

NA w/ dual post armrest styles or single post height adj. armrests

- 043SG3 K0108 Plastic Regular \$140.00
- 043SG4 K0108 Plastic Kids \$140.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	043JC20	JAY Bundled Cushion Cushion Name/Part #: _____ Width: _____ Depth: _____	Price List
--------------------------	---------	--	------------



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA