



Zippie Kidz
January 2018
EIKZ

Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: () _____ Fax: () _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email Mark For: _____	SUBMITTING FOR: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: () _____ Fax: () _____
--	---

*The HCPCS CODES provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of HCPCS CODES does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
ZIPPIE KIDZ									
<input type="checkbox"/>	013WM36	E1235	Zippie KIDZ	\$3,095					

1. FRAME STYLE									
<input type="checkbox"/>	013FC3		Standard Configuration	Std					
<input type="checkbox"/>	013FC4		Reverse Configuration <i>(Large wheels in front & casters in back. N/A w/ 20" rear wheels or 5" casters.</i>	NC	<input type="checkbox"/>	013CV1		Reverse Configuration Conversion Kit <i>(N/A w/ 20" rear wheels or 5" casters. Used to reverse configure chair w/ Large wheels in front & casters in back.</i>	\$50

2. FRAME COLOR

Standard Colors			
<input type="checkbox"/>	013C140	Aztec Gold	NC
<input type="checkbox"/>	013C2	Black <i>(N/A as base w/ Zebra)</i>	NC
<input type="checkbox"/>	013C141	Black Cherry	NC
<input type="checkbox"/>	013C10	Black Opal	NC
<input type="checkbox"/>	013C38	Blue Opal	NC
<input type="checkbox"/>	013C36	Candy Blue	NC
<input type="checkbox"/>	013C14	Candy Purple	NC
<input type="checkbox"/>	013C35	Candy Red	NC
<input type="checkbox"/>	013C142	Electric Blue	NC
<input type="checkbox"/>	013C143	Evergreen	NC
<input type="checkbox"/>	013C126	Glow	NC
<input type="checkbox"/>	013C144	Green Apple	NC
<input type="checkbox"/>	013C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	013C148	Mauve Pink	NC
<input type="checkbox"/>	013C103	Rootbeer	NC
<input type="checkbox"/>	013C114	Sparkle Silver	NC
<input type="checkbox"/>	013C145	Sunrise Orange	NC
<input type="checkbox"/>	013C113	Titanium Color Paint	NC
<input type="checkbox"/>	013C8	Yellow	NC

Standard Matte Colors			
<input type="checkbox"/>	013C116	Matte Black	NC
<input type="checkbox"/>	013C151	Matte Black Cherry	NC
<input type="checkbox"/>	013C152	Matte Purple	NC
<input type="checkbox"/>	013C150	Matte Electric Blue	NC
<input type="checkbox"/>	013C149	Matte Evergreen	NC

Custom Colors			
<input type="checkbox"/>	013C43	Mossy Oak Camouflage	\$295
<input type="checkbox"/>	013C42	Stars and Stripes	\$295
<input type="checkbox"/>	013C108	Zebra <i>(Choose base color above: White Std, N/A Black)</i>	\$295

3. FRAME WIDTH

<i>(Seat Width 1/2" Narrower)</i>									
<input type="checkbox"/>	013W1	10"	STD	<input type="checkbox"/>	013W4	13"	NC		
<input type="checkbox"/>	013W2	11"	NC	<input type="checkbox"/>	013W5	14"	NC		
<input type="checkbox"/>	013W3	12"	NC						

4. SEAT UPHOLSTERY DEPTH

Regular Frame			
<input type="checkbox"/>	013SD1	8" Seat Depth	STD
<input type="checkbox"/>	013SD2	9" Seat Depth	NC
<input type="checkbox"/>	013SD3	10" Seat Depth	NC

Long Frame			
<input type="checkbox"/>	013SD4	11" Seat Depth	STD
<input type="checkbox"/>	013SD5	12" Seat Depth	NC
<input type="checkbox"/>	013SD6	13" Seat Depth	NC
<input type="checkbox"/>	013SD7	14" Seat Depth	NC

5. SEATING

<input type="checkbox"/>	013SL1		Seat Upholstery	STD
<input type="checkbox"/>	013SL6		Omit Seat Upholstery	NC
<input type="checkbox"/>	013CU1	K0669	2" Seat Cushion	\$120

Non-Folding

<input type="checkbox"/>	013B103		Regular Height 8"-10"	STD
<input type="checkbox"/>	013B104		Tall Height 10"-12"	NC

Forward-Folding

<input type="checkbox"/>	013B101		Regular Height 8"-10"	\$60
<input type="checkbox"/>	013B102		Tall Height 10"-12"	\$60

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	013JC20	JAY Bundled Cushion		Price List
		Cushion Name or Part #		
		Width (in.):	Depth (in.):	

6. BACKREST

<input type="checkbox"/>	013SR1		Stroller Handle	\$75
<input type="checkbox"/>	013BU14		Omit Back Upholstery	-\$15
<input type="checkbox"/>	013BP8		Omit Backrest & Back Upholstery	-\$35

JAY Backs (Optional)

Ordering JAY backs on-chair, See the JAY Price List. Requires Part #

<input type="checkbox"/>	013JB02	JAY Bundled Backrest		Price List
		Backrest Name or Part #		
		Width (in.):	Height (in.):	

7. FOOTREST

<input type="checkbox"/>	013F54		Regular 8-1/2" to 10-3/4"	STD
<input type="checkbox"/>	013F46		Long 11" to 13-1/4"	NC

<input type="checkbox"/>	013LS5		Leg Strap	NC
<input type="checkbox"/>	013F38		High Mount	\$55

8. ARMRESTS (T-POST)

Adjust in 1/4" increments

<input type="checkbox"/>	013AA12		Kidz Std 4-1/2" to 5-3/4"	STD
<input type="checkbox"/>	013AM14		Foam Pad	NC
<input type="checkbox"/>	013AM1		Desk Pad	NC

<input type="checkbox"/>	013AA13		Kidz Std 5-1/2" to 6-3/4"	NC
<input type="checkbox"/>	013AM14		Foam Pad	STD
<input type="checkbox"/>	013AM1		Desk Pad	NC
<input type="checkbox"/>	013AA14		Omit Armrest	NC

9. CASTERS**Caster Wheel**

<input type="checkbox"/>	013CW1		3" Micro	NC
<input type="checkbox"/>	013CW30		3" Micro Light up	\$65
<input type="checkbox"/>	013CW5		5" Low-Profile Polyurethane	STD

Caster Options

<input type="checkbox"/>	013CF4		3/4" Ext. Fork Stem Bolt	NC
<input type="checkbox"/>	013CW12		Omit Caster Wheels	NC

Rear Wheel Size

<input type="checkbox"/>	013RW13		12" Mag (Avail. w/ 3" casters only. 12" rear wheels on reg. frame are N/A with high pull wheel locks.)	NC
<input type="checkbox"/>	013RW14		16" Mag (N/A with full poly tires)	STD
<input type="checkbox"/>	013RW3		20" Mag Wheel	NC
<input type="checkbox"/>	013RW4		20" Spoke	\$100

Rear Wheel Tire

<input type="checkbox"/>	013RT1		Pneumatic	NC
<input type="checkbox"/>	013RT2	E2211 E2213	Pneumatic Airless Insert	\$80
<input type="checkbox"/>	013RT3		Full Polyurethane	NC

10. AXLES

Axles

<input type="checkbox"/>	013AX1	Quick Release Axles	STD
<input type="checkbox"/>	013AX11	Quad Quick Release Axles <i>(Only available w/ 20" rear wheels. N/A w/ 16" rear wheel reverse configuration.)</i>	\$125

Handrim

N/A with 12" or 16" rear wheels.

<input type="checkbox"/>	013HR1	Aluminum Anodized	NC
<input type="checkbox"/>	013HR4	Plastic Coated	\$95
<input type="checkbox"/>	013HR7	Omit Handrims	NC

Wheel Lock

<input type="checkbox"/>	013WL9	High Push	STD
<input type="checkbox"/>	013WL8	High Pull <i>(12" rear wheels on reg. frame are not avail. w/ high pull wheel locks. N/A w/ 16" rear wheel reverse configuration.)</i>	NC
<input type="checkbox"/>	013WL10	Omit	NC

Extension Handle

<input type="checkbox"/>	013EH1	E0961	6" Ext. Handles <i>(Avail. w/ high-push wheel locks only.)</i>	\$75
--------------------------	--------	-------	---	------

Anti-Tip (Standard Configuration)

<input type="checkbox"/>	013AT1	E0971	Rear Anti-Tip Tubes	\$140
<input type="checkbox"/>	013AO1	E0971	Front Anti-Tip Tubes	\$150

Anti-Tip (Reverse Configuration)

<input type="checkbox"/>	013AT19	E0971	Rear Anti-Tip Tubes	\$150
<input type="checkbox"/>	013AO1	E0971	Front Anti-Tip Tubes	\$150

11. ACCESSORIES

Upper Extremity Support Tray

<input type="checkbox"/>	013A29	E0950	X-Small <i>(Suggested 10"-12" Width)</i>	\$255
<input type="checkbox"/>	013A30	E0950	Small <i>(Suggested 13"-14" Width)</i>	\$255

Miscellaneous

<input type="checkbox"/>	013A83		Tool Kit	\$45
--------------------------	--------	--	----------	------

Spoke Guard

(N/A with 12" & 16" rear wheels)

<input type="checkbox"/>	188A23	K0065	Black	\$115
<input type="checkbox"/>	188A303	K0065	Clear Zippie	\$115



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA

Customer Service: 800-333-4000 Fax: 800-300-7502 Email: orders@sunmed.com Web: www.sunrisemedical.com

©1/2018 MK-100007 REV C