



MyWay+

AN OPEN FRAME GAIT TRAINER

LECKEY

Adventure in every step!

1. **OPEN FRAME GAIT TRAINER**
supports children in an upright position, improving stepping and maximizing opportunities for exploration, interaction, and participation.
2. **CORSET-LIKE HARNESS**
for maximized trunk support and upright, hands-free walking.
3. **RE-SIZED FOR EARLY INTERVENTION**
Size 1 Frame fits users as young as 12 months old.
4. **7" INDOOR/OUTDOOR CASTERS**
with patent-pending four-way directional lock allow easy maneuverability over different surfaces. The only gait trainer that offers **SIDE-STEPPING** feature.
5. **NARROW SADDLE & DROP ATTACHMENT**
reduces abduction while addressing scissoring gait.
6. **3-WAY ADJUSTABLE HEADREST**
with adjustable laterals for additional head support.
7. **3-WAY ADJUSTABLE USER HANDLES**
for additional forearm stability.



Learn more



Video & Resources

TECHNICAL DATA	
AGE RANGE	1 to 16 years
USER SHOULDER HEIGHT	23" to 59"
USER SADDLE HEIGHT	9.5" to 32"
FRAME MAX USER WEIGHT	Size 1 55 lbs. Size 2 110 lbs. Size 3 176 lbs.
UPPER SUPPORT UNIT WEIGHT	Size 1 N/A Size 2 23 lbs. Size 3 24 lbs.
FRAME WEIGHT	Size 1 26.5 lbs. Size 2 14.5 lbs. Size 3 21 lbs.
HARNESS MAX USER WEIGHT	Size 1 to 3 110 lbs. Size 4 to 6 176 lbs.
HARNESS WAIST MEASUREMENT	18" to 35.5"

Refer to order form as specifications may vary based on configuration.



This document can be viewed in PDF format at www.SunriseMedical.com.



MK-130108 Rev A ©09.2023



800.333.4000 www.SunriseMedical.com

