



Kid Kart TLC
January 2019
EITC

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirm Via: Fax Email

Mark For: _____

SUBMITTING FOR: Quote Order

PO#: _____

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

KID KART TLC

1. FRAME

<input type="checkbox"/>	100BS2	E1232	Kid Kart TLC with Tilt <i>(Includes standard planar seat/back)</i>	\$2,995
<input type="checkbox"/>	100BS6	E1234	Kid Kart TLC with Tilt No Seating	\$2,895

Frame Options				
<input type="checkbox"/>	100KJ8	E1014	Recline Feature	\$200
<input type="checkbox"/>	100A76	E1029	Ventilator and Battery Tray	\$440

Transit System

<input type="checkbox"/>	100KJ7	K0108	Transit Option	\$300
--------------------------	--------	-------	----------------	-------

2. SIZE/SEATING

Choose Standard OR Adjustable

STANDARD				
<input type="checkbox"/>	100ST8		Standard Cushion	NC
<input type="checkbox"/>	100ST12	E2294	Wedge Cushion	\$225
<input type="checkbox"/>	100ST13	E2294	Ledge Cushion - Medium Density Contoured Foam - 7" to 13-1/2" Hip Width - 3" Lateral Hip to Thigh Supports Not avail. w/o lateral support option - Slight Medial Thigh Contour - Fixed Ischial Ledge - 2" Seat Depth Growth Kit - Removable Neoprene Cover	\$225

CUSTOM ADJUSTABLE				
<input type="checkbox"/>	100ST11	E2294	Adjustable Contour Cushion - Medium Density Contoured Foam - Adjustable Pre-Ischial Ledge - Adjustable 1" High Foam Wedges for: Lateral Hip and Thigh Support Medial Knee and Thigh Support - 3" Seat Depth Growth Kit - Removable Neoprene Cover	\$320

OR

Seat Depth Options

<input type="checkbox"/>	100SD15		Short 7" through 9"	NC
<input type="checkbox"/>	100SD17		Medium 9" through 11"	NC
<input type="checkbox"/>	100SD19		Long 11" through 13"	NC

Seat Depth Options				
<input type="checkbox"/>	100SD16		Short 8" through 10"	NC
<input type="checkbox"/>	100SD18		Long 11" through 14"	NC

Foam Options

<input type="checkbox"/>	100FM1		Conform Soft Foam	NC
<input type="checkbox"/>	100FM2		Conform Medium/Soft Foam	NC
<input type="checkbox"/>	100FM3		Conform Medium Foam	NC

Client Actual Seat Depth: _____ inches
 Set Seat Depth at: _____ inches

Client Actual Seat Depth: _____ inches
 Set Seat Depth at: _____ inches

3. WHEELS AND CASTERS

<input type="checkbox"/>	100WQ01		12-1/2" Pneumatic Rear w/ 6" Casters	NC
<input type="checkbox"/>	100WQ02	E2211 E2213	12-1/2" Pneumatic Rear w/ Airless Inserts with 6" Casters	\$75

<input type="checkbox"/>	100WQ03	E2211 E2212	14" Pneumatic Rear w/ 8" Casters	\$135
<input type="checkbox"/>	100WQ04	E2211 E2213	14" Pneumatic Rear w/ Airless Inserts with 8" Casters	\$210

4. BACKREST

2.5" "T" Back used for short cushion seat depths less than 7.5" and long cushion seat depths less than 9.5"

Choose 1" OR 2-1/2"

<input type="checkbox"/>	100BT26		"I" Back Cushion 1" Thick	NC	<input checked="" type="checkbox"/>	100BT27		"T" Back Cushion 2-1/2" Thick	\$90.00
--------------------------	---------	--	---------------------------	----	-------------------------------------	---------	--	-------------------------------	---------

5. SUPPORTS

**Priced as each, selecting these options indicates Qty 2. Hardware Bracket must be chosen with any support option as pair.*

Lateral Hip Support

Not available w/standard cushion.

<input type="checkbox"/>	100KA1/7	E0956	SM (4" x 4", 3/8" thick)*	\$100
<input type="checkbox"/>	100KA2/8	E0956	LG (5-3/8" x 4-3/4", 1/2" thick)*	\$100
<input type="checkbox"/>	100HE68	E1028	Hardware Hip Support Brackets	\$50

Lateral Thigh Support

Not available w/standard cushion.

<input type="checkbox"/>	100KB1/9	E0953	SM (4" x 4", 3/8" thick)*	\$100
<input type="checkbox"/>	100KB2/10	E0953	LG (5-3/8" x 4-3/4", 1/2" thick)*	\$100
<input type="checkbox"/>	100HE69	E1028	Hardware Lateral Thigh Brackets	\$50

Medial Support

<input type="checkbox"/>	100KG1	E0957	SM (3"H x 2-1/2"W x 4.5"L)	\$175
<input type="checkbox"/>	100KG2	E0957	LG (3-3/4"H x 4"W x 6"L)	\$175
<input type="checkbox"/>	100HE70	E1028	Hardware Medial Support Bracket	\$50

Lateral Trunk Support

Fixed Laterals (Adjustable 1"-10" widths)

<input type="checkbox"/>	100KE34		Fixed Angle Adjustable - Right	NC
<input type="checkbox"/>	100KE40		Fixed Angle Adjustable - Left	NC
<input type="checkbox"/>	100KD26	E0956	Large Straight Pad 3" x 6" - Right	\$100
<input type="checkbox"/>	100KD24	E0956	Large Straight Pad 3" x 6" - Left	\$100

Lateral Trunk Support - Cont.

Swing Away Laterals (Adjustable 1"-8" widths)

<input type="checkbox"/>	100KE39		Swing Away Bracket Adjustable-Right	\$175
<input type="checkbox"/>	100KE35		Swing Away Bracket Adjustable-Left	\$175
<input type="checkbox"/>	100KD26	E0956	Large Straight Pad 3" x 6" - Right	\$100
<input type="checkbox"/>	100KD24	E0956	Large Straight Pad 3" x 6" - Left	\$100
<input type="checkbox"/>	100KD25	E0956	Large Curved Pad 3" x 6" - Right	\$100
<input type="checkbox"/>	100KD22	E0956	Large Curved Pad 3" x 6" - Left	\$100

Anterior Trunk Support

Dimensions are seat to shoulder height.

<input type="checkbox"/>	100KI1	E0960	XX Small Butterfly Pad (< 10")	\$110
<input type="checkbox"/>	100KI2	E0960	X Small Butterfly Pad (10" - 12")	\$110
<input type="checkbox"/>	100KI3	E0960	Small Butterfly Pad (12" - 14")	\$110
<input type="checkbox"/>	100KI4	E0960	Medium Butterfly Pad (14" - 16")	\$110
<input type="checkbox"/>	100KI5	E0960	Large Butterfly Pad (>16")	\$110

Trunk Support Options

<input type="checkbox"/>	100SS1	E0960	Shoulder Strap Pads	\$25
--------------------------	--------	-------	---------------------	------

6. KEY OPTIONS

Head Support Options

Must choose Headrest Hardware from Head Support Options

<input type="checkbox"/>	100HE5	E0955	Planar	\$135.00
<input type="checkbox"/>	100HE6	E0955	Tri-Wing	\$185.00
<input type="checkbox"/>	100HE8	E0955	Micro Adj. Incl. depth brackets	\$295.00
<input type="checkbox"/>	100HQ05	E1028	Headrest hardware	\$50.00
<input type="checkbox"/>	100HQ04		Head Support Halo Strap (N/A with Planar)	\$110.00
<input type="checkbox"/>	100HQ01	E1028	Depth Adj. Brackets	\$50.00

Shoe Holders

<input type="checkbox"/>	100A68	E0952 E0951	6" Long Shoe Holder - Pair	\$120
<input type="checkbox"/>	100A69	E0952 E0951	8" Long Shoe Holder - Pair	\$120

Foot Straps

<input type="checkbox"/>	100A79	E0952	16" Foot Straps with Covers - Pair	\$35
<input type="checkbox"/>	100A80	E0952	18" Foot straps with Covers - Pair	\$35

Upper Extremity Support Tray

<input type="checkbox"/>	100A41	E0950	ABS (Light Grey)	\$125
<input type="checkbox"/>	100A39	E0950	Polycarbonate (Clear)	\$190

Pelvic Positioning Belt

Must use when ordering Anterior Trunk Supports.

<input type="checkbox"/>	100A14	E0978	1" Padded	\$45
<input type="checkbox"/>	100A15	E0978	1-1/2" Padded	\$45
<input type="checkbox"/>	100A263	E0978	1" Padded Auto Style	\$55
<input type="checkbox"/>	100A16	E0978	1-1/2" Padded Auto Style	\$55
<input type="checkbox"/>	100PP1		Set at 45°	NC
<input type="checkbox"/>	100PP2		Set at 90° (90° for Adjustable Contour Cushion only)	NC

Floor Sitter Feeder Base

<input type="checkbox"/>	100A81		Folding Fixed Height Feeder Base	\$675
--------------------------	--------	--	----------------------------------	-------

Life Support Accessories

<input type="checkbox"/>	100A74	K0105	IV Support Pole	\$145
<input type="checkbox"/>	100A711	E2208	O2 Tank Holder	\$145
<input type="checkbox"/>	100A75	K0105 E2208	IV / O2 Combination	\$220

7. COLORS

Package includes: Silver Frame and Black Upholstery

Color Package

Canopy, shell and catch all

<input type="checkbox"/>	100CL9		Purple	NC
<input type="checkbox"/>	100CL10		Soft Pink	NC
<input type="checkbox"/>	100CL11		Green	NC
<input type="checkbox"/>	100CL12		Navy Blue	NC
<input type="checkbox"/>	100CL13		Orange	NC
<input type="checkbox"/>	100CL14		Butterflies	NC
<input type="checkbox"/>	100CL15		Red Race Flags	NC

Canopy Pattern

<input type="checkbox"/>	100UF01		Cross Over	NC
<input type="checkbox"/>	100UF02		Center Stripe	NC



Sunrise Medical (US) LLC

2842 Business Park Ave. · Fresno, CA 93727 · USA