



# Quickie M6™ January 2016

EIF1

**Mark For:**

**Submitting for:**

Quote     Order

Date: \_\_\_\_\_  
 Dealer Acct #: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: (    )                      Fax: (    )  
 Confirmation Email: \_\_\_\_\_  
 Confirm Via:     Fax             Email

PO#: \_\_\_\_\_  
**ADDITIONAL SHIPPING INFORMATION**  
 Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: (    )                      Fax: (    )

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.  
 The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓ Option #    HCPCS                      Definition                      Price                      ✓ Option #    HCPCS                      Definition                      Price

**1. FRAME**

**Base Model**

<input type="checkbox"/>	115WM21	K0007	Quickie M6	\$2,945.00
--------------------------	---------	-------	------------	------------

**2. SIZE / SEATING**

**Frame Width**

<input type="checkbox"/>	115W13	E2201	22"	\$500.00
<input type="checkbox"/>	115W14	E2202	24"	\$500.00
<input type="checkbox"/>	115W15	E2202	26"	\$550.00
<input type="checkbox"/>	115W16	K0108	28"	\$550.00
<input type="checkbox"/>	115W17	K0108	30"	\$550.00

**Frame Length**

<input type="checkbox"/>	115FL9		Regular (18"-19" depth)	NC
<input type="checkbox"/>	115FL18		Long (18"-22" depth)	NC

**Seat Depth**

<input type="checkbox"/>	115SD11		18"	NC
<input type="checkbox"/>	115SD12		19"	NC
<input type="checkbox"/>	115SD13	E2203	20" (Long frame only)	\$550.00
<input type="checkbox"/>	115SD14	E2203	21" (Long frame only)	\$550.00
<input type="checkbox"/>	115SD33	E2204	22" (Long frame only)	\$900.00

**Seat Height**

<input type="checkbox"/>	115HT65		17" Includes low adjustable axle plate and shorter fork. Only available with 5" caster	\$300.00
<input type="checkbox"/>	115HT11		18.5" Only available w/ 5" casters	NC
<input type="checkbox"/>	115HT14		20"	NC

**Cushion**

<input type="checkbox"/>	115CU31	E2602	2" Cushion	\$120.00
--------------------------	---------	-------	------------	----------

**JAY Cushions (Optional)**

To order any JAY cushions on-chair, please see the JAY price list.  
 Cushion name or part number must be filled out with desired width and depth

<input type="checkbox"/>	115JC20	JAY Bundled Cushion		Price List
		Cushion Name or Part #		
		Width: _____	Depth: _____	

**3. CASTERS**

**Casters**

<input type="checkbox"/>	115CW6		5" Soft Roll	\$85.00
<input type="checkbox"/>	115CW36	E2214 E2215	8" x 2" Pneumatic	\$145.00
<input type="checkbox"/>	115CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$225.00

**Seat Height Matrix**

Caster	Front/Rear
5"	17" (Different Axle & Fork)
5"	18.5"
8"	20"

**4. REAR WHEELS**

**Rear Wheel Type**

<input type="checkbox"/>	115RW2		24" Spoke Wheel	NC
--------------------------	--------	--	-----------------	----

**Rear Wheel Tire**

<input type="checkbox"/>	115RT1	E2211 E2212	Pneumatic Tire	NC
<input type="checkbox"/>	115RT2	E2211 E2213	Pneumatic Airless Insert Tire	\$80.00

**Handrim**

<input type="checkbox"/>	115HR1		Aluminum Anodized	NC
<input type="checkbox"/>	115HR4		Plastic Coated	\$95.00
<input type="checkbox"/>	115HR7		Omit Handrim	NC

**Adjustable Axle Plate**

<input type="checkbox"/>	115AL2		Standard (for 18.5", 20" Seat Ht.)	NC
<input type="checkbox"/>	115AL28		Low (only for 17" Seat Ht.)	NC

**Axle**

<input type="checkbox"/>	115AX1		Quick Release	NC
--------------------------	--------	--	---------------	----

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
<b>5. WHEEL LOCKS</b>									
<b>Wheel Lock</b>					<b>Extension Handle</b>				
<input type="checkbox"/>	115WL1		High Mount Push Wheel Lock	NC	<input type="checkbox"/>	115EH1	E0961	6" Pair	\$75.00
<input type="checkbox"/>	115WL12		Do Not Mount Wheel Lock	NC					
<input type="checkbox"/>	115WL10		Omit Wheel Lock	NC					
<b>6. HANGERS AND FOOTREST</b>									
<b>Front Hangers</b>					<b>Extension Tubes</b>				
<input type="checkbox"/>	115H1		60° Swing Away	NC	<input type="checkbox"/>	115E1		Short (13"-15.5")	NC
<input type="checkbox"/>	115H146		70° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$95.00	<input type="checkbox"/>	115E2		Medium (15.5"-18")	NC
<input type="checkbox"/>	115H132		80° Heavy-Duty Ft. Mount Lift Off <i>Only available with 5" caster and ang adj footplates. NA w/ extension tubes</i>	\$95.00	<input type="checkbox"/>	115E5		Long (18"-20.5") NA w/ 5" caster	NC
<input type="checkbox"/>	115H41	E0990	ELR - Heavy-Duty Reinforced <i>Only available with aluminum footplates</i>	\$595.00					
<b>Lower Leg Length Range Matrix</b>									
<b>Hanger</b>		<b>Ext Tube</b>		<b>Footplate</b>					
				<b>Composite</b>	<b>Alum Angle Adj</b>	<b>Aluminum</b>			
60°				Short	11.50" - 13.25"	11.75" - 14.00"	11.25" - 13.50"		
				Medium	13.75" - 16.00"	14.00" - 16.25"	13.50" - 14.00"		
				Long	16.00" - 16.50"	16.00" - 17.75"	NA		
70°		No Ext Tube	NA	9.75" - 13.50"	NA				
80°		No Ext Tube	NA	10.75" - 14.75"	NA				
ELR				Short	NA	NA	13.00" - 17.00"		
				Medium	NA	NA	14.25" - 17.00"		
								<i>All values are ± 0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights. ELR values based on 70° hanger angle</i>	
<b>Footplates</b>									
<i>Heel loop or leg strap not included, See Footplate Options Section</i>									
<input type="checkbox"/>	115F1		Composite	NC					
<input type="checkbox"/>	115F25	K0040	Adult Angle Adj	\$185.00					
<input type="checkbox"/>	115F5		Aluminum	NC					
<b>Footplate Options</b>									
<input type="checkbox"/>	115HL1	E0951	Heel Loop	\$30.00					
<b>7. BACKREST</b>									
<b>Backrest</b>					<b>Back Height Preset</b>				
<input type="checkbox"/>	115BT34		8° w/ Push Handles <i>Adjustable (16-3/4" to 20-1/8")</i>	NC	<input type="checkbox"/>	115B91		16-3/4"	NC
					<input type="checkbox"/>	115B96		17-7/8"	NC
					<input type="checkbox"/>	115B39		19"	NC
					<input type="checkbox"/>	115B98		20-1/8"	NC
<b>8. ARMREST</b>									
<b>Armrest</b>									
<input type="checkbox"/>	115AA33	E0973	Dual Post Height Adj - Desk 10"	\$245.00					
<input type="checkbox"/>	115AA34	E0973	Dual Post Height Adj - Full 14"	\$245.00					
<b>9. KEY OPTIONS &amp; ACCESSORIES</b>									
<b>Rear Anti-Tip</b>					<b>Positioning Belt</b>				
<input type="checkbox"/>	115AT1	E0971	Rear Anti-Tip	\$140.00	<input type="checkbox"/>	115A12	E0978	Belt 2" Aircraft Buckle	\$110.00
<input type="checkbox"/>	115AT3		Omit Rear Anti-Tip	NC	<input type="checkbox"/>	115A13	E0978	Belt 2" Aircraft Buckle Padded	\$120.00
<b>Miscellaneous</b>									
<input type="checkbox"/>	115A83		Tool Kit	\$45.00	<input type="checkbox"/>	115A142	E0978	Belt 1.5" Auto Buckle	\$70.00
<b>10. COLOR</b>									
<b>Standard Frame Colors</b>					<b>Matte Frame Colors</b>				
<input type="checkbox"/>	115C140		Aztec Gold	NC	<input type="checkbox"/>	115C116		Matte Black	NC
<input type="checkbox"/>	115C2		Black <i>(Not avail as base w/Zebra)</i>	NC	<input type="checkbox"/>	115C151		Matte Black Cherry	NC
<input type="checkbox"/>	115C141		Black Cherry	NC	<input type="checkbox"/>	115C150		Matte Electric Blue	NC
<input type="checkbox"/>	115C10		Black Opal	NC	<input type="checkbox"/>	115C149		Matte Evergreen	NC
<input type="checkbox"/>	115C38		Blue Opal	NC	<input type="checkbox"/>	115C152		Matte Purple	NC
<input type="checkbox"/>	115C36		Candy Blue	NC	<b>Custom Frame Color</b>				
<input type="checkbox"/>	115C14		Candy Purple	NC	<input type="checkbox"/>	115C112		Desert Camouflage	\$295.00
<input type="checkbox"/>	115C35		Candy Red	NC	<input type="checkbox"/>	115C43		Mossy Oak Camouflage	\$295.00
<input type="checkbox"/>	115C142		Electric Blue	NC	<input type="checkbox"/>	115C42		Stars & Stripes	\$295.00
<input type="checkbox"/>	115C143		Evergreen	NC	<input type="checkbox"/>	115C108		Zebra <i>(Choose above base color) Std with White unless above base paint is chosen NA w/ Bkck, Matte Black, Bkck Opal or Titanium)</i>	\$295.00
<input type="checkbox"/>	115C126		Glow	NC					
<input type="checkbox"/>	115C144		Green Apple	NC					
<input type="checkbox"/>	115C147		Hot Sparkle Pink	NC					
<input type="checkbox"/>	115C148		Mauve Pink	NC					
<input type="checkbox"/>	115C103		Rootbeer	NC					
<input type="checkbox"/>	115C114		Sparkle Silver	NC					
<input type="checkbox"/>	115C145		Sunrise Orange	NC					
<input type="checkbox"/>	115C113		Titanium Color Paint	NC					
<input type="checkbox"/>	115C8		Yellow	NC					

