



ORDER FORM



To see QUICKIE Power Options
or Accessories, scan here or visit:
[https://www.sunrisemedical.com/
power-wheelchair-options-accessories](https://www.sunrisemedical.com/power-wheelchair-options-accessories)




Q700-UPM SEDEO® PRO ADVANCED™ EIPW25

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		State	Zip
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Seating Specifications

HCPC	K0861 / K0884
Weights Starting at: (with Batteries)	435 lbs
Power Seating	25° of CG Tilt and Power Recline with 7" mechanical sheer reduction / 8" seat elevation / 30° Anterior Tilt
Seat Type	SEDEO PRO ADVANCED
Seat Size Range	Width 16" to 20" Depth 17" to 22"
Seat-to-Floor Height Range	18" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	265 lbs. Standing
Base Size	Width 24.5", Length 36.5"
Turning Radius	19.8" (Excludes Footplates)
Drive Control	R-net PM 120A, R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24
Color	9 shroud colors and 7 color accents
Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Q700™ M Base <i>End Item #EIPW25</i>			\$
250PC335	K0861	QUICKIE Q700-UP M MPC Tilt/Recline for SEDEO PRO ADVANCED Included in Base Pricing: Q700 M Base, R-net 90amp electronics, ProGrip™ with 6-wheel independent suspension , 14" drive wheels with dual density inserts, 6" flat free casters, ISO Occupied Tie-Down brackets and 10amp charger. max speed = 6mph, weight capacity = 265lbs, STFH = 17", 3" ground clearance.	\$8,979.00
250PC336	K0884	QUICKIE Q700-UP M MPD Tilt/Recline for SEDEO PRO ADVANCED Included in Base Pricing: Q700 M Base, R-net 120amp electronics, ProGrip™ 6-wheel independent suspension , 14" drive wheels with dual density inserts, 6" flat free casters, ISO Occupied Tie-Down brackets and 10amp charger. max speed = 8mph, weight capacity = 265lbs, STFH = 17", 3" ground clearance.	\$14,297.00
User Weight			\$
250WG24	—	Up to 265 lbs.	NC
Motor Package			\$
250PF99	—	4-pole 10 kph / 6 mph (Standard for Group 3 K0861)	NC
250PF100	—	4-pole 12.5 kph / 8 mph (Standard for Group 4 K0884)	NC
250PF101	—	4-pole 12.5 kph / 8 mph (Upgrade for Group 3 K0861)	\$1,316.00
Tracking			\$
250MI01	K0108	SureTrac® Module & Motor Encoders (Improves drive tracking.)	\$1,596.00
Drive Wheels  P. 17			\$
250DW58	—	Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels  P. 17			\$
250CW55	—	Caster Wheels 6", Solid, Black	NC
Shroud Color			\$
250C119	—	Matte Black	NC
250C165	—	Candy Apple Red	NC
250C166	—	Candy Blue	NC
250C167	—	Candy Purple	NC
250C176	—	Green Apple	NC
250C178	—	Pearl Pink	NC
250C182	—	Silver	NC
250C125	—	Pearl White	NC
250C170	—	American Flag	\$590.00

BASE

OPTION	HCPC	DESCRIPTION	RETAIL
Color Accents			\$
250EZ20	—	Apple Green	NC
250EZ22	—	Black	NC
250EZ21	—	Blue	NC
250EZ23	—	Orange	NC
250EZ24	—	Red	NC
250EZ19	—	White	NC
250EZ25	—	Silver	NC
Service Decals			\$
250D10	—	Army	NC
250D11	—	Air Force	NC
250D13	—	Marines	NC
250D14	—	Navy	NC
250D15	—	Purple Heart	NC
ISO Transit Tie Downs			\$
250T050	—	ISO Transit Package	NC
LED Lights  P. 17			\$
250LI04	—	LED Lights and Indicators	\$1,373.00
250LI09	—	LED Fender Lights	NC
250LI02	—	Omit LED Lights	NC
Seat Heights			\$
<i>Measured from front of seat pan to floor.</i>			
250HT56	—	18"	NC
250HT105	—	18.5"	NC
250HT57	—	19"	NC
Batteries and Chargers			\$
Group 24 Battery Options (Recommended for max daily range)			
250AG25	E2363	GP24 Flat Top Pair w/Strap MK Gel	\$1,116.00
250AG26	E2363	GP24 Flat Top Pair w/Strap Brightway AGM	\$840.00
Chargers			
250BJ28	—	10A Dual Mode (Compact)	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Power Seating Options	P. 17	\$
-----------------------	---	----

Seat Elevator (Required)

250PI06	E2298	Seat Elevator, Powered, 8"	\$3,717.00
---------	-------	----------------------------	------------

Seat Tilt/Recline (Required)

250PA51	E1007	Posterior 0° - 25° Tilt and 180° Recline Powered, w/ 7" Shear Reduction	\$14,292.00
---------	-------	---	-------------

Seat Anterior Tilt / Standing (Must Choose One)

250TG03	E2301	Standing	\$12577.00
---------	-------	----------	------------

Back Frame Type	\$
-----------------	----

250BF09	—	Standing Seat Frame	NC
---------	---	---------------------	----

Back Cushion Type	\$
-------------------	----

250BT138	—	JAY Syncra® <i>(Choose Seating in section B)</i>	NC
250BT135	—	JAY® J3® <i>(Choose Seating in section C)</i>	NC
250BT137	—	SEDEO Aftermarket Adapter <i>(Height adjustable round tubes for traditional mounting)</i>	\$515.00
250BT133	—	OMIT <i>(Shear back plate will have 4 mounting posts for personalized mounting. No headrest mount.)</i>	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

B - JAY® Syncra™ Seating**B - JAY® Syncra™ Back** \$

250JP14	E2620	JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	\$1,116.00
250JP15	E2621	JAY® Syncra™ Back 22" Wide with 3" Contour Depth	\$1,230.00

B - JAY® Syncra™ Back Width \$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

250JQ07	—	16" Wide (Adjustment range 16 to 17" Wide)*	NC
250JQ08	—	18" Wide (Adjustment range 18 to 19" Wide)*	NC
250JQ09	—	20" Wide (Adjustment range 20 to 21" Wide)*	NC
250JQ10	—	22" Wide (Adjustment range 22 to 23" Wide)	NC

**Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.*

B - JAY® Syncra™ Back Shell Options and Height Ranges \$

JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown based on the back frame type.

JAY® Syncra™ Shell Height

250JR42	—	16" Tall Shell (22" to 24" Back Height)	NC
250JR43	—	18" Tall Shell (24" to 26" Back Height)	NC
250JR44	—	20" Tall Shell (26" to 28" Back Height)	NC
250JR45	—	22" Tall Shell (28" to 30" Back Height)	NC

JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)

250JR50	—	22" (16" Shell Height)	NC
250JR51	—	23" (16" Shell Height)	NC
250JR52	—	24" (16" Shell Height)	NC
250JR53	—	24" (18" Shell Height)	NC
250JR54	—	25" (18" Shell Height)	NC
250JR55	—	26" (18" Shell Height)	NC
250JR56	—	26" (20" Shell Height)	NC
250JR57	—	27" (20" Shell Height)	NC
250JR58	—	28" (20" Shell Height)	NC
250JR67	—	28" (22" Shell Height)	NC
250JR68	—	29" (22" Shell Height)	NC
250JR69	—	30" (22" Shell Height)	NC

B - JAY® Syncra™ Back Upholstery \$

250BU96	—	Syncra™ Premium Leatherette	NC
250BU97	—	Syncra™ Stretch	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

B- JAY J3 Back \$
B- JAY J3 Back – Contour P. 17

250JP01	E2613	J3 PA Back for 14"-20" Wide (J3PAXXXYYZZ)	\$613.00
250JP11	E2615	J3 PL Back for 14"-20" Wide (J3PLXXYYZZ)	\$715.00
250JP12	E2620	J3 PDC Back for 14"-20" Wide (J3PDCXXYYZZ)	\$858.00
250JP13	E2620	J3 PD Back for 14"-20" Wide (J3PDXXYYZZ)	\$858.00

B- JAY J3 Back – Width *Back can be ordered a different width than seat pan.* \$

250JQ01	–	14"	NC
250JQ02	–	16"	NC
250JQ03	–	18"	NC
250JQ04	–	20"	NC

B- JAY J3 Shell Options and Height Ranges \$

JAY J3 shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear.
Back height can be adjusted in the field within the ranges shown.

JAY J3 Shell Height**PRO ADVANCED**

250JR41	–	15" Tall Shell (Mid Thoracic Tall)	21"-24"	Installed Back Height	NC
250JR14	–	16.5" Tall Shell (Upper Thoracic Short)	21"-24"		NC
250JR15	–	18" Tall Shell (Upper Thoracic Med)	22"-26"		NC
250JR16	–	19.5" Tall Shell (Upper Thoracic Tall)	23"-27"		NC
250JR18	–	18" tall shell (<i>Shoulder Height</i>) for 22" wide	21"-25"		NC
250JR17	–	21" Tall Shell (<i>Shoulder Height</i>)	25"-28"		NC
250JR27	–	22.5" tall shell (<i>Shoulder Height</i>)	26"-29"		NC
250JR28	–	24" tall shell (<i>Shoulder Height</i>)	28"-30"		NC

JAY J3 Installed Back Height (From Pan to Top of Back)

250JR29	–	21"	NC
250JR19	–	22"	NC
250JR20	–	23"	NC
250JR21	–	24"	NC
250JR22	–	25"	NC
250JR23	–	26"	NC
250JR24	–	27"	NC
250JR25	–	28"	NC
250JR30	–	29"	NC
250JR31	–	30"	NC

B- JAY Anterior Thoracic Support P. 18 *Classic-Structure Pull Down Harness (4 point strap, only avail. w/JAY backrest)* \$

250K10	E0960	Medium (ATS8302)	\$174.00
250K11	E0960	Large (ATS8303)	\$174.00

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)
\$
JAY® Lateral Thoracic Support Pads

250JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$144.00
250JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$144.00
250JW36	E0956	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	\$144.00
250JW39	E0956	Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	\$144.00
250JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$144.00
250JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$144.00

JAY® Lateral Thoracic Support Pad Upholstery

250NU01	—	Premium Leatherette	NC
250NU02	—	Stretch	NC

JAY® Lateral Thoracic Swing-Away Mount

250JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$286.00
250JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$286.00

Headrests (Requires Pad, Upholstery, & Mount)

P. 18
\$
Headrest Pad

250HE36	E0955	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$184.00
250HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$234.00
250HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$234.00
250HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$245.00
250HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (<i>Adjustable from 7 to 16" Wide x 4" Tall</i>)	\$360.00

Headrest Pad Upholstery

250NU10	—	SEDEO® Upholstery	NC
250NU04	—	Syncra™ Premium Leatherette	NC
250NU05	—	Syncra™ Stretch	NC
250NU11	—	Whitmyer® Lycra	NC

Headrest Mount

250MH32	E1028	SEDEO® Multi-Position Headrest Mount with Quick Release (<i>Only available with SEDEO® Pad</i>)	\$229.00
250MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$188.00
250MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$229.00
250MH35	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$121.00

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
Push Handles			\$
250PH23	—	Push Handles <i>(Not available with Q2 holder)</i>	NC
Seat Frame Rails			\$
250SN14	—	Seat Rail Short <i>(For effective seat depth 14" to 18")</i>	NC
250SN15	—	Seat Rail Long <i>(For effective seat depth 17" to 22")</i>	NC
Seat Cushion Type			\$
250ST81	—	JAY® Syncra® Cushion (See Section E)	NC
250ST74	—	JAY® (See Section D)	NC
250ST75	—	Omit	NC
C - SEDEO PRO ADVANCED Seat Size and D - JAY® Seat Size			\$
Seat Pan Width (Seat cushion width is dependent upon bolster selection. See thigh bolster section for limitations.)			
250R53	—	16"	NC
250R54	—	17"	NC
250R55	—	18"	NC
250R56	—	19"	NC
250R57	—	20"	NC
Seat Pan Depth (Front of Seat Pan to Front of Back Cushion Effective Depth)			
250RF45	—	17"	NC
250RF46	—	18"	NC
250RF47	—	19"	NC
250RF48	—	20"	NC
250RF49	—	21" <i>(Not available with JAY® Syncra® Cushion)</i>	NC
250RF50	—	22"	NC
E - JAY® Syncra™ Cushion (JAY Your Way modifications available through Customer Service (800) 333-4000)			\$
250CU89	E2622	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology <i>(16" to 20" wide)</i>	\$618.00
JAY® Syncra™ Cushion Cover Upholstery			
250YR73	—	Syncra™ Premium Leatherette	NC
250YR74	—	Syncra™ Stretch	NC
D- JAY Seat Cushion  P. 18			\$
250CU59	E2607	JAY Ion® (Sizes 16" to 20" Wide)	\$451.00
250CU61	E2607	JAY Union® (Sizes 16" to 20" Wide)	\$510.00
250CU64	E2622	JAY Fusion® (Sizes 16" to 20" Wide)	\$504.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

JAY® Lateral Thigh Supports (Must pick pad, upholstery, and mount) \$

JAY® Lateral Thigh Support Pads

250KF16	E0953	Lateral Thigh Support Pad Standard, Right (8.25" x 3.25")	\$144.00
250KF17	E0953	Lateral Thigh Support Pad Standard, Left (8.25" x 3.25")	\$144.00
250KF18	E0953	Lateral Thigh Support Pad Large, Right (12" x 3.25")	\$144.00
250KF19	E0953	Lateral Thigh Support Pad Large, Left (12" x 3.25")	\$144.00

JAY® Lateral Thigh Support Pad Upholstery

250NU07	—	Syncra™ Premium Leatherette	NC
250NU08	—	Syncra™ Stretch	NC

JAY® Lateral Thigh Removable Mount

250KF20	E1028	Lateral Thigh Removable Mount - Right	\$286.00
250KF21	E1028	Lateral Thigh Removable Mount - Left	\$286.00

Traditional Knee and Thigh Swing-Away Support Pads \$

Lateral Knee Support Pads and Hardware (Must pick pad and mount)

250XY05	E0953	Lateral Knee Support Pad, Right (4" round)	\$201.00
250XY06	E0953	Lateral Knee Support Pad, Left (4" round)	\$201.00
250XX08	E1028	Lateral Swing-Away Bracket Knee - Right	\$263.00
250XX09	E1028	Lateral Swing-Away Bracket Knee - Left	\$263.00

Lateral Thigh Support Pads and Hardware (Must pick pad and mount)

250KF22	E0953	Lateral Thigh Support Pad Short, Right (5" H x 5" L)	\$144.00
250KF23	E0953	Lateral Thigh Support Pad Short, Left (5" H x 5" L)	\$144.00
250KF24	E0953	Lateral Thigh Support Pad Long, Right (5" H x 7" L)	\$144.00
250KF25	E0953	Lateral Thigh Support Pad Long, Left (5" H x 7" L)	\$144.00
250KF26	E1028	Lateral Thigh Removable Mount - Right	\$286.00
250KF27	E1028	Lateral Thigh Removable Mount - Left	\$286.00

Armrest Structure \$

250EG51	—	Reclining, Flip-Back, Dual Link Armrest (Height adjustable 8" to 15" including length and lateral adjustment)	NC
---------	---	---	----

Arm Pads P. 18 \$

SEDEO® Armrest Pad

250AB43	—	300 mm / 12", Integral Foam, Right	NC
250AE43	—	300 mm / 12", Integral Foam, Left	NC
250AB44	—	400 mm / 16", Integral Foam, Right	NC
250AE44	—	400 mm / 16", Integral Foam, Left	NC

Comfort Arm Pad (Foam)

250AB46	—	2" Wide Right	\$86.00
250AE46	—	2" Wide Left	\$86.00

Comfort Arm Pad (Gel)

250AB48	—	2" Wide Right	\$172.00
250AE48	—	2" Wide Left	\$172.00

Ottobock® Arm Troughs

250AE06	—	Ottobock Arm Trough, Med, Left	\$133.00
250AB06	—	Ottobock Arm Trough, Med, Right	\$133.00
250AE07	—	Ottobock Arm Trough, Large, Left	\$133.00
250AB07	—	Ottobock Arm Trough, Large, Right	\$133.00
250OL01	—	Ottobock Hand Pad, Med, Left	\$58.00
250OR01	—	Ottobock Hand Pad, Med, Right	\$58.00
250OL02	—	Ottobock Hand Pad, Large, Left	\$58.00
250OR02	—	Ottobock Hand Pad, Large, Right	\$58.00
250AE78	—	Ottobock 1 Piece Arm Trough, Left	\$186.00
250AB78	—	Ottobock 1 Piece Arm Trough, Right	\$186.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Centermount Legrest			\$
250H173	E1012	Power Centermount Standing Legrest, 90° - 0° Angle, 7.5" Extension (12.5" to 17.5" Leg Length)	\$3,717.00
Footplates			
250F127	K0040	Dual, Height and Angle Adjustable, Standard (13.5" to 14.5" wide X 10" long)	\$193.00
250F128	K0040	Dual, Height and Angle Adjustable, Large (15.5" to 16.5" wide X 12" long)	\$193.00
Calf Pad Options			
250DI34	—	Calf Pads, Pair, Standard	NC
Knee Block (Required for UP)			
250DI37	E0957	Knee Blocks: Dual, Height, Width and Angle Adjustable	\$406.00
250DI55	E1028	Knee Block Quick-Release Removable Hardware	\$138.00
Pelvic Positioning Belts  P. 18			\$
250A835	—	Positioning Belt, Aircraft Buckle 50mm x 1300mm / 2" x 50" (Included with UP)	NC
250A836	E0978	Positioning Belt, Aircraft Buckle, Padded 50mm x 1300mm / 2" x 50"	\$86
Chest Straps			\$
250ZR06	—	Chest Strap, 2-Point, Padded (Included with UP)	NC
250ZR07	—	Additional Chest Strap, 2-Point, Padded	\$286.00
Chest Bars			\$
250ZR08	—	Chest Bar, Detachable, Round Pad	\$567.00
250ZR09	—	Chest Bar, Detachable, Flat Pad, 200 mm x 100mm	\$603.00
Accessories			\$
Ventilator Option			
250VT11	E1029	Vent Tray Fixed (Pulmonetic LTV Models)	\$531.00
O2 Tank Option			
250A703	E2208	O2 Tank Holder to Backrest Styles (Not available with fixed push handles)	\$372.00
Wheelchair Tray Table			
2507A08	E0950	Tray, Swing-Away	\$400.00
Backpack			
250A1	—	Black Utility	\$62.00
Misc. Accessories			
2508A52	—	Ram Mount Phone	\$98.00
2508A53	—	Ram Mount Tablet	\$400.00
2508A54	—	Ram Mount Cup Holder	\$98.00
2508A55	—	Rear View Mirror	\$98.00
2508A59	—	Transfer Bar (Pair)	\$400.00
2508A73	—	Armrest Bag	\$98.00
2509A10	—	Utility Hooks	\$59.00

JOYSTICK / ELECTRONICS *See Last Pages For Specialty Input Devices images*

OPTION	HCPC	DESCRIPTION	RETAIL
R-net Joystick (Remote)  P. 20			\$
250J116	—	Color Joystick (2.3" Color LCD Screen Display, CJM1)	NC
250J118	—	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and Dual (iOS/Android) Bluetooth Control	\$1,206.00
250J119	—	Omit Remote Control (Required with arm mounted specialty controls and table top control)	NC
Joystick Mounting			\$
250JM31	—	Fixed on Armrest	NC
250JM32	E1028	Parallel Swing-Away	\$338.00
250JM33	E1028	Link-It Mid-Line Mount on Armrest (Swing-Away)	\$972.00
Joystick Mount (Location)			\$
250HM16	—	Right Hand Armrest	NC
250HM17	—	Left Hand Armrest	NC
Display Language Selection R-net			\$
250L1	—	English	NC
250L3	—	French	NC
250L7	—	Spanish	NC
Joystick Handle (Optional)  P. 19			\$
250JH13	E2323	Ball Handle	\$144.00
250JH15	E2323	T-Handle	\$144.00
250JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	\$187.00
250JH25	E2323	Bodypoint – U-Shaped Joystick Handle (4 inch)	\$187.00
250JH26	E2323	Bodypoint – Grooved Mushroom Joystick	\$187.00
Electronic Programmers  P. 19			\$
250PR10	—	R-net PC Programmer w/ Dongle	\$601.00
R-net Power Modules  P. 19			\$
250ER17	—	Power Module, R-net, 90A, Expandable	\$1,425.00
250EX15	E2377	Expandable 90A Controller - R-net	\$749.00
250EX16	E2313	Harness for 90A Expandable Controller	\$676.00
250ER10	—	Power Module, R-net, 120A, Expandable (Only available with Group 4)	\$1,425.00
250EX09	E2377	Expandable 120A Controller - R-net	\$749.00
250EX10	E2313	Harness for 120A Expandable Controller	\$676.00
R-net Seating Module Control			\$
250E039	E2311	UP Actuator control module, R-Net, CxSM, Thru Drive 2+ Actuators	\$4,519.00
Attendant Controls  P. 19			\$
2500A05	E2331	Dual Control Module for Attendant	\$1,487.00
2500A06	E2331	R-net Color LCD (CJSM1) Control Attendant Only	\$1,487.00
Attendant Mounting Location			\$
250DN1	—	Right Back Cane	NC
250DN2	—	Left Back Cane	NC

JOYSTICK / ELECTRONICS *See Last Pages For Specialty Input Devices images*

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Ctrl+5® Switchbox Options \$

250QG09	—	Ctrl+5 Switch Box w/ Free Assignable Functions and Switch-It™ Remote Seating™	NC
---------	---	---	----

Ctrl+5 Mounting Location

250EY09	—	Right Armrest	NC
250EY10	—	Left Armrest	NC
250EY11	—	Back Cane, Right	\$172.00
250EY12	—	Back Cane, Left	\$172.00

Ctrl+5 Switch Type P. 19 \$

250EY03	—	Push Buttons	NC
250EY17	K0108	Toggle Switches	\$343.00

R-net Programming Selection \$**Biometric Repositioning**

250PM03	K0108	Include Biometric Repositioning (Requires Memory Position #1, #2, #3 Below)	NC
---------	-------	---	----

Memory Seating

Memory Seating feature allows for up to 6 personalized seating positions that are accessed either through the joystick, Ctrl+5, joystick shortcut keys or programmable jacks on Omni2 and joystick. Choose desired presets below. Default positions are dealer adjustable.

250PM06	—	Include #1 Memory Seating (Seated: Required as return position from all other positions)	NC
250PM07	K0108	Include #2 Memory Seating (Supine)	NC
250PM08	K0108	Include #3 Memory Seating (Pressure Relief)	NC
250PM09	K0108	Include #4 Memory Seating (Forward Tilt)	NC
250PM10	K0108	Include #5 Memory Seating (Elevate) (Required)	NC
250PM11	K0108	Include #6 Memory Seating (Stand)	NC

Assignable Buttons

250PM12	—	Include Assignable Buttons	NC
---------	---	----------------------------	----

PRAR Pressure Relief Alert and Recording

250PM14	K0108	Include PRAR (Minimum Requirement Memory Position #1 and #3)	NC
---------	-------	--	----

Power Centermount Footplate Extension Option (Required for UP)

250IE03	K0108	Include 7.5" Powered Footplate Extension for Transfer Assist (Electronics programmed to allow footplates to be extended with or without legrest being raised.)	\$1,430.00
---------	-------	--	------------

SPECIALTY INPUT DEVICES *See last page for Specialty Input Device images.*

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

R-net Omni2 Options **P. 19** *For Expandable Controls, choose all that apply.* \$

2500N14	—	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth, Choose Mount Below	\$2,031.00
2500N15	—	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth, Bundled Omni2 Available w/ Switch-It Input Devices	NC

Omni2 Mount *If used w/ Standard Joystick, it must be mounted opposite side.* \$

250DS12	E1028	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Right	\$567.00
250DS13	E1028	Link-It with Parallel Swing-Away Omni2 Mount w/ Rear Mounted Interface, Left	\$567.00

Hand Control Packages **P. 20** *Requires Standard Joystick Omit / Only one hand control device can be ordered.* \$**Proportional Input Devices**

250EU14	E2312	Switch-It MicroGuide	\$3,769.00
250EU15	E2312	Switch-It MicroPilot	\$3,769.00
250EU17	E2328	Switch-It VersaGuide	\$2,076.00
250EU18	E2328	Switch-It VersaGuide EZ	\$2,453.00

Mounting Position (Required) **20**

250EM04	—	Link-It to the Backrest (Not Available w/ Omit Back Cushion Type)	\$1,087.00
250EM05	—	Link-It to the Fixed Joystick Mount	\$338.00
250EM06	E1028	Link-It to the Parallel Swing-Away Joystick Mount	\$567.00
250EM07	E1028	Link-It to the Midline Swing-Away	\$1,201.00
250EM08	—	Omit Mount	NC

Mounting Platforms (Optional) *Not available w/ Link-It to Backrest*

250MG16	E0950	SunDisk Hand Tray (Only available w/ MicroSeries Joysticks)	\$338.00
250MG15	E0950	Bullet Hand Tray	\$338.00

Mounting Location (Required)

250IM09	—	Right	NC
250IM10	—	Left	NC

Switch Options

—	—	See Switch Selections (Maximum 2)	—
---	---	-----------------------------------	---

Head Control Packages **P. 20** *Only one head control device can be ordered. Must select headrest mount on next page.* \$**Proportional Input Devices**

250EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads	\$6,537.00
250EU26	E2328	Switch-It DualPro Proportional Head Array - Spot Pads	\$6,537.00
250EU51	E2328	Switch-It Vigo Head Drive	\$5,537.00

Switched Input Devices

250EU28	E2330	Switch-It 3-Switch Head Array - Short Curved Wings	\$5,026.00
250EU30	E2330	Switch-It 3-Switch Head Array - Long Curved Wings	\$5,026.00

Switch Options

—	—	See Switch Selections (Maximum 2)	—
---	---	-----------------------------------	---

SPECIALTY INPUT DEVICES *See last page for Specialty Input Device images.*

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Sip 'N Puff \$

250ID73	E2330	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware (Requires headrest mount)	\$5,713.00
250ID74	K0108	Switch-It OptiStop with Sip 'N Puff	\$2,225.00

Headrest Mounts *(Applies to Head/Sip 'N Puff control packages.)* \$

250IJ13	E1028	Whitmyer® Linx Headrest Mount w/ Quick Release	\$188.00
250IJ14	E1028	Whitmyer Linx² Headrest Mount w/ Quick Release	\$229.00
250IJ15	E1028	Whitmyer Axys Headrest Mount w/ Quick Release	\$121.00

Chin Control Packages P. 20 \$**Proportional Input Devices**

250EU34	E2312	Switch-It MicroGuide (Switches mounted to the headrest only.)	\$3,769.00
250EU35	E2312	Switch-It MicroPilot (Switches mounted to the headrest only.)	\$3,769.00
250EU36	E2328	Switch-It VersaGuide	\$2,076.00
250EU37	E2328	Switch-It VersaGuide EZ	\$2,453.00

Mounting Position

250EM09	E1028	Link-It to Chin Harness, Small (150mm / 6")	\$1,430.00
250EM10	E1028	Link-It to Chin Harness, Large (220mm / 8.5")	\$1,430.00
250EM11	E1028	Link-It to Headrest S/A	\$1,430.00
250EM12	E1028	Link-It to Manual Chin Boom - Right Side	\$4,341.00
250EM13	E1028	Link-It to Manual Chin Boom - Left Side	\$4,341.00
250EM16	E1028	Omit Mount	NC

Input Device Mounting Location

250IM13	—	Right	NC
250IM14	—	Left	NC

Switch Options

—	—	See Switch Selections	—
---	---	-----------------------	---

Safety Devices \$

250SA06	—	Switch-It Remote Stop Switch - XLR Connector	\$1,070.00
---------	---	--	------------

Switched Driving Packages \$

250SW44	E2328	Switch-It® Cool Cube® with DB9 connector for proximity/proportional switch driving <i>(Pro Spot® proportional switches are required if proportional driving is desired.) (For non-proportional driving, single switches can be used.)</i>	\$3,969.00
---------	-------	--	------------

SPECIALTY INPUT DEVICES *See last page for Specialty Input Device images.*

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Switches P. 21

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, Cool Cube®, head controls

Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections
Switch / Quantity

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	\$
2.5" Single Switch - Black							\$223.00
2.5" Single Switch - Blue							\$223.00
2.5" Single Switch - Green							\$223.00
2.5" Single Switch - Red							\$223.00
2.5" Single Switch - Yellow							\$223.00
Egg Switch - Black							\$126.00
Egg Switch - Blue							\$126.00
Egg Switch - Green							\$126.00
Egg Switch - Red							\$126.00
Egg Switch - Yellow							\$126.00
MicroLite <i>(Includes color labels)</i>							\$188.00
Mini Button 1 - Blue <i>(Velcro)</i>							\$121.00
Mini Button 1 - Green <i>(Velcro)</i>							\$121.00
Mini Button 1 - Red <i>(Velcro)</i>							\$121.00
1" Single Switch - Black							\$315.00
1" Single Switch - Blue							\$315.00
1" Single Switch - Green							\$315.00
1" Single Switch - Orange							\$315.00
1" Single Switch - Red							\$315.00
1" Single Switch - Yellow							\$315.00

Electronic Switches

Flex Switch (Headrest mount only)							\$245.00
-----------------------------------	--	--	--	--	--	--	----------

ProSpot - Proportional Switches

ProSpot Proportional Switch - Red							\$543.00
ProSpot Proportional Switch - Green							\$543.00
ProSpot Proportional Switch - Blue							\$543.00
ProSpot Proportional Switch - Yellow							\$543.00
ProSpot Proportional Switch - Orange							\$543.00

Section B – Switch Mounting Position

Link-It to the Backrest (Up to 1 per side)							\$744.00
Headrest Mount - Gooseneck & Rigid (Up to 2 per side)							\$744.00
Link-It to the Fixed Joystick Mount (Up to 2 per side)							\$567.00
Link-It to the Parallel Swing-Away Joystick Mount (Up to 2 per side)							\$567.00
Link-It to the Seat Rail (Up to 2 per side) <i>(Not avail w/knee & thigh supports)</i>							\$451.00
Omit Mount (includes Velcro)							NC

Section C – Switch Mounting Location

Right side of chair							NC
Left side of chair							NC

Joystick Integrated Switches (Pair)

250J003	–	Proximity Switches on Flex Mount, Pair (Used with VersaGuide/VersaGuide EZ)	\$1,138.00
---------	---	---	------------

built4me

Q700-UP M with SEDEO PRO ADVANCED

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year

