



Zippie ZM-310™ Group 3

K0848, K0856
EIPW22, January 2019

Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City, State, Zip: _____
Marked For: _____	Phone: _____

The HCPCS CODES provided are based on code verification through PDAC or our interpretation of the code definitions.
The use of HCPCS CODES does not ensure coverage or payment. For coverage information, verify with the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. Zippie ZM-310™ GROUP 3 CATEGORY (Select One)

Zippie ZM-310™ BC (Group 3 K0848/K0011)			
<input type="checkbox"/> Zippie ZM-310™ BC w/ ASAP II Seat			\$7,375.00
195PC106	K0848 K0011	Base Package	\$7,375.00

Zippie ZM-310™ SC (Group 3 K0856/K0011)			
<input type="checkbox"/> Zippie ZM-310™ SC for SPOT Tilt			\$7,375.00
195PC107	K0856 K0011	Base Package	\$7,375.00

BC Seating Specifications			
HCPC Code	K0848/ K0011		
Typical Chair Weight	W/ battery	235 lbs.	W/O Battery 157 lbs.
Transit	WC19 w/optional crash tested belt		
Seat Type	ASAP II / Optional Manual tilt and/or Manual Recline		
Seat size range	Width	10" - 16"	Depth 10" - 19"
STFH range	15.5" - 21.6"		
Power Seating	9" Seat elevator		

SC Seating Specifications			
HCPC Code	K0856 / K0011		
Typical Chair Weight	W/ battery	260 lbs.	W/O Battery 182 lbs.
Transit	WC19 w/optional crash tested belt		
Seat Type	ASAP II / Optional Manual Recline		
Seat size range	Width	10" - 16"	Depth 10" - 19"
STFH range	17.6" - 21.8"		
Power Seating	50° CG tilt / optional 12" seat elevator		



Zippie ZM-310 Base Specifications

Drive Wheel Position	Midwheel Drive		
User weight Cap	165 lbs.		
Base size	Width	24"	Length 34" Turn Radius 18"
Drive Control	VR-2 90A (non-expand) R-net EL90 (Expandable)		
Driving Range	Up to 13.5 miles		
Drive Motors	4-pole 6 mph		

Drive Wheels	13" Drive Wheels w/ flat free insert (Black)		
Tires/Casters	7" Front / 6" rear (Black)		
Charger	8A Dual Mode Charger (Compact) CTE		
Battery Size	Group 22		
Shroud Colors	12 standard colors, 3 upgrade prints		
Ground Clearance	3"		

2. BASE CHAIR OPTIONS

Shroud Color				
<input type="checkbox"/>	195EZ01	-	Candy Apple Red	N/C
<input type="checkbox"/>	195EZ02	-	High Gloss Black	N/C
<input type="checkbox"/>	195EZ03	-	Candy Blue	N/C
<input type="checkbox"/>	195EZ04	-	Silver	N/C
<input type="checkbox"/>	195EZ05	-	Copper	N/C
<input type="checkbox"/>	195EZ08	-	Green Apple	N/C
<input type="checkbox"/>	195EZ09	-	Yellow	N/C
<input type="checkbox"/>	195EZ10	-	Pearl Pink	N/C
<input type="checkbox"/>	195EZ11	-	Candy Purple	N/C
<input type="checkbox"/>	195EZ12	-	Matte Black	N/C
<input type="checkbox"/>	195EZ13	-	Midnight Blue	N/C
<input type="checkbox"/>	195EZ15	-	Pearl White	N/C
<input type="checkbox"/>	195EZ17	-	Pink Camo	\$515.00
<input type="checkbox"/>	195EZ18	-	Grunge Skulls	\$515.00
<input type="checkbox"/>	195EZ16	-	American Flag	\$515.00

Drive wheel/tire package				
<input type="checkbox"/>	195DW60	-	Cast Wheel Solid Dual Density Inserts Black Knobby Tires and Black Casters	N/C

User weight				
<input type="checkbox"/>	195WG17	-	< 75 lbs.	N/C
<input type="checkbox"/>	195WG18	-	76-100 lbs.	N/C
<input type="checkbox"/>	195WG19	-	101-125 lbs.	N/C
<input type="checkbox"/>	195WG20	-	126-165 lbs.	N/C

Transit Option (Required)				
<input type="checkbox"/>	195WN03	-	WC19 Transit Package (order belt separately below)	N/C

WC19 Transit Lap Belt				
<input type="checkbox"/>	195WN02	-	WC19 Transit Lap Belt	\$210.00

Pre-Set Seat Frame Angle				
<input type="checkbox"/>	195FA10	-	0°	N/C
<input type="checkbox"/>	195FA3	-	3° N/A with +3" or +4" STFH or Lift	N/C
<input type="checkbox"/>	195FA4	-	6° N/A with +3" or +4" STFH or Lift/Tilt	N/C

Seat Height				
<i>(Measured in the front) Obtain Quote to verify est. STFH</i>				
<input type="checkbox"/>	195HT44	-	Lowest STFH	NC
<input type="checkbox"/>	195HT45	-	+1" Seat to Floor Height	NC
<input type="checkbox"/>	195HT46	-	+2" Seat to Floor Height	NC
<input type="checkbox"/>	195HT47	-	+3" STFH	NC
<input type="checkbox"/>	195HT48	-	+4" STFH	NC

Motor Package				
<input type="checkbox"/>	195PF96	-	165 lbs. 6mph 4-Pole	Std

Ventilator Options				
<i>Not available w/ Transit</i>				
<input type="checkbox"/>	195VT4	E1029	Fixed Vent Tray (Pulmonetic LTV Models)	\$465.00
<input type="checkbox"/>	195VT9	E1029	Lift off Vent Tray (Pulmonetic LTV Models)	\$465.00

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 tilt				
	0	+1"	+2"	+3"	+4"
BC ASAP II No power seating*	15.5	16.5	17.5	18.5	19.5
BC ASAP II 9" Lift only	18.4	19.4	20.4	21.4	22.4
BC ASAP II Manual Tilt	17.6	18.6	19.6	20.6	21.6
SC ASAP II Tilt only	17.6	18.6	19.6	20.6	21.6
SC ASAP II Tilt and 12" Lift combo	17.8	18.8	19.8	20.8	21.8

*For seat widths 12"-16" add .5"

3. BATTERIES and CHARGERS

Batteries				
<input type="checkbox"/>	195AG21	E2361	GP22 Flat Top Batteries Pair w/ strap	\$525.00
<input type="checkbox"/>	195AG05	-	Omit batteries (Includes harness)	N/C

Battery Charger				
<input type="checkbox"/>	195BJ19	-	8A Dual Mode Charger (Compact) CTE	N/C

4. ELECTRONICS

Drive Control Options

Non-Expandable Controllers

<input type="checkbox"/>	195ER01	Non-expandable VR2-90 controller (N/A with High Speed Motor package)	N/C
--------------------------	---------	---	-----

Expandable Controllers

<input type="checkbox"/>	195ER02	R-Net Expandable Drive Control Upgrade (90A)		\$980.00
	195EX01	E2377	Expandable Controller - R-net	\$540.00
	195EX02	E2313	Harnessing for Expandable Controller	\$440.00

VR2 Joystick (non-expandable)

<input type="checkbox"/>	195J94	-	Non Expandable VR2 Joystick	N/C
--------------------------	--------	---	-----------------------------	-----

***R-net Joystick

<input type="checkbox"/>	195J70	-	Std LED Joystick for R-net Ctrl's	N/C
<input type="checkbox"/>	195J108	-	CJSM1: Joystick Standard: w/ Stereo Jacks & Color Display	N/C
<input type="checkbox"/>	195J114	-	CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys, IR, BT, and Stereo Jack	\$1,055.00
<input type="checkbox"/>	195J95	-	Omit JS. (For use with OMNI)	N/C

*Joystick Mount (Fixed, Swing Away, Retractable, Mid)

<input type="checkbox"/>	195JM10	-	Fixed Mount	N/C
<input type="checkbox"/>	195JM11	E1028	Swing Away Mount	\$245.00
<input type="checkbox"/>	195JM15	E1028	Retractable Joystick mount w/ 45° stops	\$585.00
<input type="checkbox"/>	195JM14	E1028	Mid Line Mount	\$245.00

Joystick Mount (Location)

Height Adjustable.

<input type="checkbox"/>	195HM1	-	Right Side Mount	N/C
<input type="checkbox"/>	195HM2	-	Left Side Mount	N/C

Joystick Handle (Optional)

<input type="checkbox"/>	195JH1	E2323	Ball Handle	\$80.00
<input type="checkbox"/>	195JH2	E2323	Foam Ball Handle	\$80.00
<input type="checkbox"/>	195JH4	E2323	T-Handle	\$80.00
<input type="checkbox"/>	195JH6	E2323	Goalpost	\$80.00

Attendant Control

<input type="checkbox"/>	1950A01	E2331	Attendant Control for VR-2	\$1,025.00
<input type="checkbox"/>	1950A02	E2331	Attendant Control for R-net	\$1,025.00

Attendant Mounting location

<input type="checkbox"/>	195DN1	-	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	195DN2	-	Mounting Location: Left Back Cane	N/C

R-net OMNI and Output Options

(For expandable controls, choose all that apply)

<input type="checkbox"/>	195ON1	-	R-net OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below (required for Specialty Input Devices)	\$1,005.00
<input type="checkbox"/>	195ON2	-	R-net Output Module (ECM)	\$1,025.00

OMNI Mount

(If used with Joystick must be mounted opposite side)

<input type="checkbox"/>	195DS04	E1028	Swing Away Left Side OMNI Mount	\$245.00
<input type="checkbox"/>	195DS03	E1028	Swing Away Right Side OMNI Mount	\$245.00

Switch inputs for R-Net Joysticks or OMNI

Drive or Mode / Profile Control

<input type="checkbox"/>	195SW22	-	Stereo to 2 Mono Splitter	\$65.00
<input type="checkbox"/>	195SW1	-	Buddy Button Switch	\$80.00
<input type="checkbox"/>	195SW4	-	Micro Light Switch	\$80.00
<input type="checkbox"/>	195SW24	-	Egg Switch-Black	\$100.00
<input type="checkbox"/>	195SW21	-	Dual toggle switch and bracket	\$375.00
<input type="checkbox"/>	195SW35	-	Switch-It Mini 1 button on Mono plug	\$200.00
<input type="checkbox"/>	195SW36	-	Switch-It Mini 2 button on 2 Mono plugs	\$200.00
<input type="checkbox"/>	195SW37	-	Switch-It Mini 2 button on 1 Stereo plug	\$200.00
<input type="checkbox"/>	195SW38	-	Switch-It Mini 4 button on DB9 connector	\$200.00

Single Switch Mount Bracket

<input type="checkbox"/>	195EM01	-	Single Switch Mount Bracket	\$200.00
--------------------------	---------	---	-----------------------------	----------

Electronic Programmers

<input type="checkbox"/>	195PR2	-	VR-2 Hand Held Programmer	\$260.00
<input type="checkbox"/>	195PR5	-	R-net PC Programmer with Dongle	\$415.00

5. Specialty Driver Control Options

Requires expandable R-net controls and OMNI unless otherwise noted. N/A with Upholstered backs.

5a. Proportional Control Packages

Head Control Proportional Packages			
□ Switch-It Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)			\$5,535.00
195ID72	E2328	Switch-It Dual Pro Head Control w/ facial pad (Adult)	\$5,195.00
195IJ09	E1028	S/A Boom Mount with Egg Switch Control	\$340.00
□ Switch-It Dual Pro Proportional Head Control with Spot Pads Package			\$5,535.00
195ID70	E2328	Switch-It Dual Pro Head Control w/ Spot pads	\$5,195.00
195IJ09	E1028	S/A Boom Mount with Egg Switch Control	\$340.00

Hand Control Proportional Joystick Packages			
□ Switch-It Versa-Guide Compact Joystick Midline with Gatlin Package			\$1,895.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IJ07	E1028	S/A Gatlin Mount	\$245.00
□ Switch-It Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package			\$2,420.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$525.00
195IJ07	E1028	S/A Gatlin Mount	\$245.00
□ Switch-It Micro-Pilot Joystick with SpaceDisc Package			\$3,610.00
195ID62	E2312	Micro-Pilot Proportional Joystick	\$2,995.00
195IJ08	K0108	S/A SpaceDisc mount	\$615.00
□ Switch-It Micro-Guide Joystick with SpaceDisc Package			\$3,610.00
195ID63	E2312	Micro-Guide Proportional Joystick	\$2,995.00
195IJ08	K0108	S/A SpaceDisc mount	\$615.00

Chin Control Proportional Joystick Packages			
□ Switch-It Versa-Guide Chin Control Compact Joystick and Egg Switch Package			\$2,246.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM07	E1028	S/A Mount off backrest	\$256.00
195IJ04	E1028	S/A Boom Mount with egg switch	\$340.00
□ Switch-It Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package			\$2,431.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$525.00
195IM07	E1028	S/A Mount off backrest	\$256.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"			\$2,790.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$525.00
195IM05	K0108	Bib mount Small	\$615.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"			\$2,790.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$525.00
195IM06	K0108	Bib mount Large	\$615.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"			\$2,605.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM05	K0108	Bib mount Small	\$615.00
195IJ04	E1028	SA Boom Mount with egg switch	\$340.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"			\$2,605.00
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,650.00
195IM06	K0108	Bib mount Large	\$615.00
195IJ04	E1028	SA Boom Mount with egg switch	\$340.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"			\$3,950.00
195ID62	E2312	Micro-Pilot Proportional Joystick	\$2,995.00
195IM05	K0108	Bib mount Small	\$615.00
195IJ04	E1028	S/A Boom Mount with egg switch	\$340.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"			\$3,950.00
195ID62	E2312	Micro-Pilot Proportional Joystick	\$2,995.00
195IM06	K0108	Bib mount Large	\$615.00
195IJ04	E1028	S/A Boom Mount with egg switch	\$340.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6"			\$3,950.00
195ID63	E2312	Micro-Guide Proportional Joystick	\$2,995.00
195IM05	K0108	Bib mount Small	\$615.00
195IJ04	E1028	S/A Boom Mount with egg switch	\$340.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8"			\$3,950.00
195ID63	E2312	Micro-Guide Proportional Joystick	\$2,995.00
195IM06	K0108	Bib mount Large	\$615.00
195IJ04	E1028	S/A Boom Mount with egg switch	\$340.00

5. Specialty Driver Control Options (cont.)

5b. Switch Control

Requires Headrest and/or Headrest Mount

3 Switch Proximity Head Arrays *Choose headrest mount below*

<input type="checkbox"/> Switch-It Head Array Package 3 Switch Curved (Long) \$4,535.00			
195II13	E2330	Switch-It Proximity Head Array Long	\$4,195.00
195IJ09	E1028	S/A Boom Mount with Egg Switch Control	\$340.00
<input type="checkbox"/> Switch-It Head Array Package 3 switch curved (Stubby) \$4,535.00			
195II14	E2330	Switch-It Proximity Head Array Stubby	\$4,195.00
195IJ09	E1028	S/A Boom Mount with Egg Switch Control	\$340.00

5 Switch Proximity Head Array (No Egg Switch Needed) *Choose "Headrest Mount for Special Controls"*

<input type="checkbox"/>	195II04	E2330	Switch-It Head Array 5 switch flat	\$4,295.00
--------------------------	---------	-------	------------------------------------	------------

Multiple Switches

<input type="checkbox"/>	195SW31	-	Switch-It 3 Naked proximity switches on DB9 connector (requires OMNI)	\$2,675.00
<input type="checkbox"/>	195SW32	-	Switch-It 4 Naked proximity switches on DB9 connector (requires OMNI)	\$3,290.00
<input type="checkbox"/>	195SW33	-	Switch-It 5 Naked proximity switches on DB9 connector (requires OMNI)	\$3,705.00
<input type="checkbox"/>	195SW23	-	4-way toggle switch button DB9 Connector. Includes bracket	\$550.00

Sip-N-Puff Package *Choose headrest and HR mount below*

<input type="checkbox"/> Quickie Sip n' Puff System \$2,180.00			
195ID34	E2325	Sip-N-Puff Kit	\$1,465.00
195IH1	E2326	Sip-N-Puff Tubing and SA Hardware	\$375.00
195IJ10	E1028	S/A Boom Mount with egg switch	\$340.00

Headrest Mount for Specialty Controls

<input type="checkbox"/>	195MH12	E1028	Axys Mount & Bracket for Specialty Controls	\$230.00
<input type="checkbox"/>	195MH13	E1028	Cobra Mount & Bracket for Specialty Controls	\$290.00
<input type="checkbox"/>	195MH14	E1028	Linx Mount & Bracket for Specialty Controls	\$230.00

Headrest Pad for Specialty controls

<input type="checkbox"/>	195HE79	E0955	Plush Pad 8" for Specialty controls	\$230.00
<input type="checkbox"/>	195HE80	E0955	Plush Pad 10" for Specialty controls	\$230.00
<input type="checkbox"/>	195HE81	E0955	Plush Pad 14" for Specialty controls	\$230.00

Extra Swing Away Head Boom Egg Switch and Mount

<input type="checkbox"/>	195IJ05	E1028	Extra SA Boom Mount with egg switch #1	\$340.00
<input type="checkbox"/>	195IJ06	E1028	Extra SA Boom Mount with egg switch #2	\$340.00

Single Switches

<input type="checkbox"/>	195IB01	-	Ribbon Switch	\$100.00
<input type="checkbox"/>	195IB02	-	Buddy Button Switch	\$80.00
<input type="checkbox"/>	195IB03	-	Disc Switch	\$165.00
<input type="checkbox"/>	195IB04	-	Wobble Switch	\$515.00
<input type="checkbox"/>	195IB05	-	Micro Light Switch	\$80.00
<input type="checkbox"/>	195IB06	-	Egg Switch-Black	\$100.00
<input type="checkbox"/>	195IB07	-	Switch-It Mini 1 button on Mono plug	\$200.00

6. SEAT TYPE

<input type="checkbox"/>	195ST56	-	A.S.A.P II (Adjustable Seating and Positioning) Choose Seat Frame Size Platform	N/C
--------------------------	---------	---	---	-----

ASAP II Frame size

(each frame is width and depth adjustable within given range)

<input type="checkbox"/>	195FU16	-	Adjustable 10x10 to 13x16 (Z)	N/C
<input type="checkbox"/>	195FU17	-	Adjustable 12x12 to 16X19 (A)	N/C

7. SEAT SIZE

Seat Size Matrix. *Seat type is chosen in section 6*

Z = ASAP II Seat Frame Adjustable from 10x10 to 13x16

A = ASAP II Seat Frame Adjustable from 12x12 to 16X19

		ASAPII Seat Width							
		10	11	12	13	14	15	16	
Seat Depth	10	Z	Z	Z	Z				
	11	Z	Z	Z	Z				
	12	Z	Z	Z,A	Z,A	A	A	A	
	13	Z	Z	Z,A	Z,A	A	A	A	
	14	Z	Z	Z,A	Z,A	A	A	A	
	15	Z	Z	Z,A	Z,A	A	A	A	
	16	Z	Z	Z,A	Z,A	A	A	A	
	17				A	A	A	A	A
	18				A	A	A	A	A
	19				A	A	A	A	A

Seat Width				
<input type="checkbox"/>	195R45	-	10" Seat width	N/C
<input type="checkbox"/>	195R46	-	11" Seat width	N/C
<input type="checkbox"/>	195R21	-	12" Seat width	N/C
<input type="checkbox"/>	195R22	-	13" Seat width	N/C
<input type="checkbox"/>	195R23	-	14" Seat width	N/C
<input type="checkbox"/>	195R24	-	15" Seat width	N/C
<input type="checkbox"/>	195R5	-	16" Seat width	N/C

Seat Depth				
<input type="checkbox"/>	195RF40	-	10" Seat depth	N/C
<input type="checkbox"/>	195RF41	-	11" Seat depth	N/C
<input type="checkbox"/>	195RF20	-	12" Seat depth	N/C
<input type="checkbox"/>	195RF21	-	13" Seat depth	N/C
<input type="checkbox"/>	195RF22	-	14" Seat depth	N/C
<input type="checkbox"/>	195RF23	-	15" Seat depth	N/C
<input type="checkbox"/>	195RF5	-	16" Seat depth	N/C
<input type="checkbox"/>	195RF6	-	17" Seat depth	N/C
<input type="checkbox"/>	195RF7	-	18" Seat depth	N/C
<input type="checkbox"/>	195RF8	-	19" Seat depth	N/C

8. SEAT CUSHIONS

JAY / Zip Cushions (Optional)

To order any JAY / Zip cushions on-chair, please see the JAY/Zip price list.

Cushion name or part number must be filled out w/ desired width & depth

JAY / Zip Bundled Cushion			
<input type="checkbox"/>	195JC20	Cushion Name or Part #	Price List
		Width: _____ Depth: _____	

9. SEATING SYSTEMS

Non-Powered Seating

<input type="checkbox"/>	195PJ8	K0108	Manual Tilt 45° CG Tilt <i>Must be ordered with Stroller handle N/A with height adjustable back cane</i>	\$850.00
--------------------------	--------	-------	---	----------

Power Seating Options (SC Models)

<input type="checkbox"/>	195PJ2	E1002	SPOT 50° CG Tilt	\$4,950.00
--------------------------	--------	-------	------------------	------------

Powered Seat Elevator - (AKA: Power Lift)

<input type="checkbox"/>	195PA8	E2300	Power Seat Elevate (9")	\$2,500.00
<input type="checkbox"/>	195PA15	E2300	Power Seat Elevate for use with tilt (12") <i>(Must be ordered with SPOT tilt)</i>	\$2,800.00

Power Seating Control

Powered Seating Switch Control Method

For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics

Single Actuator Switch Control (Lift or Tilt)

<input type="checkbox"/>	195EO10	-	Switch Control 1 actuator with <u>Dual Toggle</u>	N/C
<input type="checkbox"/>	195EO11	-	Switch Control 1 actuator with <u>Buddy Button</u>	N/C
<input type="checkbox"/>	195EO19	-	Switch Control 1 actuator with Switch-It Mini 1 button	N/C
<input type="checkbox"/>	195EO20	-	Switch Control 1 actuator with Switch-It Mini 2 button	N/C

2-3 Actuator Switch Control (Lift and Tilt)

<input type="checkbox"/>	195EO21	-	Switch Control 2-3 actuator with <u>Dual Toggle</u>	\$500.00
<input type="checkbox"/>	195EO12	-	Switch Control 2-3 actuator with <u>4-way toggle</u>	\$750.00
<input type="checkbox"/>	195EO23	-	Switch Control 2-3 actuator with Switch-It Mini 2 button	\$500.00
<input type="checkbox"/>	195EO22	-	Switch Control 2-3 actuator with Switch-It Mini 4 button	\$500.00

Switch Mount

(for use with switch control powered seating)

<input type="checkbox"/>	195EY01	-	Switch Mount Right	N/C
<input type="checkbox"/>	195EY02	-	Switch Mount Left	N/C

Powered Seating Thru-drive Control Method

Thru-Drive 1 or 2 Actuator

<input type="checkbox"/>	195EO26	E2310	Thru-Drive 1 or 2 Actuator(s) For operating Tilt and/or Lift thru Joystick. Includes Ctrl+5 when ordered with Expandable Controls (R- net)	\$1,420.00
--------------------------	---------	-------	---	------------

Powered Seating Thru-drive Control Options

Only available with Expandable Controls (R-net)

Assignable Buttons

"Assignable Button" feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle

<input type="checkbox"/>	195PM01	-	Enable Assignable buttons	N/C
--------------------------	---------	---	---------------------------	-----

Ctrl+5 Switchbox Options

Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices



<input type="checkbox"/>	195QG01	-	Ctrl+5 Switchbox mounted behind joystick	N/C
<input type="checkbox"/>	195QG02	-	Ctrl+5 Switchbox mounted opposite joystick	N/C
<input type="checkbox"/>	195QG03	-	Ctrl+5 Switchbox mounted on back cane	N/C
<input type="checkbox"/>	195QG05	-	Omit Ctrl+5	(-\$500.00)

Ctrl+5 Switch Type

<input type="checkbox"/>	195EY03	-	Push buttons	N/C
<input type="checkbox"/>	195EY04	-	Toggle switches	\$250.00

10. BACKREST AND BACK OPTIONS FOR ASAP II SEAT

Backrest Frame Type

<input type="checkbox"/>	195BT8	-	Angle Adjustable	N/C
<input type="checkbox"/>	195BT98	E1225	Manual Recline 40° (not width adjust.)	\$440.00

Manual Recline Adjustment Method Option

Back Cane Height 23". Fits all seat types. N/A with Cantilever Arms

<input type="checkbox"/>	195MN01	-	Tool	N/C
<input type="checkbox"/>	195MN02	-	Knob (N/A with dual post arms)	N/C

Angle Adjustable Backrest Options

(N/A with Manual recline)

<input type="checkbox"/>	195OD8	-	Non-Folding Backrest	N/C
<input type="checkbox"/>	195OD9	-	Folding Backrest	\$115.00

Cane Angle

<input type="checkbox"/>	195BP40	-	Straight Back Post	N/C
<input type="checkbox"/>	195BP41	-	8° Bend Back Post (Avail. only on Zip or J3 Back) N/A with Manual recline	\$105.00

Cane Adjustable height

<input type="checkbox"/>	195B1	-	Fixed Low (15.5")	N/C
<input type="checkbox"/>	195B5	-	Fixed Med (17.5")	N/C
<input type="checkbox"/>	195B9	-	Fixed Tall (19.5")	N/C
<input type="checkbox"/>	195B243	-	Adjustable Low 14" - 18"	\$105.00
<input type="checkbox"/>	195B244	-	Adjustable High 20" - 26"	\$105.00

Towel Bar Depth

<input type="checkbox"/>	195QT03	-	Shallow	N/C
<input type="checkbox"/>	195QT04	-	Deep (Ext. 4" from back canes for deep backs)	N/C

Back Type (Solid)

<input type="checkbox"/>	195BU86	-	S.P.O.T. Curved back with rail for mounting laterals	N/C
<input type="checkbox"/>	195BU87	E2613	JAY SPO Posterior (Curved back with T-Foam inserts and mounting rails)	\$490.00
<input type="checkbox"/>	195BU14	-	Omit Backrest	(-\$15.00)

JAY J3 / ZIP Option

<input type="checkbox"/>	195JB01	See J3 OF	JAY J3 / Zip (select OMIT above)	See J3 OF
--------------------------	---------	-----------	----------------------------------	-----------

Back Height - Required for SPOT Curved and SPO Backs

(Limited combinations of widths) Measured from pan to top of solid back

<input type="checkbox"/>	195B118	-	12"	Available on 10" - 16" wide seats	N/C
<input type="checkbox"/>	195B132	-	14"		N/C
<input type="checkbox"/>	195B36	-	16"	Available on 12" - 16" wide seats	N/C
<input type="checkbox"/>	195B38	-	18"		N/C
<input type="checkbox"/>	195B40	-	20"		N/C


Push Handle Options

<input type="checkbox"/>	195PH14	-	Stroller Handle: Single height adjustable Bar Push Handle with 2 utility hooks. Required with Manual Tilt (N/A with adjustable height cane)	\$55.00
--------------------------	---------	---	---	---------

11. HEADREST OPTIONS

Choose One Each

Whitmyer Headrest Mount				
<input type="checkbox"/>	195MH10	E1028	Axys Mount & Bracket	\$230.00
<input type="checkbox"/>	195MH11	E1028	Cobra Mount and Bracket	\$290.00
<input type="checkbox"/>	195MH1	E1028	Linx Mount and Bracket	\$230.00



Whitmyer Headrest Pad				
<input type="checkbox"/>	195HE65	E0955	Plush Head Support Pad 8"	\$230.00
<input type="checkbox"/>	195HE66	E0955	Plush Head Support Pad 10"	\$230.00
<input type="checkbox"/>	195HE82	E0955	Adjust-a-plush pad Narrow	\$230.00
<input type="checkbox"/>	195HE83	E0955	Adjust-a-plush pad Standard	\$230.00
<input type="checkbox"/>	195HE84	E0955	Contour Cradle Small	\$230.00
<input type="checkbox"/>	195HE85	E0955	Contour Cradle Medium	\$230.00
<input type="checkbox"/>	195HE11	-	Omit headrest pad only	(-\$15.00)

12. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)				
<i>Available on all Solid Backs</i>				
<input type="checkbox"/>	LAT 3x3 S/A FLUSH MOUNTS - PR		\$805.00	
	195ZV20	E0956	LATERAL 3x3 Flush Mount - PAIR	\$260.00
	195ZS07	E0884	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/>	LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR		\$805.00	
	195ZV14	E0956	LATERAL 4X4 W/ 1" OFFSET - PAIR	\$260.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/>	LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR		\$805.00	
	195ZV15	E0956	LATERAL 4X5 W/ 1" OFFSET - PAIR	\$260.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR	\$545.00
<input type="checkbox"/>	LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR		\$805.00	
	195ZV16	E0956	LATERAL 4X6 W/ 1" OFFSET - PAIR	\$260.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR	\$545.00

Thigh Support w/ Removable Mount				
<i>Available with Cantilever Arms, N/A with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms</i>				
<input type="checkbox"/>	JAY Rectangular Thigh Support - Left		\$335.00	
	195TV3	E0956	5" x 6" Thigh Support - Left	\$105.00
	195TX1	E1028	Removable Mount - Left	\$230.00
<input type="checkbox"/>	JAY Rectangular Thigh Support - Right		\$335.00	
	195TV4	E0956	5" x 6" Thigh Support - Right	\$105.00
	195TX2	E1028	Removable Mount - Right	\$230.00

JAY Anterior Trunk Support Style <i>(Choose Size below)</i>				
<input type="checkbox"/>	195JJ12	E0960	Classic Dynamic	\$230.00
<input type="checkbox"/>	195JJ13	E0960	Classic Structured	\$230.00
<input type="checkbox"/>	195JJ14	E0960	Contoured Dynamic	\$230.00
<input type="checkbox"/>	195JJ15	E0960	Contoured Structured	\$230.00
<input type="checkbox"/>	195ZR03	E0960	Padded Chest Strap	\$230.00

JAY Anterior Trunk Support Size				
<input type="checkbox"/>	195JK01	-	Small	NC
<input type="checkbox"/>	195JK02	-	Medium	NC

JAY Pelvic Positioning Belts				
<input type="checkbox"/>	195ZZ19	-	2 Point Pelvic Belt 1' Push Button Padded	\$95.00
<input type="checkbox"/>	195ZZ20	-	2 Point Pelvic Belt 1' Side Release Padded	\$75.00
<input type="checkbox"/>	195ZZ21	-	2 Point Pelvic Belt 1.5" Push Button Padded	\$95.00
<input type="checkbox"/>	195ZZ22	-	2 Point Pelvic Belt 1.5" Side Release Padded	\$75.00
<input type="checkbox"/>	195ZZ23	-	4 Pnt V Contour Pelvic Belt 1.5" Push Btn	\$175.00

Swing Away Knee Adductor Supports				
<i>Available only with Fixed and Power Center Mount Legrests</i>				
<i>Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms</i>				
Knee Adductor Supports - LEFT (choose Pad and Arm)			\$405.00	
Pad Size				
<input type="checkbox"/>	195KB5	E0956	Knee Adductor Pad SM 3"- Left	\$175.00
S/A Arm				
<input type="checkbox"/>	195XX04	E1028	Swing Away Short - Left	\$230.00
Knee Adductor Supports - RIGHT (choose Pad and Arm)			\$405.00	
Pad Size				
<input type="checkbox"/>	195KB6	E0956	Knee Adductor Pad SM 3"- Right	\$175.00
S/A Arm				
<input type="checkbox"/>	195XX05	E1028	Swing Away Short - Right	\$230.00

Medial Thigh Support				
<input type="checkbox"/>	Small JAY Style Wedge Medial Thigh Support (Height adjustable)		\$425.00	
	195KH12	E0957	Small Wedge Pad (3" x 4")	\$175.00
	195LM3	E1028	Push Button Flip Down Bracket	\$250.00
<input type="checkbox"/>	Medium JAY Style Wedge Medial Thigh Support (Height adjustable)		\$425.00	
	195KH13	E0957	Medium Wedge Pad (4.5" x 5")	\$175.00
	195LM3	E1028	Push Button Flip Down Bracket	\$250.00
<input type="checkbox"/>	Small JAY Style OVAL Medial Thigh Support (Height adjustable)		\$425.00	
	195KH14	E0957	Small Wedge Pad (3" x 4")	\$175.00
	195LM3	E1028	Push Button Flip Down Bracket	\$250.00
<input type="checkbox"/>	Medium JAY Style Oval Medial Thigh Support (Height adjustable)		\$425.00	
	195KH15	E0957	Medium Wedge Pad (4.5" x 5")	\$175.00
	195LM3	E1028	Push Button Flip Down Bracket	\$250.00

13. ARMREST AND ARMPAD OPTIONS

Armrest Mounting for K0848 Models with ASAP II			
Dual Post Armrests			
<input type="checkbox"/>	195EG01	-	Short Dual Post Flip Back 10" Height N/C
<input type="checkbox"/>	195EG21	-	Long Dual Post Flip Back 10" Height N/C
<input type="checkbox"/>	195EG02	E0973	Shrt Dual Post Flip Back HT Adj 10"-15" \$230.00
<input type="checkbox"/>	195EG22	E0973	Lng Dual Post Flip Back HT Adj 10"-15" \$230.00
Cantilever Armrests			
<i>Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) N/A with Manual Recline Back</i>			
<input type="checkbox"/>	195EG03	E0973	Standard Cantilever Arm w/ 2 Point Locking Post, (9"-10" Ht Range) N/A w/ seat depths below 16" \$240.00
<input type="checkbox"/>	195EG04	E0973	Standard Cantilever Arm w/ 2 Point Locking Post, (11"-12" Ht Range) N/A w/ seat depths below 16" \$240.00
<input type="checkbox"/>	195EG05	E0973	STD Cantilever Arm (9"-12" Ht Range) \$240.00
<input type="checkbox"/>	195EG06	E0973	Low Cantilever Arm (7"-10" Ht Range) Choose ASAPII pad type and size below \$240.00
Single Post Armrests			
<input type="checkbox"/>	195EG07	E0973	Single Post Height Adjustable - Std \$240.00
<input type="checkbox"/>	195EG08	E0973	Single Post Height Adjustable - Low \$240.00

Armrest Mounting for K0856 SC Models with Tilt			
Cantilever Armrests			
<i>Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) N/A with Manual Recline Back</i>			
<input type="checkbox"/>	195EG12	-	STD Cantilever Arm (9"-12" Ht Range) N/C
<input type="checkbox"/>	195EG13	-	Low Cantilever Arm (7"-10" Ht Range) Choose ASAPII pad type and size below N/C
<input type="checkbox"/>	195EG10	-	Standard Cantilever Arm w/ 2 Point Locking Post, (9"-10" Ht Range) N/A w/ seat depths below 16" \$240.00
<input type="checkbox"/>	195EG11	-	Standard Cantilever Arm w/ 2 Point Locking Post, (11"-12" Ht Range) N/A w/ seat depths below 16" \$240.00
Single Post Armrests			
<input type="checkbox"/>	195EG14	-	Single Post Height Adjustable - Std \$240.00
<input type="checkbox"/>	195EG15	-	Single Post Height Adjustable - Low \$240.00
Dual Post Armrests			
<input type="checkbox"/>	195EG09	-	Shrt Dual Post Flip Back HT Adj 10"-15" \$240.00
<input type="checkbox"/>	195EG23	-	Lng Dual Post Flip Back HT Adj 10"-15" \$240.00

Arm Pad Options for ASAPII			
Pad Type / Size Options for Dual Post / Single Post / Low Cantilever			
Left Arm Pad (N/C)		Right Arm Pad (N/C)	
<input type="checkbox"/>	195AE01	-	Desk Length Waterfall (10") <input type="checkbox"/> 195AB01
<input type="checkbox"/>	195AE02	-	Full Waterfall (14") N/A w/ Low Cant <input type="checkbox"/> 195AB02
<input type="checkbox"/>	195AE03	-	Desk Length Standard (10") <input type="checkbox"/> 195AB03
<input type="checkbox"/>	195AE04	-	Full Standard (14") N/A w/ Low Cant <input type="checkbox"/> 195AB04
Pad Setting Options for Standard Cantilever Arm			
Left Arm Pad (N/C)		Right Arm Pad (N/C)	
<input type="checkbox"/>	195AE11	-	Full Length Cantilever (14") <input type="checkbox"/> 195AB11
<input type="checkbox"/>	195AE10	-	Desk Length Cantilever (10") <input type="checkbox"/> 195AB10

14. HANGERS AND FOOTRESTS FOR ASAP II

Hanger and Extension Tubes Options (see footplate compatibility in the brackets)

Swing-away Hangers

Non-Powered Swing Away Legrests

Sold as Pairs

<input type="checkbox"/>	195H2	-	70° Swing Away (B,C)	N/C
<input type="checkbox"/>	195H106	-	80° Swing In-Out (Front Mount) (B,C)	N/C
<input type="checkbox"/>	195H107	-	90° Swing In-Out (Front Mount) (B,C)	N/C
<input type="checkbox"/>	195H138	E0990	Manual ELR S/A (B,C) (Avail. with Fixed Centermount, choose below)	\$280.00

Swing Away Footplates

<input type="checkbox"/>	B - Adult Angle Adjust w/Heel Loop for use with 65° and 70° hanger (Pair)			\$220.00
	195F282	K0040	Adult Angle Adjustable Footplate	\$170.00
	195HL10	E0951	Heel Loop	\$50.00
<input type="checkbox"/>	C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)			\$205.00
	195F283	K0040	Kids Angle Adjustable Footplate	\$170.00
	195LS13	K0038	Leg Strap	\$35.00
<input type="checkbox"/>	C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)			\$205.00
	195F286	K0040	Kids Angle Adjustable Footplate	\$170.00
	195LS13	K0038	Leg Strap	\$35.00
<input type="checkbox"/>	195F279	Omit Footplate and extension tube		N/C

Omit Hangers

<input type="checkbox"/>	195H144	-	Omit Hangers	N/C
--------------------------	---------	---	---------------------	-----

Hanger Extensions (See range listed in matrix below)

NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance

<input type="checkbox"/>	195E16	-	High Mount Hanger	\$55.00
<input type="checkbox"/>	195E1	-	Short	N/C
<input type="checkbox"/>	195E2	-	Medium	N/C
<input type="checkbox"/>	195E5	-	Long	N/C
<input type="checkbox"/>	195E23	-	2" Extension for swing in/out hanger	\$55.00
<input type="checkbox"/>	195E24	-	4" Extension for swing in/out hanger	\$55.00
<input type="checkbox"/>	195E21	-	Omit Extension	N/C

Footrest Accessories

<input type="checkbox"/>	1955A01	E0952	Toe Loop (Pair) <i>Only avail on (opt B)</i>	\$50.00
<input type="checkbox"/>	195LS14	K0038	Extra Leg Strap	\$35.00
<input type="checkbox"/>	1952A01	E0951 E0952	Shoe holders Small Pair w/ padded strap	\$150.00
<input type="checkbox"/>	1952A02	E0951 E0952	Shoe holders Medium Pair w/ padded strap	\$150.00

Centermount Legrests

Center Mount Legrests

<input type="checkbox"/>	195H134	-	Fixed Center Mount Short w/ footplate	N/C
<input type="checkbox"/>	195H135	-	Fixed Center Mount Long w/ footplate	N/C
<input type="checkbox"/>	195H103	-	Angle adjustable Center Mount Short (3" - 5.25") Choose Footplate (I or J) <i>infinite and Independent leg length and angle adjustable</i>	\$315.00
<input type="checkbox"/>	195H104	-	Angle adjustable Center Mount Medium (3" - 11.25") Choose Footplate (I or J) <i>infinite and Independent leg length and angle adjustable</i>	\$315.00
<input type="checkbox"/>	195H108	-	Angle adjustable Center Mount Long (3" - 13.25") Choose Footplate (I or J) <i>infinite and Independent leg length and angle adjustable</i>	\$315.00

Centermount Footplates

<input type="checkbox"/>	I - Angle adj Centermount Footplate pair Small			\$170.00
	195F125	K0040	<i>Small Angle Adjustable Footplates 5.75" x 6"</i>	\$170.00
<input type="checkbox"/>	J - Angle adj Centermount Footplate pair Large			\$170.00
	195F126	K0040	<i>Large Angle Adjustable Footplates 5.75" x 10"</i>	\$170.00

Lower Leg Length Matrix (measured seat pan to the top of the footplate surface)

Hanger Type	High Mount	Short	Medium	Long
70° Swing Away	7" - 14"	13.5" - 16"	16" - 18.5"	18" - 20.5"
80° Swing in/out	7" - 14"	13.5" - 16"	14.5" - 19"	15" - 18"
90° Swing in/out	5.5" - 9" Footplate mounts to front of hanger.			
Manual ELR S/A		12.5" - 17"	14.5" - 19"	17" - 21"
Angle adj Cntrmt		3" - 7"	3" - 13"	
Fixed Cntrmt		12.5" - 15.5"		14.5" - 19.5"

Footplate Matrix

Hanger Type	B	C	I
70° Swing Away	X	X	
80° Swing in/out	X	X	
90° Swing in/out	X	X	
Manual ELR S/A	X	X	
Ang Adj Cntrmt			X
Fixed Cntrmt			

15. ACCESSORIES

Wheelchair Tray Table

<input type="checkbox"/>	1957A01	E0950	TRAY TABLE SM TUFFAK 12" - 14" PWR	\$230.00
<input type="checkbox"/>	1957A02	E0950	TRAY TABLE MED TUFFAK 15" - 17" PWR	\$230.00

Backpack

<input type="checkbox"/>	195A1	-	Black Utility	\$55.00
<input type="checkbox"/>	195A2	-	Black/Charcoal Kids	\$50.00



Accessories

<input type="checkbox"/>	1958A09	-	Ram Phone Mount on Round JS tube	\$50.00
<input type="checkbox"/>	1958A10	-	Ram Tablet Mount on Round JS tube	\$200.00
<input type="checkbox"/>	1958A11	-	Ram Cup holder	\$50.00
<input type="checkbox"/>	1958A12	-	Switch-It Remote Stop Switch	\$895.00

Mounting Accessories (available with All Seats)

Only one accessory can be mounted at a time

<input type="checkbox"/>	195A603	E2207	Cane/Crutch Holder	\$260.00
<input type="checkbox"/>	195A604	E2208	O2 holder	\$260.00
<input type="checkbox"/>	195A605	K0108	Walker holder	\$260.00



Sunrise Medical (US) LLC
2842 Business Park Ave.
Fresno, CA 93727 · USA