



QUICKIE P222 SE BD/SD
Solid Seat (No Power or Single Power Seating Options)
EIPW2 - Effective January 2018

Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City,State,Zip: _____
Marked For: _____	Phone: _____

*The HCPCS CODES provided are based on code verification through PDAC or our interpretation of the code definitions.
 The use of HCPCS CODES does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. QUICKIE P222 SE GROUP 4 CATEGORY

Quickie P222 SE			
<input type="checkbox"/>	116PC19	K0868 K0014	Quickie P222 SE BD (Solid Seat - No Power) \$10,076.00
<input type="checkbox"/>	116PC1	K0877 K0014	Quickie P222 SE SD (solid seat/single power seat option) \$10,076.00
<input type="checkbox"/>	116PC3	K0877 K0014	Quickie P222 SE (Base Only) \$9,383.00

Performance Package			
<input type="checkbox"/>	116PF8		8.5 (8.5 mph / 350 LB Capacity) (Max. Wt. cap. For CG Tilt 250lbs) N/C

Power CG Tilt (53° Tilt)			
<i>See electronic section for control options</i>			
<input type="checkbox"/>	116ST2	E1002	P222 SE Power CG Tilt Max. Wt. capacity for CG Tilt 250lbs. Not available with base only purchase. \$5,045.00

Tilt Weight Capacity			
<input type="checkbox"/>	116LB3		Up to 175 lb Weight Capacity N/C
<input type="checkbox"/>	116LB4		176 to 250 lb Weight Capacity N/C

LED Light Package (Rnet only)			
<input type="checkbox"/>	116LI04		LED Light Package (NA Base only) \$1,030.00

2. FRAME COLOR OPTIONS

Standard Colors			
<input type="checkbox"/>	116C140		Aztec Gold NC
<input type="checkbox"/>	116C2		Black (Not avail as base w/Zebra) NC
<input type="checkbox"/>	116C141		Black Cherry NC
<input type="checkbox"/>	116C10		Black Opal NC
<input type="checkbox"/>	116C38		Blue Opal NC
<input type="checkbox"/>	116C36		Candy Blue NC
<input type="checkbox"/>	116C14		Candy Purple NC
<input type="checkbox"/>	116C35		Candy Red NC
<input type="checkbox"/>	116C142		Electric Blue NC
<input type="checkbox"/>	116C143		Evergreen NC
<input type="checkbox"/>	116C126		Glow NC
<input type="checkbox"/>	116C144		Green Apple NC
<input type="checkbox"/>	116C147		Hot Sparkle Pink NC
<input type="checkbox"/>	116C148		Mauve Pink NC

<input type="checkbox"/>	116C103		Rootbeer NC
<input type="checkbox"/>	116C114		Sparkle Silver NC
<input type="checkbox"/>	116C145		Sunrise Orange NC
<input type="checkbox"/>	116C113		Titanium Color Paint NC
<input type="checkbox"/>	116C8		Yellow NC

Standard Matte Colors			
<input type="checkbox"/>	116C116		Matte Black NC
<input type="checkbox"/>	116C151		Matte Black Cherry NC
<input type="checkbox"/>	116C152		Matte Purple NC
<input type="checkbox"/>	116C150		Matte Electric Blue NC
<input type="checkbox"/>	116C149		Matte Evergreen NC

Custom Colors			
<input type="checkbox"/>	116C43		Mossy Oak Camouflage \$305.00
<input type="checkbox"/>	116C42		Stars and Stripes \$305.00
<input type="checkbox"/>	116C112		Desert Camo \$305.00
<input type="checkbox"/>	116C108		Zebra (Choose base color above: White Std, n/a black) \$305.00

3. BATTERIES and CHARGERS

Batteries (Gel Cells)			
<input type="checkbox"/>	116AG03	E2363	GP24 MK Gel Flat Top Batteries Pair w/strap \$975.00
<input type="checkbox"/>	116AG22	E2363	GP24 Brightway AGM Flat Top Batteries Pair w/strap \$735.00
<input type="checkbox"/>	116AG05	-	Omit batteries Includes harness & straps for standard post batteries N/C

Battery Charger			
<input type="checkbox"/>	116BJ19	-	8A Dual Mode (Compact) CTE N/C
<input type="checkbox"/>	116BJ3	-	Omit Charger -\$40.00

4. SEATING

Seat Width			
<input type="checkbox"/>	116S5	14" (not avail. w/ 19" & 20" seat depths)	N/C
<input type="checkbox"/>	116S6	15" (not avail. w/ 19" & 20" seat depths)	N/C
<input type="checkbox"/>	116S7	16" (not avail. w/19" & 20" seat depths and dual post arms)	N/C
<input type="checkbox"/>	116S8	17"	N/C
<input type="checkbox"/>	116S9	18"	N/C
<input type="checkbox"/>	116S10	19"	N/C
<input type="checkbox"/>	116S11	20"	N/C
<input type="checkbox"/>	116S27	22" (not avail. w/CG Power Tilt; Base width is 26.5" w/ 12"x4" tires)	N/C
<input type="checkbox"/>	116S28	24" (not avail. w/CG Power Tilt; Base width is 26.5" w/ 12"x4" tires)	N/C

Seat Height			
<input type="checkbox"/>	116HT10	18" (not avail. w/long ext. tubes; add 1.75" for CG Tilt)	N/C
<input type="checkbox"/>	116HT12	19" (add 1.75" for CG Tilt)	N/C
<input type="checkbox"/>	116HT14	20" (add 1.75" for CG Tilt)	N/C

Seat Base Type			
<input type="checkbox"/>	116SL2	Seat Sling Parapak (n/a w/CG Power Tilt)	N/C
<input type="checkbox"/>	116SL5	Seat Pan Aluminum	N/C
<input type="checkbox"/>	116SL6	Omit Seat Sling	-\$15.00

Seat Depth			
<input type="checkbox"/>	116SD7	14"	N/C
<input type="checkbox"/>	116SD8	15"	N/C
<input type="checkbox"/>	116SD9	16"	N/C
<input type="checkbox"/>	116SD10	17"	N/C
<input type="checkbox"/>	116SD11	18"	N/C
<input type="checkbox"/>	116SD12	19"	N/C
<input type="checkbox"/>	116SD73	20"	N/C

5. SEAT CUSHIONS

Standard Cushion			
<input type="checkbox"/>	116CU1	2" Cushion.	\$120.00

JAY Cushions (Optional)			
<i>To order any JAY cushions on-chair, please see the JAY price list.</i>			
<i>Cushion name or part number must be filled out w/ desired width & depth</i>			
<input type="checkbox"/>	116JC20	JAY Bundled Cushion Cushion Name or Part # Width: _____ Depth: _____	Price List

6. BACKS

Backrest Type			
<input type="checkbox"/>	116BT74	Folding Angle Adj. Backrest	\$115.00
<input type="checkbox"/>	116BT3	Non-Folding Angle Adj. Backrest	N/C

Back Height (Upholstery)			
<input type="checkbox"/>	116B4	Low 15.5"	N/C
<input type="checkbox"/>	116B8	Medium 17.5"	N/C
<input type="checkbox"/>	116B13	Tall 19.5"	N/C

Back Upholstery			
<input type="checkbox"/>	116BU27	Parapak Upholstery (select back height)	N/C
<input type="checkbox"/>	116BU26	K0669 Jay Solid Backrest (select back height; Decrease seat depth by 1")	N/C
<input type="checkbox"/>	116BU14	Omit Back Upholstery	-\$15.00

Jay Solid Backrest Back Height			
<input type="checkbox"/>	116B36	16"	N/C
<input type="checkbox"/>	116B37	17"	N/C
<input type="checkbox"/>	116B38	18"	N/C
<input type="checkbox"/>	116B39	19"	N/C
<input type="checkbox"/>	116B40	20"	N/C
<input type="checkbox"/>	116B41	21"	N/C
<input type="checkbox"/>	116B42	22"	N/C
<input type="checkbox"/>	116B43	23"	N/C
<input type="checkbox"/>	116B44	24"	N/C

Back Post Type			
<input type="checkbox"/>	116BP1	Straight Back Post	N/C
<input type="checkbox"/>	116BP41	8 Bend Back Post - Not avail. w/ Jay Solid backrest	\$105.00

Jay 3 Option			
<input type="checkbox"/>	116JB01	Jay J3 (select OMIT Back Uph above)	See J3 OF

7. POSTURAL SUPPORTS

Headrest (Incl. on Specified Specialty Controls)

JAY Solid Backrest Only

<input type="checkbox"/>	Headrest (8" x 4.5")		\$460.00
	116HE65	E0955 Curved Headrest 8" x 4.5"	\$230.00
	116MH1	E1028 Removable Mount	\$230.00
<input type="checkbox"/>	Headrest (9" x 4.5")		\$460.00
	116HE66	E0955 Curved Headrest 9" x 4.5"	\$230.00
	116MH1	E1028 Removable Mount	\$230.00

Pad Style

JAY Solid Backrest Only

<input type="checkbox"/>	Lateral Pads Curved		\$460.00
	116KE1	E0956 Lateral Pad Curved Pair	\$235.00
	116LM1	E1028 Swing Away Mount Pair	\$225.00
<input type="checkbox"/>	Lateral Pads Straight		\$450.00
	116KE2	E0956 Lateral Pad Straight Pair	\$225.00
	116LM1	E1028 Swing Away Mount Pair	\$225.00

Pad Size

<input type="checkbox"/>	116KD1	Lateral Support 3" x 6" (pr)	N/C
<input type="checkbox"/>	116KD2	Lateral Support 4" x 5" (pr)	N/C
<input type="checkbox"/>	116KD3	Lateral Support 4" x 6" (pr)	N/C
<input type="checkbox"/>	116KD4	Omit Lateral Support	-\$25.00

8. ARMRESTS

Armrest (Pair)

Single Post Height Adj. Armrests are N/C on Power CG Tilt

<input type="checkbox"/>	116AA131	Dual Post Flip Back Short (14-16in)	N/C
<input type="checkbox"/>	116AA132	Dual Post Flip Back Long (16-20in)	N/C
<input type="checkbox"/>	116AA51	E0973 Dual Post Hgt Adj - Short (14-16 in)	\$285.00
<input type="checkbox"/>	116AA52	E0973 Dual Post Hgt Adj - Long (16-20 in)	\$285.00
<input type="checkbox"/>	116AA8	E0973 Single Post Height Adjustable	\$270.00
<input type="checkbox"/>	116AA14	Omit Armrests and mounts	N/C

Armpads (Pair)

<input type="checkbox"/>	116AM7	Desk Length Waterfall (10") (n/a. w/ sgl. post Armrest w/Jay Solid Backrest)	N/C
<input type="checkbox"/>	116AM8	Full Length Waterfall (14") (n/a. w/ sgl. post Armrest w/Jay Solid Backrest)	N/C
<input type="checkbox"/>	116AM1	Desk Length Standard (10")	N/C
<input type="checkbox"/>	116AM2	Full Length Standard (14")	N/C

9. FOOTRESTS

Hanger (Pair)

<input type="checkbox"/>	116H12	70° Fixed Straight (only avail. w/ solid 1-piece footplate.; n/a. w/ CG Power Tilt)	N/C
<input type="checkbox"/>	116H2	70° Swing Away Hanger (Not. Avail w/ solid 1-piece footplate.)	N/C
<input type="checkbox"/>	116H13	80° Fixed Tapered (only avail. w/ solid 1-piece footplate.; n/a. w/ CG Power Tilt)	N/C
<input type="checkbox"/>	116H37	90° Swing Away (Only available with kids angle adjustable footplates)	N/C
<input type="checkbox"/>	116H16	E0990 Elevating Hanger (Req. aluminum footplate or angle adj.)	\$280.00
<input type="checkbox"/>	116H32	Omit Hangers	-\$15.00

Extension Tube

Not available with 90° hangers.

<input type="checkbox"/>	116E1	Short (13" - 15")	N/C
<input type="checkbox"/>	116E2	Medium (15" - 18") Not available with multi-position footplate.	N/C
<input type="checkbox"/>	116E5	Long (18" - 20") Not available with 70° or 80° hangers.	N/C

Footrest (Pair)

<input type="checkbox"/>	Composite w/ Heel Loops		N/C
	Not Avail. w/ fixed hangers		
	116F1	- Composite Footplate	N/C
	116HL1	- Heel Loop	\$30.00
<input type="checkbox"/>	Adult Angle Adjust. w/ heel loop		\$220.00
	Not Avail. w/ fixed hangers		
	116F282	K0040 Adult Angle Adj. Footplate	\$170.00
	116HL10	E0951 Heel Loops	\$50.00
<input type="checkbox"/>	Kids Angle Adjust. w/ heel loop		\$205.00
	Only avail. w/ 90° hangers		
	116F283	K0040 Kids Angle Adj. Footplate	\$170.00
	116LS13	E0951 Leg Straps	\$35.00
<input type="checkbox"/>	116F287	K0042 Aluminum (Only avail. with ELRs)	\$35.00
<input type="checkbox"/>	116F60	K0041 Solid One Piece not avail. w/Power CG Tilt	\$35.00
<input type="checkbox"/>	Two-piece Locking w/Heel Loop		\$220.00
	116F58	K0040 Two-piece Locking Footplate	\$170.00
	116HL10	E0951 Heel Loop	\$50.00
<input type="checkbox"/>	Multi-Positional w/Heel Loop		\$225.00
	Not Avail. w/ fixed hangers		
	116F47	K0040 Multi-Positional Footplate	\$175.00
	116HL10	E0951 Heel Loop	\$50.00
<input type="checkbox"/>	116F6	Omit Foot Plate w/Heel Loop	N/C

Footrest Accessories

<input type="checkbox"/>	116LS3	K0038 Extra Leg Strap, Black	\$35.00
--------------------------	--------	------------------------------	---------

10. ELECTRONICS

Drive Control Options

<input type="checkbox"/>	116ER04	R-Net Non-Expandable Drive Control (120A)	NC
<input type="checkbox"/>	R-Net Expandable Drive Control Upgrade (120A)		\$980.00
	116EX03	E2377 Expandable 120A Controller - R-net	\$540.00
	116EX04	E2313 Harnessing for 120A Expandable Controller	\$440.00

R-net Joystick (for expandable electronics)

<input type="checkbox"/>	116J70	-	Std LED Joystick for Expandable Ctrl's	N/C
<input type="checkbox"/>	116J108		Joystick w/Stereo Jacks & Color Display	\$210.00
<input type="checkbox"/>	116J88	-	Joystick W/Toggle, Spd Pot, Stereo Jacks & Color Display	\$410.00
<input type="checkbox"/>	116J95	-	Omit JS for use with OMNI	N/C

Joystick Mount (Fixed or Swing Away)

<input type="checkbox"/>	116JM10	-	Fixed Mount	N/C
<input type="checkbox"/>	116JM11	E1028	Swing Away Mount	\$245.00

Joystick Mount (Location)

Height Adjustable.

<input type="checkbox"/>	116HM1	-	Right Side Mount	N/C
<input type="checkbox"/>	116HM2	-	Left Side Mount	N/C

Joystick Handle (Optional)

<input type="checkbox"/>	116JH1	E2323	Ball Handle	\$80.00
<input type="checkbox"/>	116JH2	E2323	Foam Ball Handle	\$80.00
<input type="checkbox"/>	116JH4	E2323	T-Handle	\$80.00
<input type="checkbox"/>	116JH6	E2323	Goalpost	\$80.00

R-net OMNI and Output Options

(For expandable controls, choose all that apply)

<input type="checkbox"/>	116ON1	-	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$1,005.00
--------------------------	--------	---	--	------------

OMNI Mount

(If used with Joystick must be mounted opposite side)

<input type="checkbox"/>	116DS03	E1028	Swing Away Left Side OMNI Mount	\$245.00
<input type="checkbox"/>	116DS04	E1028	Swing Away Right Side OMNI Mount	\$245.00

Switch inputs for Joysticks or OMNI

Drive or Mode / Profile Control

<input type="checkbox"/>	116SW22		Stereo to 2 Mono Splitter	\$65.00
<input type="checkbox"/>	116SW6		Ribbon Switch	\$100.00
<input type="checkbox"/>	116SW1		Buddy Button Switch	\$80.00
<input type="checkbox"/>	116SW2		Disc Switch	\$165.00
<input type="checkbox"/>	116SW9		Wobble Switch	\$515.00
<input type="checkbox"/>	116SW4		Micro Light Switch	\$80.00
<input type="checkbox"/>	116SW24		Egg Switch-Black	\$100.00
<input type="checkbox"/>	116SW21		Dual toggle switch only	\$175.00
<input type="checkbox"/>	116SW23		4-way toggle switch w/5th button DB9 Connector.	\$175.00

Switch Mount Bracket

<input type="checkbox"/>	116EM01		Single Switch Mount Bracket	\$200.00
<input type="checkbox"/>	116EM02		Dual toggle switch mount bracket	\$200.00
<input type="checkbox"/>	116EM03		4-way toggle switch mount bracket.	\$200.00

Specialty Driver Control Options

Requires expandable controls and OMNI. All packages must be added to the specialty control on-chair package. Not available with Upholstered backs.

<input type="checkbox"/>	Quickie Sip n' Puff System		\$2,300.00
	116ID34	E2325 Sip N' Puff Kit	\$1,465.00
	116I01	E0955 Headrest	\$230.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
	116IH1	E2326 Sip n' Puff Tubing and Hardware	\$375.00
<input type="checkbox"/>	ASL Head Array Package (Adult)		\$4,140.00
	116I02	E2330 ASL Proximity Head Array (Adult)	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	ASL Head Array Package (Ped)		\$4,140.00
	116I03	E2330 ASL Proximity Head Array (Ped)	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	Switch-It Head Array Package (Adult)		\$4,140.00
	116I04	E2330 Switch-it Proximity Head Array	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	Switch-It Head Array Package (Pediatric)		\$4,140.00
	116I05	E2330 Switch-it Proximity Head Array	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	Quickie Proportional Chin Control Package		\$4,350.00
	116ID45	E2312 Proportional Joystick	\$4,120.00
	116IJ03	E1028 ASL Proportional S/A Joystick Mount	\$230.00

Attendant Control

<input type="checkbox"/>	116A02	E2331	Attendant Control for R-net	\$1,025.00
--------------------------	--------	-------	-----------------------------	------------

Attendant Mounting location

<input type="checkbox"/>	116DN1		Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	116DN2		Mounting Location: Left Back Cane	N/C

Electronic Programmers

<input type="checkbox"/>	116PR5		R-net PC Programmer with Dongle	\$415.00
<input type="checkbox"/>	116PR6		DTT Rnet cables included.	\$1,030.00

11. POWERED SEATING

Powered Seating Control Method (TILT)

(For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.)

<input type="checkbox"/>	116EO10		Switch Control 1 actuator w/dual toggle and mount. (For operating Tilt)	N/C
<input type="checkbox"/>	116EO11		Switch Control 1 actuator w/buddy button & mount. (For operating Tilt)	N/C

Switch Mount

(for use with switch control powered seating)

<input type="checkbox"/>	116EY01	-	Switch Mount Right	N/C
<input type="checkbox"/>	116EY02	-	Switch Mount Left	N/C

Powered Seating Thru-drive Control Method.

"Assignable Button" feature listed below (available on Rnet Systems) allows for direct "short-cut" access to various seating functions via joystick buttons. Functions accessed by depressing the button for more than 1 second. By selecting the thru-drive w/assignable button the chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

Thru-Drive 1 Actuator options.

<input type="checkbox"/>	116EO15	E2310	Thru Drive 1 Actuator w/out assignable buttons. For use with R-net (For operating Tilt thru Joystick.)	\$1,420.00
<input type="checkbox"/>	116EO13	E2310	Thru Drive 1 Actuator with assignable buttons. (For operating Tilt thru Joystick and/or Joystick. Requires R-net.)	\$1,420.00

12. WHEELS

Drive Wheel Type

<input type="checkbox"/>	116DW5		12" x 4" Pneumatic Tire	N/C
<input type="checkbox"/>	116DW6		12" x 4" Pneumatic w/Ins.	N/C
<input type="checkbox"/>	116DW27		12" x 3" Pneumatic Tire	N/C
<input type="checkbox"/>	116DW28		12" x 3" Pneumatic w/Ins.	N/C

The 12" x 3" wheel size reduces the overall width specification by approx. 2" and also alters the overall appearance of the chair.

Drive Wheel Tire

<input type="checkbox"/>	116DT1		Knobby Tire	N/C
<input type="checkbox"/>	116DT2		V Groove Tire (n/a with 12"x3" tires)	N/C

Wheel Lock

<input type="checkbox"/>	116WL26		Push to Lock	\$90.00
<input type="checkbox"/>	116WC1		Do Not Mount	N/C

Wheel Lock Extension Handle

Wheel-lock extensions handles are not available with 22" & 24" seats

<input type="checkbox"/>	116EH1	E0961	Extension Handle 6"	\$65.00
--------------------------	--------	-------	---------------------	---------

13. CASTERS

Front Caster Wheel

<input type="checkbox"/>	116CW17		8" Pneumatic	N/C
<input type="checkbox"/>	116CW20		8" Airless with Insert	N/C
<input type="checkbox"/>	116CW29		9" Pneumatic (not avail. With Suspension Fork)	\$90.00
<input type="checkbox"/>	116CW28		9" Pneumatic w/Insert (not avail. With Suspension Fork)	\$90.00

Caster Fork

Suspension Fork Options are only available with 8" casters.

<input type="checkbox"/>	116CF31		Solid Caster Fork	N/C
<input type="checkbox"/>	116CF35	E1016	Suspension Fork Med. 150-250 lbs. (Only avail. w/ 8" front casters)	\$270.00
<input type="checkbox"/>	116CF36	E1016	Suspension Fork Heavy. 251-350 lbs. (Only avail. w/ 8" front casters)	\$270.00

14. MANUAL CONVERSION KIT

<input type="checkbox"/>	116CV2	Manual Wheelchair Conversion Kit	\$1,025.00
--------------------------	--------	----------------------------------	------------

Axle

<input type="checkbox"/>	116AX1	Quick Release	N/C
<input type="checkbox"/>	116AX11	Quad Quick Release	\$95.00

Wheel Rim/Size

<input type="checkbox"/>	116RW1	24" Mag	N/C
<input type="checkbox"/>	116RW2	24" Spoke (not avail w/14" seat width)	N/C
<input type="checkbox"/>	116RW7	Omit Manual Rear Wheel	-\$55.00

Rear Wheel Tire

<input type="checkbox"/>	116RT1	E2211 E2212	Pneumatic	N/C
<input type="checkbox"/>	116RT2	E2211 E2213	Pneumatic Airless Insert	\$85.00
<input type="checkbox"/>	116RT3	E2218	Full Polyurethane	N/C
<input type="checkbox"/>	116RT4		Low Polyurethane	N/C
<input type="checkbox"/>	116RT6	E2211 E2212	High Pressure Clincher (not available with mag wheels)	\$100.00

Hand Rim

<input type="checkbox"/>	116HR1	Aluminum Anodized	N/C
<input type="checkbox"/>	116HR4	Plastic Coated	\$90.00
<input type="checkbox"/>	116HR7	Omit hand rim	-\$20.00

Projections

<input type="checkbox"/>	116HR8	8 Oblique	\$200.00
--------------------------	--------	-----------	----------

Wheel Locks

<input type="checkbox"/>	116WC2	Push to Lock Wheel Lock	\$45.00	
<input type="checkbox"/>	116WC1	Do not mount	N/C	
<input type="checkbox"/>	116ET1	E0961	6" Extension Handles	\$55.00

15. ACCESSORIES

Backpack

<input type="checkbox"/>	116A1	Black Utility Backpack	\$55.00
<input type="checkbox"/>	116A2	Black/Charcoal Kids	\$50.00

Seat Pouch

<input type="checkbox"/>	116A8	Black/Charcoal	\$50.00
--------------------------	-------	----------------	---------

Positioning Belt

<input type="checkbox"/>	116A13	E0978	2" Aircraft Buckle Padded	\$125.00
<input type="checkbox"/>	116A12	E0978	2" Aircraft Buckle	\$115.00
<input type="checkbox"/>	116A142	E0978	1.5" Auto Buckle	\$75.00

Wheel Chair Tray Table

Requires standard Armpads to mount.

<input type="checkbox"/>	116A31	E0950	Medium Tuffak (15"-17") Requires std. Armpads to mount tray	\$260.00
<input type="checkbox"/>	116A32	E0950	Large Tuffak (18"-20") Requires std. Armpads to mount tray	\$260.00
<input type="checkbox"/>	116A33	E0950	X-Large Tuffak (21"-22") Requires std. Armpads to mount tray	\$260.00

Miscellaneous

<input type="checkbox"/>	16A8105A1	Caddy	\$50.00
<input type="checkbox"/>	116A114	Dual Post Armrest Pouch	\$25.00

Additional Specifications

Weight:	115 lbs. (without batteries)
Overall Width	Using 12" x 4" Drive Tire: 24.5" using 14"-20" Seat Width 26.5" using 22" and 24" Seat Width
Overall Length:	31.5" (Without Footrests) Deduct 2" if using 12" x 3" Tire
Range:	As far as 30 miles (With Group 24 batteries depending on conditions and user environment)
Top Speed:	Up to 8.5 mph



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA