



## PERSONALIZE YOUR JAY<sup>®</sup> OR AES CUSTOM SEATING

**Vinyl** is water resistant, easy to clean and offers some anti-bacterial and anti-mildew benefits. It is recommended for planar seats and backs. Vinyl can be used on non-contact surfaces to add color to your seating system.


**Stretch** upholstery is recommended for contoured seats and backs, anti-thrust seats and similar items. Darlexx<sup>®</sup> and Dartex<sup>®</sup> allow four-way stretch that is anti-microbial, anti-fungal, waterproof and breathable. Dartex is two-sided allowing you to choose between a cloth surface or a smooth surface that is easy to wipe clean. Neoprene, sometimes referred to as “wet-suit material”, allows four-way stretch. Neoprene is a soft nylon material with rubber backing. It is waterproof but not breathable.

**Cordora** is a heavy-duty, textured nylon fabric that dries quickly if wetted.


**Nylon** is a finer weave and smoother texture than Cordora but offers the same excellent durability.

**3DX<sup>™</sup>** Microclimatic protects the skin by dissipating heat and moisture.

### Stretch




### Cordora



### Nylon

### 3DX<sup>™</sup>

### Texture Vinyl



See page 2 for other Vinyl options



## Vinyl

