



**Breezy 600**  
EIB60, January 2020

**ACCOUNT INFORMATION**

Account No: \_\_\_\_\_ Drop Ship/Ship To: \_\_\_\_\_  
 Date: \_\_\_\_\_ Name: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_ Address: \_\_\_\_\_  
 Buyer: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Phone Number: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**1. FRAME**

**BREEZY 600**

Hangers not included in base price

- 086WM20 K0004 Breezy 600 \$1,565.00

**2. SIZE / SEATING**

**Seat to Floor Height Range**

- Standard (15" - 20 1/2", see matrix) Std
- Ultra Hemi (13" - 18", see matrix) NC

**Seat Width**

- 086W5 14" Std
- 086W7 16" Std
- 086W9 18" Std
- 086W42 E2201 20" \$140.00

**Seat Depth**

- 086FD7 16" Std
- 086FD9 18" NC

**Seat Height**

- 086HT1 13" \$95.00
- 086HT2 14" \$95.00
- 086HT3 14 1/2 " \$95.00
- 086HT4 15" \$95.00
- 086HT5 15 1/2" \$95.00
- 086HT6 16" \$95.00
- 086HT7 16 1/2" \$95.00
- 086HT8 17" NC
- 086HT9 17.5" NC
- 086HT10 18" NC
- 086HT11 18.5" NC
- 086HT12 19" <sup>12</sup> NC
- 086HT13 19.5" <sup>12</sup> NC
- 086HT14 20" <sup>12</sup> NC
- 086HT16 20.5" <sup>12</sup> NC

**Seat Sling**

- 086SL3 Padded Seat Sling \$15.00
- 086SL4 Non-padded Seat Sling <sup>2</sup> Std
- 086SL6 Omit Seat Sling -\$15.00

**Cushion**

To order any JAY cushions on-chair, please see JAY price list.

- 086JC01 JAY Basic Price List
- 086JC02 JAY Easy Price List
- 086JC03 JAY Soft Combi P Price List
- 086JC04 JAY J2 Price List
- 086JC05 JAY J2 Deep Contour Price List
- 086JC06 JAY Lite Price List
- 086JC07 JAY J3 Price List
- 086CU1 K0669 2" Cushion \$150.00

### 3. CASTERS

#### CASTERS

<input type="checkbox"/>	5" Polyurethane	NC
<input type="checkbox"/>	8" Polyurethane <sup>16</sup>	Std
<input type="checkbox"/>	8" Pneumatic <sup>6</sup>	\$85.00

**Footnotes:**

- <sup>1</sup> Not available when using Standard frame.
- <sup>2</sup> For use with cushion.
- <sup>3</sup> Only available with standard hemi extension tube.
- <sup>4</sup> Not available with Ultra Hemi frame.
- <sup>5</sup> Not available with 8" pneumatic casters.
- <sup>6</sup> Not available with 20" rear wheels when using Standard frame.
- <sup>7</sup> Only available with 8" casters.
- <sup>8</sup> Not available with 14" seat width.
- <sup>9</sup> Pull-to-lock wheel locks required.
- <sup>10</sup> Not available with composite footplates.
- <sup>11</sup> Not available with articulating legrest.
- <sup>12</sup> Must use 1-1/2" stem.
- <sup>13</sup> Not available with ELR.
- <sup>14</sup> Only available with ELR
- <sup>15</sup> To order a J3 Back please select Omit Back Upholstery.  
Please note: The J3 Back requires back canes or push handles.
- <sup>16</sup> Not available with 70° hangers due to interference

#### Caster Options

<input type="checkbox"/>	1-1/2" Fork Stem (Double Bearing)	\$35.00
<input type="checkbox"/>	086CP13 Double Bearing Caster Plate	\$35.00

### 4. REAR WHEELS

#### Rear Wheel Size

<input type="checkbox"/>	086WS4	20"	NC
<input type="checkbox"/>	086WS5	22"	NC
<input type="checkbox"/>	086WS6	24"	Std

#### Rear Wheel Tires

<input type="checkbox"/>	086RT4	Low-profile Polyurethane <sup>17,21</sup>	Std
<input type="checkbox"/>	086RT29	E2218 Full-profile Polyurethane <sup>18,21</sup>	\$50.00
<input type="checkbox"/>	086RT35	E2211 Pneumatic <sup>21</sup>	\$165.00
<input type="checkbox"/>	086RT33	E2211 Pneumatic w/Airless Insert <sup>21,22</sup>	\$215.00

#### Hand Rims

<input type="checkbox"/>	086HR5	Composite <sup>21</sup>	Std
<input type="checkbox"/>	086HR2	Aluminum - painted <sup>21</sup>	NC
<input type="checkbox"/>	086HR1	Aluminum - anodized <sup>21</sup>	\$0.00
<input type="checkbox"/>	086HR4	Plastic Coated <sup>21</sup>	\$110.00
<input type="checkbox"/>	086HR7	Omit Hand Rims <sup>21</sup>	NC

#### Axle

<input type="checkbox"/>	086AL4	Threaded Offset Axles	Std
<input type="checkbox"/>	086AX43	Quick-Release Rear Axles	\$95.00
<input type="checkbox"/>	086OA2	E0958 One-Arm Drive <sup>23,24</sup>	\$850.00
<input type="checkbox"/>	086HM2	Left Mount	
<input type="checkbox"/>	086HM1	Right Mount	

#### Anti-Tip Tubes

<input type="checkbox"/>	086AT1	E0971 Rear Anti-Tip Tubes	\$150.00
--------------------------	--------	---------------------------	----------

#### Projections<sup>19, 21</sup>

<input type="checkbox"/>	086HR8	Oblique - 8	\$325.00
--------------------------	--------	-------------	----------

#### Wheel Locks

<input type="checkbox"/>	086WL1	High-Push Breezy-Style	Std
<input type="checkbox"/>	086WL2	High-Pull Breezy-Style	NC
<input type="checkbox"/>	086WL22	High-Push Quickie-Style	\$40.00
<input type="checkbox"/>	086WL21	High-Pull Quickie-Style	\$40.00

#### Wheel Lock Options

<input type="checkbox"/>	086EH1	E0961 6" Ext. Handles with Bungee	\$75.00
<input type="checkbox"/>	086GA1	E0974 Grade Aids <sup>20,25</sup>	\$150.00

**See Sunrise Seating Solutions on reverse side.**

**TOTAL CHAIR PRICE: \$ \_\_\_\_\_**

**FOOTNOTES:**

- <sup>17</sup> Not available with Projection Hand Rims.
- <sup>18</sup> "No Charge" option when used with Projection Hand Rims.
- <sup>19</sup> Not available with Low-profile polyurethane tires. Full-profile Polyurethane tires available as "No Charge" option. Requires 24" rear wheels.
- <sup>20</sup> Not available with Polyurethane tires.
- <sup>21</sup> Not available with 24" Mag SLT
- <sup>22</sup> Must order Pneumatic tire with Flat Free Insert
- <sup>23</sup> Not available with Ultra Hemi frame.
- <sup>24</sup> Not available with Quick-Release Rear Axles.
- <sup>25</sup> Grade Aids require Quickie style wheel locks

## 5. FOOTREST

Hanger Type			
<input type="checkbox"/>	086H4	60° Swing-Away Hemi <sup>1,3</sup>	\$140.00
<input type="checkbox"/>	086H131	60° Swing-Away <sup>4</sup>	\$140.00
<input type="checkbox"/>	086H2	70° Swing-Away <sup>4,5</sup>	\$140.00
<input type="checkbox"/>	086H24	K0053 Articulating Legrests <sup>4,6,7,8,9</sup>	\$370.00
<input type="checkbox"/>	086H14	E0990 Elevating Legrests <sup>4,10</sup>	\$315.00

Footplates			
<input type="checkbox"/>	086F1	E0951 Composite w/Heel Loops <sup>13</sup>	Std
<input type="checkbox"/>	086F5	Aluminum <sup>14</sup>	NC
<input type="checkbox"/>	086F3	E0951 K0040 Adult Angle-Adj. w/Heel Loops <sup>16</sup>	\$190.00

Footrest Extension Tubes			
<input type="checkbox"/>	086E1	Short (13"-15 1/2") <sup>4,11</sup>	NC
<input type="checkbox"/>	086E2	Medium (15 1/2" - 18") <sup>4,6</sup>	Std
<input type="checkbox"/>	086E44	Long (18" - 20 1/2") <sup>4,6,12</sup>	\$35.00
<input type="checkbox"/>	086E21	Omit Ext. Tubes & Footplates	NC

## 6. BACKREST / SEATING

### BACKREST (16"-19", 8° w/Push Handles)

<input type="checkbox"/>	JAY J3 Back <sup>15</sup>		
	086JB01	Please see the J3 Back order form for part #'s.	
<input type="checkbox"/>	086B12	Adjustable w/tall back canes <sup>1</sup>	\$15.00
<input type="checkbox"/>	086BU14	Omit Back Upholstery	-\$15.00

## 7. ARMRESTS

Flip Back Armrests			
<input type="checkbox"/>		Fixed Height, Desk Length	Std
<input type="checkbox"/>	086AA5	Fixed Height, Full Length	NC
<input type="checkbox"/>	086AA6	E0973 Flip Back Height-Adjustable, Desk	\$295.00
<input type="checkbox"/>		E0973 Flip Back Height-Adjustable, Full	\$295.00

## 8. KEY OPTIONS AND ACCESSORIES

Backpack & Seat Pouch			
<input type="checkbox"/>	086A1	Utility Back Pack (Black/Charcoal)	\$75.00
<input type="checkbox"/>	086A8	Seat Pouch (Black/Charcoal)	\$50.00

Wheelchair Tray Table - Tuffak			
<input type="checkbox"/>	086A31	E0950 Medium (15"-17")	\$250.00
<input type="checkbox"/>	086A32	E0950 Large (18"-20")	\$250.00

Pelvic Positioning Belt			
<input type="checkbox"/>	086A12	E0978 Aircraft Buckle (2" wide)	\$110.00
<input type="checkbox"/>	086A142	E0978 1.5" Auto Buckle	\$70.00

Miscellaneous			
<input type="checkbox"/>	086A82	Caddy	\$50.00
<input type="checkbox"/>	086A51	E2207 Crutch Holder	\$195.00
<input type="checkbox"/>	086A83	Tool Kit	\$55.00

## 9. COLOR

Frame Color			
<input type="checkbox"/>	086C6	Forest Green	Std
<input type="checkbox"/>	086C16	Silver Texture	NC
<input type="checkbox"/>	086C1	Blue	NC
<input type="checkbox"/>	086C2	Black	NC



Sunrise Medical (US) LLC  
2842 Business Park Ave. · Fresno, CA 93727 · USA

The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Standard** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
15.0"	5" lp	N/A	N/A	N/A	N/A
16.0"	5" lp	5" lp	N/A	N/A	N/A
16.5"	5" lp		N/A	N/A	N/A
17.0"	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	N/A
17.5"	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp	5" lp & 1-1/2" stem or 8" lp
18.0"	N/A	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem
18.5"	N/A	8" lp *20" full poly not recommended	8" lp	8" lp	8" lp
19.0"	N/A	N/A	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu
19.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem
20.0"	N/A	N/A	N/A	N/A	8" lp & 1-1/2" stem
20.5"	N/A	N/A	N/A	N/A	8" pneu & 1-1/2" stem

The chart below specifies which casters and rear wheels are necessary to achieve the front seat heights listed below for **Ultra Hemi** frames.

Seat Height	Rear Wheels				
	20" Low Poly	20" Full Poly/Pneu 22" Low Poly	22" Full Poly/Pneu	24" Low Poly	24" Full Poly/Pneu
13.0"	5" lp	N/A	N/A	N/A	N/A
14.0"	5" lp	5" lp	N/A	N/A	N/A
14.5"	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	5" lp & 1-1/2" stem	N/A
15.0"	8" lp	8" lp	8" lp	8" lp	N/A
15.5"	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem	8" lp or 5" lp & 1-1/2" stem
16.0"	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu	8" lp or 8" pneu
16.5"	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem	8" lp & 1-1/2" stem
17.0"	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu	8" lp & 1-1/2" stem or 8" pneu
17.5"	N/A	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu	1-1/2" stem & 8" lp or 8" pneu
18.0"	N/A	N/A	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem	8" pneu & 1-1/2" stem