



Advance Bath Chair

August 2020

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

 Confirm Via: Fax Email

 Submitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Advance Bath Chair

Weight Capacity: 159 lbs

Packages

Package includes Chair Frame, Chair Cover, Hip Belt and Chest Belt. Substitutions are not available for packages.

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> PABC1P	Size 1 Advance Bath Chair Package - Pink		\$804
<input type="checkbox"/> PABC1B	Size 1 Advance Bath Chair Package - Royal Blue		\$804
<input type="checkbox"/> PABC2P	Size 2 Advance Bath Chair Package - Pink		\$809
<input type="checkbox"/> PABC2B	Size 2 Advance Bath Chair Package - Royal Blue		\$809
<input type="checkbox"/> PABC3P	Size 3 Advance Bath Chair Package - Pink		\$863
<input type="checkbox"/> PABC3B	Size 3 Advance Bath Chair Package - Royal Blue		\$863
<input type="checkbox"/> PABC4P	Size 4 Advance Bath Chair Package - Pink		\$880
<input type="checkbox"/> PABC4B	Size 4 Advance Bath Chair Package - Royal Blue		\$880

Chair Frame

Selection Required if an Advance Bath Chair Package was not selected.

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> LAB/1/FR	Size 1 Chair Frame		\$473
<input type="checkbox"/> LAB/2/FR	Size 2 Chair Frame		\$478
<input type="checkbox"/> LAB/3/FR	Size 3 Chair Frame		\$521
<input type="checkbox"/> LAB/4/FR	Size 4 Chair Frame		\$538

Chair Cover

Selection Required. Must match chair frame size.

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> LAB/1/D-PK	Size 1 Cover - Pink		\$148
<input type="checkbox"/> LAB/1/D-RB	Size 1 Cover - Royal Blue		\$148
<input type="checkbox"/> LAB/2/D-PK	Size 2 Cover - Pink		\$148
<input type="checkbox"/> LAB/2/D-RB	Size 2 Cover - Royal Blue		\$148
<input type="checkbox"/> LAB/3/D-PK	Size 3 Cover - Pink		\$159
<input type="checkbox"/> LAB/3/D-RB	Size 3 Cover - Royal Blue		\$159
<input type="checkbox"/> LAB/4/D-PK	Size 4 Cover - Pink		\$159
<input type="checkbox"/> LAB/4/D-RB	Size 4 Cover - Royal Blue		\$159

Accessories

Abduction Belts must match chair frame size.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LAB/1/E-PK	Hip Belt - Pink		\$89
<input type="checkbox"/>	LAB/1/E-RB	Hip Belt - Royal Blue		\$89
<input type="checkbox"/>	LAB/1/A-PK	Chest Belt - Pink		\$94
<input type="checkbox"/>	LAB/1/A-RB	Chest Belt - Royal Blue		\$94
<input type="checkbox"/>	LAB/1/B-PK	Size 1 and 2 Abduction Belt - Pink		\$90
<input type="checkbox"/>	LAB/1/B-RB	Size 1 and 2 Abduction Belt - Royal Blue		\$90
<input type="checkbox"/>	LAB/2/B-PK	Size 3 and 4 Abduction Belt - Pink		\$91
<input type="checkbox"/>	LAB/2/B-RB	Size 3 and 4 Abduction Belt - Royal Blue		\$91
<input type="checkbox"/>	LAB/1/C-PK	Head Laterals - Pink		\$90
<input type="checkbox"/>	LAB/1/C-RB	Head Laterals - Royal Blue		\$90
<input type="checkbox"/>	LSF/1	Shower Trolley		\$777

Assembly

	Item Code	Description	Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee	\$120

Technical Specifications - Imperial

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	30" / 41"	37" / 49"	45" / 57"	53" / 65"
Head Support	6"	6"	8.5"	8.5"
Back Support Height	14"	16"	18"	22"
Seat Depth	9"	11"	13"	16"
Calf Support	10"	10"	13"	13"
Seat Height	2" - 17"	2" - 17"	2" - 17"	2" - 17"
Maximum Weight	159 lbs	159 lbs	159 lbs	159 lbs
Total Width	17"	17"	17"	17"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3	Size 4
Age	1 - 5 years	4 - 9 years	8 - 14 years	12 - 18 years
User Height (Max / Min)	750 mm / 1050 mm	950 mm / 1250 mm	1150 mm / 1450 mm	1350 mm / 1650 mm
Head Support	150 mm	150 mm	215 mm	215 mm
Back Support Height	350 mm	400 mm	460 mm	550 mm
Seat Depth	225 mm	275 mm	335 mm	400 mm
Calf Support	260 mm	260 mm	340 mm	340 mm
Seat Height	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm	50 mm - 425 mm
Maximum Weight	72 kgs	72 kgs	72 kgs	72 kgs
Total Width	430 mm	430 mm	430 mm	430 mm



Sunrise Medical (US) LLC 2842 Business Park Ave. • Fresno, CA 93727 • USA