



JAY[®] BASIC CUSHION

E2601 / E2602

February 2023

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY[®] Basic Cushion - E2601 / E2602

JAY BASIC CUSHION Less than 22" in Width

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	331-MJ	E2601	JAY BASIC Cushion 14" X 14"	\$105
<input type="checkbox"/>	308	E2601	JAY BASIC Cushion 14" X 16"	\$105
<input type="checkbox"/>	332-MJ	E2601	JAY BASIC Cushion 14" X 18"	\$105
<input type="checkbox"/>	356	E2601	JAY BASIC Cushion 15" X 16"	\$105
<input type="checkbox"/>	333-MJ	E2601	JAY BASIC Cushion 15" X 18"	\$105
<input type="checkbox"/>	305-MJ	E2601	JAY BASIC Cushion 16" X 16"	\$105
<input type="checkbox"/>	334-MJ	E2601	JAY BASIC Cushion 16" X 17"	\$105
<input type="checkbox"/>	306	E2601	JAY BASIC Cushion 16" X 18"	\$105
<input type="checkbox"/>	376	E2601	JAY BASIC Cushion 17" X 16"	\$105
<input type="checkbox"/>	335-MJ	E2601	JAY BASIC Cushion 17" X 17"	\$105
<input type="checkbox"/>	336	E2601	JAY BASIC Cushion 17" X 18"	\$105
<input type="checkbox"/>	300	E2601	JAY BASIC Cushion 18" X 16"	\$105
<input type="checkbox"/>	337	E2601	JAY BASIC Cushion 18" X 17"	\$105
<input type="checkbox"/>	301	E2601	JAY BASIC Cushion 18" X 18"	\$105
<input type="checkbox"/>	338	E2601	JAY BASIC Cushion 18" X 20"	\$105
<input type="checkbox"/>	339	E2601	JAY BASIC Cushion 19" X 16"	\$105
<input type="checkbox"/>	330	E2601	JAY BASIC Cushion 19" X 17"	\$105
<input type="checkbox"/>	341	E2601	JAY BASIC Cushion 19" X 18"	\$105
<input type="checkbox"/>	343	E2601	JAY BASIC Cushion 19" X 20"	\$105
<input type="checkbox"/>	344	E2601	JAY BASIC Cushion 20" X 16"	\$105
<input type="checkbox"/>	309	E2601	JAY BASIC Cushion 20" X 18"	\$105
<input type="checkbox"/>	345	E2601	JAY BASIC Cushion 20" X 20"	\$105

JAY BASIC CUSHION 22" in Width and Greater

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	328	E2602	JAY BASIC Cushion 22" X 18"	\$229
<input type="checkbox"/>	320	E2602	JAY BASIC Cushion 22" X 20"	\$229
<input type="checkbox"/>	346	E2602	JAY BASIC Cushion 22" X 22"	\$229
<input type="checkbox"/>	348	E2602	JAY BASIC Cushion 24" X 18"	\$229
<input type="checkbox"/>	340	E2602	JAY BASIC Cushion 24" X 20"	\$229
<input type="checkbox"/>	347	E2602	JAY BASIC Cushion 24" X 22"	\$229
<input type="checkbox"/>	349	E2602	JAY BASIC Cushion 24" X 24"	\$229


JAY Basic Size Matrix

		WIDTH												
		12	13	14	15	16	17	18	19	20	21	22	23	24
DEPTH	12	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	13	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	14	JYW	JYW		JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW
	15	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW
	16	JYW	JYW								JYW	JYW	JYW	JYW
	17	JYW	JYW	JYW	JYW					JYW	JYW	JYW	JYW	JYW
	18	JYW	JYW								JYW	JYW		JYW
	19	N/A	N/A	N/A	N/A	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW	JYW
	20	N/A	N/A	N/A	N/A	JYW	JYW				JYW		JYW	
	21	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	JYW	JYW	JYW	JYW	JYW
	22	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	JYW	JYW		JYW	
	23	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	JYW	JYW	JYW
24	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	JYW	JYW		

	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

JAY Basic Measuring Guide

Measurements made from completed cushion with covers.



Width

X = Width

Note: Measurement made from outside edge to outside edge of cushion.

Ex P/N: 300 X = 18"

Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower edge of cushion.

Ex P/N: 300 Y = 16"

JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.

Please refer to the "JAY YOUR WAY" product catalog for complete listing of modifications

Please visit www.sunparts.us for spare parts information.

