

All-terrain MWD power wheelchair.

EIMM

Customer				Ship To			
Account				Name			
Account #			Date	Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City	State		
Buyer				ZIP	Phone		
Mark For				Email / Fax			

Please submit orders to: magicmobility@sunmed.com

Seating Specifications


HCPK	K0848 K0856 K0861 K0850 K0858 K0862 K0849 K0857 K0851 K0859	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	17", 18", 19"

Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 90A controller
User Weight Capacity	400 lbs. 340 lbs. with Elevate or Tilt	Drive Motors	V6 4-pole 6mph
V6 Base Size: All Terrain	Width 28" Length 41.4" Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels 14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3" Length 41.5" Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	4 frame colors
		Ground Clearance	3.5"

Client Measurements

a)	Top of head to seat	in
b)	Top of shoulder to seat	in
c)	Underarm to seat	in
d)	Elbow to seat	in
e)	Top of leg to seat	in
f)	Shoulder width	in
g)	Chest width	in
i)	Hip width	in
j)	Actual sitting depth	in
m)	Foot length	in
n)	Knee to foot	in
o)	Seat width	in



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL		
SEATING OPTIONS			\$		
<input type="checkbox"/> MM-RHBSB		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC		
<input type="checkbox"/> MM-MPSSB		MPS Seat and Back - Captain style seating	NC		
<input type="checkbox"/> MM-MPSBRHBS		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC		
BASE			\$		
Frontier V6 Base		End Item Number: EIMM	\$		
<input type="checkbox"/> MM-V6NP	K0848	Magic Mobility - Frontier V6 No Power Seating Options	\$ 15,724.00		
<input type="checkbox"/> MM-V6SP	K0856	Magic Mobility - Frontier V6 Single Power Seating Options	\$ 15,724.00		
<input type="checkbox"/> MM-V6MP	K0861	Magic Mobility - Frontier V6 Multi Power Seating Options	\$ 15,724.00		
<input type="checkbox"/> MM-V6NPHD	K0850	Magic Mobility - Frontier V6 No Power Seating Options - HD	\$ 15,724.00		
<input type="checkbox"/> MM-V6SPHD	K0858	Magic Mobility - Frontier V6 Single Power Seating Options - HD	\$ 15,724.00		
<input type="checkbox"/> MM-V6MPHD	K0862	Magic Mobility - Frontier V6 Multi Power Seating Options - HD	\$ 15,724.00		
<input type="checkbox"/> MM-V6NPC	K0849	Magic Mobility - Frontier V6 No Power Seating Options - MPS (Captain)	\$ 15,724.00		
<input type="checkbox"/> MM-V6SPC	K0857	Magic Mobility - Frontier V6 Single Power Seating Options - MPS (Captain)	\$ 15,724.00		
<input type="checkbox"/> MM-V6MPC		Magic Mobility - Frontier V6 Multi Power Seating Options - MPS (Captain)	\$ 15,724.00		
<input type="checkbox"/> MM-V6NPCHD	K0851	Magic Mobility - Frontier V6 No Power Seating Options - HD MPS (Captain)	\$ 15,724.00		
<input type="checkbox"/> MM-V6SPCHD	K0859	Magic Mobility - Frontier V6 Single Power Seating Options - HD MPS (Captain)	\$ 15,724.00		
<input type="checkbox"/> MM-V6MPCHD		Magic Mobility - Frontier V6 Multi Power Seating Options - HD MPS (Captain)	\$ 15,724.00		
Frame Color			\$		
<input type="checkbox"/> MM-BKFRM	Black Frame Color	NC	<input type="checkbox"/> MM-AGFRM	Mossy Green Frame Color	NC
<input type="checkbox"/> MM-RDFRM	Red Frame Color	NC	<input type="checkbox"/> MM-SBFRM	Royal Blue Frame Color	NC
Batteries			\$		
<input type="checkbox"/> MM-BAT24G	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 932.00		
<input type="checkbox"/> MM-BAT34G	E2359	Group 34 Gel Battery (Pair) - Installed	\$ 798.00		
Battery Charger			\$		
<input type="checkbox"/> MM-BATCHG		Battery Charger	NC		
Drive Wheels			\$		
<input type="checkbox"/> MM-ATDW		All Terrain Drive Wheels and Tires	NC		
<input type="checkbox"/> MM-HPDW		Hybrid Drive Wheels and Tires - Pneumatic	NC		
<input type="checkbox"/> MM-HSDW		Hybrid Drive Wheels and Tires - Solid	\$ 137.00		
Tie Downs			\$		
<input type="checkbox"/> MM-V6TK		Crash Tested V6 Tie Down Kit - Always included	NC		
Castor Options			\$		
<input type="checkbox"/> MM-CASTP		Pneumatic Castor 2.50-4 - Non-marking	NC		
<input type="checkbox"/> MM-CASTS		Foam Filled Castor 2.50-4 - Non-marking	\$ 67.00		

BASE

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Optional Drive Tire Options			\$
<input type="checkbox"/> MM-TS2		Tire Sealant for 2 Pneumatic Drive Tires	\$ 34.00
Base Accessories			\$
<input type="checkbox"/> MM-01595		Jack for Drive Tire	\$ 180.00

SEATING: SEATING SYSTEM

No Power Seat Options \$

<input type="checkbox"/> MM-NP		No Power Seating System - 400 lb. weight capacity	\$ 3,675.00
<input type="checkbox"/> MM-NPFO		No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 4,725.00

Power Tilt \$

Utilized without Power Recline.

<input type="checkbox"/> MM-TLT	E1002	50° Power Tilt in Space - 340 lb. weight capacity	\$ 6,798.00
---------------------------------	-------	---------------------------------------------------	-------------

Power Seat Elevate \$

Only available when Power Tilt or Power Tilt and Recline is selected.

<input type="checkbox"/> MM-ELVT	E2300	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 3,938.00
----------------------------------	-------	--------------------------------------------------------	-------------

Power Recline \$

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

<input type="checkbox"/> MM-PRMPS	E1003	150° Power Recline for MPS Seating - 400 lb. weight capacity	\$ 4,836.00
<input type="checkbox"/> MM-PRRHB	E1003	150° Power Recline with Shear Reduction for REHAB Seating - 400 lb. weight capacity	\$ 4,836.00

Power Tilt and Recline \$

<input type="checkbox"/> MM-TLTPRMP	E1006	50° Tilt & 150° Power Recline - MPS Seating - 340 lb. weight capacity	\$ 10,735.00
<input type="checkbox"/> MM-TLTPRRHB	E1007	50° Tilt & 150° Power Recline with Shear Reduction - REHAB Seating - 340 lb. weight capacity	\$ 12,353.00

Power Seat Function \$

<input type="checkbox"/> MM-SFCK	E2310	Single Seat Function Control Kit - Must be selected if ordering 1 actuator function	\$ 1,573.00
<input type="checkbox"/> MM-MFCK	E2311	Multiple Seat Function Control Kit - Must be selected if ordering 2 or more actuator functions	\$ 3,203.00

Seat To Floor Height \$

V6 Seat to Floor Heights start at 17" +/- 0.25"

<input type="checkbox"/> MM-STFH		Lowest Seat To Floor Height	NC
<input type="checkbox"/> MM-STFH1		Seat To Floor Height +1"	NC
<input type="checkbox"/> MM-STFH2		Seat To Floor Height +2"	NC

Power Seat Options \$

<input type="checkbox"/> MM-CRTL5JS		Ctrl+5 - 5 Button Box (Joystick Side) - For controlling power seat functions. Includes USB port	\$ 843.00
<input type="checkbox"/> MM-CRTL5OS		Ctrl+5 - 5 Button Box (Opposite Side) - For controlling power seat functions. Includes USB port	\$ 843.00

SEATING: SEAT SIZES

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
MPS (Captain) Seat Cushion (width x depth)			\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
<input type="checkbox"/> MM-01673		16" x 16" Vinyl Seat	NC
<input type="checkbox"/> MM-01674		16" x 18" Vinyl Seat	NC
<input type="checkbox"/> MM-01675		18" x 16" Vinyl Seat	NC
<input type="checkbox"/> MM-01676		18" x 18" Vinyl Seat	NC
<input type="checkbox"/> MM-01678		18" x 20" Vinyl Seat	NC
<input type="checkbox"/> MM-01679		20" x 18" Vinyl Seat	NC
<input type="checkbox"/> MM-01680		20" x 20" Vinyl Seat	NC
<input type="checkbox"/> MM-01681		22" x 18" Vinyl Seat	\$ 411.00
<input type="checkbox"/> MM-01682		22" x 20" Vinyl Seat	\$ 411.00

REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes

\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> MM-02792STD		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> MM-02792DP		Deep Adjustable Flat Seat Pan - 20" to 21" deep	\$ 411.00

Deep Seat Pan Extension Kits

<input type="checkbox"/> MM-05690		23" Deep Seat Pan Extension Kit - Requires selection of MM-02792DP	\$ 725.00
<input type="checkbox"/> MM-05741		24" Deep Seat Pan Extension Kit - Requires selection of MM-02792DP	\$ 725.00

REHAB Seat Width

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

SEATING: SEATBACK OPTIONS - SELECT ONE

MPS Backrest Options - Comes standard with Manual Recline

\$

Breathable Cloth MPS Backrest

<input type="checkbox"/> MM-01360		Cloth Backrest - 16-20" wide and 22" tall	NC
<input type="checkbox"/> MM-01520		Cloth Tall Backrest - 16-20" wide and 26" tall	\$ 495.00
<input type="checkbox"/> MM-01361		Cloth Wide Backrest - 22-23" wide and 22" tall	\$ 298.00
<input type="checkbox"/> MM-04598		Cloth Wide and Tall Backrest - 22-23" wide and 26" tall	\$ 590.00

Wipeable Vinyl MPS Backrest

<input type="checkbox"/> MM-01360VI		Vinyl Backrest - 16-20" wide and 22" tall	\$ 270.00
<input type="checkbox"/> MM-01520VI		Vinyl Tall Backrest - 16-20" wide and 26" tall	\$ 708.00
<input type="checkbox"/> MM-01361VI		Vinyl Wide Backrest - 22-23" wide and 22" tall	\$ 359.00
<input type="checkbox"/> MM-04598VI		Vinyl Wide and Tall Backrest - 22-23" wide and 26" tall	\$ 837.00

SEATING: SEATBACK OPTIONS - SELECT ONE

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
MPS Back Accessories			\$
<input type="checkbox"/> MM-04784		Push Handle - Needed for O2 holder, backpack, and cane holder	\$ 141.00
<input type="checkbox"/> MM-02533		Accessory Mount - Needed for aftermarket headrest and/or attendant control	\$ 107.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or plat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
<input type="checkbox"/> J3PA	E2613	JAY J3 PA (Posterior) - 2.2" of lateral thoracic support (up to 20" wide)	\$ 592.00
Select Back Height:	<input type="checkbox"/> 16.5" <input type="checkbox"/> 18"	<input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24"	
<input type="checkbox"/> J3PAW	E2614	JAY J3 PA (Posterior) - 2.2" of lateral thoracic support (up to 22" wide)	\$ 810.00
Select Back Height:	<input type="checkbox"/> 18" <input type="checkbox"/> 21"		
<input type="checkbox"/> J3PL	E2615	JAY J3 PL (Posterior Lateral) - 3" of lateral thoracic support (up to 20" wide)	\$ 683.00
Select Back Height:	<input type="checkbox"/> 16.5" <input type="checkbox"/> 18"	<input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24"	
<input type="checkbox"/> J3PD	E2620	JAY J3 PD (Posterior Deep) - 6" of lateral thoracic support (up to 20" wide)	\$ 811.00
Select Back Height:	<input type="checkbox"/> 16.5" <input type="checkbox"/> 18"	<input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24"	
MPS Backrest Options for use with REHAB Back Canes			\$
<input type="checkbox"/> MM-03996		Cloth MPS Removable Backrest - Works with REHAB back canes	\$ 1,752.00
<input type="checkbox"/> MM-03996VI		Vinyl MPS Removable Backrest - Works with REHAB back canes	\$ 1,752.00
REHAB Flip Forward Backrest Option			\$
<input type="checkbox"/> MM-FLPFWD		REHAB Flip Forward Backrest Option	\$ 263.00
REHAB Seating Cushion Options			\$
<input type="checkbox"/> JION	E2607	JAY Ion	\$ 405.00
<input type="checkbox"/> JU	E2607	JAY Union	\$ 454.00
<input type="checkbox"/> JFUSION	E2622	JAY Fusion	\$ 450.00
Headrest Options - MPS Style			\$
<input type="checkbox"/> MM-02770	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 365.00
<input type="checkbox"/> MM-02770VI	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 365.00
Headrest Options - REHAB Style			\$
<input type="checkbox"/> J3PLSH10	E0955	WHITMYER PLUSH Headrest Pad - Requires J3HEDX (10" x 5")	\$ 164.00
<input type="checkbox"/> J3HEDX	E1028	WHITMYER AXYS Headrest w/ Removable Hardware - Depth adjustable 0.5" backward & 7" forward	\$ 224.00
Armrest Style Options			\$
<input type="checkbox"/> MM-ARMR		Removable Armrests (Post Style) - Attaches to seat pan. Height adjustable 2"	NC
<input type="checkbox"/> MM-ARMFR		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 391.00
<input type="checkbox"/> MM-ARMC		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 785.00
<input type="checkbox"/> MM-ARMC2		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,042.00

SEATING: SEATBACK OPTIONS - SELECT ONE

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Armrest Pad Options			\$
<input type="checkbox"/> MM-ARMPAD01		3" x 2" Desk Length Pads (Pair)	NC
<input type="checkbox"/> MM-ARMPAD02		4" x 16" Tapered Pads (Pair)	\$ 96.00
<input type="checkbox"/> MM-ARMPAD03		3" x 12" Rounded Leatherette Pad (Each)	\$ 258.00
<input type="checkbox"/> MM-ARMPAD04		3" x 14" Rounded Leatherette Pad (Each)	\$ 258.00
Lateral Support Options - REHAB Seat Only			\$
<input type="checkbox"/> J3SAMNTR	E1028	Single Swing-Away Lateral Hardware - Right	\$ 269.00
<input type="checkbox"/> J3SAMNTL	E1028	Single Swing-Away Lateral Hardware - Left	\$ 269.00
<input type="checkbox"/> J3LAT6X4R	E1028	Single 6x4 Lateral Pad - Right	\$ 128.00
<input type="checkbox"/> J3LAT6X4L	E1028	Single 6x4 Lateral Pad - Left	\$ 128.00
Positioning Options			\$
Positioning Hardware - Compatible with REHAB Seat Pan only			
<input type="checkbox"/> MM-FDHIP	E1028	Flip Down Hip Guide Hardware - Fixed to the seat pan near the back	\$ 607.00
<input type="checkbox"/> MM-RMHIP	E1028	Removable Hip Guide Hardware - Fixed to the seat pan near the back	\$ 494.00
<input type="checkbox"/> MM-RMKNE	E1028	Removable Knee Adductors Hardware - Fixed to the seat pan near the back	\$ 494.00
<input type="checkbox"/> MM-FDKNE	E1028	Flip Down Knee Adductors - Fixed to the seat pan near the back	\$ 607.00
Positioning Pad			
<input type="checkbox"/> MM-4X6PAD	E0956	4" x 6" Positioning Pad	\$ 372.00
<input type="checkbox"/> MM-4X10PAD	E0956	4" x 10" Positioning Pad	\$ 372.00
Positioning Pad Accessories			
<input type="checkbox"/> MM-4X6GEL	E0956	4" x 6" Gel Cover	\$ 270.00
Chest Harness Options			\$
<input type="checkbox"/> MM-SH120M-B2	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 151.00
<input type="checkbox"/> MM-SH120L-B2	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 151.00
<input type="checkbox"/> MM-50-CH110	E0960	Healthwares Padded Butterfly Positioning Chest Strap - 7" at chest	\$ 174.00
<input type="checkbox"/> MM-SH331L-B1	E0960	Bodypoint Stayflex Butterfly Harness - 3.5" at chest	\$ 236.00
Pelvic Positioning Options - Select One			\$
Pelvic Positioning			
<input type="checkbox"/> MM-33140	E0978	Therafin Push Button Seat Restraint - 2" x 64"	NC
<input type="checkbox"/> MM-33148	E0978	Therafin Extra Large Push Button Seat Restraint - 2" x 80"	\$ 57.00
<input type="checkbox"/> MM-ABSB	E0978	Therafin Airline Buckle Seat Restraint - 2" x 80"	\$ 146.00
<input type="checkbox"/> MM-EB205-L62	K0108	Bodypoint Evoflex Pelvic Stabilizer - Requires HW312-KIT. Sitting hip width 14-18"	\$ 107.00
<input type="checkbox"/> MM-WC19SB6K	E0978	Crash Tested Seat Belt Kit - 600mm	\$ 938.00
<input type="checkbox"/> MM-WC19SB7K	E0978	Crash Tested Seat Belt Kit - 700mm	\$ 938.00
Pelvic Positioning Options			
<input type="checkbox"/> MM-HW312-KIT		Hardware for Bodypoint Pelvic Stabilizer	\$ 79.00
<input type="checkbox"/> MM-EBE100		Extension for Bodypoint Evoflex Pelvic Stabilizer	\$ 47.00

LEGREST OPTIONS

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Fixed Legrest			\$
Fixed Center Mount Legrest			
<input type="checkbox"/> MM-040541	K0040	V6 - Center Mount Flip-Up Footplate - 1 Piece Footplate	\$ 404.00
<input type="checkbox"/> MM-040542	K0040	V6 - Center Mount Flip-Up Footplate - 2 Piece Footplate	\$ 404.00
Swing-Away Legrest			
<input type="checkbox"/> MM-03528	K0040	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 706.00
Fixed Center Mount Legrest Calf Pad Options			\$
<input type="checkbox"/> MM-05547SM	E0995	Fixed Center Mount Calf Pads - Small - With breathable fabric (13-15")	\$ 511.00
<input type="checkbox"/> MM-05547ME	E0995	Fixed Center Mount Calf Pads - Medium - With breathable fabric (16-19")	\$ 511.00
<input type="checkbox"/> MM-05547LG	E0995	Fixed Center Mount Calf Pads - Large - With breathable fabric (20-23")	\$ 511.00
Fixed Legrest Positioning Strap Options			\$
<input type="checkbox"/> MM-SP102M	K0039	Bodypoint Medium Calf Strap - Side mount only. Length 36" and Width 3.5"	\$ 107.00
<input type="checkbox"/> MM-SP102L	K0039	Bodypoint Large Calf Strap - Side mount only. Length 40" and Width 3.5"	\$ 124.00
<input type="checkbox"/> MM-04597	E0951	Heel and Side Support - Attaches to the footplate stem and footplate	\$ 584.00
<input type="checkbox"/> MM-50-SR230	E0952	Foot Strap - 1.5" wide strap	\$ 79.00
Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
<input type="checkbox"/> MM-05761	E1012	V6 & V4 FWD - Center Mount Power Elevating - 1 Piece Footplate - Includes calf pads (11" x 11.5")	\$ 3,708.00
<input type="checkbox"/> MM-06636	E1012	V6 & V4 FWD - Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads (5" x 10.5")	\$ 3,595.00
Swing-Away Power Elevating Legrest			
<input type="checkbox"/> MM-03526	E1010	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate Must be ordered with MM-SFCK or MM-MFCK depending on number of options selected	\$ 4,045.00
Power Legrest Positioning Options			\$
<input type="checkbox"/> MM-05324	K0042	2" Metal Heel Guards for 2 Piece Footplates Requires MM-06636 or MM-03526. May not fit with all seat sizes	\$ 753.00
<input type="checkbox"/> MM-827A	E0952	Foot Straps for 2 Piece Footplates (Pair) - 11-14"	\$ 79.00
<input type="checkbox"/> MM-HL	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 112.00

ELECTRONICS

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Joystick			\$
<input type="checkbox"/> MM-PG-JSMLEDL		LED R-net Joystick - No LCD display or mono jacks	NC
<input type="checkbox"/> MM-PG-JSMD		Color LCD Joystick - With mono jack ports for external switches	\$ 354.00
<input type="checkbox"/> MM-PG-JSM2IRBT		Infra-Red Color LCD Joystick with Bluetooth	\$ 528.00
Expandable Power Module			\$
Must select both items			
<input type="checkbox"/> MM-PM90	E2377	R-net 90 amp Expandable Controller	\$ 731.00
<input type="checkbox"/> MM-PG-SL	E2313	R-net Harness for Expandable Controller	\$ 511.00
Joystick Handle Type			\$
Comes standard with PG Cone-Shaped Joystick Handle			
<input type="checkbox"/> MM-PC101B	E2323	Bodypoint U-Shaped Joystick Handle - 3"	\$ 169.00
<input type="checkbox"/> MM-PC102B	E2323	Bodypoint U-Shaped Joystick Handle - 4"	\$ 169.00
<input type="checkbox"/> MM-PC110B	E2323	Bodypoint Grooved Mushroom Joystick Handle	\$ 158.00
Joystick Mount			\$
<input type="checkbox"/> MM-03069R		Fixed Joystick Mount - Right - Depth adjustable 8.25"	NC
<input type="checkbox"/> MM-03069L		Fixed Joystick Mount - Left - Depth adjustable 8.25"	NC
<input type="checkbox"/> MM-04607	E1028	Swing-Away Joystick Mount - Right - Includes latch to lock in the forward position	\$ 562.00
<input type="checkbox"/> MM-04608	E1028	Swing-Away Joystick Mount - Left - Includes latch to lock in the forward position	\$ 562.00
<input type="checkbox"/> MM-JMMIDR	E1028	Stealth Midline Joystick Mount - Right - 2 position mount (Rotation 0° - 180°)	\$ 1,360.00
<input type="checkbox"/> MM-JMMIDL	E1028	Stealth Midline Joystick Mount - Left - 2 position mount (Rotation 0° - 180°)	\$ 1,360.00
Joystick Accessories			\$
<input type="checkbox"/> MM-03330AM	E2331	Attendant Control Joystick for MPS Seating with Mounting Kit	\$ 1,882.00
<input type="checkbox"/> MM-03329AM	E2331	Attendant Control Joystick for REHAB Seating with Mounting Kit	\$ 1,882.00
<input type="checkbox"/> 139920		Tash Buddy Button - Outer edge force 3.4 oz. Center force 7 oz. (D:2.5" & H:0.7")	\$ 174.00
<input type="checkbox"/> MM-ESBLK		Egg Switch - Requires CPS9. Single function push button, scroll through modes or profiles	\$ 146.00
<input type="checkbox"/> MM-CPS0-L-HW		Swing-Away Egg Switch Hardware	\$ 449.00
Omni 2			\$
<input type="checkbox"/> MM-05436R		OMNI 2 Display - Right - For alternative drive systems	\$ 2,708.00
<input type="checkbox"/> MM-05436L		OMNI 2 Display - Left - For alternative drive systems	\$ 2,708.00
Alternative Drive Controls			\$
<input type="checkbox"/> MM-HDR-4M-S-C	E2327	Switch-It Head Array - Three switch head array	\$ 5,955.00
<input type="checkbox"/> MM-03335	E2373	Chin Control for MPS Seating System with Mount	\$ 3,764.00
<input type="checkbox"/> MM-03540	E2373	Chin Control for REHAB Seating System with Mount	\$ 3,764.00

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
ACCESSORIES			\$
<input type="checkbox"/> MM-01605		Waterproof Wheelchair Cover	\$ 146.00
<input type="checkbox"/> MM-01606		Waterproof Joystick Cover	\$ 45.00
<input type="checkbox"/> MM-01603		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 34.00
<input type="checkbox"/> MM-001B		Deluxe Backpack - Requires push handle (MM-04784) on MPS seat to attach (13.5" x 15" x 4")	\$ 96.00
<input type="checkbox"/> MM-PP0052		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 158.00
<input type="checkbox"/> MM-06806	E2208	Articulating Oxygen Holder for MPS Seating System	\$ 1,040.00
<input type="checkbox"/> MM-06807	E2208	Articulating Oxygen Holder for REHAB Seating System	\$ 1,040.00
<input type="checkbox"/> MM-05275	E2207	Cane Holder for MPS Seating System	\$ 669.00
<input type="checkbox"/> MM-05276	E2207	Cane Holder for REHAB Seating System	\$ 669.00
<input type="checkbox"/> MM-03197		Luggage Rack for V6 & V4 FWD - Attaches to base behind the chair (10.5" x 13")	\$ 332.00
<input type="checkbox"/> MM-02855		Lighting Kit for V6 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$ 1,180.00
<input type="checkbox"/> MM-03199		Mud Guards for V6 - All terrain drive tires Only (pair). Requires WH145BPB	\$ 669.00
<input type="checkbox"/> MM-04705		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 533.00
<input type="checkbox"/> MM-PP0133		Cell Phone Holder - Adjustable to fit most cell phones	\$ 292.00
<input type="checkbox"/> MM-RAM117BU		Fishing Rod Holder	\$ 292.00
<input type="checkbox"/> MM-DDT		Hand-held Programmer	\$ 1,607.00
Be Adaptive Products			
<input type="checkbox"/> MM-LM100		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 1,489.00
<input type="checkbox"/> MM-TM100		Be Adaptive Trigger Mechanism	\$ 825.00
<input type="checkbox"/> MM-FM200		Be Adaptive Limited Mobility Fishing Mount	\$ 1,489.00