



Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		State	
Buyer				ZIP		Phone	
Mark For				Email / Fax			

Please submit orders to: magicmobility@sunmed.com

Seating Specifications


HCPC	K0848 K0856 K0861 K0850 K0858 K0862	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
	K0849 K0857 K0851 K0859	Seat Type	MPS and Rehab Style
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	STF Height Range	17", 18", 19"
Seat Size Range	Width: 13-23" Depth: 13-24"		

Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 90A controller
User Weight Capacity	400 lbs. 340 lbs. with Elevate or Tilt	Drive Motors	4-pole 6mph
V6 Base Size: All Terrain	Width 28" Length 41.4" Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels 14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3" Length 41.5" Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	4 Frame Colors
		Ground Clearance	3.5"

Client Measurements

- | | | |
|----------------------------|-------|----|
| a) Top of head to seat | _____ | in |
| b) Top of shoulder to seat | _____ | in |
| c) Underarm to seat | _____ | in |
| d) Elbow to seat | _____ | in |
| e) Top of leg to seat | _____ | in |
| f) Shoulder width | _____ | in |
| g) Chest width | _____ | in |
| i) Hip width | _____ | in |
| j) Actual sitting depth | _____ | in |
| m) Foot length | _____ | in |
| n) Knee to foot | _____ | in |
| o) Seat width | _____ | in |



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL		
Seating Options			\$		
☐ 600SF4		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC		
☐ 600SF5		MPS Seat and Back - Captain style seating	NC		
☐ 600SF6		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC		
Base					
Frontier V6 Base - Transit Included		End Item Number: EIPW60	\$		
☐ 600MW67	K0848	Magic Mobility - Frontier V6 No Power Seating Options	\$ 18,587.00		
☐ 600MW68	K0856	Magic Mobility - Frontier V6 Single Power Seating Options	\$ 18,587.00		
☐ 600MW69	K0861	Magic Mobility - Frontier V6 Multi Power Seating Options	\$ 18,587.00		
☐ 600MW70	K0850	Magic Mobility - Frontier V6 No Power Seating Options - HD	\$ 18,587.00		
☐ 600MW71	K0858	Magic Mobility - Frontier V6 Single Power Seating Options - HD	\$ 18,587.00		
☐ 600MW72	K0862	Magic Mobility - Frontier V6 Multi Power Seating Options - HD	\$ 18,587.00		
☐ 600MW73	K0849	Magic Mobility - Frontier V6 No Power Seating Options - MPS (Captain)	\$ 18,587.00		
☐ 600MW74	K0857	Magic Mobility - Frontier V6 Single Power Seating Options - MPS (Captain)	\$ 18,587.00		
☐ 600MW75		Magic Mobility - Frontier V6 Multi Power Seating Options - MPS (Captain)	\$ 18,587.00		
☐ 600MW76	K0851	Magic Mobility - Frontier V6 No Power Seating Options - HD MPS (Captain)	\$ 18,587.00		
☐ 600MW77	K0859	Magic Mobility - Frontier V6 Single Power Seating Options - HD MPS (Captain)	\$ 18,587.00		
☐ 600MW78		Magic Mobility - Frontier V6 Multi Power Seating Options - HD MPS (Captain)	\$ 18,587.00		
Drive Wheels					
☐ 600DW80		All Terrain Drive Wheels and Tires	NC		
☐ 600DW81		Hybrid Drive Wheels and Tires - Pneumatic	NC		
☐ 600DW82		Hybrid Drive Wheels and Tires - Solid	\$ 163.00		
Optional Drive Tire Options					
☐ 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	\$ 40.00		
Caster Wheel Option					
☐ 600CA14	Pneumatic Caster	NC	☐ 600CA15	Solid Foam Filled Caster	\$ 80.00
Base Accessories					
☐ 6008A91		Jack for Drive Tire	\$ 213.00		
Frame Color					
☐ 600C2	Black Frame Color	NC	☐ 600C21	Mossy Green Frame Color	NC
☐ 600C3	Red Frame Color	NC	☐ 600C195	Royal Blue Frame Color	NC
Seat To Floor Height					
\$					
V6 Seat to Floor Heights start at 17" +/- 0.25"					
☐ 600HT44		Lowest Seat To Floor Height	NC		
☐ 600HT45		Seat To Floor Height +1"	NC		
☐ 600HT46		Seat To Floor Height +2"	NC		
Batteries					
☐ 600AG25	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 1,138.00		
☐ 600AG27	E2359	Group 34 Gel Battery (Pair) - Installed	\$ 1,097.00		
Battery Charger					
☐ 600BJ28		10A Battery Charger	NC		
☐ 600BJ38		12A Battery Charger	\$ 495.00		

Seating: Seat Functions

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
No Power Seat Options			\$
☐ 600ST86		No Power Seating System - 400 lb. weight capacity	\$ 4,345.00
☐ 600ST87		No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 5,586.00
Power Tilt			\$

Utilized without Power Recline.

☐ 600PA64	E1002	50° Power Tilt in Space - 340 lb. weight capacity	\$ 8,458.00
Power Seat Elevate			\$

Only available when Power Tilt or Power Tilt and Recline is selected.

☐ 600PI16	E2298	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 4,750.00
Power Recline			\$

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

☐ 600PA65	E1003	150° Power Recline Only - 300 lb. weight capacity	\$ 7,225.00
Power Tilt and Recline			\$
☐ 600PA66	E1006	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 13,709.00

Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)			\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
☐ 600SI24		16" x 18" Vinyl Seat	NC
☐ 600SI25		18" x 18" Vinyl Seat	NC
☐ 600SI26		18" x 20" Vinyl Seat	NC
☐ 600SI27		20" x 18" Vinyl Seat	NC
☐ 600SI28		20" x 20" Vinyl Seat	NC
☐ 600SI29		22" x 20" Vinyl Seat	\$ 486.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

☐ 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 486.00
☐ 600PM15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 858.00

REHAB Seat Width

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options			\$
-------------------------------	--	--	----

Note cushion size will be determined by Rehab Seat size selected above

☐ 600CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 460.00
☐ 600CU60	E2608	JAY® Ion® (Size 22" Wide)	\$ 578.00
☐ 600CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$ 520.00
☐ 600CU62	E2608	JAY® Union® (Size 22" Wide)	\$ 635.00
☐ 600CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 514.00
☐ 600CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$ 613.00

Seating: Seat Sizes			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)			\$
Lateral Thigh Support Pads			
☐ 600KF22	E0953	Lateral Thigh Support Pad Short - Right (5x5)	\$ 147.00
☐ 600KF23	E0953	Lateral Thigh Support Pad Short - Left (5x5)	\$ 147.00
Lateral Thigh Support Mounts			
☐ 600KF26	E1028	Lateral Thigh Removable Mount - Right	\$ 292.00
☐ 600KF27	E1028	Lateral Thigh Removable Mount - Left	\$ 292.00
Seating: Seatback Options			
MPS Backrest Options - Comes standard with Manual Recline			\$
Backrest Width			
☐ 600BN46		Standard Backrest Width - 16-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22-23" Wide	\$ 353.00
Backrest Height			
☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	\$ 586.00
Backrest Material			
☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	\$ 319.00
MPS Back Accessories			\$
☐ 600PH33		Push Handle - Needed for O2 holder	\$ 167.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 127.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
☐ 600JP01	E2313	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 625.00
☐ 600JP02	E2614	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 928.00
☐ 600JP11	E2615	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 729.00
☐ 600JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 875.00
JAY® J3® Back Shell Options and Height Ranges			
☐ 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15		18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR17		21" tall shell (Shoulder Height)	NC
☐ 600JR27		22.5" tall shell (Shoulder Height)	NC
☐ 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			\$
☐ 600D23		REHAB Flip Forward Backrest Option	\$ 311.00

Seating: Seatback Options			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
☐ 600JW35	E0956	Lateral Thoracic (L x H) 4" x 5" Pad - Right	\$ 147.00
☐ 600JW38	E0956	Lateral Thoracic (L x H) 4" x 5" Pad - Left	\$ 147.00
☐ 600JW37	E0956	Lateral Thoracic (L x H) 6" x 4" Pad - Right	\$ 147.00
☐ 600JW40	E0956	Lateral Thoracic (L x H) 6" x 4" Pad - Left	\$ 147.00
JAY® Lateral Thoracic Swing-Away Mount			
☐ 600JW41	E1034	Lateral Thoracic Swing Away Bracket - Right	\$ 292.00
☐ 600JW42	E1034	Lateral Thoracic Swing Away Bracket - Left	\$ 292.00
JAY® Anterior Thoracic Support			\$
☐ 600ZR10	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 158.00
☐ 600ZR11	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 158.00
☐ 600K10	E0960	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 177.00
☐ 600K11	E0960	JAY® Anterior Thoracic Support - Large ATS8303	\$ 177.00
Pelvic Positioning Options - Select One			\$
Pelvic Positioning			
☐ 600A848		Positioning Belt, Auto Style Buckle, 2" x 60"	NC
☐ 600A849	E0978	Positioning Belt, Aircraft Buckle, Padded 2" x 50"	\$ 88.00
☐ 600A850	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 203.00
☐ 600A851	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 224.00
Headrest Options - MPS Style			\$
☐ 600NU12	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 432.00
☐ 600NU13	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 432.00
Headrest Options - REHAB Style			\$
Headrest Pad			
☐ 600HE39	E0955	Whitmyer® Plush Headrest Pad, 10"	\$ 239.00
☐ 600HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 367.00
Headrest Pad			
☐ 600MH33	E1033	Whitmyer® Linx Headrest Mount with Quick Release	\$ 192.00
☐ 600MH34	E1033	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 234.00
Armrest Style Options			\$
☐ 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 462.00
☐ 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 928.00
☐ 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,232.00
Armrest Pad Options			\$
☐ 600AB81		10" Waterfall Armrest Pad - Right	NC
☐ 600AE81		10" Waterfall Armrest Pad - Left	NC
☐ 600AB82		14" Waterfall Armrest Pad - Right	NC
☐ 600AE82		14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
☐ 600H210	K0040	Center Mount (14.5" - 17" Leg Length) and Flip-Up Footplate (11" x 10")	\$ 477.00
Center Mount Fixed Legrest - Accessories			
☐ 600DI59	E0995	Fixed Center Mount Calf Pads	\$ 604.00
☐ 600LS16	E0952	Foot Strap - 1.5" wide strap	\$ 93.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
☐ 600H213	E1012	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 4,384.00
Center Mount Power Elevating Legrest Extension			
☐ 600E42		Center Mount Power Elevating Extension - Standard (15" - 17" Leg Length)	NC
☐ 600E43		Center Mount Power Elevating Extension - Short (10" - 15" Leg Length)	NC
Center Mount Power Elevating Legrest Footplates			
☐ 600F117		Center Mount Power Elevating Footplate - Standard (5" x 12" each)	NC
☐ 600F118		Center Mount Power Elevating Footplate - Short (4.5" x 7" each)	NC
Swing-Away Legrest			\$
Swing-Away Legrest			
☐ 600H211	K0040	Swing-Away Legrest (12" - 17" Leg Length) - flip-up footplate (6" x 8.5")	\$ 835.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
☐ 600H214	E1010	Swing-Away Power Elevating Legrest (12" - 17" Leg Length), flip-up height adj footplate (6" x 8.5")	\$ 4,783.00
S/A Power Legrest Positioning Options			
☐ 600HL15	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 132.00

Accessories			\$
☐ 6008A93		Waterproof Wheelchair Cover	\$ 173.00
☐ 6008A94		Waterproof Joystick Cover	\$ 54.00
☐ 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 40.00
☐ 600A221		Deluxe Backpack	\$ 113.00
☐ 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 187.00
☐ 600A719	E2208	Articulating Oxygen Holder	\$ 1,230.00
☐ 600A612	E2207	Cane Holder - Includes push handles	\$ 791.00
☐ 6009A15		Luggage Rack for V6 - Attaches to base behind the chair (10.5" x 13")	\$ 393.00
☐ 6009A17		Lighting Kit for V6 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$ 1,395.00
☐ 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 630.00
☐ 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 345.00
☐ 600A614		Fishing Rod Holder	\$ 345.00

Be Adaptive Products			
☐ 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 1,760.00
☐ 6008A99		Be Adaptive Trigger Mechanism	\$ 975.00
☐ 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 1,760.00

Electronics									
✓ OPTION		HCPCS		DESCRIPTION				RETAIL	
Joystick									\$
☐	600J121			Standard LED Joystick (w/o Screen Display)				NC	
☐	600J116			Color Joystick (2.3" Color LCD Screen Display, CJSM1)				NC	
☐	600J118			High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display				\$ 1,230.00	
Joystick Mount									\$
☐	600JM34			Fixed on Armrest - Depth adjustable 8.25"				NC	
☐	600JM45	E1032		Swing-Away Joystick Mount - Includes latch to lock in the forward position				\$ 665.00	
☐	600JM46	E1032		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)				\$ 1,608.00	
Joystick Mount - Location									\$
☐	600HM16			Right Hand Armrest				NC	
☐	600HM17			Left Hand Armrest				NC	
Display Language Selection R-net									\$
☐	600L1	English	NC	☐	600L3	French	NC	☐	600L7 Spanish NC
Joystick Handle Type									\$
Comes standard with PG Cone-Shaped Joystick Handle									
☐	600JH13	E2323		Ball Handle				\$ 147.00	
☐	600JH26	E2323		Bodypoint Grooved Mushroom Joystick Handle				\$ 191.00	
☐	600JH24	E2323		Bodypoint – Rubber Dome Joystick Handle				\$ 191.00	
☐	600JH25	E2323		Bodypoint – U-Shaped Joystick Handle (4 inch)				\$ 191.00	
Expandable Power Module									\$
Must select both items									
☐	600EX15	E2377		Expandable 90A Controller - R-net				\$ 786.00	
☐	600EX16	E2313		Harness for 90A Expandable Controller				\$ 710.00	
Power Seat Function									\$
☐	600EO37	E2310		Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function				\$ 2,306.00	
☐	600EO38	E2311		Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions				\$ 4,174.00	
Power Seat Control Option									\$
☐	600QG11			Ctrl+5 Switch Box with Voice Activation				\$ 1,459.00	
☐	600QG12			Ctrl+5 Switch Box with Voice Activation (Multi-power only)				NC	
Ctrl+5 Switch Type									\$
☐	600EY22			Push Button				NC	
☐	600EY23	K0108		Toggle Switches				\$ 350.00	
Power Seat Control Option - Ctrl+5 Mounting location									\$
☐	600EY19	Right Armrest		NC	☐	600EY20	Left Armrest		NC
Attendant Controls									\$
☐	6000A05	E2331		Dual Control Module for Attendant - Includes accessory mount and push handle				\$ 1,561.00	
Attendant Control Mounting Location									\$
☐	600DN5	Right Back Cane		NC	☐	600DN6	Left Back Cane		NC
R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth									\$
☐	600ON14			R-net Omni2 Only Input Control Module				\$ 2,072.00	
☐	600ON15			R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices				NC	
R-net Omni2 -Mounting Location									\$
☐	600DS23	Right Hand Armrest		NC	☐	600DS24	Left Hand Armrest		NC
Head Control Packages									\$
☐	600EU25	E2328		Switch-It DualPro® Proportional Head Array - Facial Pads				\$ 6,668.00	
☐	600EU26	E2328		Switch-It DualPro® Proportional Head Array - Spot Pads				\$ 6,668.00	