

MWD power wheelchair



Customer				Ship To			
Account				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		State	
Buyer				ZIP	Phone		
Mark For				Email / Fax			

Please submit orders to: magicmobility@sunmed.com

Seating Specifications


HPCP	K0848 K0856 K0861 K0850 K0858 K0862 K0849 K0857 K0851 K0859	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	17", 18", 19"

Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 90A controller
User Weight Capacity	400 lbs. 340 lbs. with Elevate or Tilt	Drive Motors	4-pole 6mph
V6 Base Size: All Terrain	Width 28" Length 41.4" Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels 14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3" Length 41.5" Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	4 Frame Colors
		Ground Clearance	3.5"

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL		
Seating Options			\$		
<input type="checkbox"/> 600SF4		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC		
<input type="checkbox"/> 600SF5		MPS Seat and Back - Captain style seating	NC		
<input type="checkbox"/> 600SF6		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC		
Base					
Frontier V6 Base		End Item Number: EIPW60	\$		
<input type="checkbox"/> 600MW67	K0848	Magic Mobility - Frontier V6 No Power Seating Options	\$ 17,702.00		
<input type="checkbox"/> 600MW68	K0856	Magic Mobility - Frontier V6 Single Power Seating Options	\$ 17,702.00		
<input type="checkbox"/> 600MW69	K0861	Magic Mobility - Frontier V6 Multi Power Seating Options	\$ 17,702.00		
<input type="checkbox"/> 600MW70	K0850	Magic Mobility - Frontier V6 No Power Seating Options - HD	\$ 17,702.00		
<input type="checkbox"/> 600MW71	K0858	Magic Mobility - Frontier V6 Single Power Seating Options - HD	\$ 17,702.00		
<input type="checkbox"/> 600MW72	K0862	Magic Mobility - Frontier V6 Multi Power Seating Options - HD	\$ 17,702.00		
<input type="checkbox"/> 600MW73	K0849	Magic Mobility - Frontier V6 No Power Seating Options - MPS (Captain)	\$ 17,702.00		
<input type="checkbox"/> 600MW74	K0857	Magic Mobility - Frontier V6 Single Power Seating Options - MPS (Captain)	\$ 17,702.00		
<input type="checkbox"/> 600MW75		Magic Mobility - Frontier V6 Multi Power Seating Options - MPS (Captain)	\$ 17,702.00		
<input type="checkbox"/> 600MW76	K0851	Magic Mobility - Frontier V6 No Power Seating Options - HD MPS (Captain)	\$ 17,702.00		
<input type="checkbox"/> 600MW77	K0859	Magic Mobility - Frontier V6 Single Power Seating Options - HD MPS (Captain)	\$ 17,702.00		
<input type="checkbox"/> 600MW78		Magic Mobility - Frontier V6 Multi Power Seating Options - HD MPS (Captain)	\$ 17,702.00		
Drive Wheels			\$		
<input type="checkbox"/> 600DW80		All Terrain Drive Wheels and Tires	NC		
<input type="checkbox"/> 600DW81		Hybrid Drive Wheels and Tires - Pneumatic	NC		
<input type="checkbox"/> 600DW82		Hybrid Drive Wheels and Tires - Solid	\$ 155.00		
Optional Drive Tire Options			\$		
<input type="checkbox"/> 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	\$ 38.00		
Caster Wheel Option			\$		
<input type="checkbox"/> 600CA14	Pneumatic Caster	NC	<input type="checkbox"/> 600CA15	Solid Foam Filled Caster	\$ 76.00
Base Accessories			\$		
<input type="checkbox"/> 6008A91		Jack for Drive Tire	\$ 203.00		
Frame Color			\$		
<input type="checkbox"/> 600C2	Black Frame Color	NC	<input type="checkbox"/> 600C21	Mossy Green Frame Color	NC
<input type="checkbox"/> 600C3	Red Frame Color	NC	<input type="checkbox"/> 600C195	Royal Blue Frame Color	NC
Seat To Floor Height			\$		
V6 Seat to Floor Heights start at 17" +/- 0.25"					
<input type="checkbox"/> 600HT44		Lowest Seat To Floor Height	NC		
<input type="checkbox"/> 600HT45		Seat To Floor Height +1"	NC		
<input type="checkbox"/> 600HT46		Seat To Floor Height +2"	NC		
Batteries			\$		
<input type="checkbox"/> 600AG25	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 1,116.00		
<input type="checkbox"/> 600AG27	E2359	Group 34 Gel Battery (Pair) - Installed	\$ 1,075.00		
Battery Charger			\$		
<input type="checkbox"/> 600BJ28		10A Battery Charger	NC		
Power Retractable Docking Pin			\$		
<input type="checkbox"/> 6009A14		Retractable Docking Pin	\$ 3,071.00		

Seating: Seat Functions

<input checked="" type="checkbox"/> OPTION	HCPCS	DESCRIPTION	RETAIL
No Power Seat Options			\$
<input type="checkbox"/> 600ST86		No Power Seating System - 400 lb. weight capacity	\$ 4,138.00
<input type="checkbox"/> 600ST86		No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 5,320.00
Power Tilt			\$

Utilized without Power Recline.

<input type="checkbox"/> 600PA64	E1002	50° Power Tilt in Space - 340 lb. weight capacity	\$ 8,292.00
Power Seat Elevate			\$

Only available when Power Tilt or Power Tilt and Recline is selected.

<input type="checkbox"/> 600PI16	E2300	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 4,750.00
Power Recline			\$

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

<input type="checkbox"/> 600PA65	E1003	150° Power Recline Only - 300 lb. weight capacity	\$ 6,750.00
Power Tilt and Recline			\$
<input type="checkbox"/> 600PA66	E1006	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 12,750.00

Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)			\$
---	--	--	----

MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.

<input type="checkbox"/> 600SI24		16" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI25		18" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI26		18" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI27		20" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI28		20" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI29		22" x 20" Vinyl Seat	\$ 463.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 463.00
<input type="checkbox"/> 600PM15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 817.00

REHAB Seat Width

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options			\$
--------------------------------------	--	--	----

Note cushion size will be determined by Rehab Seat size selected above

<input type="checkbox"/> 600CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 451.00
<input type="checkbox"/> 600CU60	E2608	JAY® Ion® (Size 22" Wide)	\$ 567.00
<input type="checkbox"/> 600CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$ 510.00
<input type="checkbox"/> 600CU62	E2608	JAY® Union® (Size 22" Wide)	\$ 623.00
<input type="checkbox"/> 600CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 504.00
<input type="checkbox"/> 600CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$ 601.00

Lateral Thigh Support Pads and Hardware (must pick pad and mount) \$

Lateral Thigh Support Pads

<input type="checkbox"/> 600KF22	E0953	Lateral Thigh Support Pad Short - Right (5"x5")	\$ 144.00
<input type="checkbox"/> 600KF23	E0953	Lateral Thigh Support Pad Short - Left (5"x5")	\$ 144.00

Lateral Thigh Support Mounts

<input type="checkbox"/> 600KF26	E1028	Lateral Thigh Removable Mount - Right	\$ 286.00
<input type="checkbox"/> 600KF27	E1028	Lateral Thigh Removable Mount - Left	\$ 286.00

Seating: Seatback Options

MPS Backrest Options - Comes standard with Manual Recline \$

Backrest Width

<input type="checkbox"/> 600BN46		Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47		Wide Backrest Width - 22-23" Wide	\$ 336.00

Backrest Height

<input type="checkbox"/> 600B80		Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81		Tall Backrest Height - 26" Tall	\$ 558.00

Backrest Material

<input type="checkbox"/> 600BU106		Cloth Backrest	NC
<input type="checkbox"/> 600BU107		Vinyl Backrest	\$ 304.00

MPS Back Accessories \$

<input type="checkbox"/> 600PH33		Push Handle - Needed for O2 holder	\$ 159.00
<input type="checkbox"/> 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 121.00

REHAB Seating Backrest Options \$

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

JAY J3 Back Contour

<input type="checkbox"/> 600JP01	E2313	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 613.00
<input type="checkbox"/> 600JP02	E2614	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 910.00
<input type="checkbox"/> 600JP11	E2615	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 715.00
<input type="checkbox"/> 600JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 858.00

JAY® J3® Back Shell Options and Height Ranges

<input type="checkbox"/> 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15		18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR17		21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27		22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28		24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option \$

<input type="checkbox"/> 600D23		REHAB Flip Forward Backrest Option	\$ 296.00
---------------------------------	--	------------------------------------	-----------

Seating: Seatback Options

<input checked="" type="checkbox"/> OPTION	HPCPS	DESCRIPTION	RETAIL
--	-------	-------------	--------

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only \$

JAY® Lateral Thoracic Support Pads

<input type="checkbox"/> 600JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 144.00
<input type="checkbox"/> 600JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 144.00
<input type="checkbox"/> 600JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 144.00
<input type="checkbox"/> 600JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 144.00

JAY® Lateral Thoracic Swing-Away Mount

<input type="checkbox"/> 600JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$ 286.00
<input type="checkbox"/> 600JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$ 286.00

JAY® Anterior Thoracic Support				\$
<input type="checkbox"/>	600ZR10	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 155.00
<input type="checkbox"/>	600ZR11	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 155.00
<input type="checkbox"/>	600K10	E0960	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 174.00
<input type="checkbox"/>	600K11	E0960	JAY® Anterior Thoracic Support - Large ATS8303	\$ 174.00
Pelvic Positioning Options - Select One				\$
Pelvic Positioning				
<input type="checkbox"/>	600A848		Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
<input type="checkbox"/>	600A849	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 86.00
<input type="checkbox"/>	600A850	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 199.00
<input type="checkbox"/>	600A851	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 220.00
Headrest Options - MPS Style				\$
<input type="checkbox"/>	600NU12	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 411.00
<input type="checkbox"/>	600NU13	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 411.00
Headrest Options - REHAB Style				\$
Headrest Pad				
<input type="checkbox"/>	600HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$ 234.00
<input type="checkbox"/>	600HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 360.00
Headrest Pad				
<input type="checkbox"/>	600MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$ 188.00
<input type="checkbox"/>	600MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 229.00
Armrest Style Options				\$
<input type="checkbox"/>	600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 440.00
<input type="checkbox"/>	600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 884.00
<input type="checkbox"/>	600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,173.00
Armrest Pad Options				\$
<input type="checkbox"/>	600AB80		4" x 16" Tapered Pads - Right	\$ 54.00
<input type="checkbox"/>	600AE80		4" x 16" Tapered Pads - Left	\$ 54.00
<input type="checkbox"/>	600AB81		10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE81		10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/>	600AB82		14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE82		14" Waterfall Armrest Pad - Left	NC
Legrest Options				
<input checked="" type="checkbox"/>	OPTION	HCPCS	DESCRIPTION	RETAIL
Center Mount Fixed Legrest				\$
Center Mount Fixed Legrest				
<input type="checkbox"/>	600H210	K0040	Center Mount Flip-Up Footplate	\$ 454.00
Center Mount Fixed Legrest - Calfpads				
<input type="checkbox"/>	600DI59	E0995	Fixed Center Mount Calf Pads	\$ 575.00
Center Mount Power Elevating Legrest				\$
Center Mount Power Elevating Legrest				
<input type="checkbox"/>	600H213	E1012	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 4,175.00
Center Mount Power Elevating Legrest Extension				
<input type="checkbox"/>	600E42		Center Mount Power Elevating Extension - Standard	NC
<input type="checkbox"/>	600E43		Center Mount Power Elevating Extension - Short	NC
Center Mount Power Elevating Legrest Footplates				
<input type="checkbox"/>	600F117		Center Mount Power Elevating Footplate - Standard	NC
<input type="checkbox"/>	600F118		Center Mount Power Elevating Footplate - Short	NC

Swing-Away Legrest	\$
--------------------	----

Swing-Away Legrest

<input type="checkbox"/> 600H211	K0040	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$	795.00
----------------------------------	-------	--	----	--------

S/A Non-Power Positioning Options

<input type="checkbox"/> 600LS15	E0951	Heel and Side Support - Attaches to the footplate stem and footplate	\$	657.00
<input type="checkbox"/> 600LS16	E0952	Foot Strap - 1.5" wide strap	\$	89.00

Swing-Away Power Elevating Legrest	\$
------------------------------------	----

Swing-Away Power Elevating Legrest

<input type="checkbox"/> 600H214	E1010	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$	4,555.00
----------------------------------	-------	--	----	----------

S/A Power Legrest Positioning Options

<input type="checkbox"/> 600HL15	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$	126.00
----------------------------------	-------	--	----	--------

Accessories	\$
-------------	----

<input type="checkbox"/> 6008A93		Waterproof Wheelchair Cover	\$	165.00
<input type="checkbox"/> 6008A94		Waterproof Joystick Cover	\$	51.00
<input type="checkbox"/> 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$	38.00
<input type="checkbox"/> 600A221		Deluxe Backpack	\$	108.00
<input type="checkbox"/> 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$	178.00
<input type="checkbox"/> 600A719	E2208	Articulating Oxygen Holder	\$	1,171.00
<input type="checkbox"/> 600A612	E2207	Cane Holder - Includes push handles	\$	753.00
<input type="checkbox"/> 6009A15		Luggage Rack for V6 - Attaches to base behind the chair (10.5" x 13")	\$	374.00
<input type="checkbox"/> 6009A17		Lighting Kit for V6 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$	1,329.00
<input type="checkbox"/> 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$	600.00
<input type="checkbox"/> 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$	329.00
<input type="checkbox"/> 600A614		Fishing Rod Holder	\$	329.00

Be Adaptive Products

<input type="checkbox"/> 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$	1,676.00
<input type="checkbox"/> 6008A99		Be Adaptive Trigger Mechanism	\$	929.00
<input type="checkbox"/> 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$	1,676.00

Electronics	
-------------	--

<input checked="" type="checkbox"/> OPTION	HCPCS	DESCRIPTION		RETAIL
--	-------	-------------	--	--------

Joystick	\$
----------	----

<input type="checkbox"/> 600J121		Standard LED Joystick (w/o Screen Display)		NC
<input type="checkbox"/> 600J116		Color Joystick (2.3" Color LCD Screen Display, CJSM1)		NC
<input type="checkbox"/> 600J118		High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$	1,206.00

Joystick Mount	\$
----------------	----

<input type="checkbox"/> 600JM34		Fixed on Armrest - Depth adjustable 8.25"		NC
<input type="checkbox"/> 600JM45	E1028	Swing-Away Joystick Mount - Includes latch to lock in the forward position	\$	633.00
<input type="checkbox"/> 600JM46	E1028	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)	\$	1,531.00

Joystick Mount - Location	\$
---------------------------	----

<input type="checkbox"/> 600HM16		Right Hand Armrest		NC
<input type="checkbox"/> 600HM17		Left Hand Armrest		NC

Display Language Selection R-net	\$
----------------------------------	----

<input type="checkbox"/> 600L1	English	NC	<input type="checkbox"/> 600L3	French	NC	<input type="checkbox"/> 600L7	Spanish		NC
--------------------------------	---------	----	--------------------------------	--------	----	--------------------------------	---------	--	----

Joystick Handle Type	\$
----------------------	----

Comes standard with PG Cone-Shaped Joystick Handle

<input type="checkbox"/> 600JH13	E2323	Ball Handle	\$	144.00
<input type="checkbox"/> 600JH26	E2323	Bodypoint Grooved Mushroom Joystick Handle	\$	187.00
<input type="checkbox"/> 600JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	\$	187.00
<input type="checkbox"/> 600JH25	E2323	Bodypoint – U-Shaped Joystick Handle (4 inch)	\$	187.00

Expandable Power Module					\$	
Must select both items						
<input type="checkbox"/>	600EX15	E2377	Expandable 90A Controller - R-net		\$ 749.00	
<input type="checkbox"/>	600EX16	E2313	Harness for 90A Expandable Controller		\$ 676.00	
Power Seat Function					\$	
<input type="checkbox"/>	600EO37	E2310	THRU DRIVE CNTRL 1 ACTUATOR- Must be selected if ordering 1 actuator function		\$ 2,196.00	
<input type="checkbox"/>	600EO38	E2311	THRU DRIVE CNTRL 2+ ACTUATORS- Must be selected if ordering 2 or more actuator functions		\$ 3,975.00	
Power Seat Control Option					\$	
<input type="checkbox"/>	600QG11		Ctrl+5 Switch Box with Free Assignable Functions		\$ 1,430.00	
<input type="checkbox"/>	600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)		NC	
Ctrl+5 Switch Type					\$	
<input type="checkbox"/>	600EY22		Push Button		NC	
<input type="checkbox"/>	600EY23	K0108	Toggle Switches		\$ 343.00	
Power Seat Control Option - Mounting location					\$	
<input type="checkbox"/>	600EY19	CTRL+5 MNT ARMREST RIGHT	NC	<input type="checkbox"/>	600EY20 CTRL+5 MNT ARMREST LEFT	NC
Attendant Controls					\$	
<input type="checkbox"/>	6000A05	E2331	Dual Control Module for Attendant - Includes accessory mount and push handle		\$ 1,487.00	
Attendant Control Mounting Location					\$	
<input type="checkbox"/>	600DN5	Right Back Cane	NC	<input type="checkbox"/>	600DN6 Left Back Cane	NC
R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth					\$	
<input type="checkbox"/>	600ON14		R-net Omni-2 Only Input Control Module		\$ 2,031.00	
<input type="checkbox"/>	600ON15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices		NC	
R-net Omni-2 -Mounting Location					\$	
<input type="checkbox"/>	600DS23	Right Hand Armrest	NC	<input type="checkbox"/>	600DS24 Left Hand Armrest	NC
Head Control Packages					\$	
<input type="checkbox"/>	600EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads		\$ 6,537.00	
<input type="checkbox"/>	600EU26	E2328	Switch-It DualPro® Proportional Head Array - Spot Pads		\$ 6,537.00	