



Totstander

February 2023

Mark For: _____

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax EmailSubmitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Totstander
HCPCS: E0638
Totstander

Selection required. Includes Adjustable Chest and Knee Supports, Lap Belt and Sandals.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LFS/T-MW	Totstander - Meadow		\$1,244

Technical Specifications - Imperial

Specification	Size 1
Age (approximate)	1 - 5 years
Maximum User Weight	48 lbs
Chest Support Height (Min / Max)	18.5" / 26"
Base	18" x 18"

Technical Specifications - Metric

Specification	Size 1
Age (approximate)	1 - 5 years
Maximum User Weight	22 kg
Chest Support Height (Min / Max)	470 mm / 660.5 mm
Base	450 mm x 450 mm



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA