



# Everyday Activity Seat January 2022

**Mark For:**

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via:  Fax  Email

**Submitting for:**  Quote  Order

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**Everyday Activity Seat****HCPCS: E1399**

*Size 1 Weight Capacity: 66 lbs | Size 2 Weight Capacity: 110 lbs | Size 3 Weight Capacity: 220 lbs*

**Base Frame**

*Selection Required.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-6102	Size 1 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	E1399	\$2,901
<input type="checkbox"/>	151-6202	Size 2 Seat Base, Backrest and Hi-Low Chassis with Gas Spring	E1399	\$3,011
<input type="checkbox"/>	151-6302	Size 3 Seat Base, Backrest and Hi-Low Hydraulic Assist Chassis	E1399	\$5,581

**Chair Cover**

*Selection Required. Must match base frame size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-6104-03VIN	Size 1 Vinyl Seat and Back Covers - Blue		\$348
<input type="checkbox"/>	151-6206-03VIN	Size 2 Vinyl Seat and Back Covers - Blue		\$393
<input type="checkbox"/>	151-6304-03VIN	Size 3 Vinyl Seat and Back Covers - Blue		\$653

**Lateral Support**

*Laterals must match base frame size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-6106-07VIN	Size 1 Rigid Chest Laterals - Black Vinyl (Pair)		\$365
<input type="checkbox"/>	151-6210-07VIN	Size 2 Rigid Chest Laterals - Black Vinyl (Pair)		\$365
<input type="checkbox"/>	151-6306-07VIN	Size 3 Rigid Chest Laterals - Black Vinyl (Pair)		\$498
<input type="checkbox"/>	151-6108-07VIN	Size 1 Swing-Away Chest Laterals - Black Vinyl (Pair)		\$454
<input type="checkbox"/>	151-6212-07VIN	Size 2 Swing-Away Chest Laterals - Black Vinyl (Pair)		\$459
<input type="checkbox"/>	151-6308-07VIN	Size 3 Swing-Away Chest Laterals - Black Vinyl (Pair)		\$602
<input type="checkbox"/>	151-6218-07VIN	Size 1 and 2 Seat Protraction Pads - Black Vinyl (Pair)		\$442
<input type="checkbox"/>	151-6314-07VIN	Size 3 Seat Protraction Pads - Black Vinyl (Pair)		\$476
<input type="checkbox"/>	137-808-05	Chest Support - Size 1 & 2 <i>Requires rigid or swing-away laterals</i>		\$174
<input type="checkbox"/>	133-625	Chest Support - Size 3 <i>Requires rigid or swing-away laterals</i>		\$174

## Upper Extremity Positioning

*Tray requires Armrests. Selections must match base frame size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-6117	Size 1 Armrests (Pair)		\$379
<input type="checkbox"/>	151-6226	Size 2 Armrests (Pair)		\$379
<input type="checkbox"/>	151-3747	Size 3 Armrests (Pair)		\$533
<input type="checkbox"/>	151-1644	Size 1 Tray - Plastic Gray		\$379
<input type="checkbox"/>	151-2644	Size 2 Tray - Plastic Gray		\$379
<input type="checkbox"/>	151-3767	Size 3 Tray - Plastic Black		\$474

## Head Supports

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-6208-07VIN	Size 1 and 2 Flat Headrest Support and Cushion - Black Vinyl <i>Must match base frame size</i>		\$354
<input type="checkbox"/>	151-6209-07VIN	Size 1 and 2 Laterals for Flat Headrest - Black Vinyl (Pair) <i>Must match base frame size. Includes width adjustable lateral pads &amp; covers. Requires Flat Headrest Support and Cushion.</i>		\$299
<input type="checkbox"/>	151-3751-05	Size 3 Contoured Headrest - Black Neoprene <i>Must match base frame size</i>		\$461
<input type="checkbox"/>	151-1756	Size 1 and 2 Contoured Headrest Adapter Bracket		\$79
<input type="checkbox"/>	120-859	Size 1 and 2 Contoured Headrest Support Hardware		\$358
<input type="checkbox"/>	137-800-05	Size 1 and 2 Contoured Headrest Cushion - Black		\$164

## Foot Positioning

*Footplate selection required with Sandals*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-1745	Size 1 One Piece Footplate - Metal		\$379
<input type="checkbox"/>	151-2745	Size 2 One Piece Footplate - Metal		\$379
<input type="checkbox"/>	151-3745	Size 3 One Piece Footplate - Metal		\$644
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)		\$254
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)		\$277
<input type="checkbox"/>	152-3600	Large Sandals w/ Straps (Pair)		\$299
<input type="checkbox"/>	151-1760	Size 1 and 2 Sandal Raiser 1" (Each) <i>Requires Sandal Selection, must match Sandal size.</i>	Quantity: (1-6)	\$56
<input type="checkbox"/>	151-3760	Size 3 Sandal Raiser 1" (Each) <i>Requires Sandal Selection, must match Sandal size.</i>	Quantity: (1-6)	\$64

## Pelvic Fixation

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	151-1660	Size 1 Pelvic Belt 4-Point		\$139
<input type="checkbox"/>	151-2660	Size 2 Pelvic Belt 4-Point		\$149
<input type="checkbox"/>	151-3660	Size 3 Pelvic Belt 4-Point		\$160
<input type="checkbox"/>	151-6115-07VIN	Size 1 Side Pads - Black Vinyl (Pair)		\$492
<input type="checkbox"/>	151-6222-07VIN	Size 2 Side Pads - Black Vinyl (Pair)		\$586
<input type="checkbox"/>	151-6317-07VIN	Size 3 Side Pads - Black Vinyl (Pair)		\$751
<input type="checkbox"/>	151-6224-07VIN	Size 2 Adduction Side Pad - Black Vinyl		\$282
<input type="checkbox"/>	151-6320-07VIN	Size 3 Adduction Side Pad - Black Vinyl		\$282
<input type="checkbox"/>	151-6120-07VIN	Size 1 Wedge Shape Pommel - Black Vinyl		\$199
<input type="checkbox"/>	151-6229-07VIN	Size 2 Wedge Shape Pommel - Black Vinyl		\$238
<input type="checkbox"/>	151-6334-07VIN	Size 3 Wedge Shape Pommel - Black Vinyl		\$248

## Assembly

	Item Code	Description	Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee	\$126

**Technical Specifications - Imperial**

Specification	Size 2	Size 3
Max Load	110 lbs	220 lbs
Weight (Seat and Hi-Low Chassis)	75 lbs	108 lbs
Seat Height (Hi-Low Chassis)	12.6" - 23.6"	17.7" - 23.6"
Seat Width	10.2" - 15"	13.4" - 19.7"
Seat Depth	9" - 15"	13" - 18.9"
Lateral Distance (Chest)	6.7" - 9.4"	9.4" - 14.6"
Back Height	15.7" - 20.5"	17.7" - 24.8"
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	7.9" - 11.8"	7" - 12.7"
Lower Leg Length	9.8" - 17.7"	9.8" - 17.7"
Lower Leg Angle	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone

**Technical Specifications - Metric**

Specification	Size 2	Size 3
Max Load	50 kgs	100 kgs
Weight (Seat and Hi-Low Chassis)	34 kgs	49 kgs
Seat Height (Hi-Low Chassis)	320 mm - 600 mm	450 mm - 600 mm
Seat Width	260 mm - 380 mm	340 mm - 500 mm
Seat Depth	230 mm - 380 mm	330 mm - 480 mm
Lateral Distance (Chest)	170 mm - 240 mm	240 mm - 370 mm
Back Height	400 mm - 520 mm	450 mm - 630 mm
Back Angle	30° Supine / 10° Prone	30° Supine / 10° Prone
Armrest Height	200 mm - 300 mm	180 mm - 325 mm
Lower Leg Length	250 mm - 450 mm	250 mm - 450 mm
Lower Leg Angle	-10° to +50°	-10° to +50°
Footplate Angle	-15° to +10°	-15° to +10°
Tilt-In-Space	30° Supine / 10° Prone	30° Supine / 10° Prone



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA