



ORDER FORM

Pronestander

Prone for success with this angle adjustable stander.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Pronestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
----------	-------------	-------	--------

Pronestander			\$
--------------	--	--	----

Includes Flexible Laterals, PU Knee Supports, Footplate, Tray, Tray Cover and Activity Bowl.

Must select pelvic band and sandal size below.

<input type="checkbox"/>	LPS/1-MW	Size 1 Pronestander - Meadow	E0638	\$3,241
<input type="checkbox"/>	LPS/2-MB	Size 2 Pronestander - Mondeo Blue	E0638	\$3,543
<input type="checkbox"/>	LPS/3-MB	Size 3 Pronestander - Mondeo Blue	E0638	\$4,149

Pelvic Band			\$	
-------------	--	--	----	--

<input type="checkbox"/>	LPS/1/B-MW	Size 1 Pelvic Band - Meadow		\$194
<input type="checkbox"/>	LPS/2/B-MB	Size 2 Pelvic Band - Mondeo Blue		\$213
<input type="checkbox"/>	LPS/3/B-MB	Size 3 Pelvic Band - Mondeo Blue		\$231

Additional Knee Supports			\$	
--------------------------	--	--	----	--

<input type="checkbox"/>	LPS/1/K-MW	Size 1 Upholstered Knee Cups - Meadow [Pair]		\$132
<input type="checkbox"/>	LPS/2/K-MB	Size 2 Upholstered Knee Cups - Mondeo Blue [Pair]		\$132
<input type="checkbox"/>	LPS/3/K-MB	Size 3 Upholstered Knee Cups - Mondeo Blue [Pair]		\$132
<input type="checkbox"/>	LPS/4/K-MB	Size 4 Upholstered Knee Cups - Mondeo Blue [Pair]		\$132

Foot Positioning			\$	
------------------	--	--	----	--

<input type="checkbox"/>	152-1600	Small Sandals with Straps [Pair]		\$273
<input type="checkbox"/>	152-2600	Medium Sandals with Straps [Pair]		\$298
<input type="checkbox"/>	152-3600	Large Sandals with Straps [Pair]		\$321
<input type="checkbox"/>	117-787	Sandal Raiser 1" [Each] Requires Sandal Selection.	QTY	\$95

Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55" / 73"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Plate Height	18.5" - 27"	27" - 38"	38" - 53"
Base	32" x 24"	42" x 26"	50" x 28"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	800 mm / 1150 mm	1100 mm / 1450 mm	1400 mm / 1850 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Plate Height	470 mm - 680 mm	680 mm - 980 mm	970 mm - 1140 mm
Base	820 mm x 600 mm	1070 mm x 650 mm	1270 mm x 700 mm
Activity Tray	400 mm x 525 mm	450 mm x 575 mm	500 mm x 600 mm

Sunrise Medical (US) LLC

2842 Business Park Avenue
 Fresno, CA 93727 USA
 Tel.: 1-800-333-4000
 Fax: 1-800-300-7502
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129853 Rev. E ©01/2024



800.333.4000

www.SunriseMedical.com