



# JAY® Protector

January 2022

**Mark For:** \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Dealer Acct #: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Confirm Via:  Fax  Email

**Submitting for:**  Quote  Order  
 PO#: \_\_\_\_\_

### ADDITIONAL SHIPPING INFORMATION

Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

## JAY® PROTECTOR

- 140 JAY Protector - Adult**  
 K0108 10" - 18" Seat Widths \$279
- 143-MJ JAY Protector - Junior**  
 K0108 8" - 12" Seat Widths \$253

### Specifications

	Adult	Junior
Width	15"	11"
Depth	9"	7"
Height	2"	1.5"
Accommodates Seat Widths	10" - 18"	8" - 12"
Insert	JAY Flow™ Fluid Pad	JAY Flow™ Fluid Pad
Outer Cover	Moisture Resistant	Moisture Resistant
Product Weight	2 lbs	1 lbs
Weight Capacity	250 lbs	250 lbs



The JAY® Protector is a small, portable wheelchair seating accessory that allows clients to experience the benefits of JAY on a limited basis when sitting or moving outside of their wheelchair.



Sunrise Medical (US) LLC  
 2842 Business Park Ave. · Fresno, CA 93727 · USA