



# Quickie GP/GPV

January 2021

EIGP70

### Mark For:

Date: \_\_\_\_\_  
 Dealer Acct #: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Confirm Via:  Fax  Email

### Submitting for:

Quote  Order

PO#: \_\_\_\_\_

### ADDITIONAL SHIPPING INFORMATION

Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

## 1. FRAME

### Base Model

For sizes 14" or less please use HCPCS E1235

<input type="checkbox"/>	024WM23	K0005	Quickie GP-70°	\$1,953.00
<input type="checkbox"/>	024WM23	K0005	Quickie GPV-70° (NA 12"-13")	\$1,953.00

### Frame Angle

<input type="checkbox"/>	024FA6	60° (N/A 18 1/4" and 19 1/4" frame length)	NC
<input type="checkbox"/>	024FA7	70°	NC
<input type="checkbox"/>	024FA1	75° (N/A with rigid footrest)	NC
<input type="checkbox"/>	024FA8	80° (N/A with rigid footrest)	NC

## 2. SIZE / SEATING

### Seat Width

Seat Width 1/2" Narrower

<input type="checkbox"/>	024W3		12"	NC
<input type="checkbox"/>	024W4		13"	NC
<input type="checkbox"/>	024W5		14"	NC
<input type="checkbox"/>	024W6		15"	NC
<input type="checkbox"/>	024W7		16"	NC
<input type="checkbox"/>	024W8		17"	NC
<input type="checkbox"/>	024W9		18"	NC
<input type="checkbox"/>	024W41		19"	\$128.00
<input type="checkbox"/>	024W11	E2201	20"	\$510.00
<input type="checkbox"/>	024W12	E2201	21"	\$510.00
<input type="checkbox"/>	024W13	E2201	22"	\$510.00

### Seat Depth

21 & 22" seat depths only available with 16" or wider chairs

<input type="checkbox"/>	024SD5		12"	NC
<input type="checkbox"/>	024SD6		13"	NC
<input type="checkbox"/>	024SD7		14"	NC
<input type="checkbox"/>	024SD8		15"	NC
<input type="checkbox"/>	024SD9		16"	NC
<input type="checkbox"/>	024SD10		17"	NC
<input type="checkbox"/>	024SD11		18"	NC
<input type="checkbox"/>	024SD72		19"	\$77.00
<input type="checkbox"/>	024SD13	E2203	20"	\$561.00
<input type="checkbox"/>	024SD14	E2203	21"	\$561.00
<input type="checkbox"/>	024SD33	E2204	22"	\$918.00

### Frame Length

<input type="checkbox"/>	024FL30	18-1/4" (Avail. on 12" Sling Depth)	NC
<input type="checkbox"/>	024FL31	19-1/4" (Avail. on 13" Sling Depth)	NC
<input type="checkbox"/>	024FL32	20-1/4" (Avail. on 14", 15" Sling Depth)	NC
<input type="checkbox"/>	024FL33	21-1/4" (Avail. on 14", 15", 16" Sling Depth)	NC
<input type="checkbox"/>	024FL34	22-1/4" (Avail. on 16", 17" Sling Depth)	NC
<input type="checkbox"/>	024FL35	23-1/4" (Avail. on 17", 18" Sling Depth)	NC
<input type="checkbox"/>	024FL36	24-1/4" (Avail. on 18" Sling Depth)	NC
<input type="checkbox"/>	024FL37	25-1/4" (Avail. on 19" Sling Depth)	NC
<input type="checkbox"/>	024FL38	26-1/4" (Avail. on 20" Sling Depth)	NC
<input type="checkbox"/>	024FL39	27-1/4" (Avail. on 21" Sling Depth)	NC
<input type="checkbox"/>	024FL40	28-1/4" (Avail. on 22" Sling Depth)	NC

### Frame Height

See Appendix 11 for height chart. Short frame N/A w/Xtender

<input type="checkbox"/>	024FH2	Short (front 9.25"-rear 8.25") N/A with Xtender Power Assist	\$112.00
<input type="checkbox"/>	024FH5	Standard (front 10" - rear 9")	NC
<input type="checkbox"/>	024FH7	Tall (front 11" - rear 10")	\$112.00

### Cushion

<input type="checkbox"/>	024CU9	Reinforced Seat Sling	\$31.00
<input type="checkbox"/>	024CU1	K0669 2" Cushion	\$153.00
<input type="checkbox"/>	024CU2	Omit Cushion	NC
<input type="checkbox"/>	024CU4	Omit Cushion and Seat Sling	-\$15.00

### JAY Bundled Cushions - Optional

To order any JAY cushions on-chair, please see JAY Price List. Cushion name or part number must be filled out with desired width & depth.

<input type="checkbox"/>	024JC20	JAY Bundled Cushion Cushion Name or Part #: _____ Width: _____ Depth: _____	Price List
--------------------------	---------	---	------------

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
<b>3. CASTERS</b>									

Caster Wheels				
<input type="checkbox"/>	024CW1		3" Micro	NC
<input type="checkbox"/>	024CW3		4" Low Profile Polyurethane	NC
<input type="checkbox"/>	024CW31		4" Micro Lighted	\$71.00
<input type="checkbox"/>	024CW6	E2219	5" Soft Roll Semi-Pneumatic	\$92.00
<input type="checkbox"/>	024CW5		5" Low Profile Polyurethane	NC
<input type="checkbox"/>	024CW7	E2219	5" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$87.00
<input type="checkbox"/>	024CW33	E2214 E2215	6" Pneumatic	\$77.00
<input type="checkbox"/>	024CW34		6" Polyurethane	NC
<input type="checkbox"/>	024CW35	E2219	6" x 1.5" Semi-Pneumatic <i>only avail. For use in bottom fork holes</i>	\$87.00
<input type="checkbox"/>	024CW21	E2214 E2215	8" Pneumatic <i>Only available on 60° frame option</i>	\$87.00

Caster Wheels - Continued				
<input type="checkbox"/>	024CW32		8" Polyurethane <i>Only available on 60° frame option</i>	NC
<input type="checkbox"/>	024CW36	E2214 E2215	8" x 2" Pneumatic <i>Only available on 60° frame option</i>	\$148.00
<input type="checkbox"/>	024CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert <i>Only available on 60° frame option</i>	\$230.00

Caster Fork				
<input type="checkbox"/>	024CS5		3/4" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	024CS9		1-1/2" Longer Fork Stem Bolt	NC
<input type="checkbox"/>	024CS3	K0073	Caster Pin Locks	\$199.00

### 4. REAR WHEELS

Xtender Power Assist				
<i>RSH 15"-19" with 22" wheels, RSH 16"-20" with 24" wheels. Spare batteries must match set up. Anti-Tips included. Dual Drive mode available w/B4ME</i>				
<input type="checkbox"/>	024XT6	E0986	22" Xtender Power Assist JWX-2 <b>NICKEL</b>	\$8,206.00
<input type="checkbox"/>	024XT7	E0986	22" Xtender Power Assist JWX-2 <b>LITHIUM</b>	\$9,022.00
<input type="checkbox"/>	024XT4	E0986	24" Xtender Power Assist JWX-2 <b>NICKEL</b>	\$8,206.00
<input type="checkbox"/>	024XT5	E0986	24" Xtender Power Assist JWX-2 <b>LITHIUM</b>	\$9,022.00
<i>Please select one of the following handrim options</i>				
<input type="checkbox"/>	024XH5		Plastic Coated	NC
<input type="checkbox"/>	024XH6		Stainless Steel	NC

Rear Wheel Tire				
<input type="checkbox"/>	024RT1		Pneumatic	NC
<input type="checkbox"/>	024RT2	E2211 E2213	Pneumatic Airless Insert <i>N/A with 26" wheels &amp; Xtender</i>	\$87.00
<input type="checkbox"/>	024RT3		Full Poly <i>N/A with 26" wheels</i>	NC
<input type="checkbox"/>	024RT6	E2211 E2212	High-Pressure Clincher	\$97.00

Hand Rim				
<input type="checkbox"/>	024HR1		Aluminum Anodized	NC
<input type="checkbox"/>	024HR4		Plastic Coated <i>-N/A on all rims</i>	\$112.00
<input type="checkbox"/>	024HR7		Omit Hand Rim	NC

Xtender - Other Options				
<input type="checkbox"/>	024XA10		Spare NiMH Power Pack - NICKEL	\$1,423.00
<input type="checkbox"/>	024XA11		Spare Li-ion Power Pack - LITHIUM	\$1,831.00

Projection Hand Rim				
<i>Projections N/A w/ plastic coated handrims. N/A w/ Spoke guards</i>				
<input type="checkbox"/>	024HR8		8 Oblique (24"/26")	\$332.00

Rear Wheel Size				
<i>Spinergy N/A w/ Xtender, spoke guards. Contact Built 4 Me for Spinergy Color Options</i>				
<input type="checkbox"/>	024RW8		22" Mag	NC
<input type="checkbox"/>	024RW9		22" Spoke	\$204.00
<input type="checkbox"/>	024RW1		24" Mag	NC
<input type="checkbox"/>	024RW2		24" Spoke	NC
<input type="checkbox"/>	024RW11		26" Spoke <i>(N/A with Xtender Power Assist)</i>	\$204.00
<input type="checkbox"/>	024RW42		24" Spinergy (540) SPOX Black	\$862.00
<input type="checkbox"/>	024RR25		Omit Rear Wheel and Axle Intended Size for Set Up :	-\$55.00

Axle				
<input type="checkbox"/>	024AX1		Quick Release Steel	NC
<input type="checkbox"/>	024AX11		Quad Quick Release	\$128.00

Axle Plate				
<input type="checkbox"/>	024AL2		Standard Adjustable	NC
<input type="checkbox"/>	024AL13		Double Slotted Adjustable	\$61.00
<input type="checkbox"/>	024AL6	E0959	Amputee Adjustable <i>N/A with padded swing-away armrests or stroller handles</i>	\$179.00

### 5. WHEEL LOCKS

Wheel Lock				
<input type="checkbox"/>	024WL1		High Mount Push	NC
<input type="checkbox"/>	024WL2		High Mount Pull	NC
<input type="checkbox"/>	024WL5		High Mount Scissor <i>N/A w/ extension handles, grade aides</i>	NC
<input type="checkbox"/>	024WL12		Do Not Mount Wheel Locks	NC
<input type="checkbox"/>	024WL10		Omit	NC

Extension Handle				
<input type="checkbox"/>	024EH1	E0961	6" Pair	\$77.00

Grade Aid				
<input type="checkbox"/>	024GA1	E0974	Grade Aid <i>N/A w/ poly, high-pressure clincher tires, or Xtender</i>	\$153.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
<b>6. FOOTREST</b>									

Footplates				
<i>Heel loop or leg strap not included, See Footplate Options Section</i>				
<b>1. GP FOOTREST OPTIONS</b>				
<input type="checkbox"/>	024F56	K0040	Rigid N/A w/ Ext. Tubes; N/A with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	K0040	Platform N/A w/ Ext. Tubes	\$97.00
<input type="checkbox"/>	024F1	K0040	Composite Footplate	NC
<input type="checkbox"/>	024F3	K0040	Adult Angle Adj.	\$189.00
<input type="checkbox"/>	024F48	K0040	Platform Flip Up	\$194.00
<b>2. GPV FOOTREST OPTIONS</b>				
<input type="checkbox"/>	024F56	K0040	Rigid N/A w/ Ext. Tubes; N/A with 75° or 80° frame angle	NC
<input type="checkbox"/>	024F50	K0040	Platform N/A w/ Ext. Tubes	\$97.00
<input type="checkbox"/>	024F48	K0040	Platform Flip Up	\$194.00

Footrest Extension Tubes				
<input type="checkbox"/>	024E1		Short (14" to 16-1/2")	NC
<input type="checkbox"/>	024E2		Medium (16-1/2" to 19")	NC
<input type="checkbox"/>	024E5		Long (19" to 21-1/2")	NC
<input type="checkbox"/>	024E16		High Mount Only avail. with adult angle adjustable footplates	\$56.00

Footplate Options (Ext Tube & Direct Mount)				
<input type="checkbox"/>	024HL1	E0951	Heel Loop	\$41.00
<input type="checkbox"/>	024HL5	K0038	Leg Strap	\$41.00
<input type="checkbox"/>	024HL4		Omit Heel Loop	NC
<input type="checkbox"/>	024HL6		Omit Leg Strap	NC
<input type="checkbox"/>	024F7		Omit Footplate & Ext Tube	-\$15.00
<input type="checkbox"/>	024IG2		Impact Guard Neoprene	\$92.00

**7. BACKREST**

Back Type				
<input type="checkbox"/>	<b>024BT3</b>		<b>Non-Folding Backrest (Fixed)</b> <i>N/A with push handles or adj. upholstery</i>	
<input type="checkbox"/>	024BB1		8" Back Height	NC
<input type="checkbox"/>	024BB2		9" Back Height	NC
<input type="checkbox"/>	024BB3		10" Back Height	NC
<input type="checkbox"/>	024BB4		11" Back Height	NC
<input type="checkbox"/>	024BB5		12" Back Height	NC
<input type="checkbox"/>	024BB6		13" Back Height	NC
<input type="checkbox"/>	024BB7		14" Back Height	NC
<input type="checkbox"/>	024BB8		15" Back Height	NC
<input type="checkbox"/>	024BB9		16" Back Height	NC
<input type="checkbox"/>	<b>024BT1</b>		<b>Folding Backrest (Adjustable)</b>	
<input type="checkbox"/>	024BA1		11" - 14" Back Height	\$66.00
<input type="checkbox"/>	024BA2		12" - 16" Back Height	\$66.00
<input type="checkbox"/>	024BA3		14" - 18" Back Height	\$66.00
<input type="checkbox"/>	024BA4		16" - 20" Back Height	\$66.00

Back Upholstery				
<input type="checkbox"/>	024BU1		Standard	NC
<input type="checkbox"/>	024BU8		Adjustable	\$179.00
<input type="checkbox"/>	024BU16		Reinforced (N/A w/ adj. upholstery)	\$31.00
<input type="checkbox"/>	024BU14		Omit Back Upholstery	-\$15.00

Backrest Options				
<input type="checkbox"/>	024BT33		Folding w/ Integral Push Handles	NC
<input type="checkbox"/>	024A85		Bolt-On Push Handle (Folding only)	\$117.00
<input type="checkbox"/>	024SR2		Regular Stroller Handle <i>N/A with amputee axle plate</i>	\$138.00
<input type="checkbox"/>	024SR3		Tall Stroller Handle <i>N/A with amputee axle plate</i>	\$138.00

JAY Back (Optional)				
<input type="checkbox"/>	024JB01		J3 Back (See J3 Order Form)	Price List

**8. ARMREST / SIDEGUARDS**

Single Post				
<input type="checkbox"/>	<b>024AA9</b>	<b>E0973</b>	<b>Single Post Ht. Adjust. Low</b> <i>Must be removed to fold back down</i>	<b>\$291.00</b>
<input type="checkbox"/>	024AM1		Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2		Full Length - Classic (14")	NC
<input type="checkbox"/>	<b>024AA8</b>	<b>E0973</b>	<b>Single Post Ht. Adjust. Std</b> <i>Must be removed to fold back down</i>	<b>\$291.00</b>
<input type="checkbox"/>	024AM1		Desk Length - Classic (10")	NC
<input type="checkbox"/>	024AM2		Full Length - Classic (14")	NC

Padded Swing Away Armrest				
<i>All Padded Swing-Away N/A with amputee axle plate or Xtender Power Assist</i>				
<input type="checkbox"/>	024AA19	E0973	Padded S-A Medium	\$230.00
<input type="checkbox"/>	024AA16	E0973	Padded S-A Adult	\$230.00

Side Guard				
<input type="checkbox"/>	024SG3	K0108	Plastic Regular	\$148.00
<input type="checkbox"/>	024SG4	K0108	Plastic Kids	\$148.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
<b>9. KEY OPTIONS AND ACCESSORIES</b>									

Rear Anti-Tip				
<input type="checkbox"/>	024AT1	E0971	Rear Anti-Tip	\$153.00

Backpack				
<input type="checkbox"/>	024A1		Quickie Backpack Black	\$77.00
<input type="checkbox"/>	024A2		Quickie Kids Backpack Black	\$61.00

Seat Pouch				
<input type="checkbox"/>	024A8		Quickie Seat Pouch Black	\$51.00

Positioning Belts				
<input type="checkbox"/>	024A142	E0978	Belt 1.5" Auto Buckle w/Clamps	\$71.00
<input type="checkbox"/>	024A12	E0978	Belt 2" Aircraft Buckle w/Clamps	\$112.00
<input type="checkbox"/>	024A13	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$122.00

Spoke Guards				
<i>N/A w/ Xtender Power Assist; Only available w/ Mag and Spoke wheels</i>				
<input type="checkbox"/>	024A23	K0065	Black Spoke Guards	\$117.00
<input type="checkbox"/>	024A24	K0065	Clear Spoke Guards	\$117.00

Crutch Holder				
<input type="checkbox"/>	024A53	E2207	Left Side	\$199.00
<input type="checkbox"/>	024A52	E2207	Right Side	\$199.00

Extra Leg Strap				
<input type="checkbox"/>	024A55	K0038	Extra Leg Strap Black	\$31.00

Miscellaneous				
<input type="checkbox"/>	024A83		Tool Kit	\$56.00
<input type="checkbox"/>	024A999		Quickie Water Bottle & Holder	\$56.00

Tray Table				
<input type="checkbox"/>	024A31	E0950	Tray Table	\$260.00

<b>10. COLOR</b>									
------------------	--	--	--	--	--	--	--	--	--

Standard Frame Colors				
<input type="checkbox"/>	024C140		Aztec Gold	NC
<input type="checkbox"/>	024C2		Black <i>(Not avail as base w/Zebra)</i>	NC
<input type="checkbox"/>	024C141		Black Cherry	NC
<input type="checkbox"/>	024C10		Black Opal	NC
<input type="checkbox"/>	024C38		Blue Opal	NC
<input type="checkbox"/>	024C36		Candy Blue	NC
<input type="checkbox"/>	024C14		Candy Purple	NC
<input type="checkbox"/>	024C35		Candy Red	NC
<input type="checkbox"/>	024C142		Electric Blue	NC
<input type="checkbox"/>	024C143		Evergreen	NC
<input type="checkbox"/>	024C126		Glow	NC
<input type="checkbox"/>	024C144		Green Apple	NC
<input type="checkbox"/>	024C147		Hot Sparkle Pink	NC
<input type="checkbox"/>	024C148		Mauve Pink	NC
<input type="checkbox"/>	024C103		Rootbeer	NC
<input type="checkbox"/>	024C114		Sparkle Silver	NC
<input type="checkbox"/>	024C145		Sunrise Orange	NC
<input type="checkbox"/>	024C113		Titanium Color Paint	NC
<input type="checkbox"/>	024C8		Yellow	NC

Matte Frame Colors				
<input type="checkbox"/>	024C116		Matte Black	NC
<input type="checkbox"/>	024C151		Matte Black Cherry	NC
<input type="checkbox"/>	024C150		Matte Electric Blue	NC
<input type="checkbox"/>	024C149		Matte Evergreen	NC
<input type="checkbox"/>	024C152		Matte Purple	NC

Personalized Frame Color				
<input type="checkbox"/>	024C112		Desert Camouflage	\$403.00
<input type="checkbox"/>	024C43		Mossy Oak Camouflage	\$403.00
<input type="checkbox"/>	024C42		Stars & Stripes	\$403.00
<input type="checkbox"/>	024C108		Zebra <i>(Choose above base color) Std with White unless above base paint is chosen N/A w/ Blck, Matte Black, Blck Opal or Titanium.</i>	\$403.00

11. APPENDIX

**GP Series Rear STF Chart**

\*All measurements given are an estimate and have been rounded to the nearest .5"

SHORT Frame (9 1/4" Front / 8 1/4" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
			22"	Mag / Spoke	Full Poly	14.5	15	15.5
22"	Mag / Spoke	Pneum / Insert	14.5	15	16	16	17	17.5
24"	Mag / Spoke	Full Poly	15	15.5	16.5	16.5	17	17.5
24"	Spoke / Spinerger	Clincher	15	16	16.5	16.5	17	18
24"	Mag / Spoke / Spinerger	Pneum / Insert	15.5	16	17	17	17.5	18
26"	Spoke / Spinerger	Clincher	16	17	17.5	17.5	18	19
26"	Spoke / Spinerger	Pneumatic	16.5	17	17.5	18	18.5	19

TALL Frame (11" Front / 10" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
			22"	Mag / Spoke	Full Poly	16	16.5	17.5
22"	Mag / Spoke	Pneum / Insert	16.5	17	17.5	18	18.5	19
24"	Mag / Spoke	Full Poly	17	17.5	18	18.5	18.5	19.5
24"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	17	17.5	18.5	19	19	20
26"	Spoke	Clincher	18	18.5	19	19.5	20	20.5
26"	Spoke	Pneumatic	18	19	19.5	19.5	20	21

STANDARD Frame (10" Front / 9" Rear)			Rear Axle Plate Hole Position					
Wheel	Rim	Tire	Standard			Inverted		
			Top	Middle	Bottom	Top	Middle	Bottom
			22"	Mag / Spoke	Full Poly	15	15.5	16.5
22"	Mag / Spoke	Pneum / Insert	15.5	16	17	17	17.5	18
24"	Mag / Spoke	Full Poly	16	16.5	17	17	18	18.5
24"	Spoke / Spinerger	Clincher	16	16.5	17	17.5	18	18.5
24"	Mag / Spoke / Spinerger	Pneum / Insert	16	16.5	17.5	17.5	18	18.5
26"	Spoke / Spinerger	Clincher	17	17.5	18	18.5	19	19.5
26"	Spoke / Spinerger	Pneumatic	17	18	18.5	18.5	19.5	20


Maximum Dump 2 3/4"

\*\*\* Custom Frame Height Through Configurator - REAR 6 1/2 - 13" (Max frame squeeze = 3 1/2") Calculate Accordingly

**GP Series Front STF Chart**

\*All measurements given are an estimate and have been rounded to the nearest .5"

Frame Dimensions



SHORT Frame (9 1/4" Front / 8 1/4" Rear)				3"	4" Poly 4" Lighted	5" Poly 5" x 1.5 5" Soft Roll	6" Poly 6" x 1.5	6" Pneu	8" Poly	8" Pneu	8x2
4 1/4"	STD	Bottom		15.50	16.00	16.50	17.00	17.00	x	x	x
4 1/4"	3/4"	Bottom		16.50	17.00	17.00	17.50	18.00	x	x	x
4 1/4"	1 1/2"	Bottom		17.00	17.50	18.00	18.50	18.50	x	x	x
5 1/4"	STD	Top		x	x	16.50	17.00	17.00	x	x	x
5 1/4"	STD	Bottom		16.50	17.00	17.5	18.00	18.00	19.00	20.5	x
5 1/4"	3/4"	Top		x	x	17.50	18.00	18.00	x	x	x
5 1/4"	3/4"	Bottom		17.50	18.00	18.50	18.50	19.00	19.50	20.00	x
5 1/4"	1 1/2"	Top		x	x	18.00	18.50	18.50	x	x	x
5 1/4"	1 1/2"	Bottom		18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6"	STD	Top		x	x	17.50	18.00	18.00	18.50	19.00	x
6"	STD	Bottom		17.50	18.00	18.50	19.00	19.00	19.50	20.00	x
6"	3/4"	Top		x	x	18.00	18.50	18.50	19.50	19.50	x
6"	3/4"	Bottom		18.00	18.50	19.00	19.50	19.50	20.50	20.50	x
6"	1 1/2"	Top		x	x	19.00	19.50	19.50	20.00	20.50	x
6"	1 1/2"	Bottom		19.00	19.50	20.00	20.50	20.50	21.00	21.50	x
7"	STD	Top		x	x	18.50	19.00	19.00	19.50	20.00	x
7"	STD	Bottom		18.50	19.00	19.50	20.00	20.00	20.50	21.00	x
7"	3/4"	Top		x	x	19.00	19.50	19.50	20.50	20.50	x
7"	3/4"	Bottom		19.00	19.50	20.00	20.50	20.50	21.50	21.50	x
7"	1 1/2"	Top		x	x	20.00	20.50	20.50	21.00	21.50	x
7"	1 1/2"	Bottom		20.00	20.50	21.00	21.00	21.50	22.00	22.50	x
8x2	STD	Top		x	x	x	x	x	x	x	19.00
8x2	STD	Bottom		x	x	x	x	x	x	x	20.00
8x2	3/4"	Top		x	x	x	x	x	x	x	19.50
8x2	3/4"	Bottom		x	x	x	x	x	x	x	20.50
8x2	1 1/2"	Top		x	x	x	x	x	x	x	20.50
8x2	1 1/2"	Bottom		x	x	x	x	x	x	x	21.00



Sunrise Medical (US) LLC

2842 Business Park Ave. · Fresno, CA 93727 · USA

STANDARD Frame (10" Front / 9" Rear) ADD 3/4" to each measurement above

TALL Frame (11" Front / 10" Rear) ADD 1 3/4" to each measurement above

\*\*\* Custom Frame Height through Configurator - Front 9 - 13" - Calculate Accordingly (Minimum Front STF = 15 1/2" / Max = 26")