



# Quickie Pulse™ 6 BC Group 3

K0848, K0011  
EIPW22

Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City, State, Zip: _____
Marked For: _____	Phone: _____

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

## 1. QUICKIE Pulse™ 6 GROUP 3 CATEGORY (Select One)

<b>QUICKIE Pulse™ 6 BC (Group 3 K0848, K0011)</b>			
<input type="checkbox"/>	Quickie Pulse™ 6 BC w/ASAP II Seat (Fig 1.1)		<b>\$6,661.00</b>
<input checked="" type="checkbox"/>	195PC74	K0848 / K0011	Base Package \$6,661.00



Fig 1.1 Pulse™ 6 BC

BC Seating Specifications			
<b>HCPC Code</b>	K0848 / K0011		
<b>Typical Chair Weight</b>	<b>W/battery</b>	235lbs	<b>W/O Battery</b> 157lbs
<b>Transit</b>	ISO occupied approved		
<b>Seat Type</b>	ASAP II		
<b>Seat size range</b>	<b>Width</b>	12"-22"	<b>Depth</b> 12"-22"
<b>STFH range</b>	<b>ASAP II</b>	16.2" - 20.2"	
<b>Power Seating</b>	9" Seat elevator		

QUICKIE Pulse™ 6 Base Specifications					
<b>Drive Wheel Position</b>	Midwheel Drive			<b>Drive Wheels</b>	13" Drive Wheels w/ flat free insert (Black)
<b>User weight Cap</b>	300 lb			<b>Tires/Casters</b>	7" Front / 6" rear (Black)
<b>Base size</b>	<b>Width</b>	24"	<b>Length</b> 34"	<b>Charger</b>	8A Dual Mode Charger (Compact) CTE
			<b>Turn Radius</b> 18"	<b>Battery Size</b>	Group 22
<b>Drive Control</b>	VR-2 90A (non-expand)		R-net EL90 (Expandable)	<b>Shroud Colors</b>	12 standard colors, 3 upgrade prints.
<b>Driving Range</b>	Up to 13.5 miles			<b>Ground Clearance</b>	3"
<b>Drive Motors</b>	4-pole 6 mph				

## 2. BASE CHAIR OPTIONS

Shroud Color				
<input type="checkbox"/>	195EZ01	-	Candy Apple Red	N/C
<input type="checkbox"/>	195EZ02	-	High Gloss Black	N/C
<input type="checkbox"/>	195EZ03	-	Candy Blue	N/C
<input type="checkbox"/>	195EZ04	-	Silver	N/C
<input type="checkbox"/>	195EZ05	-	Copper	N/C
<input type="checkbox"/>	195EZ08	-	Green Apple	N/C
<input type="checkbox"/>	195EZ09	-	Yellow	N/C
<input type="checkbox"/>	195EZ10	-	Pearl Pink	N/C
<input type="checkbox"/>	195EZ11	-	Candy Purple	N/C
<input type="checkbox"/>	195EZ12	-	Matte Black	N/C
<input type="checkbox"/>	195EZ13	-	Midnight Blue	N/C
<input type="checkbox"/>	195EZ15	-	Pearl White	N/C
<input type="checkbox"/>	195EZ06	-	Digital Camo.	\$525.00
<input type="checkbox"/>	195EZ07	-	Woodland Camo.	\$525.00
<input type="checkbox"/>	195EZ14	-	Carbon Fiber Print	\$525.00
<input type="checkbox"/>	195EZ17	-	Pink Camo	\$525.00
<input type="checkbox"/>	195EZ18	-	Grunge Skulls	\$525.00
<input type="checkbox"/>	195EZ16	-	American Flag	\$525.00

Drive wheel/tire package				
<input type="checkbox"/>	195DW60	-	Cast Wheel Solid Dual Density Inserts <b>Black</b> Knobby Tires and <b>Black</b> Casters	Std

User weight				
<input type="checkbox"/>	195WG06	-	< 150lb.	N/C
<input type="checkbox"/>	195WG13	-	151-200lb	N/C
<input type="checkbox"/>	195WG14	-	201-250lb	N/C
<input type="checkbox"/>	195WG03	-	251-300lb	N/C

Ventilator Options				
<i>Not available w/ Transit or Power Legs</i>				
<input type="checkbox"/>	195VT4	E1029	Fixed Vent Tray (Pulmonetic LTV Models)	\$474.00

Motor Package				
<input type="checkbox"/>	195PF88	-	300lbs 6mph 4-Pole	Std

Transit Option				
<input type="checkbox"/>	195TO40	-	Omit Transit	NC
<input type="checkbox"/>	195TO25	-	ISO Transit Package	\$316.00

Pre-Set Seat Frame Angle				
<input type="checkbox"/>	195FA10	-	0°	N/C
<input type="checkbox"/>	195FA3	-	3° NA with +3" or +4" STFH or lift	N/C
<input type="checkbox"/>	195FA4	-	6° NA with +3" or +4" STFH or Lift	N/C

Seat Height				
<i>(Measured in the front) - Obtain Quote to verify est. STFH</i>				
<input type="checkbox"/>	195HT44		Position "0"	NC
<input type="checkbox"/>	195HT45		+1" Seat to Floor Height	NC
<input type="checkbox"/>	195HT46		+2" Seat to Floor Height	NC
<input type="checkbox"/>	195HT47		+3" Seat to Floor Height	NC
<input type="checkbox"/>	195HT48		+4" Seat to Floor Height	NC

STFH Matrix					
Seating Type and Power Option Combinations	Seat-to-Floor Height at 0 Pre-tilt				
	Low	+1"	+2"	+3"	+4"
No power seating.	16"	17"	18"	19"	20"
Lift only.	18.4"	19.4"			

Adding seat angle will increase STFH by 1" for each incremental change.

## 3. BATTERIES and CHARGERS

Batteries				
<input type="checkbox"/>	195AG21	E2361	GP22 Flat Top Batteries Pair w/strap	\$536.00
<input type="checkbox"/>	195AG05	-	Omit batteries (Includes harness)	N/C

Battery Charger				
<input type="checkbox"/>	195BJ19	-	8A Dual Mode Charger (Compact) CTE	N/C

## 4. ELECTRONICS

Drive Control Options			
Non-Expandable Controllers			
<input type="checkbox"/>	195ER01	Non-expandable VR2-90 controller (NA with High Speed Motor package)	N/C
Expandable Controllers			
<input type="checkbox"/>	195ER02	R-Net Expandable Drive Control Upgrade (90A)	\$1,000.00
<input type="checkbox"/>	195EX01	E2377 Expandable Controller - R-net	\$551.00
<input type="checkbox"/>	195EX02	E2313 Harnessing for Expandable Controller	\$449.00

VR2 Joystick (non-expandable)			
<input type="checkbox"/>	195J94	- Non Expandable VR2 Joystick	N/C

R-net Joystick			
<input type="checkbox"/>	195J70	- Std LED Joystick for R-net Ctrl's	N/C
<input type="checkbox"/>	195J108	- CJSM1: Joystick Standard: w/Stereo Jacks & Color Display	N/C
<input type="checkbox"/>	195J114	- CJSM2: Joystick Advanced: 3.5" Color Display, Paddle Switches, Shortcut Keys, IR, BT, and Stereo Jack.	\$1,076.00
<input type="checkbox"/>	195J95	- Omit JS. For use with Omni2	N/C

Joystick Mount (Fixed, Swing Away, Retractable, Mid)			
<input type="checkbox"/>	195JM10	- Fixed Mount	N/C
<input type="checkbox"/>	195JM11	E1028 Swing Away Mount	\$250.00
<input type="checkbox"/>	195JM15	E1028 Retractable Joystick mount w/45° stops	\$597.00
<input type="checkbox"/>	195JM14	E1028 Mid Line Mount	\$250.00

Joystick Mount (Location)			
<i>Height Adjustable.</i>			
<input type="checkbox"/>	195HM1	- Right Side Mount	N/C
<input type="checkbox"/>	195HM2	- Left Side Mount	N/C

Joystick Handle (Optional)			
<input type="checkbox"/>	195JH1	E2323 Ball Handle	\$82.00
<input type="checkbox"/>	195JH2	E2323 Foam Ball Handle	\$82.00
<input type="checkbox"/>	195JH4	E2323 T-Handle	\$82.00
<input type="checkbox"/>	195JH6	E2323 Goalpost	\$82.00

Attendant Control			
<input type="checkbox"/>	1950A01	E2331 Attendant Control for VR-2	\$1,046.00
<input type="checkbox"/>	1950A02	E2331 Attendant Control for R-net	\$1,046.00

Attendant Mounting location			
<input type="checkbox"/>	195DN1	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	195DN2	Mounting Location: Left Back Cane	N/C

R-net Omni2 and Output Options			
<i>(For expandable controls, choose all that apply)</i>			
<input type="checkbox"/>	195ON16	- R-net Omni2 Input Control Module and 3.5" LCD display with built in Infrared, Bluetooth and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$1,076.00
<input type="checkbox"/>	195ON2	E2351 R-net Output Module (ECM)	\$1,046.00

Omni2 Mount			
<i>(If used with Joystick must be mounted opposite side)</i>			
<input type="checkbox"/>	195DS17	E1028 Swing Away Left Side Omni2 Mount w/ Rear Mounted Interface	\$250.00
<input type="checkbox"/>	195DS16	E1028 Swing Away Right Side Omni2 Mount w/ Rear Mounted Interface	\$250.00

Switch inputs for R-Net Joysticks or Omni2			
<i>Drive or Mode / Profile Control</i>			
<input type="checkbox"/>	195SW22	Stereo to 2 Mono Splitter	\$66.00
<input type="checkbox"/>	195SW1	Buddy Button Switch	\$82.00
<input type="checkbox"/>	195SW4	Micro Light Switch	\$82.00
<input type="checkbox"/>	195SW24	Egg Switch-Black	\$102.00
<input type="checkbox"/>	195SW21	Dual toggle switch and bracket	\$383.00
<input type="checkbox"/>	195SW35	Switch-It Mini 1 button on Mono plug	\$204.00
<input type="checkbox"/>	195SW36	Switch-It Mini 2 button on 2 Mono plugs	\$204.00
<input type="checkbox"/>	195SW37	Switch-It Mini 2 button on 1 Stereo plug	\$204.00
<input type="checkbox"/>	195SW38	Switch-It Mini 4 button on DB9 connector	\$204.00

Single Switch Mount Bracket			
<input type="checkbox"/>	195EM01	Single Switch Mount Bracket	\$204.00

Electronic Programmers			
<input type="checkbox"/>	195PR2	VR-2 Hand Held Programmer	\$265.00
<input type="checkbox"/>	195PR5	R-net PC Programmer with Dongle	\$423.00

**5. Specialty Driver Control Options**

*Requires expandable R-net controls and Omni2 unless otherwise noted. Not available with Upholstered backs.*

**5a. Proportional Control Packages**

<b>Head Control Proportional Packages</b>			
□ Switch-It Dual Pro Proportional Head Control with Whitmyer and Facial Pads (Adult)			<b>\$5,646.00</b>
195ID72	E2328	Switch-It Dual Pro Head Control w/facial pad (Adult)	\$5,299.00
195J09	E1028	S/A Boom Mount with Egg Switch Control	\$347.00
□ Switch-It Dual Pro Proportional Head Control with Spot Pads Package			<b>\$5,646.00</b>
195ID70	E2328	Switch-It Dual Pro Head Control w/Spot pads	\$5,299.00
195J09	E1028	S/A Boom Mount with Egg Switch Control	\$347.00

<b>Hand Control Proportional Joystick Packages</b>			
□ Switch-It Versa-Guide Compact Joystick Midline with Gatlin Package			<b>\$1,933.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195J07	E1028	S/A Gatlin Mount	\$250.00
□ Switch-It Versa-Guide Compact Joystick Midline and prox flex switches with Gatlin Package			<b>\$2,469.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$536.00
195J07	E1028	S/A Gatlin Mount	\$250.00
□ Switch-It Micro-Pilot Joystick with SpaceDisc Package			<b>\$3,682.00</b>
195ID62	E2312	Micro-Pilot Proportional Joystick	\$3,055.00
195J08	K0108	S/A SpaceDisc mount	\$627.00
□ Switch-It Micro-Guide Joystick with SpaceDisc Package			<b>\$3,682.00</b>
195ID63	E2312	Micro-Guide Proportional Joystick	\$3,055.00
195J08	K0108	S/A SpaceDisc mount	\$627.00

<b>Chin Control Proportional Joystick Packages</b>			
□ Switch-It Versa-Guide Chin Control Compact Joystick and Egg Switch Package			<b>\$2,291.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM07	E1028	S/A Mount off backrest	\$261.00
195J04	E1028	S/A Boom Mount with egg switch	\$347.00
□ Switch-It Versa-Guide Chin Control Compact Joystick and Proximity Flex Switches Package			<b>\$2,480.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$536.00
195IM07	E1028	S/A Mount off backrest	\$261.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex switches Package Std 6"			<b>\$2,846.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$536.00
195IM05	K0108	Bib mount Small	\$627.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Proximity Flex Switches Package Large 8"			<b>\$2,846.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM08	K0108	Dual Flex Proximity Switches (requires no force)	\$536.00
195IM06	K0108	Bib mount Large	\$627.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Std 6"			<b>\$2,657.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM05	K0108	Bib mount Small	\$627.00
195J04	E1028	SA Boom Mount with egg switch	\$347.00
□ Switch-It Versa-Guide Bib Chin Control Compact Joystick with Egg Switch Package Large 8"			<b>\$2,657.00</b>
195ID71	E2373	Switch-It Versa-Guide Compact Proportional Joystick	\$1,683.00
195IM06	K0108	Bib mount Large	\$627.00
195J04	E1028	S/A Boom Mount with egg switch	\$347.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Std 6"			<b>\$4,029.00</b>
195ID62	E2312	Micro-Pilot Proportional Joystick	\$3,055.00
195IM05	K0108	Bib mount Small	\$627.00
195J04	E1028	S/A Boom Mount with egg switch	\$347.00
□ Switch-It Micro-Pilot Bib Chin Control with Egg Switch Package Large 8"			<b>\$4,029.00</b>
195ID62	E2312	Micro-Pilot Proportional Joystick	\$3,055.00
195IM06	K0108	Bib mount Large	\$627.00
195J04	E1028	S/A Boom Mount with egg switch	\$347.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Std 6"			<b>\$4,029.00</b>
195ID63	E2312	Micro-Guide Proportional Joystick	\$3,055.00
195IM05	K0108	Bib mount Small	\$627.00
195J04	E1028	S/A Boom Mount with egg switch	\$347.00
□ Switch-It Micro-Guide Bib Chin Control with Egg Switch Package Large 8"			<b>\$4,029.00</b>
195ID63	E2312	Micro-Guide Proportional Joystick	\$3,055.00
195IM06	K0108	Bib mount Large	\$627.00
195J04	E1028	S/A Boom Mount with egg switch	\$347.00

**5. Specialty Driver Control Options (cont.)**

**5b. Switch Control**

*Requires Headrest and/or Headrest Mount*

**3 Switch Proximity Head Arrays** *Choose headrest mount below*

<input type="checkbox"/> <b>Switch-It Head Array Package 3 Switch Curved (Long)</b> <b>\$4,626.00</b>			
195II13	E2330	Switch-It Proximity Head Array Long	\$4,279.00
195IJ09	E1028	S/A Boom Mount with Egg Switch Control	\$347.00
<input type="checkbox"/> <b>Switch-It Head Array Package 3 switch curved (Stubby)</b> <b>\$4,626.00</b>			
195II14	E2330	Switch-It Proximity Head Array Stubby	\$4,279.00
195IJ09	E1028	S/A Boom Mount with Egg Switch Control	\$347.00

**5 Switch Proximity Head Array (No Egg Switch Needed)**

*Choose "Headrest Mount for Special Controls"*

<input type="checkbox"/>	195II04	E2330	Switch-It Head Array 5 switch flat	\$4,381.00
--------------------------	---------	-------	------------------------------------	------------

**Multiple Switches**

<input type="checkbox"/>	195SW31		Switch-It 3 Naked proximity switches on DB9 connector. (requires Omni2)	\$2,729.00
<input type="checkbox"/>	195SW32		Switch-It 4 Naked proximity switches on DB9 connector. (requires Omni2)	\$3,356.00
<input type="checkbox"/>	195SW33		Switch-It 5 Naked proximity switches on DB9 connector. (requires Omni2)	\$3,779.00
<input type="checkbox"/>	195SW23		4-way toggle switch button DB9 Connector. Includes bracket	\$561.00

**Sip N Puff Package** *Choose headrest and HR mount below*

<input type="checkbox"/> <b>Quickie Sip n' Puff System</b> <b>\$2,224.00</b>			
195ID34	E2325	Sip N' Puff Kit	\$1,494.00
195IH1	E2326	Sip n' Puff Tubing and SA Hardware	\$383.00
195IJ10	E1028	S/A Boom Mount with egg switch	\$347.00

**Headrest Mount for Specialty Controls**

<input type="checkbox"/>	195MH12	E1028	Axys Mount & Bracket for Specialty Controls	\$235.00
<input type="checkbox"/>	195MH13	E1028	Cobra Mount & Bracket for Specialty Controls	\$296.00
<input type="checkbox"/>	195MH14	E1028	Linx Mount & Bracket for Specialty Controls	\$235.00

**Headrest Pad for Specialty controls**

<input type="checkbox"/>	195HE79	E0955	Plush Pad 8" for Specialty controls	\$235.00
<input type="checkbox"/>	195HE80	E0955	Plush Pad 10" for Specialty controls	\$235.00
<input type="checkbox"/>	195HE81	E0955	Plush Pad 14" for Specialty controls	\$235.00

**Extra Swing Away Head Boom Egg Switch and Mount**

	195IJ05	E1028	Extra SA Boom Mount with egg switch #1	\$347.00
	195IJ06	E1029	Extra SA Boom Mount with egg switch #2	\$347.00

**Single Switches**

<input type="checkbox"/>	195IB01		Ribbon Switch	\$102.00
<input type="checkbox"/>	195IB02		Buddy Button Switch	\$82.00
<input type="checkbox"/>	195IB03		Disc Switch	\$168.00
<input type="checkbox"/>	195IB04		Wobble Switch	\$525.00
<input type="checkbox"/>	195IB05		Micro Light Switch	\$82.00
<input type="checkbox"/>	195IB06		Egg Switch-Black	\$102.00
<input type="checkbox"/>	195IB07		Switch-It Mini 1 button on Mono plug	\$204.00

## 6. SEAT TYPE

<input type="checkbox"/>	195ST56	-	Quickie A.S.A.P II (Adjustable Seating and Positioning) Choose Seat Frame Size Platform.	N/C
--------------------------	---------	---	--	-----

### Frame size

*(each frame is width and depth adjustable within given range.)*

<input type="checkbox"/>	195FU17	-	Adjustable 12x12 to 16X19 ( A )	N/C
<input type="checkbox"/>	195FU18	-	Adjustable 16x12 to 20x20 ( B )	N/C
<input type="checkbox"/>	195FU19	-	Adjustable 18x16 to 22x22 ( C )	N/C

## 7. SEAT SIZE

### Seat Size Matrix. *Seat type is chosen in section 5.*

A = Seat Frame Adjustable from 12x12 to 16X19

B = Seat Frame Adjustable from 16X12 to 20X20

C = Seat Frame Adjustable from 18X16 to 22X22

		Seat Width										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	22"	
Seat Depth	12"	A	A	A	A	A,B	B	B	B	B		
	13"	A	A	A	A	A,B	B	B	B	B		
	14"	A	A	A	A	A,B	B	B	B	B		
	15"	A	A	A	A	A,B	B	B	B	B		
	16"	A	A	A	A	A,B	B	B,C	B,C	B,C	C	
	17"	A	A	A	A	A,B	B	B,C	B,C	B,C	C	
	18"	A	A	A	A	A,B	B	B,C	B,C	B,C	C	
	19"	A	A	A	A	A,B	B	B,C	B,C	B,C	C	
	20"					B	B	B,C	B,C	B,C	C	
	21"							C	C	C	C	
	22"							C	C	C	C	

Seat Width					
<input type="checkbox"/>	195R21		12"	Seat width	N/C
<input type="checkbox"/>	195R22		13"	Seat width	N/C
<input type="checkbox"/>	195R23		14"	Seat width	N/C
<input type="checkbox"/>	195R24		15"	Seat width	N/C
<input type="checkbox"/>	195R5		16"	Seat width	N/C
<input type="checkbox"/>	195R6		17"	Seat width	N/C
<input type="checkbox"/>	195R7		18"	Seat width	N/C
<input type="checkbox"/>	195R8		19"	Seat width	N/C
<input type="checkbox"/>	195R9		20"	Seat width	N/C
<input type="checkbox"/>	195R25		22"	Seat width	N/C

Seat Depth					
<input type="checkbox"/>	195RF20		12"	Seat depth	N/C
<input type="checkbox"/>	195RF21		13"	Seat depth	N/C
<input type="checkbox"/>	195RF22		14"	Seat depth	N/C
<input type="checkbox"/>	195RF23		15"	Seat depth	N/C
<input type="checkbox"/>	195RF5		16"	Seat depth	N/C
<input type="checkbox"/>	195RF6		17"	Seat depth	N/C
<input type="checkbox"/>	195RF7		18"	Seat depth	N/C
<input type="checkbox"/>	195RF8		19"	Seat depth	N/C
<input type="checkbox"/>	195RF24		20"	Seat depth	N/C
<input type="checkbox"/>	195RF25		21"	Seat depth	N/C
<input type="checkbox"/>	195RF26		22"	Seat depth	N/C

## 8. SEAT CUSHIONS

### JAY Cushion Bundling

<input type="checkbox"/>	195JC14	-	JAY Cushion On Chair Order Choose Cushion in section 17c.	NC
--------------------------	---------	---	--	----

## 9. POWERED SEATING

Powered Seat Elevator- 9" Travel			
<input type="checkbox"/>	195PA8	E2300	Power Seat Elevate (adds 3" to STFH) \$2,550.00

### Power Seating Control

Powered Seating <u>Switch</u> Control Method			
<i>(For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.)</i>			

Single Actuator Switch Control (Lift)			
<input type="checkbox"/>	195EO10	-	Switch Control 1 actuator with <u>Dual Toggle</u> N/C
<input type="checkbox"/>	195EO11	-	Switch Control 1 actuator with <u>Buddy Button</u> N/C
<input type="checkbox"/>	195EO19	-	Switch Control 1 actuator with Switch-It Mini 1 button N/C
<input type="checkbox"/>	195EO20	-	Switch Control 1 actuator with Switch-It Mini 2 button N/C

2-3 Actuator Switch Control (Legs only)			
<input type="checkbox"/>	195EO21	-	Switch Control 2-3 actuator with <u>Dual Toggle</u> \$510.00
<input type="checkbox"/>	195EO12	-	Switch Control 2-3 actuator with <u>4-way toggle</u> \$765.00
<input type="checkbox"/>	195EO23	-	Switch Control 2-3 actuator with Switch-It Mini 2 button \$510.00
<input type="checkbox"/>	195EO22	-	Switch Control 2-3 actuator with Switch-It Mini 4 button \$510.00

Switch Mount			
<i>(for use with switch control powered seating)</i>			
<input type="checkbox"/>	195EY01	-	Switch Mount Right N/C
<input type="checkbox"/>	195EY02	-	Switch Mount Left N/C

Powered Seating <u>Thru-drive</u> Control Method.			
Thru-Drive 1 Actuator			
<input type="checkbox"/>	195EO28	E2310	Thru Drive 1 Actuator (For operating Lift thru Joystick.) Includes Ctrl+5 when ordered with R-net. \$1,958.00
Thru-Drive 2+ Actuator (Requires R-net)			
<input type="checkbox"/>	195EO29	E2311	Thru Drive 2+ Actuators: (For operating power legs) Includes Ctrl+5 \$3,478.00

Powered Seating Thru-drive Control Options			
Assignable Buttons (Requires R-net)			

"Assignable Button" feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

<input type="checkbox"/>	195PM01	-	Enable Assignable buttons N/C
--------------------------	---------	---	-------------------------------

Ctrl+5 Switchbox Options (Requires R-net)			
---	--	--	--



Ctrl+5 allows direct control to 5 power seating functions (10 total switches). Switches can also be used with assignable buttons for greater flexibility. Box includes 3-Amp USB port for charging portable devices

<input type="checkbox"/>	195QG01	-	Ctrl+5 Switchbox mounted behind joystick N/C
<input type="checkbox"/>	195QG02	-	Ctrl+5 Switchbox mounted opposite joystick N/C
<input type="checkbox"/>	195QG03	-	Ctrl+5 Switchbox mounted on back cane N/C
<input type="checkbox"/>	195QG05	-	Omit Ctrl+5 -\$500.00

Ctrl+5 Switch Type			
<input type="checkbox"/>	195EY03	-	Push buttons N/C
<input type="checkbox"/>	195EY04	-	Toggle switches \$255.00

## 10. BACKREST AND BACK OPTIONS

Backrest Frame Type			
<input type="checkbox"/>	195BT8	-	Angle Adjustable N/C
<input type="checkbox"/>	195BT98	E1225	Manual Recline 40° (not width adjust.) \$449.00

Manual Recline Adjustment Method Option			
<i>Back Height 23". Fits all seat types. NA with Cantilever Arms</i>			
<input type="checkbox"/>	195MN01		Tool N/C
<input type="checkbox"/>	195MN02		Knob (NA with dual post arms.) N/C

Angle Adjustable Backrest Options			
<i>(NA with Manual recline)</i>			
<input type="checkbox"/>	195OD8		Non-Folding Backrest N/C
<input type="checkbox"/>	195OD9		Folding Backrest \$117.00

Back Cane Angle			
<input type="checkbox"/>	195BP40		Straight Back Post N/C
<input type="checkbox"/>	195BP41		8° Bend Back Post (Avail. only w/ Uph or J3 Back) NA with Manual recline \$107.00

Cane Adjustable height			
<i>(Upholstery Height Measured pan to top of upholstery.)</i>			
<input type="checkbox"/>	195B1		Fixed Low (15.5") N/C
<input type="checkbox"/>	195B5		Fixed Med (17.5") N/C
<input type="checkbox"/>	195B9		Fixed Tall (19.5") N/C
<input type="checkbox"/>	195B243		Adjustable Low 14" - 18" \$107.00
<input type="checkbox"/>	195B244		Adjustable High 20" - 26" \$107.00

Towel Bar Depth			
<input type="checkbox"/>	195QT03		Shallow (NA with Upholstery backs.) N/C
<input type="checkbox"/>	195QT04		Deep (Ext. 4" from back canes for deep backs) N/C

Back Type (Upholstery or Solid)			
<input type="checkbox"/>	195BU1	-	Upholstery (Not avail Adj cane) N/C
<input type="checkbox"/>	195BU70	E2611	Tension adj Upholstery Ballistic Nylon. w/1" width adjustment. \$332.00
<input type="checkbox"/>	195BU71	E2611	Tension adj Upholstery 3DX Vented. w/1" width adjustment. \$347.00
<input type="checkbox"/>	195BU76	-	S.P.O.T. Curved N/C
<input type="checkbox"/>	195BU77	E2613	JAY SPO Posterior \$500.00
<input type="checkbox"/>	195BU14	-	Omit Upholstery -\$15.00

JAY J3 Option (select Omit above)			
<input type="checkbox"/>	195JB08	-	JAY J3 On Chair Order Choose back, headrest, laterals in section 17a,b,c NC




Back Height - Required for SPOT Curved and SPO Backs.			
<i>(Limited combinations of widths) Measured from pan to top of solid back.</i>			
<input type="checkbox"/>	195B132	-	14" N/C
<input type="checkbox"/>	195B36	-	16" N/C
<input type="checkbox"/>	195B38	-	18" N/C
<input type="checkbox"/>	195B40	-	20" N/C
<input type="checkbox"/>	195B42	-	22" N/C
<input type="checkbox"/>	195B44	-	24" N/C
<input type="checkbox"/>	195B140	-	26" N/C

## 11. HEADREST OPTIONS

Choose One Each.

Whitmyer Headrest Mount			
<input type="checkbox"/>	195MH10	E1028	Axys Mount & Bracket \$235.00
<input type="checkbox"/>	195MH11	E1028	Cobra Mount and Bracket \$296.00
<input type="checkbox"/>	195MH1	E1028	Linx Mount and Bracket \$235.00

Whitmyer Headrest Pad			
<input type="checkbox"/>	195HE65	E0955	Plush Head Support Pad 8" \$235.00
<input type="checkbox"/>	195HE66	E0955	Plush Head Support Pad 10" \$235.00
<input type="checkbox"/>	195HE78	E0955	Plush Head Support Pad 14" \$235.00
<input type="checkbox"/>	195HE82	E0955	Adjusta Plush Head Support Pad Narrow \$235.00
<input type="checkbox"/>	195HE83	E0955	Adjusta Plush Head Support Pad std \$235.00



## 12. POSTURAL SUPPORTS

LATERAL THORACIC SUPPORT (Swing Away)			
<i>Available on all Solid Backs.</i>			
<input type="checkbox"/>	<b>LAT 4X4 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$821.00</b>
	195ZV14	E0956	LATERAL 4X4 W/ 1" OFFSET - PAIR \$265.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR \$556.00
<input type="checkbox"/>	<b>LAT 4X5 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$821.00</b>
	195ZV15	E0956	LATERAL 4X5 W/ 1" OFFSET - PAIR \$265.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR \$556.00
<input type="checkbox"/>	<b>LAT 4X6 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$821.00</b>
	195ZV16	E0956	LATERAL 4X6 W/ 1" OFFSET - PAIR \$265.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR \$556.00
<input type="checkbox"/>	<b>LAT 6X4 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$821.00</b>
	195ZV17	E0956	LATERAL 6X4 W/ 1" OFFSET - PAIR \$265.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR \$556.00
<input type="checkbox"/>	<b>LAT 6X5 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$821.00</b>
	195ZV18	E0956	LATERAL 6X5 W/ 1" OFFSET - PAIR \$265.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR \$556.00
<input type="checkbox"/>	<b>LAT 6X6 S/A MOUNTS W/ 1" OFFSET - PR</b>		<b>\$821.00</b>
	195ZV19	E0956	LATERAL 6X6 W/ 1" OFFSET - PAIR \$265.00
	195ZS07	E1028	LATERAL SWING AWAY MOUNT - PAIR \$556.00

Swing Away Knee Adductor Supports			
<i>Available only with Fixed and Power Center Mount Legrests.</i>			
<i>Not available with Dual Post Arm Rests or 2-point Locking Cantilever Arms</i>			
<b>Knee Adductor Supports - LEFT (choose Pad and Arm)</b>			<b>\$414.00</b>
<b>Pad Size</b>			
<input type="checkbox"/>	195KB5	E0956	Knee Adductor Pad SM- Left \$179.00
<input type="checkbox"/>	195KB7	E0956	Knee Adductor Pad LG- Left \$179.00
<b>S/A Arm Size</b>			
<input type="checkbox"/>	195XX04	E1028	Swing Away Short - Left \$235.00
<input type="checkbox"/>	195XX06	E1028	Swing Away Long - Left \$235.00
<b>Knee Adductor Supports - RIGHT (choose Pad and Arm)</b>			<b>\$414.00</b>
<b>Pad Size</b>			
<input type="checkbox"/>	195KB6	E0956	Knee Adductor Pad SM- Right \$179.00
<input type="checkbox"/>	195KB8	E0956	Knee Adductor Pad LG- Right \$179.00
<b>S/A Arm Size</b>			
<input type="checkbox"/>	195XX05	E1028	Swing Away Short - Right \$235.00
<input type="checkbox"/>	195XX07	E1028	Swing Away Long - Right \$235.00

CHEST STRAP			
<input type="checkbox"/>	195ZR03	E0960	Padded Chest Strap \$235.00

Long Thigh Support w/Removable Mount			
<i>Available with Cantilever Arms, NA with Manual Recline Back, Dual Post, Single Post or Cantilever w/ 2 pt Locking arms.</i>			
<input type="checkbox"/>	<b>SPOT Long Thigh Support - Left</b>		<b>\$342.00</b>
	195TV2	E0956	Long Thigh Support - Left \$107.00
	195TX1	E1028	Removable Mount - Left \$235.00
<input type="checkbox"/>	<b>SPOT Long Thigh Support - Right</b>		<b>\$342.00</b>
	195TV1	E0956	Long Thigh Support - Right \$107.00
	195TX2	E1028	Removable Mount - Right \$235.00

## 13. ARMREST AND ARMPAD OPTIONS

Armrest Mounting			
Dual Post Armrests			
<input type="checkbox"/>	195EG01		Short Dual Post Flip Back 10" Height. N/C
<input type="checkbox"/>	195EG21		Long Dual Post Flip Back 10" Height. N/C
<input type="checkbox"/>	195EG02	E0973	Shrt Dual Post Flip Back HT Adj 10"-15" \$235.00
<input type="checkbox"/>	195EG22	E0973	Lng Dual Post Flip Back HT Adj 10"-15" \$235.00
Cantilever Armrests			
<i>(Adjustable length arm pad included (choose Pad setting below, except for Low Arm Option) NA with Manual Recline Back</i>			
<input type="checkbox"/>	195EG03	E0973	<b>Standard Cantilever Arm</b> w/ 2 Point Locking Post. (9"-10" Ht Range.) NA w/ seat depths below 16" \$245.00
<input type="checkbox"/>	195EG04	E0973	<b>Standard Cantilever Arm</b> w/ 2 Point Locking Post. (11"-12" Ht Range.) NA w/ seat depths below 16" \$245.00
<input type="checkbox"/>	195EG05	E0973	<b>STD Cantilever Arm</b> 9"-12" Ht Range \$245.00
<input type="checkbox"/>	195EG06	E0973	<b>Low Cantilever Arm</b> 7"-10" Ht Range (Choose pad type and size below) \$245.00
Single Post Armrests			
<input type="checkbox"/>	195EG07	E0973	Single Post Height Adjustable - Std. \$245.00
<input type="checkbox"/>	195EG08	E0973	Single Post Height Adjustable - Low \$245.00

Arm Pad Options			
<b>Pad Type / Size Options for Dual Post / Single Post / Low Cantilever</b>			
<b>Left Arm Pad (N/C)</b>		<b>Right Arm Pad (N/C)</b>	
<input type="checkbox"/>	195AE01	-	Desk Length Waterfall (10") <input type="checkbox"/> 195AB01
<input type="checkbox"/>	195AE02	-	Full Waterfall (14") NA w/Low Cant <input type="checkbox"/> 195AB02
<input type="checkbox"/>	195AE03	-	Desk Length Standard (10") <input type="checkbox"/> 195AB03
<input type="checkbox"/>	195AE04	-	Full Standard (14") NA w/Low Cant <input type="checkbox"/> 195AB04
<b>Pad Setting Options for Standard Cantilever Arm</b>			
<b>Left Arm Pad (N/C)</b>		<b>Right Arm Pad (N/C)</b>	
<input type="checkbox"/>	195AE11	-	Full Length Cantilever (14") <input type="checkbox"/> 195AB11
<input type="checkbox"/>	195AE10	-	Desk Length Cantilever (10") <input type="checkbox"/> 195AB10



## 14. HANGERS AND FOOTRESTS

### Swing-away Hangers

#### Non-Powered Swing Away Legrests

*Sold as Pairs*

<input type="checkbox"/>	195H2		70° Swing Away (B)	N/C
<input type="checkbox"/>	195H6		90° Swing Away (B) (5.5-9")	N/C
<input type="checkbox"/>	195H138	E0990	Manual ELR S/A (B) (Avail. with Fixed Centermount, choose below)	\$286.00

#### Power Elevating Legrests (2 Actuators)

<input type="checkbox"/>	195H126	E1010	Swing Away Power ELR w/ Articulation Pair (B footplates) (N/A for 12" or 13" wide)	\$2,315.00
--------------------------	---------	-------	--	------------

#### Swing Away Footplates

<input type="checkbox"/>			<b>B - Adult Angle Adjust w/Heel Loop for use with 70° hanger and ELR (Pair)</b>	<b>\$224.00</b>
	195F282	K0040	Adult Angle Adjustable Footplate	\$173.00
	195HL10	E0951	Heel Loop	\$51.00
<input type="checkbox"/>			<b>B - Adult Angle Adjust w/Heel Loop for use with 90° hanger (Pair)</b>	<b>\$224.00</b>
	195F285	K0040	Adult Angle Adjustable Footplate	\$173.00
	195HL10	E0951	Heel Loop	\$51.00
<input type="checkbox"/>			<b>C - Kids Angle Adjust w/Leg Strap for use with 70° hanger (Pair)</b>	<b>\$209.00</b>
	195F283	K0040	Kids Angle Adjustable Footplate	\$173.00
	195LS13	K0038	Leg Strap	\$36.00
<input type="checkbox"/>			<b>C - Kids Angle Adjust w/Leg Strap for use with 80° and 90° hanger(Pair)</b>	<b>\$209.00</b>
	195F286	K0040	Kids Angle Adjustable Footplate	\$173.00
	195LS13	K0038	Leg Strap	\$36.00
<input type="checkbox"/>	195F279		Omit Footplate and extension tube	N/C

#### Omit Hangers.

<input type="checkbox"/>	195H144		Omit Hangers	N/C
--------------------------	---------	--	--------------	-----

#### Hanger Extensions. (See range listed in matrix below)

*NOTE: Not all extension tube will be available; depending on the seat height, some extension tube options may not offer enough ground clearance.*

<input type="checkbox"/>	195E16	K0037	High Mount Hanger	\$56.00
<input type="checkbox"/>	195E1		Short	N/C
<input type="checkbox"/>	195E2		Medium	N/C
<input type="checkbox"/>	195E5		Long	N/C

#### Footrest Accessories

<input type="checkbox"/>	1955A01	E0952	Toe Loop (Pair) Only avail on (opt B).	\$51.00
<input type="checkbox"/>	195LS14	K0038	Extra Leg Strap	\$36.00

### Centermount Legrests

#### Non-Powered Center Mount Legrests

##### Fixed Center Mount Legrests (includes footplate)

**Plus options below**

<input type="checkbox"/>	195H134	-	Fixed Center Mount Short	N/C
<input type="checkbox"/>	195H135	-	Fixed Center Mount Long	N/C

#### Angle Adjustable Center Mount Legrests

<input type="checkbox"/>	195H104	-	Angle adjustable Center Mount Medium 3" - 11.25" Choose Footplate (I or J) infinite and Independent leg length and angle adjustable	\$321.00
<input type="checkbox"/>	195H108	-	Angle adjustable Center Mount Long 3" - 13.25" Choose Footplate (I or J) infinite and Independent leg length and angle adjustable	\$321.00

#### Power Centermount (2 Actuators)

<input type="checkbox"/>	195H127	K0108	Power Centermount ELR w/7.5" Power Articulation (F) Centermount footplates w/ programming options below	\$2,366.00
--------------------------	---------	-------	---	------------

#### Power Centermount Footplate Extension Option (For Thru-Drive Operation)

<input type="checkbox"/>	195IE01	K0108	<b>Enable</b> 7.5" Powered footplate extension for transfer assist. (Electronics programmed to allow footplates to be extended <b>with or without</b> legrest being raised)	\$1,275.00
<input type="checkbox"/>	195IE02	-	<b>Disable</b> Footplate extension for transfers. (Electronics programmed to <b>ONLY</b> extend footplates when legs are raised.)	N/C

#### Centermount Footplates

<input type="checkbox"/>			<b>F - Centermount Dual Footplate (Pair)</b> (Individually Height (3"-19") and Angle Adjustable)	<b>\$173.00</b>
	195F315	K0040	CMT Angle Adjustable Footplates	\$173.00
	195HL14	E0951	Padded and Adjustable Heel Block	N/C
<input type="checkbox"/>			<b>I - Angle adj Centermount Footplate pair Small</b>	<b>\$173.00</b>
	195F125	K0040	Small Angle Adjustable Footplates 5.75" x 6"	\$173.00
<input type="checkbox"/>			<b>J - Angle adj Centermount Footplate pair Large</b>	<b>\$173.00</b>
	195F126	K0040	Large Angle Adjustable Footplates 5.75" x 10"	\$173.00

Lower Leg Length Matrix (measured seat pan to the top of the footplate surface)						Footplate Matrix			
Hanger Type	Standard	High Mount	Short	Medium	Long	B	F	I	
70° Swing Away		7"-14"	13.5"- 16"	16"-18.5"	18"-20.5"	X			
90° Swing Away		5.5-9" Footplate mounts to front of hanger.						X	
Manual ELR S/A			12.5"-17"	14.5"-19"	17"-21"	X			
Power ELR S/A			12"-14"	14"-17"	16.5"-19"	X			
Fixed Centermount			12.5"-15.5"		14.5"-19.5"	Single Footplate			
Angle Adj Centermount				3"-13"	3"-15"			X	
Power Centermount	3"-19"						X		

## 15. ACCESSORIES

#### Wheelchair Tray Table

<input type="checkbox"/>	1957A01	E0950	TRAY TABLE SM TUFFAK 12"-14" PWR	\$235.00
<input type="checkbox"/>	1957A02	E0950	TRAY TABLE MED TUFFAK 15"-17" PWR	\$235.00
<input type="checkbox"/>	1957A03	E0950	TRAY TABLE LRG TUFFAK 18"-20" PWR	\$235.00

#### Backpack

<input type="checkbox"/>	195A1		Black Utility	\$56.00
--------------------------	-------	--	---------------	---------

#### Accessories

<input type="checkbox"/>	1958A09		Ram Phone Mount on Round JS tube	\$51.00
<input type="checkbox"/>	1958A10		Ram Tablet Mount on Round JS tube	\$204.00
<input type="checkbox"/>	1958A11		Ram Cup holder	\$51.00
<input type="checkbox"/>	1958A12		Switch-It Remote Stop Switch	\$913.00

#### Positioning Belt

<input type="checkbox"/>	195EJ01		Auto Style Buckle	N/C
<input type="checkbox"/>	195EJ02	E0978	Belt 2" Aircraft Buckle	\$41.00
<input type="checkbox"/>	195EJ03	E0978	Belt 2" Aircraft Buckle Padded	\$71.00

#### Mounting Accessories

*Only one accessory can be mounted at a time.*

<input type="checkbox"/>	195A603	E2207	Cane/Crutcher Holder	\$265.00
<input type="checkbox"/>	195A604	E2208	O2 holder	\$265.00
<input type="checkbox"/>	195A605	K0108	Walker holder	\$265.00



# JAY<sup>®</sup> Signature Series Backrest and Accessories

Note: Not all sizes and configurations shown below, for full offering visit [www.sunrisemedical.com](http://www.sunrisemedical.com)

## 17a. JAY BACKREST

Seat width must can be equal to or 2" less than selection in section 7.


### JAY J3 Back PA (2.2" Contour Depth)

UT Shells- Heavy Duty TH 2 Point Hardware

<input type="checkbox"/>	195JT01	E2613	J3PAUTS16TH	16" Wide J3 Back w/ 16.5" Shell Height	\$592.00
<input type="checkbox"/>	195JT02	E2613	J3PAUTS18TH	18" Wide J3 Back w/ 16.5" Shell Height	\$592.00
<input type="checkbox"/>	195JT03	E2613	J3PAUTS20TH	20" Wide J3 Back w/ 16.5" Shell Height	\$592.00
<input type="checkbox"/>	195JT04	E2613	J3PAUTM16TH	16" Wide J3 Back w/ 18" Shell Height	\$592.00
<input type="checkbox"/>	195JT05	E2613	J3PAUTM18TH	18" Wide J3 Back w/ 18" Shell Height	\$592.00
<input type="checkbox"/>	195JT06	E2613	J3PAUTM20TH	20" Wide J3 Back w/ 18" Shell Height	\$592.00
<input type="checkbox"/>	195JT07	E2614	J3PAUTM22TH	22" Wide J3 Back w/ 18" Shell Height	\$811.00
<input type="checkbox"/>	195JT08	E2613	J3PAUTT16TH	16" Wide J3 Back w/ 19.5" Shell Height	\$592.00
<input type="checkbox"/>	195JT09	E2613	J3PAUTT18TH	18" Wide J3 Back w/ 19.5" Shell Height	\$592.00
<input type="checkbox"/>	195JT10	E2613	J3PAUTT20TH	20" Wide J3 Back w/ 19.5" Shell Height	\$592.00

SH Shells- TH 2 Point Hardware

<input type="checkbox"/>	195JT11	E2613	J3PASHS18TH	18" Wide J3 Back w/ 21" Shell Height	\$592.00
<input type="checkbox"/>	195JT12	E2613	J3PASHS20TH	20" Wide J3 Back w/ 21" Shell Height	\$592.00
<input type="checkbox"/>	195JT13	E2614	J3PASHS22TH	22" Wide J3 Back w/ 21" Shell Height	\$811.00




### JAY J3 Back PD (6" Contour Depth)

UT Shells- Heavy Duty TH 2 Point Hardware

<input type="checkbox"/>	195JT14	E2620	J3PDUTM16TH	16" Wide J3 Back w/ 18" Shell Height	\$765.00
<input type="checkbox"/>	195JT15	E2620	J3PDUTM18TH	18" Wide J3 Back w/ 18" Shell Height	\$765.00
<input type="checkbox"/>	195JT16	E2620	J3PDUTM20TH	20" Wide J3 Back w/ 18" Shell Height	\$765.00

SH Shells- Heavy Duty TH 2 Point Hardware

<input type="checkbox"/>	195JT17	E2620	J3PDSHS16TH	16" Wide J3 Back w/ 21" Shell Height	\$765.00
<input type="checkbox"/>	195JT18	E2620	J3PDSHS18TH	18" Wide J3 Back w/ 21" Shell Height	\$765.00
<input type="checkbox"/>	195JT19	E2620	J3PDSHS20TH	20" Wide J3 Back w/ 21" Shell Height	\$765.00




### JAY J3 Back PL (3" Contour Depth)

UT Shells- Heavy Duty TH 2 Point Hardware

<input type="checkbox"/>	195JT20	E2613	J3PLUTS16TH	16" Wide J3 Back w/ 16.5" Shell Height	\$683.00
<input type="checkbox"/>	195JT21	E2613	J3PLUTS18TH	18" Wide J3 Back w/ 16.5" Shell Height	\$683.00
<input type="checkbox"/>	195JT22	E2613	J3PLUTS20TH	20" Wide J3 Back w/ 16.5" Shell Height	\$683.00
<input type="checkbox"/>	195JT23	E2613	J3PLUTM16TH	16" Wide J3 Back w/ 18" Shell Height	\$683.00
<input type="checkbox"/>	195JT24	E2613	J3PLUTM18TH	18" Wide J3 Back w/ 18" Shell Height	\$683.00
<input type="checkbox"/>	195JT25	E2613	J3PLUTM20TH	20" Wide J3 Back w/ 18" Shell Height	\$683.00
<input type="checkbox"/>	195JT26	E2613	J3PLUTT16TH	16" Wide J3 Back w/ 19.5" Shell Height	\$683.00
<input type="checkbox"/>	195JT27	E2613	J3PLUTT18TH	18" Wide J3 Back w/ 19.5" Shell Height	\$683.00
<input type="checkbox"/>	195JT28	E2613	J3PLUTT20TH	20" Wide J3 Back w/ 19.5" Shell Height	\$683.00

SH Shells- TH 2 Point Hardware

<input type="checkbox"/>	195JT29	E2613	J3PLSHS18TH	18" Wide J3 Back w/ 21" Shell Height	\$683.00
<input type="checkbox"/>	195JT30	E2613	J3PLSHS20TH	20" Wide J3 Back w/ 21" Shell Height	\$683.00



### JAY Back

To order any **JAY Back not listed above**, please see the JAY price list. Cushion name or part number must be filled out w/ desired width and depth

<input type="checkbox"/>	195JB01	JAY Bundled Back	Back Name or Part #	Width:	Height:	Price List
--------------------------	---------	------------------	---------------------	--------	---------	------------

## 17b. JAY BACKREST ACCESSORIES


### JAY J3 Back Headrest

*AXYS Mount and Pad (includes removable mount)*

<input type="checkbox"/>	195HD01	E0955/E1028	J3HEDX10	AXYS Headrest with 10" Plush Pad	\$389.00
<input type="checkbox"/>	195HD02	E0955/E1028	J3HEDX14	AXYS Headrest with 14" Plush Pad	\$402.00
<input type="checkbox"/>	195HD05	E0955/E1028	J3HEDXAAPS	AXYS Headrest with Narrow Adjust-a-Plush	\$490.00
<input type="checkbox"/>	195HD06	E0955/E1028	J3HEDXAAPN	AXYS Headrest with Standard Adjust-a-Plush	\$490.00

*LINX Mount and Pad (includes removable mount)*

<input type="checkbox"/>	195HD03	E0955/E1028	J3HEDL10	LINX Headrest with 10" Plush Pad	\$450.00
<input type="checkbox"/>	195HD04	E0955/E1028	J3HEDL14	LINX Headrest with 14" Plush Pad	\$463.00
<input type="checkbox"/>	195HD07	E0955/E1028	J3HEDLAAPS	LINX Headrest with Narrow Adjust-a-Plush	\$551.00
<input type="checkbox"/>	195HD08	E0955/E1028	J3HEDLAAPN	LINX Headrest with Standard Adjust-a-Plush	\$551.00



### JAY J3 Back Laterals

*Swing Away Laterals With Pad - Pair*

<input type="checkbox"/>	195JW02	E0956/E1028	J3LAT4X4SA	SwingAway Lat with 4"x4" Pad (PR)	\$794.00
<input type="checkbox"/>	195JW03	E0956/E1028	J3LAT4X5SA	SwingAway Lat with 4"x5" Pad (PR)	\$794.00
<input type="checkbox"/>	195JW04	E0956/E1028	J3LAT4X6SA	SwingAway Lat with 4"x6" Pad (PR)	\$794.00

## 17c. JAY CUSHIONS

*Seat width must match selection in section 7. Seat depth can be +/- 2" of selection in section 7*

### JAY Fusion Cushion

<input type="checkbox"/>	195JU01	E2622	JFUSION1616	JAY Fusion 16" x 16"	\$449.00
<input type="checkbox"/>	195JU02	E2622	JFUSION1618	JAY Fusion 16" x 18"	\$449.00
<input type="checkbox"/>	195JU03	E2622	JFUSION1816	JAY Fusion 18" x 16"	\$449.00
<input type="checkbox"/>	195JU04	E2622	JFUSION1818	JAY Fusion 18" x 18"	\$449.00
<input type="checkbox"/>	195JU05	E2622	JFUSION1820	JAY Fusion 18" x 20"	\$449.00
<input type="checkbox"/>	195JU06	E2622	JFUSION2018	JAY Fusion 20" x 18"	\$449.00
<input type="checkbox"/>	195JU07	E2622	JFUSION2020	JAY Fusion 20" x 20"	\$449.00

### JAY X2 Cushion

<input type="checkbox"/>	195JU22	E2622	JX21616	JAY X2 16" x 16"	\$434.00
<input type="checkbox"/>	195JU23	E2622	JX21618	JAY X2 16" x 18"	\$434.00
<input type="checkbox"/>	195JU24	E2622	JX21816	JAY X2 18" x 16"	\$434.00
<input type="checkbox"/>	195JU25	E2622	JX21818	JAY X2 18" x 18"	\$434.00
<input type="checkbox"/>	195JU26	E2622	JX21820	JAY X2 18" x 20"	\$434.00
<input type="checkbox"/>	195JU27	E2622	JX22020	JAY X2 20" x 20"	\$434.00

### JAY Union Cushion

<input type="checkbox"/>	195JU15	E2607	JU1616	JAY UNION 16" x 16"	\$454.00
<input type="checkbox"/>	195JU16	E2607	JU1618	JAY UNION 16" x 18"	\$454.00
<input type="checkbox"/>	195JU17	E2607	JU1816	JAY UNION 18" x 16"	\$454.00
<input type="checkbox"/>	195JU18	E2607	JU1818	JAY UNION 18" x 18"	\$454.00
<input type="checkbox"/>	195JU19	E2607	JU1820	JAY UNION 18" x 20"	\$454.00
<input type="checkbox"/>	195JU20	E2607	JU2018	JAY UNION 20" x 18"	\$454.00
<input type="checkbox"/>	195JU21	E2607	JU2020	JAY UNION 20" x 20"	\$454.00

### JAY Ion Cushion

<input type="checkbox"/>	195JU08	E2607	JION1616	JAY ION 16" x 16"	\$403.00
<input type="checkbox"/>	195JU09	E2607	JION1618	JAY ION 16" x 18"	\$403.00
<input type="checkbox"/>	195JU10	E2607	JION1816	JAY ION 18" x 16"	\$403.00
<input type="checkbox"/>	195JU11	E2607	JION1818	JAY ION 18" x 18"	\$403.00
<input type="checkbox"/>	195JU12	E2607	JION1820	JAY ION 18" x 20"	\$403.00
<input type="checkbox"/>	195JU13	E2607	JION2018	JAY ION 20" x 18"	\$403.00
<input type="checkbox"/>	195JU14	E2607	JION2020	JAY ION 20" x 20"	\$403.00

### JAY Cushion

*To order any **JAY cushions not listed above**, please see the JAY price list. Cushion name or part number must be filled out w/ desired width and depth*

<input type="checkbox"/>	195JC20	JAY Bundled Cushion	Cushion Name or Part #	<input style="width: 100px;" type="text"/>	Width:	<input style="width: 50px;" type="text"/>	Depth:	<input style="width: 50px;" type="text"/>	Price List
--------------------------	---------	---------------------	------------------------	--	--------	---	--------	---	------------

## 16. Personalized Options

Looking for more options? Contact our "built 4 me" department.

Below are common upgrades that utilize Ready Made Quotes. These options can be quickly quoted to add a personal touch to your Pulse power chair.

- Lights and Indicators

*Some Restrictions Apply*



To Inquire about any unique personalizing of Quickie products, contact:  
 Phone: (800) 333-4000, Option #3 or Fax: (800) 333-9011



**Sunrise Medical (US) LLC**  
 2842 Business Park Ave.  
 Fresno, CA 93727 · USA