



ORDER FORM

To see QUICKIE Power Options
or Accessories, scan here or visit:
[https://www.sunrisemedical.com/
power-wheelchair-options-accessories](https://www.sunrisemedical.com/power-wheelchair-options-accessories)



Q700M


with SEDEO® LITE™

EIPW25

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		State	Zip
Mark For		Phone	

Client Measurements

a) Top of head to seat		in
b) Top of shoulder to seat		in
c) Underarm to seat		in
d) Elbow to seat		in
e) Top of leg to seat		in
f) Shoulder width		in
g) Chest width		in
i) Hip width		in
j) Actual sitting depth		in
m) Foot length		in
n) Knee to foot		in
o) Seat width		in



Notes

Chair Specifications

HCPC	K0848 / K0850 / K0856 / K0858 / K0868 / K0877
Weights Starting at: (with Batteries)	NPC: 294 lbs. SPC: 318 lbs.
Power Seating	50° of CG Tilt 12" seat elevator
Seat Type	SEDEO® LITE™
Seat Size Range	Width 15" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs. 350 lbs. HD
Base Size	Width 25", Length 36.5"
Minimum Turning Radius	19.8" (excludes footplates)
Drive Control	R-net PM 120A / R-net PM 90A
Driving Range	Up to 18 miles (w/ Group 24 Batteries)
Drive Motors	4-pole 6 mph, 5 mph HD, 8 mph
Drive Wheels	14" w/flat-free inserts (black)
Casters	6" front / 6" rear (black)
Charger	10A Dual Mode (Compact)
Battery Size	Group 24 / Group 34
Color	9 shroud colors and 7 color accents
Maximum Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Mid Wheel Drive Base End Item Number EIPW25 (LITE™ Seat)			\$
☐ 250PC267	K0848	Quickie Q700 M NPC No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$9,159.00
☐ 250PC268	K0856	Quickie Q700 M SPC Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	\$9,160.00
Group 3 Heavy Duty Mid Wheel Drive Base End Item Number EIPW25 (LITE™ Seat)			\$
☐ 250PC269	K0850	Quickie Q700 M NPC HD No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$11,784.00
☐ 250PC270	K0858	Quickie Q700 M SPC HD Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 90amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 5mph, Weight Capacity = 350lbs., STFH = 16.5".	\$11,784.00
Group 4 Standard Duty Mid Wheel Drive Base End Item Number EIPW25 (LITE™ Seat)			\$
☐ 250PC271	K0868	Quickie Q700 M NPD No Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$14,583.00
☐ 250PC272	K0877	Quickie Q700 M SPD Single Power Option for SEDEO® LITE™ Seat Included in Base Pricing : Q700 M Base, R-net 120amp Electronics for No Power Seat Functions, ProGrip™ Six-Wheel Independent Suspension, 14" Drive Wheels with Dual Density Inserts, 6" Flat Free Casters, ISO Occupied Tie-Down Brackets and 10amp Charger. Max Speed = 8mph, Weight Capacity = 300lbs., STFH = 16.5".	\$14,583.00

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
User Weight			\$
☐ 250WG23	–	Up to 300 lbs.	NC
☐ 250WG15	–	301 to 350 lbs.	NC
Motor Package			\$
☐ 250PF98	–	4-pole 8 kph / 5 mph <i>(Standard for Group 3 HD K0850, K0858)</i>	NC
☐ 250PF99	–	4-pole 10 kph / 6 mph <i>(Standard for Group 3 K0848, K0856)</i>	NC
☐ 250PF100	–	4-pole 12.5 kph / 8 mph <i>(Standard for Group 4 K0868, K0877)</i>	NC
☐ 250PF101	–	4-pole 12.5 kph / 8 mph <i>(Upgrade for Group 3 K0848, K0856)</i>	\$1,342.00
Tracking			\$
☐ 250MI01	K0108	SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i>	\$1,628.00
Drive Wheels			\$
☐ 250DW58	–	Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			\$
☐ 250CW55	–	Caster Wheels 6", Solid, Black	NC
Shroud Color			\$
☐ 250C119	–	Matte Black	NC
☐ 250C125	–	Pearl White	NC
☐ 250C165	–	Candy Apple Red	NC
☐ 250C166	–	Candy Blue	NC
☐ 250C167	–	Candy Purple	NC
☐ 250C176	–	Green Apple	NC
☐ 250C178	–	Pearl Pink	NC
☐ 250C182	–	Silver	NC
☐ 250C170	–	American Flag	\$602.00
Color Accents			\$
<i>Note: Color Accents and Shroud Colors may not match.</i>			
☐ 250EZ20	–	Apple Green	NC
☐ 250EZ22	–	Black	NC
☐ 250EZ21	–	Blue	NC
☐ 250EZ23	–	Orange	NC
☐ 250EZ24	–	Red	NC
☐ 250EZ19	–	White	NC
☐ 250EZ25	–	Silver	NC

BASE				
OPTION	HCPC	DESCRIPTION	RETAIL	
ISO Transit Tie Downs				\$
☐ 250T050	–	ISO Transit Package	NC	
LED Lights				\$
☐ 250LI08	–	LED Fender Lights	\$566.00	
☐ 250LI02	–	Omit LED Lights	NC	
Seat Heights				\$

Note: Available seat-to-floor heights are dependent on 2 options: **1- Battery Size** and **2- Powered Seating**.
Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.

		Seat Height	Power Seating	No Power Seating	
☐ 250HT103	–	16.5"	Group 34	Group 34	NC
☐ 250HT55	–	17"	Group 34	GP24 or GP34	NC
☐ 250HT104	–	17.5"	Group 34	GP24 or GP34	NC
☐ 250HT56	–	18"	GP24 or GP34	GP24 or GP34	NC
☐ 250HT105	–	18.5"	GP24 or GP34	GP24 or GP34	NC
☐ 250HT57	–	19"	GP24 or GP34	GP24 or GP34	NC

Batteries and Chargers

Group 24 Battery Options (Recommended for max daily range)

☐ 250AG25	E2363	GP24 Flat Top Pair w/Strap MK Gel	\$1,138.00
☐ 250AG26	E2363	GP24 Flat Top Pair w/Strap AGM	\$857.00

Group 34 Battery Options (Required for lower seat-to-floor heights)

☐ 250AG27	E2359	GP34 Flat Top Pair w/Straps MK Gel	\$1,097.00
☐ 250AG29	E2359	GP34 Flat Top Pair w/Straps AGM	\$828.00

Omit Batteries

☐ 250AG24	–	Battery Harnesses Only (Omit Batteries)	NC
----------------------------------	---	---	----

Chargers

☐ 250BJ28	–	10A Dual Mode (Compact)	NC
----------------------------------	---	-------------------------	----

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating Options			\$
Seat Elevator			
☐ 250PI08	E2298	C-Me® Seat Elevator, Powered, 300 mm / 12" (Not available with Group 4) <i>C-Me speed up to 3 mph.</i>	\$3,791.00
Seat Tilt			
☐ 250PA35	E1002	Posterior 50° Tilt	\$8,458.00
☐ 250PA36	E1002	Posterior 50° HD Tilt	\$9,626.00
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)			\$
☐ 250FA10	—	0° (Adds 0" to STFH)	NC
☐ 250FA3	—	3° (Adds 1" to STFH)	NC
SEDEO® LITE™ Backrest Frame (Angle Adjustable 86° - 106°)			\$
Back Frame Width			
☐ 250BF21	—	Width Adjustable, 385 - 485 mm / 15 - 19"	NC
☐ 250BF22	—	Width Adjustable, 460 - 510 mm / 18 - 20"	NC
Back Cane Height			
☐ 250B263	—	Standard, 18" - 20"	NC
☐ 250B264	—	Tall, 21"-23"	NC
☐ 250PH30	—	Height Adjustable, Short w/Trailing Push Handles, 20" to 22"	NC
☐ 250PH31	—	Height Adjustable, Tall w/Trailing Push Handles, 23" to 26"	NC
Backrest Frame Option			\$
☐ 250BF23	—	Flip Down, Quick Release	\$147.00
☐ 250BF24	—	Removeable	NC
Back Cushion Type			\$
☐ 250BT148	—	SEDEO® LITE™ Sling Back with Tension Adjustable Straps (Choose Seating in Section A)	NC
☐ 250BT135	—	JAY® J3® (Choose Seating in Section B)	NC
☐ 250BT133	—	OMIT (For personalized mounting. No headrest mount)	NC

A - SEDEO® LITE™ Sling Back Dimensions

Installed Back Height - (From Pan to Top of Back)			
☐ 250JR71	—	18"	NC
☐ 250JR72	—	19"	NC
☐ 250JR73	—	20"	NC
☐ 250JR74	—	21"	NC
☐ 250JR75	—	22"	NC
☐ 250JR76	—	23"	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
B - JAY® J3® Back			
B - JAY® J3® Back - Contour			\$
☐ 250JP01	E2613	J3® PA Back for 16"-20" Wide	\$625.00
☐ 250JP11	E2615	J3® PL Back for 16"-20" Wide	\$729.00
☐ 250JP13	E2620	J3® PD Back for 16"-20" Wide <i>(18" and 21" shell height only)</i>	\$875.00
B - JAY® J3® Back Mounting Hardware			\$
☐ 250JP20	—	TH	NC
B - JAY® J3® Back Width			\$
<i>Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.</i>			
☐ 250JQ02	—	16"	NC
☐ 250JQ03	—	18"	NC
☐ 250JQ04	—	20"	NC
B - JAY® J3® Back Shell Options and Height Ranges			\$
JAY® J3® Shell Height			
☐ 250JR14	—	16.5" tall shell <i>(Upper Thoracic Short)</i>	NC
☐ 250JR15	—	18" tall shell <i>(Upper Thoracic Med)</i>	NC
☐ 250JR16	—	19.5" tall shell <i>(Upper Thoracic Tall)</i>	NC
☐ 250JR17	—	21" tall shell <i>(Shoulder Height)</i>	NC
☐ 250JR27	—	22.5" tall shell <i>(Shoulder Height)</i>	NC
☐ 250JR28	—	24" tall shell <i>(Shoulder Height)</i>	NC
JAY® J3® Installed Back Height (From Pan to Top of Back)			
☐ 250JR29	—	21"	NC
☐ 250JR19	—	22"	NC
☐ 250JR20	—	23"	NC
☐ 250JR21	—	24"	NC
☐ 250JR22	—	25"	NC
☐ 250JR23	—	26"	NC
☐ 250JR24	—	27"	NC
☐ 250JR25	—	28"	NC
☐ 250JR30	—	29"	NC
☐ 250JR31	—	30"	NC

POWERED SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

B - JAY® Anterior Thoracic Support

\$

Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)

☐	250K10	E0960	Medium (ATS8302)	\$177.00
☐	250K11	E0960	Large (ATS8303)	\$177.00

B - JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery, & Mount)

\$

JAY® Lateral Thoracic Support Pads

☐	250JW35	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Right	\$147.00
☐	250JW38	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4"x 5" Pad - Left	\$147.00
☐	250JW36	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$147.00
☐	250JW39	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$147.00
☐	250JW37	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$147.00
☐	250JW40	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$147.00

JAY® Lateral Thoracic Support Pad Upholstery

☐	250NU01	—	Syncra™ Premium Leatherette	NC
☐	250NU02	—	Syncra™ Stretch	NC

JAY® Lateral Thoracic Swing-Away Mount

☐	250JW41	E1034	Lateral Thoracic Swing-Away Bracket - Right	\$292.00
☐	250JW42	E1034	Lateral Thoracic Swing-Away Bracket - Left	\$292.00

Headrests (Requires Pad, Upholstery, & Mount)

\$

Headrest Pad

☐	250HE36	E0955	SEDEO® Pad Small With Cover (LxH) 240 x 150 mm / 10" x 6" <i>Can only be used with Multi Position Headrest Mount with Quick Release and SEDEO® Upholstery</i>	\$188.00
☐	250HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$239.00
☐	250HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$239.00
☐	250HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$250.00
☐	250HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (<i>Adjustable from 7 to 16" Wide x 4" Tall</i>)	\$367.00
☐	250HE48	—	OMIT Headrest Pad	NC

Headrest Pad Upholstery

☐	250NU10	—	SEDEO® Upholstery	NC
☐	250NU11	—	Whitmyer® Lycra	NC

Headrest Mount

☐	250MH32	E1033	SEDEO® Multi-Position Headrest Mount with Quick Release (<i>Only available with SEDEO® Pad</i>)	\$234.00
☐	250MH33	E1033	Whitmyer® Linx Headrest Mount with Quick Release	\$192.00
☐	250MH34	E1033	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$234.00
☐	250MH35	E1033	Whitmyer® AXYS Headrest Mount with Quick Release	\$123.00
☐	250MH17	—	OMIT Headrest Mount	NC

POWERED SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Seat Cushion Type			\$
☐ 250ST80	–	JAY® (See Section C)	NC
☐ 250ST75	–	Omit	NC
B - Seat Size			\$
Seat Pan Width			
☐ 250R78	–	15"	NC
☐ 250R53	–	16"	NC
☐ 250R54	–	17"	NC
☐ 250R55	–	18"	NC
☐ 250R56	–	19"	NC
☐ 250R57	–	20"	NC
Seat Pan Depth (Front of seat pan to front of back cushion – effective depth)			
☐ 250RF51	–	16"	NC
☐ 250RF45	–	17"	NC
☐ 250RF46	–	18"	NC
☐ 250RF47	–	19"	NC
☐ 250RF48	–	20"	NC
☐ 250RF49	–	21"	NC
☐ 250RF50	–	22"	NC
C - JAY® Cushions			
C - JAY® Seat Cushion			\$
☐ 250CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$460.00
☐ 250CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$520.00
☐ 250CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$514.00
☐ 250CU104	E2622	JAY® Fusion® with Cryo® Technology (Sizes 16" - 21" Wide)	\$630.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

Knee and Thigh Swing-Away Support Pads

\$

Lateral Knee Support Pads and Hardware *(must pick pad and mount)*

☐	250XY05	E0953	Lateral Knee Support Pad, Right (4" round)	\$205.00
☐	250XY06	E0953	Lateral Knee Support Pad, Left (4" round)	\$205.00
☐	250XX08	E1028	Lateral Swing-Away Bracket Knee - Right	\$268.00
☐	250XX09	E1028	Lateral Swing-Away Bracket Knee - Left	\$268.00

Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*

☐	250KF22	E0953	Lateral Thigh Support Pad Short, Right (5"x5")	\$147.00
☐	250KF23	E0953	Lateral Thigh Support Pad Short, Left (5"x5")	\$147.00
☐	250KF24	E0953	Lateral Thigh Support Pad Long, Right (5"x7")	\$147.00
☐	250KF25	E0953	Lateral Thigh Support Pad Long, Left (5"x7")	\$147.00
☐	250KF26	E1028	Lateral Thigh Removable Mount - Right	\$292.00
☐	250KF27	E1028	Lateral Thigh Removable Mount - Left	\$292.00

Armrest Structure

\$

☐	250EG36	—	Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 13" x .625")	\$292.00
☐	250EG34	—	Single Post, SEDEO LITE, Removable (N/A with thigh support), (Height Adjustable 10" to 15")	NC
☐	250EG50	—	Dual Post, Flip Back, Removable (Height Adjustable)	\$225.00

Arm Pads

\$

SEDEO® Armrest Pad

☐	250AB62	—	300 mm / 12", Integral Foam, Right	\$58.00
☐	250AE62	—	300 mm / 12", Integral Foam, Left	\$58.00
☐	250AB63	—	400 mm / 16", Integral Foam, Right	\$58.00
☐	250AE63	—	400 mm / 16", Integral Foam, Left	\$58.00
☐	250AB61	—	SEDEO LITE Pad, Right	NC
☐	250AE61	—	SEDEO LITE Pad, Left	NC

Comfort Arm Pad (Foam)

☐	250AB46	—	2" Wide Right	\$88.00
☐	250AE46	—	2" Wide Left	\$88.00

Comfort Arm Pad (Gel)

☐	250AB48	—	2" Wide Right	\$175.00
☐	250AE48	—	2" Wide Left	\$175.00

Ottobock® Arm Troughs

☐	250AE06	—	Ottobock Arm Trough, Med (12"), Left	\$136.00
☐	250AB06	—	Ottobock Arm Trough, Med (12"), Right	\$136.00
☐	250AE07	—	Ottobock Arm Trough, Large (14"), Left	\$136.00
☐	250AB07	—	Ottobock Arm Trough, Large (14"), Right	\$136.00
☐	250OL01	—	Ottobock Hand Pad, Med (7"), Left	\$59.00
☐	250OR01	—	Ottobock Hand Pad, Med (7"), Right	\$59.00
☐	250OL02	—	Ottobock Hand Pad, Large (9"), Left	\$59.00
☐	250OR02	—	Ottobock Hand Pad, Large (9"), Right	\$59.00
☐	250AE78	—	Ottobock One-Piece Arm Trough (18"), Left	\$190.00
☐	250AB78	—	Ottobock One-Piece Arm Trough (18"), Right	\$190.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Elbow Support			\$
☐ 2503A02	–	Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>	\$179.00
☐ 2503A03	–	Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>	\$179.00
☐ 2503A06	–	Elbow Support Pad, RH 3" x 3"	\$149.00
☐ 2503A07	–	Elbow Support Pad, LH 3" x 3"	\$149.00
☐ 2503A04	E1028	Swing-Away Elbow Support Hardware, RH	\$245.00
☐ 2503A05	E1028	Swing-Away Elbow Support Hardware, LH	\$245.00

Fixed Centermount Legrest			\$
☐ 250H183	–	SEDEO® Fixed Centermount Legrest <i>(3" to 17.5" Leg Length)</i>	\$292.00
☐ 250H184	–	Fixed Centermount Legrest, Short <i>(Includes Calf Pads & Single Footplate) (11.5" to 14.5" Leg Length)</i>	NC
☐ 250H185	–	Fixed Centermount Legrest, Long <i>(Includes Calf Pads & Single Footplate) (13.5" to 18.5" Leg Length)</i>	NC
☐ 250H136	–	Angle Adjustable Centermount Legrest <i>(11.5" to 18.5" Leg Length)</i>	\$525.00

Power Centermount Legrest			\$
☐ 250H175	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Mechanical Articulation) (5" to 17.5" Leg Length)</i>	\$2,706.00
☐ 250H109	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i>	\$3,791.00

Power Centermount Legrest Accessories			\$
--	--	--	-----------

SEDEO Centermount Footplates

☐ 250F127	K0040	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$197.00
☐ 250F128	K0040	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$197.00

Footplate Accessories (Available with Power Centermounts)

☐ 250DI68	K0108	Extra Large Centermount Footplate <i>(Single Footplate Overlay 16"W x 12"D)</i>	\$199.00
☐ 250DI69	K0108	Footplate Connector Link <i>(Attaches footplates together)</i>	\$97.00
☐ 250DI57	K0108	Gel Footpads <i>(Pair, LH and RH)</i>	\$403.00

Angle Adjustable Centermount Footplates

☐ 250F45	K0040	Single / Dual Footplate	\$197.00
☐ 250F44	K0040	Single Platform Footplate	\$197.00

SEDEO Centermount Calf Pads

☐ 250DI34	–	Calf Pads, Pair, Standard Integral Foam	NC
☐ 250DI35	–	Calf Pads, Pair, Padded	\$408.00

Angle Adjustable Centermount Calf Pads

☐ 250DI49	–	Calf Pads, Pair, Standard Integral Foam	NC
☐ 250DI50	–	Extra Large Single Calf Pad	\$197.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
SEDEO® PRO™ Swing-Away Legrest (Pair)			\$
☐ 250H177	–	Legrests, Swing-Away, 70°, SEDEO® PRO™ (14" to 21" Leg Length)	NC
☐ 250H178	E0990	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)	\$408.00
☐ 250H179	E1010	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ (Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)	\$2,508.00
☐ 250DI76	–	Omit Legrests	NC
Swing-Away Legrest Footplates			
☐ 250F369	K0040	Angle Adjustable Flip-Up Footplates	\$197.00
Swing-Away Legrest Calf Supports			
☐ 250DI42	–	Calf Support, Depth Adjustable, Padded (6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)	\$460.00
Swing-Away Legrest Options			
☐ 250DI44	E0951	Heel Loops	\$58.00
☐ 250KG4	–	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	\$234.00
SEDEO® Lower Extremity Swing-Away Support			
☐ 250VD01	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	\$460.00
☐ 250VD02	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	\$460.00
Pelvic Positioning Belts			
☐ 250A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐ 250A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
☐ 250A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$88.00
☐ 250A845	K0108	Bodypoint – Evoflex Pelvic Stabilizer Belt	\$203.00
☐ 250A846	K0108	Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	\$224.00
Accessories			\$
Ventilator Option			
☐ 250VT12	E1029	Vent Tray Fixed (Pulmonetic LTV Models)	\$542.00
O2 Tank Option			
☐ 250A718	E2208	O2 Tank Holder Mounted to Backrest (D and E Styles)	\$379.00
Wheelchair Tray Table			
☐ 2507A09	E0950	Tray, Swing-Away	\$408.00
Backpack			
☐ 250A1	–	Black Utility	\$63.00
Medial Thigh Support			
☐ 2508A80	E0957	Medial Thigh Support	\$197.00
☐ 2508A81	E1028	Pommel Swing-Away Hardware	\$292.00
Misc. Accessories			
☐ 2508A59	–	Transfer Bar (Pair)	\$408.00
☐ 2508A52	–	Ram Mount Phone	\$100.00
☐ 2508A53	–	Ram Mount Tablet	\$408.00
☐ 2508A54	–	Ram Mount Cup Holder	\$100.00
☐ 2508A58	–	Rear View Mirror	\$100.00
☐ 250A611	E2207	Crutch Holder	\$147.00
☐ 2508A73	–	Bag for Armrest	\$100.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
R-net Joystick <i>(Remote)</i>			\$
☐ 250J121	–	Standard LED Joystick <i>(w/o Screen Display)</i>	NC
☐ 250J116	–	Color Joystick <i>(2.3" Color LCD Screen Display, CJS1M1)</i>	NC
☐ 250J118	–	High Resolution Color Joystick,(CJS1M2+BT) <i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,230.00
☐ 250J119	–	Omit Remote Control <i>(Required with arm mounted specialty controls)</i>	NC
Joystick Mounting			\$
☐ 250JM34	–	Fixed on Armrest	NC
☐ 250JM35	E1032	Parallel Swing-Away	\$345.00
☐ 250JM36	E1032	Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i>	\$991.00
Joystick Mount (Location)			\$
☐ 250HM16	–	Right Hand	NC
☐ 250HM17	–	Left Hand	NC
Display Language Selection R-net			\$
☐ 250L1	–	English	NC
☐ 250L3	–	French	NC
☐ 250L7	–	Spanish	NC
Joystick Handle (Optional)			\$
☐ 250JH13	E2323	Ball Handle	\$147.00
☐ 250JH15	E2323	T-Handle	\$147.00
☐ 250JH24	E2323	Bodypoint - Rubber Dome Joystick Handle	\$191.00
☐ 250JH25	E2323	Bodypoint - U-Shaped Joystick Handle <i>(4 inch)</i>	\$191.00
☐ 250JH26	E2323	Bodypoint - Grooved Mushroom Joystick Handle	\$191.00
Electronic Programmers			\$
☐ 250PR13	–	R-net PC Wireless Programmer	\$613.00
R-net Power Modules			\$
☐ 250ER11	–	Power Module, R-net, 120A, Non-Expandable <i>(Only available with Group 4)</i>	NC
☐ 250ER10	–	Power Module, R-net, 120A, Expandable <i>(Only available with Group 4)</i> <i>Required for specialty controls or 3+ actuators</i>	\$1,454.00
250EX09	E2377	Expandable 120A Controller - R-net	\$764.00
250EX10	E2313	Harness for 120A Expandable Controller	\$690.00
☐ 250ER18	–	Power Module, R-net, 90A, Non-Expandable	NC
☐ 250ER17	–	Power Module, R-net, 90A, Expandable <i>Required for specialty controls or 3+ actuators</i>	\$1,454.00
250EX15	E2377	Expandable 90A Controller - R-net	\$764.00
250EX16	E2313	Harness for 90A Expandable Controller	\$690.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
Seating Control			\$
Powered Seating Control			
☐ 250E040	–	Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
☐ 250E037	E2310	Thru-Drive Control for 1 Actuator	\$2,240.00
☐ 250E038	E2311	Thru-Drive Control for 2+ Actuators	\$4,055.00
Switch Control <i>(Available with Omit Thru-Drive only) (1 actuator for Tilt or Elevate only)</i>			
☐ 250E033	–	Dual Toggle for 1 Actuator	NC
☐ 250E034	–	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
☐ 250E041	–	Ctrl+5 Switch Box for 2 Actuators or More	NC
Additional Power Seating Control <i>(Available with Thru-Drive only)</i>			
☐ 250QG11	–	Ctrl+5 Switch Box with Voice Activation	\$1,459.00
Dual Toggle Mounting Location			
☐ 250TN1	–	Toggle Mount, Right	NC
☐ 250TN2	–	Toggle Mount, Left	NC
Ctrl+5 Mounting Location <i>(Will swing away if mounted in same location as joystick)</i>			
☐ 250EY19	–	Right Armrest	NC
☐ 250EY20	–	Left Armrest	NC
☐ 250EY21	–	Back Cane, Right	\$175.00
☐ 250EY24	–	Back Cane, Left	\$175.00
Ctrl+5 Switch Type			
☐ 250EY22	–	Push Buttons	NC
☐ 250EY23	K0108	Toggle Switches	\$350.00
Attendant Controls			\$
☐ 2500A05	E2331	Dual Control Module for Attendant	\$1,517.00
☐ 2500A06	E2331	R-net Color LCD (CJSM1), Attendant Control Only	\$1,517.00
Attendant Mounting Location			\$
☐ 250DN5	–	Right Back Cane	NC
☐ 250DN6	–	Left Back Cane	NC
R-net Programming			\$
Assignable Buttons			
<i>Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.</i>			
☐ 250PM12	–	Include Assignable Buttons	NC
Power Centermount Footplate Extension Option			
☐ 250IE03	K0108	Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.</i>	\$1,459.00



Q700 M with SEDEO® LITE™		\$
Base Price		
Total of additional equipment		
Packaging and transport		
Total Price		

Signature

Send E-Mail