



BREAKTHROUGH IN COMFORT

EXCEPTIONAL COMFORT

One of the first ever wheelchair backs designed to help manage an individual's microclimate. Internal thermal technology gives the user the ability to warm up or cool down for maximum comfort.

PROVEN DESIGN

JAY J3 wheelchair backs are designed with contoured back shells, soft foam, and adjustable components. The JAY HV (heating/ventilation) technology is seamlessly incorporated into this proven backrest.

LEADING TECHNOLOGY

Based on customer feedback, the HV Back design uses JAY's industry-leading seating technology in partnership with a leader in automotive thermal management.



REDEFINING BACK TECHNOLOGY

HV

Superior Clinical Seating

EXCEPTIONAL COMFORT

The patented internal thermal technology of the JAY J3 HV Back provides exceptional comfort.

- Evaporative cooling to dissipate moisture and pull it away from the body
- Heat system for warmth, equipped with automatic shut-off

JAY Cover Technology

The J3 HV Back cover optimizes airflow through the use of 3DX™ spacer fabric and specially designed vents.

SAFE IN ALL ENVIRONMENTS

This product safely works in all environments and includes:

- The technology comes from our partnership with a company offering years of experience in the automotive industry
- JAY has then combined that technology with rigorous wheelchair industry testing standards

EASY TO USE

Battery Powered

- Easy access battery pack concealed on the back shell
- Rechargeable battery provides 8 hours of continued use from a full charge
- Stand-alone battery charger

Simple Control

- Simple on/off switch controls the heat or fan setting
- Red light indicates heat, blue light indicates fan



J3 HV BACK SPECS

HCPCS Code

HV Insert: NOT CODED
J3 PA Back: E2613
J3 PL Back: E2615
J3 PD Back: E2620

Contour Styles (depth)

PA: 2"
PL: 3"
PD: 6"

Back Widths

16", 18", 20"

Back Heights

Mid Thoracic Med: 13.5"
Shoulder Height Std: 21"

Hardware Options

SH
TH
4 Point

Cover Material

3DX™ Spacer Fabric
and Lycra

JAY Your Way

Customizations Available



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-129437 Rev. B ©06.2017

800.333.4000 www.SunriseMedical.com

