



# Quickie Sports & Recreation

Each of us possesses a spirit which is revealed in our STRENGTH, our ENERGY, our PASSION, our ABILITY, our CHOICE, our LIBERTY and our STYLE. What makes us all alike is the presence of a spirit. What makes us all unique is how we reveal it. Quickie Sports and Recreation is dedicated to providing innovative products to help you reveal your spirit.



# Quickie Sports & Recreation



	C H A I R S			
<b>Model</b>	<b>Ti All Court</b>	<b>All Court</b>	<b>Ti Match Point</b>	<b>Match Point</b>
<b>Material</b>	Titanium	Aluminum	Titanium	Aluminum
<b>Product Weight</b>	21.5 lbs	24 lbs	20 lbs	22.5 lbs
<b>Width</b>	12"-20"	12"-20"	12"-20"	12"-20"
<b>Depth/Length</b>	10"-20"	10"-20"	10"-18"	10"-18"
<b>Height</b>	16"-21"	16"-21"	16"-21"	16"-21"
<b>Weight Capacity</b>	250 lbs	250 lbs	250 lbs	250 lbs
<b>Camber Options</b>	16° or 20°	16° or 20°	16° or 20°	16° or 20°
<b>Front Seat Height</b>	16"-21"	16"-21"	16"-21"	16"-21"
<b>Rear Seat Height</b>	13"-21"	13"-21"	13"-21"	13"-21"
<b># of Color Options*</b>	20	20	20	20
<b>Anti-Tip Style</b>	Single (std), Double	Single (std), Double	Single	Single
<b>Support</b>	Anti-tip brace	Anti-tip brace	Anti-tip brace	Anti-tip brace
<b>Chest Straps</b>	2", 4"	2", 4"	2", 4"	2", 4"
<b>Wheel Type</b>	Spinergy, Heat Treated, Cross Spoke	Spinergy, Heat Treated, Cross Spoke	Spinergy, Heat Treated, Cross Spoke	Spinergy, Heat Treated, Cross Spoke
<b>Wheel Size</b>	24", 26"	24", 26"	24", 26"	24", 26"
<b>Key Option</b>	Wing w/ standup	Wing w/ standup	Extra wide 'footprint'	Extra wide 'footprint'
<b>Backrest</b>	Folding, Non-Folding	Folding, Non-Folding	Folding, Non-Folding	Folding, Non-Folding

\* Refer to the Quickie Color Selection Brochure or website for color options.



	S H A R K
<b>Product Weight</b>	33 lbs
<b>Seat Width</b>	14.8", 16.4", 17.9"
<b>Depth/Length</b>	40.2"
<b>Height</b>	8.6"
<b>Frame Width</b>	15", 16.5", 18"
<b>Weight Capacity</b>	265 lbs
<b>Rear angle adjust</b>	30°-90°, 7 Basic, 3 fine settings
<b>Crank Length/Height</b>	6.89"
<b>Gear Speeds</b>	7", 27"
<b>Gears</b>	Sram S7™, Sram™5.0
<b>Brake</b>	Shimano-Tiagra™
<b>Back</b>	Folding, adjustable

## Shark Color Selecton

