



# JAY<sup>®</sup> Zip<sup>™</sup> Cushion

E2607

January 2022

**Mark For:**

Date: \_\_\_\_\_  
 Dealer Acct #: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Confirm Via:  Fax  Email

**Submitting for:**

Quote  Order

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**JAY ZIP CUSHION - E2607**

**JAY Zip Cushion**

<input type="checkbox"/>	JZXXYY	E2607	JAY Zip Cushion 8-16 wide	\$437.00
--------------------------	--------	-------	---------------------------	----------

Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)  
 Includes: Ionic+™ fabric outer cover, inner cover, and foam base

**JAY Zip**

2 Outer Covers included with each Cushion!




**JAY Zip Size Matrix**

		WIDTH												
		8	9	10	11	12	13	14	15	16				
DEPTH	8					N/A	N/A	N/A	N/A	N/A				
	9						N/A	N/A	N/A	N/A				
	10							N/A	N/A	N/A				
	11								N/A	N/A				
	12									N/A	N/A			
	13											N/A		
	14												N/A	
	15	N/A												
	16	N/A	N/A											
17	N/A	N/A	N/A											
18	N/A	N/A	N/A	N/A										

	=	Standard Offering
	=	Offered Via Jay Your Way
	=	Not available

**JAY Zip Measuring Guide**

Measurements made from completed cushion with covers.



Bottom View

**Width**

X = Width

Note: Measurement made from bottom outside edge to outside edge of cushion.  
 Ex: JZ1214 X = 12"

**Depth**

Y = Depth

Note: Measurement made from bottom rear edge to front lower edge of cushion.  
 Ex: JZ1214 Y = 14"

**Outer Cover**

The JAY Zip outer cover is designed to promote comfort, convenience and breathability by the use of "Ionic+™" material and 3DX Spacer Fabric.



**Inner Cover**

The water resistant inner cover, combined with the Aquaguard™ Zipper and anti wicking thread provides the ultimate protection against odor.



**Note: To order via JAY Your Way add a "M" to the front of the part number and follow the prompts. For additional modifications see Page 2 or visit [www.sunrisemedical.com](http://www.sunrisemedical.com) to view even more!**

Please visit [www.sunparts.us](http://www.sunparts.us) for spare parts information.

JAY Zip Cushion									
✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JZ0808	E2607	JAY ZIP 8" X 8"	\$437.00	<input type="checkbox"/>	JZ1212	E2607	JAY ZIP 12" X 12"	\$437.00
<input type="checkbox"/>	JZ0809	E2607	JAY ZIP 8" X 9"	\$437.00	<input type="checkbox"/>	JZ1213	E2607	JAY ZIP 12" X 13"	\$437.00
<input type="checkbox"/>	JZ0810	E2607	JAY ZIP 8" X 10"	\$437.00	<input type="checkbox"/>	JZ1214	E2607	JAY ZIP 12" X 14"	\$437.00
<input type="checkbox"/>	JZ0811	E2607	JAY ZIP 8" X 11"	\$437.00	<input type="checkbox"/>	JZ1215	E2607	JAY ZIP 12" X 15"	\$437.00
<input type="checkbox"/>	JZ0812	E2607	JAY ZIP 8" X 12"	\$437.00	<input type="checkbox"/>	JZ1216	E2607	JAY ZIP 12" X 16"	\$437.00
<input type="checkbox"/>	JZ0813	E2607	JAY ZIP 8" X 13"	\$437.00	<input type="checkbox"/>	JZ1217	E2607	JAY ZIP 12" X 17"	\$437.00
<input type="checkbox"/>	JZ0814	E2607	JAY ZIP 8" X 14"	\$437.00	<input type="checkbox"/>	JZ1218	E2607	JAY ZIP 12" X 18"	\$437.00
<input type="checkbox"/>	JZ0908	E2607	JAY ZIP 9" X 8"	\$437.00	<input type="checkbox"/>	JZ1310	E2607	JAY ZIP 13" X 10"	\$437.00
<input type="checkbox"/>	JZ0909	E2607	JAY ZIP 9" X 9"	\$437.00	<input type="checkbox"/>	JZ1311	E2607	JAY ZIP 13" X 11"	\$437.00
<input type="checkbox"/>	JZ0910	E2607	JAY ZIP 9" X 10"	\$437.00	<input type="checkbox"/>	JZ1312	E2607	JAY ZIP 13" X 12"	\$437.00
<input type="checkbox"/>	JZ0911	E2607	JAY ZIP 9" X 11"	\$437.00	<input type="checkbox"/>	JZ1313	E2607	JAY ZIP 13" X 13"	\$437.00
<input type="checkbox"/>	JZ0912	E2607	JAY ZIP 9" X 12"	\$437.00	<input type="checkbox"/>	JZ1314	E2607	JAY ZIP 13" X 14"	\$437.00
<input type="checkbox"/>	JZ0913	E2607	JAY ZIP 9" X 13"	\$437.00	<input type="checkbox"/>	JZ1315	E2607	JAY ZIP 13" X 15"	\$437.00
<input type="checkbox"/>	JZ0914	E2607	JAY ZIP 9" X 14"	\$437.00	<input type="checkbox"/>	JZ1316	E2607	JAY ZIP 13" X 16"	\$437.00
<input type="checkbox"/>	JZ0915	E2607	JAY ZIP 9" X 15"	\$437.00	<input type="checkbox"/>	JZ1317	E2607	JAY ZIP 13" X 17"	\$437.00
<input type="checkbox"/>	JZ1008	E2607	JAY ZIP 10" X 8"	\$437.00	<input type="checkbox"/>	JZ1318	E2607	JAY ZIP 13" X 18"	\$437.00
<input type="checkbox"/>	JZ1009	E2607	JAY ZIP 10" X 9"	\$437.00	<input type="checkbox"/>	JZ1411	E2607	JAY ZIP 14" X 11"	\$437.00
<input type="checkbox"/>	JZ1010	E2607	JAY ZIP 10" X 10"	\$437.00	<input type="checkbox"/>	JZ1412	E2607	JAY ZIP 14" X 12"	\$437.00
<input type="checkbox"/>	JZ1011	E2607	JAY ZIP 10" X 11"	\$437.00	<input type="checkbox"/>	JZ1413	E2607	JAY ZIP 14" X 13"	\$437.00
<input type="checkbox"/>	JZ1012	E2607	JAY ZIP 10" X 12"	\$437.00	<input type="checkbox"/>	JZ1414	E2607	JAY ZIP 14" X 14"	\$437.00
<input type="checkbox"/>	JZ1013	E2607	JAY ZIP 10" X 13"	\$437.00	<input type="checkbox"/>	JZ1415	E2607	JAY ZIP 14" X 15"	\$437.00
<input type="checkbox"/>	JZ1014	E2607	JAY ZIP 10" X 14"	\$437.00	<input type="checkbox"/>	JZ1416	E2607	JAY ZIP 14" X 16"	\$437.00
<input type="checkbox"/>	JZ1015	E2607	JAY ZIP 10" X 15"	\$437.00	<input type="checkbox"/>	JZ1417	E2607	JAY ZIP 14" X 17"	\$437.00
<input type="checkbox"/>	JZ1016	E2607	JAY ZIP 10" X 16"	\$437.00	<input type="checkbox"/>	JZ1418	E2607	JAY ZIP 14" X 18"	\$437.00
<input type="checkbox"/>	JZ1108	E2607	JAY ZIP 11" X 8"	\$437.00	<input type="checkbox"/>	JZ1512	E2607	JAY ZIP 15" X 12"	\$437.00
<input type="checkbox"/>	JZ1109	E2607	JAY ZIP 11" X 9"	\$437.00	<input type="checkbox"/>	JZ1513	E2607	JAY ZIP 15" X 13"	\$437.00
<input type="checkbox"/>	JZ1110	E2607	JAY ZIP 11" X 10"	\$437.00	<input type="checkbox"/>	JZ1514	E2607	JAY ZIP 15" X 14"	\$437.00
<input type="checkbox"/>	JZ1111	E2607	JAY ZIP 11" X 11"	\$437.00	<input type="checkbox"/>	JZ1515	E2607	JAY ZIP 15" X 15"	\$437.00
<input type="checkbox"/>	JZ1112	E2607	JAY ZIP 11" X 12"	\$437.00	<input type="checkbox"/>	JZ1516	E2607	JAY ZIP 15" X 16"	\$437.00
<input type="checkbox"/>	JZ1113	E2607	JAY ZIP 11" X 13"	\$437.00	<input type="checkbox"/>	JZ1517	E2607	JAY ZIP 15" X 17"	\$437.00
<input type="checkbox"/>	JZ1114	E2607	JAY ZIP 11" X 14"	\$437.00	<input type="checkbox"/>	JZ1518	E2607	JAY ZIP 15" X 18"	\$437.00
<input type="checkbox"/>	JZ1115	E2607	JAY ZIP 11" X 15"	\$437.00	<input type="checkbox"/>	JZ1613	E2607	JAY ZIP 16" X 13"	\$437.00
<input type="checkbox"/>	JZ1116	E2607	JAY ZIP 11" X 16"	\$437.00	<input type="checkbox"/>	JZ1614	E2607	JAY ZIP 16" X 14"	\$437.00
<input type="checkbox"/>	JZ1117	E2607	JAY ZIP 11" X 17"	\$437.00	<input type="checkbox"/>	JZ1615	E2607	JAY ZIP 16" X 15"	\$437.00
<input type="checkbox"/>	JZ1209	E2607	JAY ZIP 12" X 9"	\$437.00	<input type="checkbox"/>	JZ1616	E2607	JAY ZIP 16" X 16"	\$437.00
<input type="checkbox"/>	JZ1210	E2607	JAY ZIP 12" X 10"	\$437.00	<input type="checkbox"/>	JZ1617	E2607	JAY ZIP 16" X 17"	\$437.00
<input type="checkbox"/>	JZ1211	E2607	JAY ZIP 12" X 11"	\$437.00	<input type="checkbox"/>	JZ1618	E2607	JAY ZIP 16" X 18"	\$437.00

**JAY Your Way - Modifications**

REDUCED HEIGHT			
<input type="checkbox"/>	RH	Remove 1/2" of foam from bottom	NC

NOTCH FOR BACK POST			
<input type="checkbox"/>	NP	Notch for Back Post	NC
Specify distance from back of cushion:			
<input type="checkbox"/>	1"	<input type="checkbox"/>	2"
<input type="checkbox"/>	3"		

NOTCH FOR SEAT RAILS			
<input type="checkbox"/>	NR	Notch for Seat Rails	NC

COVER OPTIONS			
<input type="checkbox"/>	CVR	Cover Options	NC
<input type="checkbox"/>	Dartex	<input type="checkbox"/>	Incontinent
<input type="checkbox"/>	Microclimatic		

ANTERIOR OR POSTERIOR SLOPE			
<input type="checkbox"/>	WGP	Posterior Slope - Wedge	NC
Wedge glued underneath cushion (Cover sewn to match)			
Select desired height increase at front of cushion			
<input type="checkbox"/>	.5"	<input type="checkbox"/>	1"
<input type="checkbox"/>	1.5"		
<input type="checkbox"/>	WGA	Anterior Slope - Wedge	NC
Wedge glued underneath cushion (Cover sewn to match)			
Select desired height increase at rear of cushion			
<input type="checkbox"/>	.5"	<input type="checkbox"/>	1"
<input type="checkbox"/>	1.5"		

**Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.**

**Notes:**

