

Magic Mobility Frontier V6

March 2022



Order Form

| | |
|-------------|--|
| Account: | |
| Date: | |
| PO#: | |
| Buyer Name: | |
| Marked For: | |

| | |
|-------------------|--|
| Ship To: | |
| Name: | |
| Address: | |
| City, State, Zip: | |
| Phone: | |


Please submit order to: magicmobility@sunmed.com

Sunrise Medical (US) LLC
2842 Business Park
Fresno, CA 93727
1-800-333-4000



Client Measurements

| | |
|----------------------------|------|
| a) Top of head to seat | Inch |
| b) Top of shoulder to seat | Inch |
| c) Underarm to seat | Inch |
| d) Elbow to seat | Inch |
| e) Top of leg to seat | Inch |
| f) Shoulder width | Inch |
| g) Chest width | Inch |
| i) Hip width | Inch |
| j) Actual sitting depth | Inch |
| m) Foot length | Inch |
| n) Knee to foot | Inch |
| o) Seat Width | Inch |



Notes:

Seating Specifications

| | | | |
|---------------------------------------|---|-------------------------|--------------------------------|
| HCPC | K0848 / K0856 / K0861 K0850 / K0858 / K0862 K0849 / K0857 / K0851 / K0859 | Seat Type | MPS and Rehab Style |
| Typical WC Weight w/ Batteries | Starting at 320 lbs. Weights vary depending on options | Seat Size Range | Width 13"-23" Depth 13"-24" |
| | | STF Height Range | 17", 18", 19" |
| Power Seating | 50 degrees of CG Tilt / Power Recline / 12" seat elevator | | |

Drive Base Specifications

| | | | |
|---------------------------------|---|-------------------------|---|
| Drive Wheel Position | Mid wheel Drive | Drive Wheels | 14" All Terrain Drive Wheels 14" Hybrid Drive Wheels |
| User Weight Capacity | 400 lbs. 340 lbs. with Elevate or Tilt | Tires/Casters | Pneumatic or Solid inserts |
| | | Charger | Dual Mode |
| V6 Base Size All Terrain | Width 28", Length 41.5", Turning Radius 23.5" | | |
| V6 Base Size Hybrid | Width 27.3", Length 41.5", Turning Radius 23.5" | | |
| Drive Control | R-Net PM 90A controller | Battery Size | Group 24 or 34 |
| Drive Motors | V6 4-pole 6mph | Color Accents | 5 frame colors |
| | | Ground Clearance | 3.5" |

| OPTION | HCPC | Description | Retail |
|---|-------------|---|---------------|
| Seating Options | | | |
| <input type="checkbox"/> | MM-RHBSB | REHAB Seat & Back (Back canes and Adjustable Flat Seat Pan) | N/C |
| <input type="checkbox"/> | MM-MPSSB | MPS Seat & Back (Captain Style Seating) | N/C |
| <input type="checkbox"/> | MM-MPSBRHBS | MPS Back & REHAB Seat (Captain Style Backrest and Adjustable Flat Seat Pan) | N/C |
| Base | | | |
| OPTION | HCPC | Description | Retail |
| Frontier V6 Base End Item Number: EIMM | | | |
| <input type="checkbox"/> | MM-V6NP | K0848 Magic Mobility - Frontier V6 No Power Seating Options | \$14,975.00 |
| <input type="checkbox"/> | MM-V6SP | K0856 Magic Mobility - Frontier V6 Single Power Seating Options | \$14,975.00 |
| <input type="checkbox"/> | MM-V6MP | K0861 Magic Mobility - Frontier V6 Multi Power Seating Options | \$14,975.00 |
| <input type="checkbox"/> | MM-V6NPHD | K0850 Magic Mobility - Frontier V6 No Power Seating Options - HD | \$14,975.00 |
| <input type="checkbox"/> | MM-V6SPHD | K0858 Magic Mobility - Frontier V6 Single Power Seating Options - HD | \$14,975.00 |
| <input type="checkbox"/> | MM-V6MPHD | K0862 Magic Mobility - Frontier V6 Multi Power Seating Options - HD | \$14,975.00 |
| <input type="checkbox"/> | MM-V6NPC | K0849 Magic Mobility - Frontier V6 No Power Seating Options - MPS (Captain) | \$14,975.00 |
| <input type="checkbox"/> | MM-V6SPC | K0857 Magic Mobility - Frontier V6 Single Power Seating Options - MPS (Captain) | \$14,975.00 |
| <input type="checkbox"/> | MM-V6MPC | Magic Mobility - Frontier V6 Multi Power Seating Options - MPS (Captain) | \$14,975.00 |
| <input type="checkbox"/> | MM-V6NPCHD | K0851 Magic Mobility - Frontier V6 No Power Seating Options - HD MPS (Captain) | \$14,975.00 |
| <input type="checkbox"/> | MM-V6SPCHD | K0859 Magic Mobility - Frontier V6 Single Power Seating Options - HD MPS (Captain) | \$14,975.00 |
| <input type="checkbox"/> | MM-V6MPCHD | Magic Mobility - Frontier V6 Multi Power Seating Options - HD MPS (Captain) | \$14,975.00 |
| Frame Color | | | |
| <input type="checkbox"/> | MM-BKFRM | Black Frame Color | N/C |
| <input type="checkbox"/> | MM-RDFRM | Red Frame Color | N/C |
| <input type="checkbox"/> | MM-AGFRM | Mossy Green Frame Color | N/C |
| <input type="checkbox"/> | MM-SBFRM | Royal Blue Frame Color | N/C |
| Batteries | | | |
| <input type="checkbox"/> | MM-BAT24G | E2363 Group 24 Gel Battery installed - Pair | \$888.00 |
| <input type="checkbox"/> | MM-BAT34G | E2359 Group 34 Gel Battery installed - Pair | \$760.00 |
| Battery Charger | | | |
| <input type="checkbox"/> | MM-BATCHG | Battery Charger | N/C |
| Drive Wheels | | | |
| <input type="checkbox"/> | MM-ATDW | All Terrain Drive Wheels and Tires | N/C |
| <input type="checkbox"/> | MM-HPDW | Hybrid Drive Wheels and Tires - Pneumatic | N/C |
| <input type="checkbox"/> | MM-HSDW | Hybrid Drive Wheels and Tires - Solid | \$130.00 |
| Tie Downs | | | |
| <input type="checkbox"/> | MM-V6TK | Crash Tested V6 Tie Down Kit - Always Included | N/C |
| Castor Options | | | |
| <input type="checkbox"/> | MM-CASTP | Pneumatic Castor 2.50-4 - non-marking | N/C |
| <input type="checkbox"/> | MM-CASTS | Foam Filled Castor 2.50-4 - non-marking | \$64.00 |
| Optional Drive Tire Options | | | |
| <input type="checkbox"/> | MM-TS2 | Tire Sealant for 2 Pneumatic Drive Tires | \$32.00 |
| Base Accessories | | | |
| <input type="checkbox"/> | MM-01595 | Jack for Drive Tire | \$171.00 |

Seating

| OPTION | HCPC | Description | Retail | | | | | | | | | | | | | | | | | | | | |
|---|-------------------|--|-------------|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|--------------------------|-----|
| Seating System | | | | | | | | | | | | | | | | | | | | | | | |
| No Power Seat Options | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-NP | No Power Seating System (400 lb. Weight Capacity) | \$3,500.00 | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-NPFO | No Power Seating System with Power Functionality (400 lb. Weight Capacity) | \$4,500.00 | | | | | | | | | | | | | | | | | | | | |
| <i>Note: required for lights, power recline, power legrest, steering lock, or future power seat functions.</i> | | | | | | | | | | | | | | | | | | | | | | | |
| Power Tilt - Utilized without Power Recline | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-TLT E1002 | 50° Power Tilt in Space - 340 lb. Weight Capacity | \$6,474.00 | | | | | | | | | | | | | | | | | | | | |
| Power Seat Elevate - Only available when Power Tilt or Power Tilt & Recline is selected | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-ELVT E2300 | 12" Adjustable Seat Elevator - 340 lb. Weight Capacity | \$3,750.00 | | | | | | | | | | | | | | | | | | | | |
| Power Recline - Utilized without Power Tilt - Requires No Power Seating with Power Functionality | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-PRMPS E1003 | 150° Power Recline for MPS Seating - 400 lb. Weight Capacity | \$4,606.00 | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-PRRHB E1003 | 150° Power Recline w/ Shear Reduction for REHAB Seating - 400 lb. Weight Capacity | \$4,606.00 | | | | | | | | | | | | | | | | | | | | |
| Power Tilt and Recline | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-TLTPRMPS E1006 | 50° Tilt & 150° Power Recline - MPS Seating - 340 lb. Weight Capacity | \$10,224.00 | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-TLTPRRHB E1007 | 50° Tilt & 150° PR w/ Shear Reduction - REHAB Seating - 340 lb. Weight Capacity | \$11,765.00 | | | | | | | | | | | | | | | | | | | | |
| Power Seat Function | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-SFCK E2310 | Single Seat Function Control Kit | \$1,498.00 | | | | | | | | | | | | | | | | | | | | |
| <i>Note: MUST be selected if ordering 1 actuator functions</i> | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-MFCK E2311 | Multiple Seat Function Control Kit | \$3,050.00 | | | | | | | | | | | | | | | | | | | | |
| <i>Note: MUST be selected if ordering 2 or more actuator functions</i> | | | | | | | | | | | | | | | | | | | | | | | |
| Seat To Floor Height | | | | | | | | | | | | | | | | | | | | | | | |
| <i>V6 Seat to Floor Heights start at 17" +/- 0.25"</i> | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-STFH | Lowest Seat To Floor Height | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-STFH1 | Seat To Floor Height +1" | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-STFH2 | Seat To Floor Height +2" | N/C | | | | | | | | | | | | | | | | | | | | |
| Power Seat Options | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-CRTL5JS | Ctrl+5 - 5 Button Box for controlling power seat functions - Includes USB port - Joystick Side | \$803.00 | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-CRTL5OS | Ctrl+5 - 5 Button Box for controlling power seat functions - Includes USB port - Opposite Side | \$803.00 | | | | | | | | | | | | | | | | | | | | |
| Seat Sizes (width x depth) | | | | | | | | | | | | | | | | | | | | | | | |
| MPS (Captain) Seat Cushion (width x depth) | | | | | | | | | | | | | | | | | | | | | | | |
| <i>Note: MPS Backrest can be paired with the REHAB flat pan - If Flat pan is desired omit the MPS cushion and select the desired REHAB flat pan</i> | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01673 | 16" x 16" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01674 | 16" x 18" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01675 | 18" x 16" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01676 | 18" x 18" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01678 | 18" x 20" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01679 | 20" x 18" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01680 | 20" x 20" Vinyl Seat | N/C | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01681 | 22" x 18" Vinyl Seat | \$391.00 | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-01682 | 22" x 20" Vinyl Seat | \$391.00 | | | | | | | | | | | | | | | | | | | | |
| REHAB Flat Pan Seat Sizes - Adjustable in width & depth for most sizes | | | | | | | | | | | | | | | | | | | | | | | |
| REHAB Seat Pan | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-02792STD | Standard Adjustable Flat Seat Pan - 13" to 19" deep | N/C | | | | | | | | | | | | | | | | | | | | |
| <i>Note: Standard seat pan available 13"-23" wide</i> | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-02792DP | Deep Adjustable Flat Seat Pan - 20" to 21" deep | \$391.00 | | | | | | | | | | | | | | | | | | | | |
| <i>Note: Deep seat pan available 13"-23" wide</i> | | | | | | | | | | | | | | | | | | | | | | | |
| Deep Seat Pan Extension Kits | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-05690 | 23" Deep Seat Pan Extension Kit - Requires selection of MM-02792DP | \$690.00 | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | MM-05741 | 24" Deep Seat Pan Extension Kit - Requires selection of MM-02792DP | \$690.00 | | | | | | | | | | | | | | | | | | | | |
| REHAB Seat Width | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 13" | <input type="checkbox"/> | 14" | <input type="checkbox"/> | 15" | <input type="checkbox"/> | 16" | <input type="checkbox"/> | 17" | <input type="checkbox"/> | 18" | <input type="checkbox"/> | 19" | <input type="checkbox"/> | 20" | <input type="checkbox"/> | 21" | <input type="checkbox"/> | 22" | <input type="checkbox"/> | 23" | | |
| REHAB Seat Depth | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 13" | <input type="checkbox"/> | 14" | <input type="checkbox"/> | 15" | <input type="checkbox"/> | 16" | <input type="checkbox"/> | 17" | <input type="checkbox"/> | 18" | <input type="checkbox"/> | 19" | <input type="checkbox"/> | 20" | <input type="checkbox"/> | 21" | <input type="checkbox"/> | 22" | <input type="checkbox"/> | 23" | <input type="checkbox"/> | 24" |

Seatback Options - select one

MPS Backrest Options - Comes standard with Manual Recline

Breathable Cloth MPS Backrest

| | | | |
|--------------------------|----------|--|----------|
| <input type="checkbox"/> | MM-01360 | Cloth Backrest - 16"-20" Wide & 22" Tall | N/C |
| <input type="checkbox"/> | MM-01520 | Cloth Tall Backrest - 16"-20" Wide & 26" Tall | \$471.00 |
| <input type="checkbox"/> | MM-01361 | Cloth Wide Backrest - 22"-23" Wide & 22" Tall | \$284.00 |
| <input type="checkbox"/> | MM-04598 | Cloth Wide & Tall Backrest - 22"-23" Wide & 26" Tall | \$562.00 |

Wipeable Vinyl MPS Backrest

| | | | |
|--------------------------|------------|--|----------|
| <input type="checkbox"/> | MM-01360VI | Vinyl Backrest - 16"-20" Wide & 22" Tall | \$257.00 |
| <input type="checkbox"/> | MM-01520VI | Vinyl Tall Backrest - 16"-20" Wide & 26" Tall | \$674.00 |
| <input type="checkbox"/> | MM-01361VI | Vinyl Wide Backrest - 22"-23" Wide & 22" Tall | \$342.00 |
| <input type="checkbox"/> | MM-04598VI | Vinyl Tall & Wide Backrest - 22"-23" Wide & 26" Tall | \$797.00 |

MPS Back Accessories

| | | | |
|--------------------------|----------|--|----------|
| <input type="checkbox"/> | MM-04784 | Push Handle - Needed for O2 holder, Backpack, Cane holder | \$134.00 |
| <input type="checkbox"/> | MM-02533 | Accessory Mount - Needed for aftermarket headrest and/or attendant control | \$102.00 |

REHAB Seating Backrest Options

Backrest will have the width of the MPS seat or flat pan. Please choose your back height

Backrest will come with 2 point hardware with a 500 lb. capacity

| | | | | |
|--------------------------|--------------------|---|---|----------|
| <input type="checkbox"/> | J3PA | E2613 | JAY J3 PA - Posterior - 2.2" of lateral thoracic support - goes to 20" wide | \$592.00 |
| | Select Back Height | <input type="checkbox"/> 16.5" <input type="checkbox"/> 18" <input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24" | | |
| <input type="checkbox"/> | J3PAW | E2614 | JAY J3 PA - Posterior - 2.2" of lateral thoracic support - goes to 22" wide | \$810.00 |
| | Select Back Height | <input type="checkbox"/> 18" <input type="checkbox"/> 21" | | |
| <input type="checkbox"/> | J3PL | E2615 | JAY J3 PL - Posterior Lateral - 3" of lateral thoracic support - goes to 20" wide | \$683.00 |
| | Select Back Height | <input type="checkbox"/> 16.5" <input type="checkbox"/> 18" <input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24" | | |
| <input type="checkbox"/> | J3PD | E2620 | JAY J3 PD - Posterior Deep - 6" of lateral thoracic support - goes to 20" wide | \$811.00 |
| | Select Back Height | <input type="checkbox"/> 16.5" <input type="checkbox"/> 18" <input type="checkbox"/> 19.5" <input type="checkbox"/> 21" <input type="checkbox"/> 22.5" <input type="checkbox"/> 24" | | |

MPS Backrest Options For use with REHAB Back Canes

| | | | |
|--------------------------|------------|--|------------|
| <input type="checkbox"/> | MM-03996 | Cloth MPS Removable Backrest (Works with REHAB back canes) | \$1,669.00 |
| <input type="checkbox"/> | MM-03996VI | Vinyl MPS Removable Backrest (Works with REHAB back canes) | \$1,669.00 |

REHAB Flip Forward Backrest Option

| | | | |
|--------------------------|-----------|------------------------------------|----------|
| <input type="checkbox"/> | MM-FLPFWD | REHAB Flip Forward Backrest Option | \$250.00 |
|--------------------------|-----------|------------------------------------|----------|

REHAB Seating Cushion Options

| | | | | |
|--------------------------|---------|-------|------------|----------|
| <input type="checkbox"/> | JION | E2607 | JAY Ion | \$405.00 |
| <input type="checkbox"/> | JU | E2607 | JAY Union | \$454.00 |
| <input type="checkbox"/> | JFUSION | E2622 | JAY Fusion | \$450.00 |

Headrest Options - MPS Style

| | | | | |
|--------------------------|------------|-------|--|----------|
| <input type="checkbox"/> | MM-02770 | E0955 | Cloth Headrest - 12" x 4.5" x 6.25" - MPS seat back only | \$348.00 |
| <input type="checkbox"/> | MM-02770VI | E0955 | Vinyl Headrest - 12" x 4.5" x 6.25" - MPS seat back only | \$348.00 |

Headrest Options - REHAB Style

| | | | | |
|--------------------------|----------|-------|---|----------|
| <input type="checkbox"/> | J3PLSH10 | E0955 | WHITMYER PLUSH Headrest Pad - 10" x 5" - requires J3HEDX | \$164.00 |
| <input type="checkbox"/> | J3HEDX | E1028 | WHITMYER AXYS Headrest w/ Removable Hardware - depth adj .5" behind backrest and 7" forward | \$224.00 |

Armrest Style Options

| | | | |
|--------------------------|----------|---|----------|
| <input type="checkbox"/> | MM-ARMR | Removable Armrests - Post style (Attaches to the Seat Pan - Height Adjustable 2") | N/C |
| <input type="checkbox"/> | MM-ARMFR | Flip-Up Removable Armrests - Post style (Attaches to the base-height adj 3" & angle adj 90° - 155°) | \$372.00 |
| <input type="checkbox"/> | MM-ARMC | Cantilever Armrests - Reclining Flip Back (Attaches to the back - Height Adjustable 4" & Angle adj 90°) | \$748.00 |
| <input type="checkbox"/> | MM-ARMC2 | Cantilever 2 Point Quick Release Post (Attaches front of armrest to seat frame) | \$992.00 |

Armrest Pad Options

| | | | |
|--------------------------|-------------|----------------------------------|----------|
| <input type="checkbox"/> | MM-ARMPAD01 | 3" x 12" Desk Length Pads | N/C |
| <input type="checkbox"/> | MM-ARMPAD02 | 4" x 16" Tapered Pads | \$91.00 |
| <input type="checkbox"/> | MM-ARMPAD03 | 3" x 12" Rounded Leatherette Pad | \$246.00 |
| <input type="checkbox"/> | MM-ARMPAD04 | 3" x 14" Rounded Leatherette Pad | \$246.00 |

Lateral Support Options - REHAB Seat Only

| | | | | |
|--------------------------|-----------|-------|--|----------|
| <input type="checkbox"/> | J3SAMNTR | E1028 | Single Swing away Lateral Hardware - Right | \$269.00 |
| <input type="checkbox"/> | J3SAMNTL | E1028 | Single Swing away Lateral Hardware - Left | \$269.00 |
| <input type="checkbox"/> | J3LAT6X4R | E0956 | Single 6x4 Lateral Pad - Right | \$128.00 |
| <input type="checkbox"/> | J3LAT6X4L | E0956 | Single 6x4 Lateral Pad - Left | \$128.00 |

| Positioning Options | | | |
|---|--------------|-------|--|
| Positioning Hardware - Compatible with REHAB Seat pan only | | | |
| <input type="checkbox"/> | MM-FDHIP | E1028 | Flip Down Hip Guide Hardware - fixed to the seat pan near the back \$578.00 |
| <input type="checkbox"/> | MM-RMHIP | E1028 | Removable Hip Guide Hardware - fixed to the seat pan near the back \$470.00 |
| <input type="checkbox"/> | MM-RMKNE | E1028 | Removable Knee Adductors Hardware - fixed to the seat pan near the front \$470.00 |
| <input type="checkbox"/> | MM-FDKNE | E1028 | Flip Down Knee Adductors - fixed to the seat pan near the front \$578.00 |
| Positioning Pad | | | |
| <input type="checkbox"/> | MM-4X6PAD | E0956 | 4" x 6" Positioning Pad \$354.00 |
| <input type="checkbox"/> | MM-4X10PAD | E0956 | 4" x 10" Positioning Pad \$354.00 |
| Positioning Pad Accessories | | | |
| <input type="checkbox"/> | MM-4X6GEL | E0956 | 4" x 6" Gel Cover \$257.00 |
| Chest Harness Options | | | |
| <input type="checkbox"/> | MM-SH120M-B2 | E0960 | Bodypoint Monoflex Medium Chest Strap - pad 3.5" x 18.5" total length 49.5" \$144.00 |
| <input type="checkbox"/> | MM-SH120L-B2 | E0960 | Bodypoint Monoflex Large Chest Strap - pad 3.5" x 21.5" total length 52.5" \$144.00 |
| <input type="checkbox"/> | MM-50-CH110 | E0960 | Healthwares Padded Butterfly Positioning Chest Strap - 7" at Chest \$166.00 |
| <input type="checkbox"/> | MM-SH331L-B1 | E0960 | Bodypoint Stayflex Butterfly Harness - 3.5" at Chest \$225.00 |
| Pelvic Positioning Options select one | | | |
| Pelvic Positioning | | | |
| <input type="checkbox"/> | MM-33140 | E0978 | Therafin Push Button Seat Restraint - 2" x 64" N/C |
| <input type="checkbox"/> | MM-33148 | E0978 | Therafin Extra Large Push Button Seat Restraint - 2" x 80" \$54.00 |
| <input type="checkbox"/> | MM-ABSB | E0978 | Therafin Airline Buckle Seat Restraint - 2" x 80" \$139.00 |
| <input type="checkbox"/> | MM-EB205-L62 | K0108 | Bodypoint Evoflex Pelvic Stabilizer - sitting hip width 14" - 18" requires Hw312-KIT \$102.00 |
| <input type="checkbox"/> | MM-WC19SB6K | E0978 | Crash Tested Seat Belt Kit - 600mm \$893.00 |
| <input type="checkbox"/> | MM-WC19SB7K | E0978 | Crash Tested Seat Belt Kit - 700mm \$893.00 |
| Pelvic Positioning Options | | | |
| <input type="checkbox"/> | MM-HW312-KIT | | Hardware for Bodypoint Pelvic Stabilizer \$75.00 |
| <input type="checkbox"/> | MM-EBE100 | | Extension for Bodypoint Evoflex Pelvic Stabilizer \$45.00 |
| Legrest Options | | | |
| Fixed Legrest | | | |
| Fixed Center Mount Legrest | | | |
| <input type="checkbox"/> | MM-040541 | K0040 | V6 MWD - Center Mount Flip-Up Footplate - 1 piece Footplate \$385.00 |
| <input type="checkbox"/> | MM-040542 | K0040 | V6 MWD - Center Mount Flip-Up Footplate - 2 piece Footplate \$385.00 |
| Swing Away Legrest | | | |
| <input type="checkbox"/> | MM-03528 | K0040 | Side Mount Swing Away Leg Rest - w/ 2 piece flip-up footplate 8.5" x 6", 13" height adj \$672.00 |
| Fixed Center Mount Legrest Calf Pad Options | | | |
| <input type="checkbox"/> | MM-05547SM | E0995 | Fixed Center Mount Calf Pads - w/ breathable fabric - Small (13"-15") \$487.00 |
| <input type="checkbox"/> | MM-05547ME | E0995 | Fixed Center Mount Calf Pads - w/ breathable fabric - Medium (16"-19") \$487.00 |
| <input type="checkbox"/> | MM-05547LG | E0995 | Fixed Center Mount Calf Pads - w/ breathable fabric - Large (20"-23") \$487.00 |
| Fixed Legrest Positioning Strap Options | | | |
| <input type="checkbox"/> | MM-SP102M | K0039 | Bodypoint Medium Calf Strap - side mount only length 36" & width 3.5" \$102.00 |
| <input type="checkbox"/> | MM-SP102L | K0039 | Bodypoint Large Calf Strap - side mount only length 40" & width 3.5" \$118.00 |
| <input type="checkbox"/> | MM-04597 | E0951 | Heel & Side Support - Attaches to the footplate stem and footplate \$556.00 |
| <input type="checkbox"/> | MM-50-SR230 | E0952 | Foot Strap - 1.5" wide strap \$75.00 |
| Power Elevating Legrest | | | |
| Center Mount Power Elevating Legrest | | | |
| <input type="checkbox"/> | MM-05761 | E1012 | Center Mount Power Elevating - 1 pc footplate 11" x 11.5" - Includes Calf Pads \$3,531.00 |
| <input type="checkbox"/> | MM-06636 | E1012 | Center Mount Power Elevating - 2 pc footplate 5" x 10.5" - Includes Calf Pads \$3,424.00 |
| Swing Away Power Elevating Legrest | | | |
| <input type="checkbox"/> | MM-03526 | E1010 | Side Mount Power Elevating Legrest* - Swing away legrest w/ Flip Up 2" height adj footplate \$3,852.00 |
| *Must be ordered with upgraded electronics MM-SFCK or MM-MFCK depending on number of options selected | | | |
| Power Legrest Positioning Options | | | |
| <input type="checkbox"/> | MM-05324 | K0042 | 2" Metal Heel Guards for 2 Piece Footplates \$717.00 |
| <i>Note: May not fit with all seat sizes - requires MM-06636 or MM-03526</i> | | | |
| <input type="checkbox"/> | MM-827A | E0952 | Foot Straps for 2 Piece Footplates - 11"-14" (pair) \$75.00 |
| <input type="checkbox"/> | MM-HL | E0951 | Heel Loops - Attaches to the footplate - length 11.5" - 13" & width 1" \$107.00 |

Electronics

| OPTION | HCPC | Description | Retail |
|--|-------|--|------------|
| Joystick | | | |
| <input type="checkbox"/> MM-PG-JSMLEDL | | LED R-net Joystick - No LCD display or mono jacks | N/C |
| <input type="checkbox"/> MM-PG-JSMD | | Color LCD Joystick - with mono jack ports for external switches | \$337.00 |
| <input type="checkbox"/> MM-PG-JSM2IRBT | | Infra-Red Color LCD Joystick with Bluetooth | \$503.00 |
| Expandable Power Module - Must Select both items | | | |
| <input type="checkbox"/> MM-PM90 | E2377 | R-net 90 amp Expandable Controller | \$696.00 |
| <input type="checkbox"/> MM-PG-SL | E2313 | R-net Harness for Expandable Controller | \$487.00 |
| Joystick Handle Type - (Comes standard with PG Cone-Shaped Joystick Handle) | | | |
| <input type="checkbox"/> MM-PC101B | E2323 | Bodypoint U-Shaped Joystick Handle - 3" | \$161.00 |
| <input type="checkbox"/> MM-PC102B | E2323 | Bodypoint U-Shaped Joystick Handle - 4" | \$161.00 |
| <input type="checkbox"/> MM-PC110B | E2323 | Bodypoint Grooved Mushroom Joystick Handle | \$150.00 |
| Joystick Mount | | | |
| <input type="checkbox"/> MM-03069R | | Fixed Joystick Mount - Depth Adj 8.25" - Right | N/C |
| <input type="checkbox"/> MM-03069L | | Fixed Joystick Mount - Depth Adj 8.25" - Left | N/C |
| <input type="checkbox"/> MM-04607 | E1028 | Swing Away Joystick Mount - Right <i>(Includes latch to lock in the forward position)</i> | \$535.00 |
| <input type="checkbox"/> MM-04608 | E1028 | Swing Away Joystick Mount - Left <i>(Includes latch to lock in the forward position)</i> | \$535.00 |
| <input type="checkbox"/> MM-JMMIDR | E1028 | Midline S/A Joystick Mount - Rotation 0° - 180° - Right | \$1,295.00 |
| <input type="checkbox"/> MM-JMMIDL | E1028 | Midline S/A Joystick Mount - Rotation 0° - 180° - Left | \$1,295.00 |
| Joystick Accessories | | | |
| <input type="checkbox"/> MM-03330AM | E2331 | Attendant Control Joystick for MPS Seating with mounting kit | \$1,792.00 |
| <input type="checkbox"/> MM-03329AM | E2331 | Attendant Control Joystick for REHAB Seating with mounting kit | \$1,792.00 |
| <input type="checkbox"/> 139920 | | Tash Buddy Button - Diameter 2.5" & Height .7" Outer edge force 3.4 oz. Center force 7 oz. | \$166.00 |
| <input type="checkbox"/> MM-ESBLK | | Egg Switch - single function push button that can scroll through modes or profiles - requires CPS9 | \$139.00 |
| <input type="checkbox"/> MM-CPS9-L-HW | | Swing Away Egg Switch Hardware - Stealth | \$428.00 |
| Omni 2 | | | |
| <input type="checkbox"/> MM-05436R | | OMNI 2 Display - for alternative drive systems - Right | \$2,579.00 |
| <input type="checkbox"/> MM-05436L | | OMNI 2 Display - for alternative drive systems - Left | \$2,579.00 |
| Alternative Drive Controls | | | |
| <input type="checkbox"/> MM-HDR-4M-S-C | E2327 | Switch-It Head Array - three switch head array | \$5,671.00 |
| <input type="checkbox"/> MM-03335 | E2373 | Chin Control for MPS Seating System w/ mount | \$3,585.00 |
| <input type="checkbox"/> MM-03540 | E2373 | Chin Control for REHAB Seating System w/ mount | \$3,585.00 |

Accessories

| OPTION | HCPC | Description | Retail |
|--|-------------|---|------------|
| Accessories | | | |
| <input type="checkbox"/> | MM-01605 | Waterproof Wheelchair Cover | \$139.00 |
| <input type="checkbox"/> | MM-01606 | Waterproof Joystick Cover | \$43.00 |
| <input type="checkbox"/> | MM-01603 | Medical Necessity Bag - Velcro's to armrest - 10.5" x 6" x 1.25" | \$32.00 |
| <input type="checkbox"/> | MM-001B | Deluxe Backpack - 13.5" x 15" x 4" | \$91.00 |
| <i>Note: Requires Push handle (MM-04784) on MPS seat to attach backrest to chair</i> | | | |
| <input type="checkbox"/> | MM-PP0052 | Cup Holder - Attaches to the armrest opposite joystick - self leveling | \$150.00 |
| <input type="checkbox"/> | MM-06806 | E2208 Articulating Oxygen Holder for MPS Seating System | \$990.00 |
| <input type="checkbox"/> | MM-06807 | E2208 Articulating Oxygen Holder for Rehab Seating System | \$990.00 |
| <input type="checkbox"/> | MM-05275 | E2207 Cane Holder for MPS Seating System | \$637.00 |
| <input type="checkbox"/> | MM-05276 | E2207 Cane Holder for REHAB Seating System | \$637.00 |
| <input type="checkbox"/> | MM-03197 | Luggage Rack for V6 - Attaches to base behind the chair 10.5" x 13" | \$316.00 |
| <input type="checkbox"/> | MM-02855 | Lighting Kit for V6 - 2 headlights, 2 red marker lights w/ turn signals-NA with Flip Down Knee Supports | \$1,124.00 |
| <input type="checkbox"/> | MM-03199 | Mud Guards Mud Guards - All Terrain DRIVE TIRES only - (pair) requires WH145BPB | \$637.00 |
| <input type="checkbox"/> | MM-04705 | Charge Mate - 12 Volt/Car outlet and USB port/cell phone charger, with a max draw of 6 amps | \$508.00 |
| <input type="checkbox"/> | MM-PP0133 | Cell Phone Holder - Adj. to fit most cell phones | \$278.00 |
| <input type="checkbox"/> | MM-RAM117BU | Fishing Rod Holder | \$278.00 |
| <input type="checkbox"/> | MM-DDT | Hand-held Programmer | \$1,530.00 |
| <input type="checkbox"/> | MM-LM100 | Be Adaptive Limited Mobility Shooting Mount - gun or crossbow | \$1,418.00 |
| <input type="checkbox"/> | MM-TM100 | Be Adaptive Trigger Mechanism | \$786.00 |
| <input type="checkbox"/> | MM-FM200 | Be Adaptive Limited Mobility Fishing Mount | \$1,418.00 |