FIT COMES FIRST




Your favorite wheelchair back—now in over 130 PD configurations!


More Sizes. More Styles. More backs for the perfect fit.



J3[®] Posterior Deep (PD) Back Your favorite back can now fit more clients.

The 6" deep PD backrest provides midline postural support and positioning while the curved 7" tall built-in contoured laterals offer exceptional thoracic stability. Achieve optimal back placement with diverse J3 back hardware options. Easy-to-order customized options, such as multi-colored shells and covers, available via JAY Your Way.





JAY TECHNOLOGY IN THE J3 BACK

















Specifications:

HCPCS Code F2620

Available Heights

LT, MT, UT and SH (8" to 24" tall)

Widths

12" to 20"

Cover Material

3DX™ Spacer Fabric and Lycra® stretch

Hardware Options

Contact: 2-point Angle: +/- 20° Width: 2" Depth: 2.25" Weight Capacity: 300 lbs

Contact: 2-point Angle: +/- 20° Width: 2" Depth: 2.25" Weight Capacity: 300 lbs

Contact: 2-point Angle: +/- 10° Width: 2" Depth: 1.75" Weight Capacity: 300 lbs

Contact: 4-point Angle: +/-10° Width: 2" Depth: 1.38" Weight Capacity: 500 lbs



