

SERIAL NO: \_\_\_\_\_



**Customer Contact Details**

Customer Name: \_\_\_\_\_

Advisor: \_\_\_\_\_

Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Ship to Customer Address:

Shipping Name: \_\_\_\_\_

Shipping Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

**Additional Information**

Notes: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Current Wheelchair: \_\_\_\_\_

Customer Weight: \_\_\_\_\_

Customer Height: \_\_\_\_\_

Customer Disability: \_\_\_\_\_

**Other**

VAT Exemption Declaration (UK Only):  
 I declare that I am chronically sick or disabled by reason of:  
 \_\_\_\_\_  
 and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.  
 Signature: \_\_\_\_\_ Date: \_\_\_\_\_

RGK Contact Authorisation:  
 From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.  
 I would like to be contacted by RGK:   
 Signature: \_\_\_\_\_ Date: \_\_\_\_\_

Guide:

Base product

Required option

Optional Extra



Critical Information



Information Needed



Non certified  
Sub4 weight  
option


**PRODUCT STYLE**

**SBTX2010 OCTANE SUB4**  
\$8,649.00

XA070 Open design

**Max User Weight 110kg (242lbs)**

**PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM**



XW009 QUOTE ONLY (NO CAD FOR APPROVAL)  
\$0.00

XW011 ORDER (NO CAD FOR APPROVAL)  
\$0.00

XW012 QUOTE CAD drawing for approval  
\$405.00 (2 customer amendments ONLY)

XW010 ORDER CAD drawing for approval  
\$250.00 (2 customer amendments ONLY)

**KEY MEASUREMENTS**

XY010 Seat width (*Inside of sideguards*)


 mm

(Max - 440mm)  
(Min - 250mm)

XY020 Seat length


 mm

(Max - 440mm)  
(Min - 250mm)

XY030 Extra seat length (*Provides additional leg support & hand grip*)


 mm

(Max - 440mm total seat)

XY050 Front frame height


 mm

(Max - 520mm)  
(Min - 380mm)

XY060 Rear frame height


 mm

(Max - 500mm)  
(Min - 360mm)

XY070 Frame length


 mm

(Max - 625mm)  
(Min = XY020 + XY030 + 140mm)

XY072 Castor position


 mm

Optimise measurements


XA020 Ergonomic seat  
\$77.00



Specify length

 mm

(Standard = 120mm)




**FRONT FRAME AND FOOTREST**

Front Frame Options

XY080 Seat taper design

**i** Specify width  mm


(When combined with XY091, dimension is inbetween the centre of the main bends)



XY091 Parallel taper design

**i** Specify width  mm

(When combined with XY080, dimensions must be smaller by 20mm)



Footrest Design Choose One)


XY100 Round design

Castor Wheel Size	Max Height Achievable
4"	110mm



**4**  XY105 Square design

Castor Wheel Size	Max Height Achievable
4"	130mm



Footrest Measurements

XY115 Front footrest height

**!**  mm

XY120 Rear footrest height


**i**  mm

Std = 20 lower than front height

XY125 Footrest length


**i**  mm

Maximum length 140mm (Certified weight)  
(160mm and over requires heavy duty upgrade XB150)



**4**  XB200 High level footrest (Non-certified Sub4 option)  
\$135.00

Indicate front footrest height:  
 mm



Footrest Style

XB020 Tubular Titanium footrest  
(Height adjustable)

Footrest Cover

XB045 Lightweight carbon fibre cover  
(With weight saving cut-out)

Footplate width Sub4 certified  
(Std if no option selected)

**4**  No Cut-Out  
(Non Certified Sub4 option)

**4**  Footplate width Standard (optimised plate width)  
(Non Certified Sub4 option)

### CASTOR WHEELS, FORKS & WHEEL LOCKS

#### Castor Wheels

XC135 Aluminium centre 4"x1.25"  
Choose colour below:

- |  |   |
|--|---|
| <input type="checkbox"/> Silver hub & black tyre (Std) | <input type="checkbox"/> Green hub & black tyre |
| <input type="checkbox"/> Black hub & black tyre        | <input type="checkbox"/> Red hub & black tyre   |
| <input type="checkbox"/> Purple hub & black tyre       | <input type="checkbox"/> Blue hub & black tyre  |

XC140 Aluminium centre 5"x1.4" (Non certified Sub4 option)  
\$122.00 (Silver hub and black tyre only)


XC145 Aluminium centre 6"x1.4" (Non certified Sub4 option)  
\$122.00 (Silver hub and black tyre only)

#### Castor Forks

XC192 Anti-Flutter castor fork  
Note: Colour follows accent

#### Castor Pin Width

XY330 Castor pin width  
 Optimise measurements  
(Minimum width for clear castor rotation and maximum wheelbase length)



#### Wheel Lock Options

- SBTXD036 Compact scissor with lightweight fixing
- XD020 Underseat scissor (Non certified Sub4 option)
- XD050 Push to lock (Non certified Sub4 option)
- XD010 Omit wheel locks

- SBTXD085 Eagle Brakes (Indicate colour)  
\$95.00 (Non Certified Sub4 option)
- SBTXD086 Eagle Brakes with L/W fixing (Indicate colour)  
\$108.00 (Only with fixed centre of gravity)  
(Non Certified Sub4 option)

BLACK     PURPLE     RED     BLUE

Black Std if no option selected

### SIDEGUARDS

#### Sideguard

XE055 Screw on lightweight carbon fibre

#### Sideguard with Fender


XE065 Screw on lightweight carbon fibre w/fender

#### Sideguard Measurements

Individual above tyre (Follow tyre - max. 25mm)

mm

Level with tyre (Follow tyre)





#### Sideguard Measurements

10mm above tyre (Standard)

Individual above tyre (Min = 10mm)

mm

Centre of tyre




### BACKREST

#### Backrest Measurements



XY230 Backrest height

Specify height


 mm


XG010 Fixed Backrest (non-folding)



SBTXG020  
\$318.00

Folding backrest with 2 position snap lock and angle adjust (Non certified Sub4 option)

#### Ergonomic Options



XY360 Backrest angle




- 5  
  -3  
  0  
  +3  
  +5

Individual

 Deg

(Max -10/+20 Deg)




XG100 Ergonomic backrest

\$77.00



Vertical height

 mm


Rearward position

 mm


XG110 Inward taper  
\$77.00 (inside of backrest posts)




Specify width or per side

 mm


XG120 Outward flare  
\$77.00 (inside of backrest posts)



Specify width or per side

 mm


### CAMBER BAR

#### Specify Camber Angle



XY390 Camber angle



- 1  
  2  
  3

#### Centre of Gravity



SBTXH100 Single post fixed centre of gravity  
(25-125mm)


 mm


SBTXH130  
\$353.00

Square back adjustable centre of gravity  
(25-125mm) (Non Certified Sub4 option)


 mm


SBTXH140  
\$353.00

Round back adjustable centre of gravity  
(60-125mm) (Non Certified Sub4 option)


 mm

### PUSH HANDLES



XG070  
\$149.00

Screw in push handles  
(Non Certified Sub4 option)



XG080  
\$229.00

Flip down push handles  
(Non Certified Sub4 option)

### ANTI-TIP



XJ100  
\$422.00

Single sport anti tip  
(Non Certified Sub4 option)

**FRAME FINISH OPTIONS**

Frame and Footrest (Choose One)

SBTXN030 Bead Blast frame



XN010 Powder coat frame

SBTXN050 Polished frame

Paint code

Accent Colour (Choose One)

**Note:** Accent colour includes, forks & castor caps.



Black accent

Silver accent

XN250  
\$165.00

Individual Accent Colour

Red accent

Blue accent

Purple accent

Carbon accent  
(Castor caps only, black for rest)

---

**UPHOLSTERY**

Seat Upholstery

XP045 Sub4 Seat sling strap adjustable



XR050 Carbon fibre seat plate  
\$324.00  
(Non certified Sub4 option)

Backrest Upholstery

XP035 Airtech Sub4 backrest strap adjustable



XP050 Airtech backrest strap adjustable  
(Non certified Sub4 option)



XP070 Backrest upgrade to leather corners (Black)  
\$115.00  
(only with Airtech upholstery XP050)  
(Non certified Sub4 option)



XP110 Backrest with zip pocket  
\$85.00  
(only with Airtech upholstery XP050)  
(Non certified Sub4 option)



XP170 Individual backrest embroidery  
\$131.00  
Please complete page 9  
(only with Airtech upholstery XP050)

Backrest Upholstery Stitching Options

Silver

Black

Blue

Individual:

Red

Orange

Purple

---

RGK Logo Colour

White (Std)

Black

Blue

Individual:

Red

Orange

Purple

---

REAR WHEELS (* non certified weight options)	
Rear Wheel Size	
<input type="checkbox"/> XY420    24" Wheels	<input type="checkbox"/> XY430    25" Wheels
Rear Wheel Gap	
<input checked="" type="checkbox"/> XY460    Rear wheel gap <input type="checkbox"/> 15 <input type="checkbox"/> 20 <input type="checkbox"/> 25	
Rear Wheel Options	
<input type="checkbox"/> XK075    Spinergy <b>BLXL</b> 18 blade spoke <i>N/A IN SILVER HUB</i>	<input type="checkbox"/> SBTXK150    Spinergy <b>CLX</b> 18 carbon blade spoke \$674.00 <i>N/A IN SILVER HUB</i>
<input checked="" type="checkbox"/> XK226    Titanium 100mm quick release axle pins	<input type="checkbox"/> SBTXK155    Spinergy <b>CX6</b> 18 carbon blade spoke \$803.00 <i>N/A IN SILVER HUB</i>
<b>Spinergy Colour Options:</b>	
<input type="checkbox"/> Black spokes <input type="checkbox"/> White spokes <input type="checkbox"/> Yellow spokes <input type="checkbox"/> Red spokes <input type="checkbox"/> Blue spokes <input type="checkbox"/> Pink spokes <input type="checkbox"/> Purple spokes <input type="checkbox"/> Orange spokes <input type="checkbox"/> Green spokes	<input type="checkbox"/> XK180    Individual spoke design \$169.00 <i>(Specify colour)</i>  Inside: _____  Outside: _____
<b>Spinergy Anodise Options</b>	
<input type="checkbox"/> XK190    Black hub <i>(std option)</i>	<input type="checkbox"/> XK200    Silver hub \$88.00 <i>(n/a with BLXL/CLX)</i>
<input type="checkbox"/> XK202    Purple hub/rims \$263.00 <i>(LX, XLX &amp; BLXL only)</i>	<input type="checkbox"/> XK201    Red hub/rims \$263.00 <i>(LX, XLX &amp; BLXL only)</i>
<input type="checkbox"/> XK203    Blue hub/rims \$263.00 <i>(LX, XLX &amp; BLXL only)</i>	<input type="checkbox"/> XK204    Gold hub/rims \$263.00 <i>(LX, XLX &amp; BLXL only)</i>
<input type="checkbox"/> XK205    Individual hub/rims \$304.00    Select colour: _____	
Supply Without Wheels	
<input type="checkbox"/> SBTXK040    Supply without Wheels -\$455.00	<input style="width: 150px; height: 20px;" type="text" value="Wheel Type:"/> <input style="width: 150px; height: 20px;" type="text" value="Tyre Type:"/>
Handrim Position / Options	
<input type="checkbox"/> XY480    Narrow set handrim	<input type="checkbox"/> XY485    Wide set handrim
<input type="checkbox"/> XL090    Remove excess mounting tabs \$97.00 <i>(recommended for small wheel gaps)</i>	
Handrim Styles	
<input type="checkbox"/> XL010    Aluminium hard anodised 6 tab <i>(Standard if no option selected)</i>	<input type="checkbox"/> XL040    Titanium 6 tab \$231.00
<input type="checkbox"/> XL160    Elipse 6 tab \$324.00 <i>(replacement for Surge)</i>	<input type="checkbox"/> XL035    Natural Fit LT \$290.00
<input type="checkbox"/> XL020    Tetra grip 6 tab \$272.00	<input type="checkbox"/> XL150    Carbolife Gekko <input type="checkbox"/> Black <input type="checkbox"/> Silver \$573.00
<input type="checkbox"/> XL100    Para grip 6 tab \$272.00	<input type="checkbox"/> XL170    Carbolife Curve L (Black) \$614.00
Tyre Options	
<input type="checkbox"/> XM060    Schwalbe One (Lightest option) \$91.00	<input type="checkbox"/> XM130    Kenda black \$34.00
<input type="checkbox"/> XM010    Kenda <input type="checkbox"/> Grey <input type="checkbox"/> Blue <input type="checkbox"/> Red	<input type="checkbox"/> XM050    Kenda K885 mountain bike tyre (24",25" only) \$26.00 <i>(24",25" only)</i>
<input type="checkbox"/> XM020    Schwalbe Marathon Plus \$66.00	
<input type="checkbox"/> XM040    Shox solid airless \$107.00	
Inner Tube Valve Options	
<input type="checkbox"/> XY500    Schraeder inner tube valve <i>(Standard if no option selected)</i>	<input type="checkbox"/> XY490    Presta inner tube valve

\* Indicates items that are not part of the certified weight of the Tiga Sub4

**4 CUSHION (\*\* non certified weight options)**

Cushion foam options

- XP080** Cushion foam  
\$80.00  
*(select density and size below)*
- Soft (CMS040)                       Medium (CMS050)
- 50mm     75mm     100mm

Cushion cover upgrades

- XP085** Nylon cushion cover  
*(Black)*
- XP090** Airlite cushion cover  
\$62.00  
*(Black)*
- XP100** Leather cushion cover  
\$204.00  
*(Black)*

Jay Products - **PLEASE PROVIDE JAY ORDER FORM**


- JAY Cushion**  
*(Price dependant on size / option)*
- JAY Backrest**  
*(Price dependant on size / option)*
- JAY Whitmyer Headrest**  
*(Price dependant on size / option)*

**4 FRAME OPTIONS (\*\* non certified weight options)**

- XR060** Welded calf strap mounting loops  
\$97.00  
*(To locate calf strap on front frame)*

Specify height from floor

- Optimise measurements*



**4 ACCESSORIES (\*\* non certified weight options)**

- XR080** Frontwheel assembly  
\$700.00  
*(Only with XB010, TGAXB080, XB150)*
- XP120** Frame protectors Neoprene zip fasten (Black)  
\$69.00
- XP420** Frame protectors Nylon velcro fasten (Black)  
\$46.00
- XR010** Seat pouch velcro attach *(Black)*  
\$36.00
- XR200** Rucksack with RGK logo (Black)  
\$61.00
- XR020** Underseat large storage bag (Black)  
\$85.00
- XR090** Tool kit
- XS065** Sub4 Calf Strap

\*\* indicates options that, when chosen, will not enable RGK to guarantee a Sub4kg product weight



**ADDITIONAL NOTES**

**INDIVIDUAL BACKREST EMBROIDERY**

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



**BRACE BAR**

**SPEC NOTES:**