



Squiggles Plus Stander

February 2023

Mark For: Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ ST: _____ ZIP: _____ Dealer Phone: _____ Fax: _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	Submitting for: <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____ ADDITIONAL SHIPPING INFORMATION Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ ST: _____ ZIP: _____ Ship To Phone: _____ Fax: _____
------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Squiggles + Stander	HCPCS: E0641
----------------------------	---------------------

Weight Capacity: 48.4 lbs

Packages

Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	PSPSTB	Squiggles Plus Stander Package - Blue	E0641	\$6,611
<input type="checkbox"/>	PSPSTP	Squiggles Plus Stander Package - Purple	E0641	\$6,611
<input type="checkbox"/>	PSPSTR	Squiggles Plus Stander Package - Red	E0641	\$6,611
<input type="checkbox"/>	PSPSTRV	Squiggles Plus Stander Package - Red Vinyl	E0641	\$6,611

Standard Model

Selection Required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplates and tray attachment.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	184-600	Standard Model with Support Shell		\$3,100
<input type="checkbox"/>	184-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning		\$1,516

Positioning Pads

Selection Required. Includes chest and hip pad covers and knee support covers.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	184-1800-02	Squiggles Plus Stander Positioning Pads - Blue		\$689
<input type="checkbox"/>	184-1800-04	Squiggles Plus Stander Positioning Pads - Purple		\$689
<input type="checkbox"/>	184-1800-09	Squiggles Plus Stander Positioning Pads - Red		\$689
<input type="checkbox"/>	184-1800-09VIN	Squiggles Plus Stander Positioning Pads - Red Vinyl		\$689

Accessories

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	118-619	Tray - Gray	E0950	\$410
<input type="checkbox"/>	188-1854-05	Padded Tray Insert - Black	K0108	\$332
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray	K0108	\$298
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$267
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$291
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Each) <i>Requires Sandal Selection</i>	K0108	\$93

Quantity: (1-6)

Headrest

Cushion selection required with hardware

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	E1028	\$231
<input type="checkbox"/>	184-1812-02	Flat Headrest Cushions - Blue	E0955	\$122
<input type="checkbox"/>	184-1812-04	Flat Headrest Cushions - Purple	E0955	\$122
<input type="checkbox"/>	184-1812-09	Flat Headrest Cushions - Red	E0955	\$122
<input type="checkbox"/>	184-1812-09VIN	Flat Headrest Cushions - Red Vinyl (Pair)	E0955	\$122
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>	E1028	\$212
<input type="checkbox"/>	184-1814-02	Flat Headrest and Lateral Cushions - Blue	E0956	\$213
<input type="checkbox"/>	184-1814-04	Flat Headrest and Lateral Cushions - Purple	E0956	\$213
<input type="checkbox"/>	184-1814-09	Flat Headrest and Lateral Cushions - Red	E0956	\$213
<input type="checkbox"/>	184-1814-09VIN	Flat Headrest and Lateral Cushions - Red Vinyl (Pair)	E0956	\$213

Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Distance Between Midline of Knees	5.5" - 8.4"	140 mm - 215 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.5" - 31.5"	480 mm - 800 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 mm L x 390 mm W x 390 mm H
Pivot Base Weight	14.3 lbs	6.5 kgs
Pivot Base Dimensions (Unfolded)	37.4" L x 23.2" W x 20.3" H	950 mm L x 590 mm W x 515 mm H
Pivot Base Dimensions (Folded)	26" L x 23.2" W x 20.3" H	660 mm L x 590 mm W x 515 mm H
Pivot Base Angle Range	Prone 90-70° / Supine 90-160°	Prone 90-70° / Supine 90-160°



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA