

magic
Mobility™



Because the world's not flat™



Extreme
X8™

OFF-ROAD
AND OFF
THE MAP



The Magic Mobility Extreme X8 provides the freedom you've always wanted. The X8 is your partner to almost any destination. It can handle sand, snow and every terrain you can throw at it. With unrivaled off-road performance and reliability, the X8 is perfect for outdoor adventurers, explorers and rural workers. If the great outdoors is where you want to be, the Extreme X8 is waiting to take you there.

The right wheelchair for:*

- Extreme off-road adventures
- Everyday off-road country living, working and rural productivity
- Active outer urban living



"WE HAVE A SIX ACRE PROPERTY AND THE EXTREME X8 ALLOWS ME TO GET TO ALL THE UNEVEN TERRAINS AND SEE PARTS OF THE FARM THAT I HAVE NOT BEEN ABLE TO ACCESS IN FOUR YEARS! THANK YOU FOR GIVING ME MY FREEDOM AND INDEPENDENCE."

Michael Trost
Maryland USA

Frontier
V6™

ULTIMATE
OUTDOOR
ADVENTURE



Our most popular powered wheelchair worldwide, the Magic Mobility Frontier V6 All Terrain is ideal for the outdoors, yet its size makes it a capable indoor performer. If performance and versatility is what you are after, the Frontier V6 All Terrain will impress those who prefer trails over sidewalks. Choose all-terrain or hybrid wheels.

The right wheelchair for:*

- Adventurous off-road recreation
- Everyday off-road country living, working and rural productivity
- Urban outdoor living and working
- Roomy indoor spaces



"THIS CHAIR HAS MET AND EXCEEDED MY EXPECTATIONS. I'VE GONE MANY MILES AND HOURS ... IN THE TIME I'VE HAD IT. IT HAS NEVER LET ME DOWN. IT'S ALWAYS GOT ME HOME, THAT'S IMPORTANT. I HAVE CONFIDENCE WHEN I LEAVE ON MY DAILY EXCURSIONS."

Gary Hays
Texas, USA



A powered wheelchair for people who know what they want. The V4 RWD is designed for people wanting to get around the street or off-road, as well as indoors. Choose all-terrain or hybrid wheels to ensure you get the best powerchair for your lifestyle.

The right wheelchair for:*

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces

“FINALLY A WHEELCHAIR THAT GETS ME TO THE TOP OF MY DRIVEWAY. SOMETIMES IT’S THE LITTLE THINGS LIKE BEING ABLE TO COLLECT MY OWN MAIL THAT MAKE A DIFFERENCE”

Sue
United Kingdom



The V4 FWD offers a great solution to some of life’s trickiest corners. The V4 FWD is designed for people wanting to get around their homes with ease, using desks and tables, or driving frequently. Choose all-terrain or hybrid wheels to ensure you get the best powerchair for your lifestyle.

The right wheelchair for:*

- Adventurous off-road recreation
- Country living and working
- Urban outdoor living and working
- Roomy indoor spaces

“THE RIGHT WHEELCHAIR EMPOWERS PEOPLE TO PARTICIPATE IN ECONOMIC AND CIVIC LIFE.”

Peta Hooke
Melbourne, VIC



**"WE'VE BEEN
CREATING
MADE-TO-ORDER
WHEELCHAIRS
FOR OVER 25
YEARS, SO THAT
YOU CAN LIVE
THE LIFE YOU
CHOOSE"**

Ashley Daff
R&D Director



ATTENTION TO DETAIL

Working together, we can figure out how to make your wheelchair uniquely yours. We know that small touches make a big difference.



ENGINEERED TO THE LAST BOLT

If we think we can make it better, we will. Our in-house engineering team is hands on, working side by side with customers and therapists. The entire Magic team works continuously to improve your entire experience, from choosing a Magic Mobility wheelchair to its performance in the places you want to go.

BUILT TO LAST

We never compromise on quality. Your powerchair is your mobility, therefore reliability is paramount.



PROUDLY AUSTRALIAN DESIGNED AND ASSEMBLED IN THE USA

We love that our product is designed in Australia. We love meeting our customers, hearing their ever changing needs and taking it back to the drawing board.

SPECIFICATIONS	EXTREME X8	FRONTIER V6	FRONTIER V4 RWD	FRONTIER V4 FWD
Max User Weight	400 lbs., 340 lbs. with Elevate or Tilt			
Seat-to-Floor Height	18.5", 19.5", 20.5"	17", 18", 19"		
Length x Width	40.5" x 28"	All Terrain 41.5" x 28" Hybrid 41.5" x 27.3"	All Terrain 39.3" x 28" Hybrid 39.3" x 25"	All Terrain 36.8" x 28" Hybrid 36.8" x 25"
Drive Wheel Position	4 Wheel Drive	Mid Wheel Drive	Rear Wheel Drive	Front Wheel Drive
Max Speed	6 mph			
Battery Size	Group 24	Group 24 or 34		
Overall Product Weight	Starting at 320 lbs.			
Drive Wheels	14" All Terrain Drive Wheels	14" All Terrain or 14" Hybrid Drive Wheels		
Caster Wheels	N/A	Pneumatic or Solid Inserts		
Power Seat Options	50° of CG Tilt/Power Recline/12" Seat Elevator/Power Elevating Legrest			
Colors	Midnight Black, Ocean Blue, Fire Engine Red, Clover Green, Tangerine and Camouflage	Black, Red, Mossy Green, Royal Blue and Yellow		

PUNCTURE RESISTANT TIRE LINERS

Puncture resistant tire liners provide peace of mind for those concerned about leaks and punctures when going off-road. An Australian rubber manufacturer developed an aramid fiber reinforced rubber compound for puncture resistance. The liners are included on all tires.



FULL LIGHTING PACKAGES

Mounted on flexible rubber brackets that resist bumps, our lights use an extremely bright and highly efficient 24V LED lighting system that runs directly from the wheelchair batteries.



RAM® MOUNT- PHONE HOLDER

The X-Grip® provides simple mounting and a secure hold for phones and tablets of all sizes. In extreme environments, you can ensure added peace-of-mind.



ACCESSORY CHARGERS

Our accessory chargers can be fitted to all Magic Mobility power wheelchairs which use the R-Net system. Two chargers are available:

- 1 x 12V & 2 x USB ports
- 2 x USB ports



See more online.



DESIGNED
IN AUSTRALIA



ASSEMBLED
IN THE USA



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-129945 Rev.A ©09.2021

800.333.4000 www.SunriseMedical.com/Magic

