

SERIAL NO: _____



Customer Contact Details

Customer Name: _____

RGK Genius: _____

Address: _____

Phone: _____

Email: _____

Ship to Customer Address:

Shipping Name: _____

Shipping Address: _____

Phone: _____

Email: _____

Additional Information

Notes: _____

Current Wheelchair: _____

Customer Weight: _____

Customer Height: _____

Customer Disability: _____

Other

VAT Exemption Declaration (UK Only):

I declare that I am chronically sick or disabled by reason of:

and that I am receiving from RGK Wheelchairs Ltd, the wheelchair specified herein and I claim relief from Value Added Tax under group 12 of schedule 8 to Value Added Tax Act 1994.

Signature: _____ Date: _____

RGK Contact Authorisation:

From time to time, RGK would like to send you exclusive information, news and offers on everything RGK, by email or post. We will not pass any of your details to any third party.

I would like to be contacted by RGK:

Signature: _____ Date: _____

Guide:



Base product



Required option



Optional Extra



Critical Information



Information Needed

PRODUCT STYLE


ALAXZ010 ALLSTAR (FOR PLAYERS AGED 16 AND BELOW)
\$4,339.00

ALAXA088 Basketball (with wing)

ALAXA089 Tennis (no wing)

Max User Weight 100kg (220lbs)

PLEASE NOTE: ALL RGK CHAIRS ARE HAND BUILT AND HAVE A MANUFACTURING MEASUREMENT TOLERANCE OF +/- 5MM



CAD APPROVAL

CAD approved by: Customer RGK Genius

XW009 QUOTE ONLY (NO CAD FOR APPROVAL)

XW011 ORDER (NO CAD FOR APPROVAL) **\$0.00**

XW012 QUOTE CAD drawing for approval
\$483.00 (2 customer amendments ONLY)

XW010 ORDER CAD drawing for approval
\$278.00 (2 customer amendments ONLY)

KEY MEASUREMENTS

XY010 Seat width *(Inside of sideguards)*



mm

in

(Max - 500mm)
(Min - 250mm)

XY020 Seat length



mm

in

(Max - 500mm)
(Min - 250mm)

XY030 Extra seat length *(Provides additional leg support & hand grip)*



mm

in

XY050 Front frame height



mm

in

(Max - 630mm)
(Min - 360mm)

XY060 Rear frame height



mm

in

(Max - 630mm)
(Min - 360mm)

XY070 Frame length



mm

in

(Min XY020 + XY030 + 110mm)

XA020 Ergonomic seat
\$96.00







Specify length

mm

in

(Standard = 120mm)



FRONT FRAME AND SEAT

Front Frame Options


XY085 Seat taper design
\$96.00



Specify width

mm

in



Seat Options

XR040 Aluminium seat plate
\$152.00
(Riveted to frame)



(Std XY020 - 20mm)

XR050 Carbon fibre seat plate
\$407.00
(Riveted to frame)



(Std XY020 - 20mm)

FOOTREST

Fixed Footrest

XW012 0 Fixed position footrest
\$254.00


Height and Position

Front footrest height

mm

in

(min = 35mm)




Rear footrest height

mm

in


(min = 35mm)



Footrest position
(front of footplate to castor pot)

mm

in




Footrest Size

Footrest width

mm

in

(Std = 200mm)



Footrest length

mm

in

(Std = 250mm)



Footrest Plate

XB180 Aluminium footplate
(2.5mm)

XB190 Carbon fibre footplate
\$265.00

XB100 ABS moulded foot tray
\$148.00

Specify length

mm

in

Max = 380mm (Full)
Std = 330mm

Adjustable Footrest

ALAXB130 Clamp on adjustable position footrest
(100mm range)

ALAXB085 Adjustable position footrest
(50mm range)


Height and Position

Front footrest height

mm

in

(min = 40mm)




Rear footrest height

mm

in


(min = 40mm)



Footrest position
(front of footplate to castor pot)

mm

in




Footrest Size

Footrest width

mm

in

(Std = 200mm)




Footrest length

mm

in

(Std = 250mm)



Footrest Plate

XB180 Aluminium footplate
(5mm)

XB100 ABS moulded foot tray
\$148.00

Specify length

mm

in

Max = 380mm (Full)
Std = 330mm

Castor Pin Width



XY330 Castor pin width




mm

in

MAX 320mm

Optimise measurements



CAMBER BAR, WING AND ANTI-TIP

Centre of Gravity


ALAXY510 Fixed centre of gravity




Specify measurement

mm

in



ALAXY515 Backrest COG adjustment
(+/- 20mm)
(Not with welded sides)




Specify Camber Angle

XY390 Camber angle




16 18 20 22

ALAXH150 Straight axle design
\$202.00




Wing Options

ALAXA110 Round design basketball wing




ALAXA140 No Wing



ALAXA120 Straight sided design basketball wing
\$201.00



ALAXA125 Arrow design basketball wing
\$201.00



Anti-tip options

ALAXJ040 Single fixed anti tip

Specify position

Basketball (5mm inside rear wheel)

Tennis (100mm outside rear wheel)

ALAXJ060 Double fixed anti tip
\$285.00
(for Basketball)

Specify width



mm in

Optimise measurements
(Wide as possible - std if not indicated)

SIDEGUARDS

Sideguard Style

XE030 Screw on aluminium
\$202.00


XE050 Screw on carbon fibre
\$452.00

ALAXE230 Welded sport
\$431.00
(Not with backrest COG Adj - ALAXY510)

XE310 High design upgrade
\$202.00
(only with ALAXE230)



Sideguard Design (Choose One)



Straight design



Round design




Follow tyre

Sideguard Measurements


Sideguard length

mm in



Sideguard height

mm in



Sideguard height (front)

mm in



Follow Tyre Design

25mm above tyre

Individual above tyre

mm in



Level with tyre



BACKREST


Backrest Height

XY230 Fixed height backrest



Specify height


mm in



ALAXY245 Adjustable height backrest
(+/- 25mm)

Specify height

mm in




Ergonomic Options

XY360 Backrest angle

-5 -3 0 +3 +5


Individual
(Max -10/+10 Deg)



XG100 Ergonomic backrest
\$96.00
(not with ALAXY245)


Vertical height

mm in



Rearward position


mm in



XG110 Inward taper
\$96.00
(not with ALAXY245)

Specify width or per side


mm in



XG120 Outward flare
\$96.00
(not with ALAXY245)

Specify width or per side

mm in



CASTORS AND REAR WHEELS

Castor Options


	XC160 Sport 61mm castor (Black)	\$109.00	Bones Reds bearings
	ALAXC230 Extruded aluminium castor fork		

Rear Wheel Size

XY420	24" Wheels	XY430	25" Wheels	XY440	26" Wheels	XY450	700C Wheels
-------	------------	-------	------------	-------	------------	-------	-------------

Rear Wheel Gap

	XY460	Rear wheel gap
		15 20 25
		Individual mm in



Rear Wheel Options

XK030	Equalizer 36/48 spoke (Black hub, black spokes)	XK100	Spinergy XSLX 24 spoke (Tennis)
		\$905.00	
XK080	Spinergy SLX 24 spoke	XK101	Spinergy XSLX 24 spoke + caps (Basketball)
\$902.00		\$943.00	
Spinergy Colour Options:			
Black spokes	White spokes	Yellow spokes	XK180
			Individual spoke design (Specify colour)
Red spokes	Blue spokes	Pink spokes	
			Inside:
Purple spokes	Orange spokes	Green spokes	Outside:
Spinergy Anodise Options			
XK190	Black hub (std option)	XK200	Silver hub
		\$111.00	
XK201	Red hub/rims	XK205	Individual hub/rims Select colour:
\$329.00		\$329.00	\$381.00
XK202	Purple hub/rims	XK203	Blue hub/rims
\$329.00		\$329.00	
XK204	Gold hub/rims		
\$329.00			

Supply Without Wheels

XK040	Supply without Wheels
-\$222.00	



Handrim Position / Options

XY480	Narrow set handrim	XY485	Wide set handrim
XL090	Remove excess mounting tabs (recommended for small wheel gaps)	XL080	Super narrow set handrim (only with XL040 and includes XL090)
\$122.00		\$233.00	

Handrim Styles

XL010	Aluminium hard anodised 6 tab	XL040	Titanium 6 tab (Not with 700C)
		\$290.00	
XL070	Aluminium hard anodised 12 tab (24" not available)	XL050	Stainless steel 6 tab
\$122.00		\$152.00	

Tyre Options

XM010	Kenda	Grey	Blue	Red	XM080	Tufo MS3 tubeless (Blue only)
					\$95.00	
XM075	Kenda Kaliente	Red	Blue		XM095	Kenda Wide (Grey Only)
\$43.00					\$43.00	
XM150	IRC Exeracer Pro (Red only)					
\$100.00						

Spare Wheel

ALAXK500	Spare wheel Spinergy (Built up + axle pin)	XK523	Spare wheel upcharge - Individual Hubs/Rims
\$1,018.00		\$147.00	
XK520	Spare wheel upcharge - Titanium Rims	XK524	Spare wheel upcharge - Individual Spoke Colour
\$237.00		\$107.00	
XK522	Spare wheel upcharge - Remove Excess Tabs		
\$63.00			

Inner Tube Valve Options

XY500	Schraeder inner tube valve	XY490	Presta inner tube valve (Standard if no option selected)
-------	----------------------------	-------	---



Frame & Component Finish

XN200	Individual paint colour (State Colour/RAL code)
\$197.00	

	Frame Comp
--	---------------

Note: Accent colour includes, forks, castor caps, decals.



Seat Uphostery / Plate

Backrest Upholstery

Backrest Stitching Options

RGK Logo ColourModel Logo Colour

Backrest Accessories

ADD DETAILS ON PAGE 9

CUSHION									
Cushion Foam Options									
XP080 \$96.00	Cushion foam (select density and size below)								
	Soft (CMS040)		Medium (CMS050)						
	50mm	75mm	100mm	Other	mm	in			
XP340 \$90.00	Hard Plastazote cushion upgrade (Indicate max height for Basketball)								
	580mm	630mm							
Cushion Cover Options									
XP085	Nylon cushion cover (Black)				XP380 \$102.00	Neoprene cushion cover (Black)			
XP090 \$74.00	Air Lite cushion cover (Black)								
STRAPS									
Velcro straps					Ratchet Straps				
XS010 \$52.00	100mm adjustable velcro strap				XS020 \$93.00	Sport ratchet strap (Click strap) (11" x 2" pad)			
XS050 \$137.00	100mm Pro adjustable velcro strap				XS021 \$93.00	Sport ratchet strap (Single)			
XS057 \$137.00	50mm Pro knee strap (measurement from tube to tube over knees)				XS025 \$243.00	Pelvi Loc 2DA ratchet strap (For thighs and knees)			
					<i>Select size:</i> Medium (12") Large (14")				
					<i>Select size:</i> S (8") M (12") L (14") XL (17")				
					XS030 \$260.00 Pelvi Loc 3DA ratchet strap (For hips)				
					XS035 \$155.00 M2 ratchet strap (For thighs and knees)				
					XS036 \$155.00 M2 ratchet strap (Single) (For Feet)				
					8" 10" 12"				
					8" 10" 12"				
Velcro strap position									
Strap 1 Specify strap length mm in Specify distance from back post mm in									
Strap 1 Specify strap length mm in Specify distance from back post mm in									
Ratchet strap position									
Strap 1 Specify distance from back post mm in Strap 2 Specify distance from back post mm in									
XS040 \$97.00 Fixed mounting points (one set per click strap) x1 x2 x3									
XS045 \$97.00 Clamped mounting points (one set per click strap) x1 x2 x3									
FRAME OPTIONS									
XR060 \$122.00	Calf strap mounting loops				XR065 \$122.00	Strap location loops			
ACCESSORIES									
XR030 \$116.00	Triple wheel bag (Black)				XS070	Foot strap (select if required)			
XP130 \$104.00	Frame protectors full length				XS060	Calf strap (select if required)			
XA150 \$74.00	Plastic tube fitted to offensive wing				XP400	Back bar protector			
					XR090	Tool kit			

ADDITIONAL NOTES

INDIVIDUAL BACKREST EMBROIDERY

PLEASE PROVIDE JPEG IMAGE AT TIME OF ORDER AND NOTE POSITION BELOW



BRACE BAR

SPEC NOTES:

OFF SPEC REQUEST*

(please specify below the page number(s) and details to which the off spec request applies.

OFF SPEC REQUEST

APPROVED *(for internal RGK office only)*

X

RGK Internal Sign-off

Off spec requests require a feasibility study which may increase order price and lead time. Any dimension/ option/ frame request not listed on the order form will be treated as an off spec request.*

Estimated MSRP*