



Kids Quickie 2 SE - Without Seating

January 2017

EIZ10

Mark For:

Submitting for:

Quote

Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: Fax Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions. The use of **HCPCS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

KIDS QUICKIE 2 SE W/O SEATING

K0005 coding still applies for EIZ10 with certain sizes

<input type="checkbox"/>	142WM08	E1238	Kids Quickie 2 SE W/O Seating	\$2,525.00
--------------------------	---------	-------	-------------------------------	------------

FRAMES

Standard Frame - 250 lbs. Maximum User Weight

Frame Type

Uses Rear Side Frame Depth Adj. & Depth Adj. Back

<input type="checkbox"/>	142FT54	Cross Brace with Seat Rails	NC
<input type="checkbox"/>	142FT55	Cross Brace w/o Seat Rails	NC

Transit System

Incl. Stabilizer Bar. N/A with Swing-Away Stroller Handles and Extension with Cross Brace with Seat Rails.

<input type="checkbox"/>	142TO1	Transit System	\$315.00
--------------------------	--------	----------------	----------

Standard Colors

<input type="checkbox"/>	142C140	Aztec Gold	NC
<input type="checkbox"/>	142C2	Black (Not avail as base w/Zebra)	NC
<input type="checkbox"/>	142C141	Black Cherry	NC
<input type="checkbox"/>	142C10	Black Opal	NC
<input type="checkbox"/>	142C38	Blue Opal	NC
<input type="checkbox"/>	142C36	Blue Purple	NC
<input type="checkbox"/>	142C14	Candy Blue	NC
<input type="checkbox"/>	142C35	Candy Red	NC
<input type="checkbox"/>	142C142	Electric Blue	NC
<input type="checkbox"/>	142C143	Evergreen	NC
<input type="checkbox"/>	142C126	Glow	NC
<input type="checkbox"/>	142C144	Green Apple	NC
<input type="checkbox"/>	142C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	142C148	Mauve Pink	NC
<input type="checkbox"/>	142C103	Rootbeer	NC
<input type="checkbox"/>	142C114	Sparkle Silver	NC
<input type="checkbox"/>	142C145	Sunrise Orange	NC
<input type="checkbox"/>	142C113	Titanium Color Paint	NC
<input type="checkbox"/>	142C8	Yellow	NC

Standard Matte Colors

<input type="checkbox"/>	142C116	Matte Black	NC
<input type="checkbox"/>	142C151	Matte Black Cherry	NC
<input type="checkbox"/>	142C152	Matte Purple	NC
<input type="checkbox"/>	142C150	Matte Electric Blue	NC
<input type="checkbox"/>	142C149	Matte Evergreen	NC

Custom Frame Color

<input type="checkbox"/>	142C43	Mossy Oak Camouflage	\$295.00
<input type="checkbox"/>	142C42	Stars & Stripes	\$295.00
<input type="checkbox"/>	142C108	Zebra (Choose above base color) Std with White unless above base paint is chosen. NA w/ Black	\$295.00

Frame Width

<input type="checkbox"/>	142W2	11" (By Upholstery) (Avail. w/Cross Brace w/ Seat Rails Only)	NC
<input type="checkbox"/>	142W3	12"	NC
<input type="checkbox"/>	142W4	13"	NC
<input type="checkbox"/>	142W5	14"	NC
<input type="checkbox"/>	142W6	15"	\$250.00
<input type="checkbox"/>	142W7	16"	\$250.00
<input type="checkbox"/>	142W8	17" (Avail. w/Cross Brace w/ Seat Rails Only)	\$250.00
<input type="checkbox"/>	142W9	18" (Avail. w/Cross Brace w/ Seat Rails Only)	\$250.00

FRAME DEPTH AND BACK POST LOCATION

Choose 1 or 2

1. Frame Depth and Back Post Location W/O Seat Rail			
<input type="checkbox"/>	142FL1	Kids Frame (10"-16" Depth)	NC
<input type="checkbox"/>	142FD3	10" Frame Depth	NC
<input type="checkbox"/>	142FD4	11" Frame Depth	NC
<input type="checkbox"/>	142FD5	13" Frame Depth	NC
<input type="checkbox"/>	142FD6	15" Frame Depth	NC
<input type="checkbox"/>	142FD7	16" Frame Depth	NC
<input type="checkbox"/>	142FL9	Regular Frame (12"-18" Depth)	NC
<input type="checkbox"/>	142FD3	12" Frame Depth	NC
<input type="checkbox"/>	142FD4	13" Frame Depth	NC
<input type="checkbox"/>	142FD6	15" Frame Depth	NC
<input type="checkbox"/>	142FD8	17" Frame Depth	NC
<input type="checkbox"/>	142FD9	18" Frame Depth	NC
<input type="checkbox"/>		Long Frame (14"-20" Depth)	NC
<input type="checkbox"/>	142FD5	14" Frame Depth	NC
<input type="checkbox"/>	142FD6	15" Frame Depth	NC
<input type="checkbox"/>	142FD8	17" Frame Depth	NC
<input type="checkbox"/>	142FD10	19" Frame Depth	NC
<input type="checkbox"/>	142FD11	E2203 20" Frame Depth	NC

CHOOSE 1 OR 2

2. Frame Depth and Back Post Location With Seat Rail			
Frame Length			
<input type="checkbox"/>	142FL1	Kids	STD
<input type="checkbox"/>	142FL9	Regular	NC
<input type="checkbox"/>	142FL18	Long	NC
Seat Rail Length (Sling Not Included)			
<input type="checkbox"/>	142SN05	10"	NC
<input type="checkbox"/>	142SN06	11"	NC
<input type="checkbox"/>	142SN07	12"	NC
<input type="checkbox"/>	142SN08	13"	NC
<input type="checkbox"/>	142SN09	14"	NC
<input type="checkbox"/>	142SN10	15"	\$250.00
<input type="checkbox"/>	142SN11	16"	\$250.00
<input type="checkbox"/>	142SN12	17"	\$250.00
<input type="checkbox"/>	142SN13	18"	\$250.00
Back Post Location			
<input type="checkbox"/>	142BL14	1" Offset From Seat Rail	NC
<input type="checkbox"/>	142BL15	2" Offset From Seat Rail	NC

DEPTH ADJUSTABLE BACK and JAY

Stroller Handle (Height)			
<input type="checkbox"/>	142BP27	Depth Adj. Back Post Standard	STD
<input type="checkbox"/>	142BP28	Depth Adj. Back Post Medium	NC
<input type="checkbox"/>	142BP29	Depth Adj. Back Post Tall	NC

Stroller Handle Options			
<input type="checkbox"/>	142BP33	Swing Away Stroller Handles	\$100.00

JAY Backs (Optional)			
Ordering JAYbacks on-chair, please see the JAY price list. Requires part #			
<input type="checkbox"/>	142JB02	JAY Bundled Backrest Backrest Name or Part #	Price List
		Width: Height:	

JAY Cushions (Optional)		
To order any JAY cushions on-chair, please see the JAY price list. Cushion name or part number must be filled out with desired width and depth		
<input type="checkbox"/>	142JC20	JAY Bundled Cushion Cushion Name or Part #
		Width: Depth:
		Price List

ARMRESTS

<input type="checkbox"/>	142AA3	Flip Up	STD
<input type="checkbox"/>	142AA1	E0973 Ht. Adj. Locking Flip Up	\$255.00
<input type="checkbox"/>	142AA7	E0973 Ht/Length Adj. Locking Flip Up	\$265.00
<input type="checkbox"/>	142AM1	Desk Length (10")	NC
<input type="checkbox"/>	142AM2	Full Length (14")	NC
<input type="checkbox"/>	142AM9	Foam Short	NC
<input type="checkbox"/>	142AM10	Foam Long	NC

<input type="checkbox"/>	142AA8	E0973 Ht Adj. Adult	\$285.00
<input type="checkbox"/>	142AM1	Desk Length (10")	NC
<input type="checkbox"/>	142AM2	Full Length (14")	NC
<input type="checkbox"/>	142AA9	E0973 Ht/Length Adj. Kids	\$285.00
<input type="checkbox"/>	142AM1	Desk Length (10")	NC
<input type="checkbox"/>	142AM2	Full Length (14")	NC
<input type="checkbox"/>	142AA14	Omit Armrest	NC

AXLE AND WHEELS

Axle Plate

<input type="checkbox"/>	142AL2		Standard Adjustable	NC
<input type="checkbox"/>	142AL71	E0959	Amputee Adjustable	\$85.00

Axle

<input type="checkbox"/>	142AX1		Quick Release	STD
<input type="checkbox"/>	142AX11		Quad Quick Release	\$125.00

Anti-Tip

<input type="checkbox"/>	142AT1	E0971	Rear Anti-Tip	\$140.00
--------------------------	--------	-------	---------------	----------

Rear Wheel Type

<input type="checkbox"/>	142RW3		20" Mag	NC
<input type="checkbox"/>	142RW74		20" Lite Spoke	NC
<input type="checkbox"/>	142RW8		22" Mag	NC
<input type="checkbox"/>	142RW75		22" Lite Spoke	NC
<input type="checkbox"/>	142RW1		24" Mag	NC
<input type="checkbox"/>	142RW54		24" Spoke	\$100.00

Rear Wheel Tire

<input type="checkbox"/>	142RT1	E2211 E2212	Pneumatic (Requires the selection of Rear Wheel Tire/Insert)	NC
<input type="checkbox"/>	142RT29	E2218	Full Polyurethane (avail w/Mag only)	\$50.00

Rear Wheel Tire Tube/Insert

Requires a selection of Pneumatic

<input type="checkbox"/>	142RZ2	E2213	Pneumatic Airless Insert	\$80.00
<input type="checkbox"/>	142RT1	E2211 E2212	Pneumatic Tube	NC

Hand rim

<input type="checkbox"/>	142HR1		Aluminum Anodized	NC
<input type="checkbox"/>	142HR4		Plastic Coated	\$95.00
<input type="checkbox"/>	142HR7		Omit	NC

Wheel Lock

<input type="checkbox"/>	142WL1		High Mount Push	NC
<input type="checkbox"/>	142WL2		High Mount Pull	NC

Extension Handle

<input type="checkbox"/>	142EH1	E0961	6" Pair	\$75.00
--------------------------	--------	-------	---------	---------

Grade Aid

Requires a selection of Rear Wheel Tire/Insert

<input type="checkbox"/>	142GA1	E0974	Grade Aid	\$150.00
--------------------------	--------	-------	-----------	----------

ACCESSORIES

Backpack

<input type="checkbox"/>	142A1		Utility Pack - Black	\$55.00
<input type="checkbox"/>	142A2		Kids Black Backpack - Zippie	\$50.00

Seat Pouch

<input type="checkbox"/>	142A8		Black	\$45.00
--------------------------	-------	--	-------	---------

Miscellaneous

<input type="checkbox"/>	142A49	E2207	Crutch Holder	\$195.00
<input type="checkbox"/>	142A83		Tool Kit	\$45.00
<input type="checkbox"/>	142A54	K0038	Extra Leg Strap	\$30.00

Side Guard

Not available with height adj. flip back armrests

<input type="checkbox"/>	142SG4	K0108	Plastic Kids	\$125.00
<input type="checkbox"/>	142SG3	K0108	Plastic Regular	\$125.00

Positioning Belt

<input type="checkbox"/>	142A141	E0978	Belt 1.5" Auto Buckle w/Clamps	\$70.00
<input type="checkbox"/>	142A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$110.00
<input type="checkbox"/>	142A118	E0978	Belt 2" Aircraft Buckle w/Clamps	\$120.00
<input type="checkbox"/>	142A112	E0978	1" Dual Pull Side Squeeze Pad w /Clmp	\$85.00
<input type="checkbox"/>	142A814	E0978	1" Comfort Fit APB5116	\$120.00

Spoke Guard

<input type="checkbox"/>	142A24	K0065	Clear	\$115.00
<input type="checkbox"/>	142A23	K0065	Black	\$115.00

Upper Extremity Support Tray

<input type="checkbox"/>	142A29	E0950	X-Small (Suggested 10"-12")	\$250.00
<input type="checkbox"/>	142A30	E0950	Small (Suggested 13"-14")	\$250.00
<input type="checkbox"/>	142A31	E0950	Medium (Suggested 15"-17")	\$250.00
<input type="checkbox"/>	142A32	E0950	Large (Suggested 18"-20")	\$250.00

	Seat Sling to Top of Footplate Range		Seat Sling to Top of Footplate Range
60° Swing-Away - Short Ext.	13 - 16 3/4	70° "V" Swing-Away - Short Ext.	14 3/4 - 17 1/2
60° Swing-Away - Medium Ext.	15 1/2 - 19 1/4	70° "V" Swing-Away - Medium Ext.	16 3/4 - 19 1/2
60° Swing-Away - Long Ext.	18 - 21 3/4	70° "V" Swing-Away - Long Ext.	18 3/4 - 21 1/2
60° Platform Flip-Up - Short	12 3/8 - 16 3/8	70° "V" Platform Flip-Up - Short	13 3/8 - 16 1/8
60° Platform Flip-Up - Medium	14 1/2 - 18 1/2	70° "V" Platform Flip-Up - Medium	15 1/2 - 18 1/4
60° Platform Flip-Up - Long	17 - 21	70° "V" Platform Flip-Up - Long	18 - 20 3/4
70° Swing-Away - Short Ext.	13 1/2 - 16 1/4		
70° Swing-Away - Medium Ext.	15 - 18 3/4		
70° Swing-Away - Long Ext.	17 1/2 - 21 1/4		
70° Platform Flip-Up - Short	13 1/2 - 16 3/8		
70° Platform Flip-Up - Medium	15 1/4 - 18 1/4		
70° Platform Flip-Up - Long	17 3/4 - 20 3/4		

* All adjustment ranges may vary due to seat height restrictions. Recommended 2 1/2" footrest clearance to ground surface.

ANGLE ADJUSTABLE BACKREST CHART

Back Tube Height	LOW	MEDIUM	TALL
Pivot Regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot Regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot Regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot Tall set at 7"	16" 17" 18"	19" 20" 21"	N/A

Adjustable Pivot Point Preset Measurement from Seat Sling

If back height option is not chosen, middle option will be preset