



JAY[®] X2 Cushion

E2622

January 2022

Mark For: _____
 Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for: Quote Order
 PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY X2 Cushion - E2622

JAY X2 CUSHION Less than 22" in Width (300lb weight capacity)

✓	Option #	HCPCS	Definition	Actual Width	Actual Depth	Price
<input type="checkbox"/>	JX21414	E2622	JAY X2 Cushion 14" X 14"	13.75"	14"	\$461.00
<input type="checkbox"/>	JX21416	E2622	JAY X2 Cushion 14" X 16"	13.75"	16"	\$461.00
<input type="checkbox"/>	JX21515	E2622	JAY X2 Cushion 15" X 15"	15.25"	15"	\$461.00
<input type="checkbox"/>	JX21516	E2622	JAY X2 Cushion 15" X 16"	15.25"	16"	\$461.00
<input type="checkbox"/>	JX21616	E2622	JAY X2 Cushion 16" X 16"	15.75"	16"	\$461.00
<input type="checkbox"/>	JX21617	E2622	JAY X2 Cushion 16" X 17"	15.75"	17"	\$461.00
<input type="checkbox"/>	JX21618	E2622	JAY X2 Cushion 16" X 18"	15.75"	18"	\$461.00
<input type="checkbox"/>	JX21620	E2622	JAY X2 Cushion 16" X 20"	15.75"	20"	\$461.00
<input type="checkbox"/>	JX21716	E2622	JAY X2 Cushion 17" X 16"	17.25"	16"	\$461.00
<input type="checkbox"/>	JX21717	E2622	JAY X2 Cushion 17" X 17"	17.25"	17"	\$461.00
<input type="checkbox"/>	JX21816	E2622	JAY X2 Cushion 18" X 16"	17.75"	16"	\$461.00
<input type="checkbox"/>	JX21818	E2622	JAY X2 Cushion 18" X 18"	17.75"	18"	\$461.00
<input type="checkbox"/>	JX21820	E2622	JAY X2 Cushion 18" X 20"	17.75"	20.25"	\$461.00
<input type="checkbox"/>	JX22020	E2622	JAY X2 Cushion 20" X 20"	19.75"	20"	\$461.00

JAY X2 Size Matrix

		WIDTH										
		12	13	14	15	16	17	18	19	20	21	22
DEPTH	12	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	13	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	14	N/A	JYW		JYW	JYW	JYW	JYW	N/A	N/A	N/A	N/A
	15	N/A	JYW	JYW		JYW	JYW	JYW	N/A	N/A	N/A	N/A
	16	N/A	JYW						JYW	JYW	N/A	N/A
	17	N/A	N/A	N/A	JYW				JYW	JYW	JYW	N/A
	18	N/A	N/A	N/A	JYW		JYW		JYW	JYW	N/A	N/A
	19	N/A	N/A	N/A	JYW	JYW	JYW	JYW	JYW	JYW	N/A	N/A
	20	N/A	N/A	N/A	JYW				JYW		N/A	N/A
	21	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
	22	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A

<input type="checkbox"/>	=	Standard Offering
<input type="checkbox"/> JYW	=	Offered Via Jay Your Way
<input type="checkbox"/> N/A	=	Not available


JAY X2 CUSHION with Extra Cover - Less than 22" in Width (300lb weight capacity)

	Option #	HCPCS	Definition	Actual Width	Actual Depth	Price
<input type="checkbox"/>	JX21414-2	E2622	JAY X2 Cushion 14" X 14"	13.75"	14"	\$568.00
<input type="checkbox"/>	JX21416-2	E2622	JAY X2 Cushion 14" X 16"	13.75"	16"	\$568.00
<input type="checkbox"/>	JX21515-2	E2622	JAY X2 Cushion 15" X 15"	15.25"	15"	\$568.00
<input type="checkbox"/>	JX21516-2	E2622	JAY X2 Cushion 15" X 16"	15.25"	16"	\$568.00
<input type="checkbox"/>	JX21616-2	E2622	JAY X2 Cushion 16" X 16"	15.75"	16"	\$568.00
<input type="checkbox"/>	JX21617-2	E2622	JAY X2 Cushion 16" X 17"	15.75"	17"	\$568.00
<input type="checkbox"/>	JX21618-2	E2622	JAY X2 Cushion 16" X 18"	15.75"	18"	\$568.00
<input type="checkbox"/>	JX21620-2	E2622	JAY X2 Cushion 16" X 20"	15.75"	20"	\$568.00
<input type="checkbox"/>	JX21716-2	E2622	JAY X2 Cushion 17" X 16"	17.25"	16"	\$568.00
<input type="checkbox"/>	JX21717-2	E2622	JAY X2 Cushion 17" X 17"	17.25"	17"	\$568.00
<input type="checkbox"/>	JX21816-2	E2622	JAY X2 Cushion 18" X 16"	17.75"	16"	\$568.00
<input type="checkbox"/>	JX21818-2	E2622	JAY X2 Cushion 18" X 18"	17.75"	18"	\$568.00
<input type="checkbox"/>	JX21820-2	E2622	JAY X2 Cushion 18" X 20"	17.75"	20.25"	\$568.00
<input type="checkbox"/>	JX22020-2	E2622	JAY X2 Cushion 20" X 20"	19.75"	20"	\$568.00

Please visit www.sunparts.us for spare parts information.

JAY X2 Measuring Guide

Measurements made from completed cushion with covers.



Width
 X = Width
 Note: Measurement made from rear bottom edge to front edge of the cushion.
 Ex: JX21618 X=16"

Depth
 Y = Depth
 Note: Measurement made from front edge of rear cutout and front lower of cushion.
 Ex: JX21618 Y=18"

JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.

Active Fluid Bladder

AFB Active Fluid Bladder \$135.00
 Notes: 40% Larger fluid bladder that provides fluid support to the outside edges of the cushion.

Fluid Volume Change

OF Fluid Volume Overfill NC
 5% 10% 15% 20% 25% 30% 35% 40%

UF Fluid Volume Underfill NC
 -5% -10% -15% -20%

Obliquity Fluid Volume Change

OOFR Obliquity Overfill - Right Side NC
 5% 10% 15% 20% 25% 30% 35% 40%

OOFL Obliquity Overfill - Left Side NC
 5% 10% 15% 20% 25% 30% 35% 40%

OUFRR Obliquity Underfill - Right Side NC
 -5% -10% -15% -20%

OUFLL Obliquity Underfill - Left Side NC
 -5% -10% -15% -20%

Cover Options

CVR Cover Options NC
 X-static w/ Stretch Incontinent Microclimatic

Size Change

RL Reduce Length NC
 -1" -2" -3" -4"

RW Reduce Width NC
 -1"

Positioning Components

APC Add Positioning Components NC
 Notes: Add 1" tall Lateral, Medial, and Pelvic Positioning Pieces.

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:

