

ORDER FORM




**Q200R** with Rehab Seat

**EIPW35**

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		State	Zip
Mark For		Phone	

## Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



## Notes

## Chair Specifications

HCPC	K0848 / K0856
Weights starting at: (with Batteries)	213 lbs
Power Seating	22° Tilt
Seat Type	Rehab Seat w/ Fold Down Back and 12° Angle Adjustment
Seat Size Range	Width 15" to 21" Depth 16" to 20"
Seat-to-Floor Height Range	17" to 19"

## Mid-Wheel Drive Base Specifications

Drive Wheel Position	Rear-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 23", Length 42"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 15.5 miles
Drive Motors	2-pole, 6 mph / 10 kph
Drive Wheels	12" w/ Flat Free Inserts 14" w/ Pneumatic Inserts 14" w/ Flat Free Inserts
Front Casters	9" w/ Pneumatic Inserts 9" w/ Flat Free Inserts
Charger	8A Dual Mode (Compact), 10A
Battery Size	Group 22
Color	7 color accents
Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Rear Wheel Drive Base <i>End Item Number EIPW35 (Rehab Seat)</i>			\$
350PC302	K0848	Quickie Q200 R NPC No Power Option Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,793.00
350PC303	K0856	Quickie Q200 R SPC Single Power Option Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,793.00
User Weight			\$
350WG23	–	Up to 300 lbs.	NC
Motor Package			\$
350PF107	–	2-pole 6 mph / 10 kph	NC
Drive Wheels			\$
350DW77	–	Drive Wheel, 12", Solid, Knobby, Black	NC
350DW78	–	Drive Wheel, 14", Pneumatic, Knobby, Black	\$279.00
350DW79	–	Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels			\$
350CW121	–	Caster Wheels 9", Pneumatic, Black	\$279.00
350CW122	–	Caster Wheels 9", Solid, Black	NC
Color Accents			\$
350EZ20	–	Apple Green	NC
350EZ22	–	Black	NC
350EZ21	–	Blue	NC
350EZ23	–	Orange	NC
350EZ24	–	Red	NC
350EZ19	–	White	NC
350EZ25	–	Silver	NC
ISO Transit Tie Downs			\$
350T050	–	ISO Transit Package	NC
LED Lights			\$
350LI04	–	LED Lights and Indicators	\$1,341.00
350LI02	–	Omit LED Lights	NC
350LI08	–	LED Fender Lights	\$542.00

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Seat Heights <i>(Seat height is measured from front of seat pan to floor with 0° frame angle.)</i>			\$
350HT55	–	Seat Height 17" (N/A with Tilt)	NC
350HT56	–	Seat Height 18"	NC
350HT57	–	Seat Height 19" (N/A with Tilt)	NC
Batteries and Chargers			\$
<b>Group 22 Battery</b>			
350AG31	E2361	GP22 Flat Top Pair w/Straps <i>(N/A w/ Tilt)</i>	\$587.00
350AG35	E2361	GP22/50Ah Flat Top Pair w/Straps	\$652.00
<b>Chargers</b>			
350BJ30	–	10A Dual Mode (Compact)	\$106.00
350BJ19	–	8A Dual Mode (Compact)	NC
SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating Options			\$
<b>Seat Tilt</b>			
350PA56	E1002	Posterior 22° Tilt	\$7,688.00
Pre Set Seat Frame Angle <i>(No Power Seating Options, Increases Seat to Floor Height)</i>			\$
350FA10	–	0° <i>(Adds 0" to STFH)</i>	NC
350FA3	–	3° <i>(Adds 1" to STFH)</i>	NC
350FA4	–	6° <i>(Adds 2" to STFH) (N/A w/ 19" STFH)</i>	NC
Backrest Frame Type			\$
350BF27	–	Tool Angle Adjustable Back, 90° to 102°, Stroller Handle	NC
350BF30	–	Tool Angle Adjustable Back, 90° to 102°, Trailing Push Handles	NC
Backrest Frame Option			\$
350BF31	–	Fold Down, Quick Release Backrest	NC
Back Cane Height			\$
350B266	–	Tall, 23"	NC
Back Cushion Type			\$
350BT148	–	Lite Sling Back, with Tension Adjustable Straps	NC
350BT135	–	JAY® J3® (Section A)	NC
350BT133	–	Omit Back Cushion	NC

## SEATING

OPTION	HCPC	DESCRIPTION	RETAIL
<b>A - JAY® J3 Back</b>			<b>\$</b>
<b>JAY® J3 Back Contour</b>			
350JP01	E2613	J3 PA Back for 14" to 20" Wide	\$599.00
350JP11	E2615	J3 PL Back for 14" to 20" Wide	\$699.00
350JP13	E2620	J3 PD Back for 14" to 20" Wide <i>(18" and 21" shell height only)</i>	\$838.00
<b>JAY® J3 Shell Height</b>			
350JR14	–	16.5" Tall Shell <i>(Upper Thoracic Short)</i>	NC
350JR15	–	18" Tall Shell <i>(Upper Thoracic Med)</i>	NC
350JR16	–	19.5" Tall Shell <i>(Upper Thoracic Tall)</i>	NC
350JR17	–	21" Tall Shell <i>(Shoulder Height)</i>	NC
350JR27	–	22.5" Tall Shell <i>(Shoulder Height)</i>	NC
350JR28	–	24" Tall Shell <i>(Shoulder Height)</i>	NC
<b>JAY® Anterior Thoracic Support</b>			
350K10	E0960	Medium <i>(ATS8302)</i>	\$170.00
350K11	E0960	Large <i>(ATS8303)</i>	\$170.00
<b>JAY® Lateral Thoracic Swing-Away Supports</b>			
350JW35	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Right	\$141.00
350JW38	E0956	Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Left	\$141.00
350JW36	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$141.00
350JW39	E0956	Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$141.00
350JW37	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$141.00
350JW40	E0956	Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$141.00
<b>Lateral Thoracic Swing-Away Hardware</b>			
350JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$279.00
350JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$279.00
<b>Seat Type</b>			<b>\$</b>
350ST80	–	JAY® <i>(Comes with Solid Seat Pan)</i> (See Section B)	NC
350ST34	E2607	SEDEO Cushion <i>(Sizes 16", 18", &amp; 20" Wide)</i> <i>(Comes with Solid Seat Pan)</i>	\$441.00
350ST85	–	Solid Seat Pan	NC
350ST84	–	Seat Sling	NC

SEATING			
OPTION	HCPC	DESCRIPTION	RETAIL
<b>Headrests (Requires Pad, Upholstery and Mount)</b>			\$
<b>Headrest Pad</b>			
350HE36	E0955	SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" <i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i>	\$180.00
350HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	\$229.00
350HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$229.00
350HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	\$239.00
350HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium <i>(Adjustable from 7" to 16" Wide x 4" Tall)</i>	\$352.00
350HE48	–	OMIT Headrest Pad	NC
<b>Headrest Mount</b>			
350MH32	E1028	SEDEO Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO Pad)</i>	\$224.00
350MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$184.00
350MH34	E1028	Whitmyer® Linx <sup>2</sup> Headrest Mount with Quick Release	\$224.00
350MH35	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$118.00
350MH17	–	OMIT Headrest Mount	NC
<b>Seat Size</b>			\$
<b>Seat Pan Width (Seat width is determined by inside of arm pads.)</b>			
350R78	–	15" <i>(16" seat pan width)</i>	NC
350R53	–	16" <i>(16" seat pan width)</i>	NC
350R54	–	17" <i>(16" seat pan width)</i>	NC
350R55	–	18" <i>(19" seat pan width)</i>	NC
350R56	–	19" <i>(19" seat pan width)</i>	NC
350R57	–	20" <i>(19" seat pan width)</i>	NC
350R58	–	21" <i>(19" seat pan width)</i>	NC
<b>Seat Pan Depth (Front of Seat Pan to Back Structure.)</b>			
350RF51	–	16"	NC
350RF45	–	17"	NC
350RF46	–	18"	NC
350RF47	–	19"	NC
350RF48	–	20"	NC
<b>B - JAY® Seat Cushion</b>			\$
350CU59	E2607	JAY® Ion® <i>(Sizes 16" - 21" Wide)</i>	\$441.00
350CU61	E2607	JAY® Union® <i>(Sizes 16" - 21" Wide)</i>	\$498.00
350CU64	E2622	JAY® Fusion® <i>(Sizes 16" - 21" Wide)</i>	\$492.00
350CU104	E2622	JAY® Fusion® with Cryo® Technology	\$603.00

**POSITIONING SUPPORTS & ACCESSORIES**

OPTION	HCPC	DESCRIPTION	RETAIL
<b>Armrest Structure</b>			<b>\$</b>
350EG42	–	Single Post Removable <i>(Only available with Single Power, SPC)</i>	NC
350EG43	E0973	Single Post Removable	\$168.00
<b>Side Guards (Pair)</b>			<b>\$</b>
350SG49	E0953	Side Guards <i>(Depth Adjustable)</i>	\$196.00
<b>Arm Pads</b>			<b>\$</b>
350AM15	–	2" x 14" Armpad	NC
<b>Swing-Away Legrest (Pair)</b>			<b>\$</b>
350H177	–	SEDEO PRO Legrests, Swing Away, 70° <i>(14" to 21" Leg Length)</i>	NC
350H178	E0990	SEDEO PRO Legrests, Manual Elevating, Swing Away <i>(Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)</i>	\$391.00
350H179	E1010	SEDEO PRO Legrests, Power Elevating, Swing Away <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$2,403.00
350H180	–	SEDEO PRO Hanger Receivers Only	NC
<b>Swing Away Legrest Calf Supports</b>			
350DI42	–	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$441.00
<b>Footplates</b>			
350F370	K0040	Angle Adjustable Flip Up Footplates	\$189.00
<b>Accessories</b>			
350DI44	E0951	Heel Loops	\$56.00
350KG4	–	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$224.00
<b>SEDEO Lower Extremity Swing Away Support</b>			
350VD01	E1020	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$441.00
350VD02	E1020	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$441.00
<b>Pelvic Positioning Belts</b>			
350A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
350A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
350A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$84.00
<b>Accessories</b>			<b>\$</b>
<b>Backpack</b>			
350A1	–	Black Utility	\$61.00
<b>Medial Thigh Support (N/A with Sling Seat)</b>			
3508A80	E0957	Medial Thigh Support	\$189.00
3508A81	E2207	Pommel Swing-Away Hardware	\$279.00
<b>Misc. Accessories</b>			
3508A52	–	Ram Mount Phone	\$96.00
3508A53	–	Ram Mount Tablet	\$391.00
3508A54	–	Ram Mount Cup Holder	\$96.00
350A611	E2207	Crutch Holder	\$141.00
3508A58	–	Rear View Mirror	\$96.00
3509A10	–	Utility Hooks	\$58.00

**JOYSTICK / ELECTRONICS**

OPTION	HCPC	DESCRIPTION	RETAIL
<b>VR-2 Joystick</b>			<b>\$</b>
350J127	–	VR-2 Joystick	NC
<b>R-net Joystick (Remote)</b>			<b>\$</b>
350J121	–	Standard LED Joystick (w/o Screen Display)	NC
350J138	–	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	\$279.00
350J118	–	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,178.00
350J119	–	Omit Remote Control	NC
<b>Joystick Mounting</b>			<b>\$</b>
350JM34	–	Fixed on Armrest	NC
350JM35	E1028	Parallel Swing-Away	\$330.00
<b>Joystick Mount (Location)</b>			<b>\$</b>
350HM16	–	Right Hand Armrest	NC
350HM17	–	Left Hand Armrest	NC
<b>Display Language Selection R-net</b>			<b>\$</b>
350L1	–	English	NC
350L3	–	French	NC
350L7	–	Spanish	NC
<b>Joystick Handle (Optional) N/A with HD Joystick</b>			<b>\$</b>
350JH13	E2323	Ball Handle	\$141.00
350JH15	E2323	T-Handle	\$141.00
<b>Electronic Programmings</b>			<b>\$</b>
350PR10	–	R-net PC Programmer w/ Dongle	\$587.00
<b>R-net Power Modules</b>			<b>\$</b>
350ER12	–	Power Module, VR-2 90A, Non-Expandable	NC
350ER13	–	Power Module, R-net, 90A, Expandable	\$1,335.00
350EX15	E2377	Expandable 90A Controller - R-net	\$704.00
350EX16	E2313	Harness for 90A Expandable Controller	\$631.00



**JOYSTICK / ELECTRONICS**

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

**Seating Control**

\$

**Powered Seating Control**

350E040	–	Omit Thru-Drive <i>(Requires Switch Control and Standard LED Joystick Requires Ctrl+5) (N/A w/ VR-2 electronics)</i>	NC
350E037	E2310	Thru-drive Control for 1 Actuator	\$2,146.00
350E038	E2311	Thru-drive Control for 2+ Actuators	\$3,884.00

**Switch Control** *(Available with Omit Thru-Drive only) (Requires use of Omni2 or Joystick Jacks)*

350E033	–	Dual Toggle for 1 Actuator	NC
350E034	–	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
350E035	–	Switch-It Mini 1 button for 1 Actuator	NC
350E036	–	Switch-It Mini 2 button for 1 Actuator	NC
350E041	–	Ctrl+5 Switch Box for 2 Actuators or more	NC

**Additional Power Seating Control**

350QG11	–	Ctrl+5 Switch Box with Free Assignable Functions	\$1,397.00
---------	---	--	------------

**Dual Toggle Mounting Location**

350TN1	–	Toggle Mount, Right	NC
350TN2	–	Toggle Mount, Left	NC

**Ctrl+5 Mounting Location**

350EY19	–	Right Armrest	NC
350EY20	–	Left Armrest	NC
350EY21	–	Back Cane, Right	\$168.00
350EY24	–	Back Cane, Left	\$168.00

**Ctrl+5 Switch Type**

350EY22	–	Push Buttons	NC
350EY23	K0108	Toggle Switches	\$335.00

**Attendant Controls**

\$

350A05	E2331	Dual Control Module for Attendant	\$1,453.00
350A06	E2331	R-net Color LCD (CJSM1), Attendant Control Only	\$1,453.00

**Attendant Mounting Location**

\$

350DN5	–	Right Back Cane	NC
350DN6	–	Left Back Cane	NC

**R-net Programming**

\$

**Assignable Buttons**

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

350PM12	–	Include Assignable Buttons	NC
---------	---	----------------------------	----

**JOYSTICK / ELECTRONICS**

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

**Specialty Input Controls**

\$

**R-net Omni2**

350ON14	–	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$1,984.00
350ON18	–	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	\$433.00

**Omni2 Mount (If used w/ Standard Joystick, it must be mounted opposite side)**

350DS14	E1028	Link-it with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Right	\$554.00
350DS15	E1028	Link-it with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Left	\$554.00

**Proportional Hand Input Devices**

350EU14	E2312	Switch-It® MicroGuide	\$3,682.00
350EU15	E2312	Switch-It® MicroPilot	\$3,682.00

**Hand Input Device Mounting Position (Required)**

350EM05	–	Link-It to the Fixed Joystick Mount	\$330.00
350EM06	E1028	Link-It to the Parallel Swing- Away Joystick Mount	\$554.00
350EM08	–	Omit Mount	NC

**Hand Input Device Mounting Platforms (Optional)**

350MG16	E0950	SunDisk Hand Tray <i>(only available with MicroSeries joysticks)</i>	\$330.00
350MG15	E0950	Bullet Hand Tray	\$330.00

**Hand Input Device Mounting Location (Required)**

350IM09	–	Right	NC
350IM10	–	Left	NC

**Proportional Head Array Input Devices**

350EU25	E2328	Switch-It® DualPro Proportional Head Array - Facial Pads	\$6,387.00
350EU26	E2328	Switch-It® DualPro Proportional Head Array - Spot Pads	\$6,387.00

**Head Array Mounting (Required)**

350IJ13	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$184.00
350IJ14	E1028	Whitmyer® Linx <sup>2</sup> Headrest Mount with Quick Release	\$224.00
350IJ15	E1028	Whitmyer® AXYS Headrest Mount with Quick Release	\$118.00



Q200 R with Rehab Seat

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year

