



Horizon Stander

January 2022

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax Email**Submitting for:****Quote****Order**

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: _____

Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Horizon Stander**HCPCS: E0641****Frame**

Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals, Standard Knees and Basic Footplate.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/178BBEXP	Size 2 Horizon Support Frame - Blue	E1399	\$2,670
<input type="checkbox"/>	LHZ/2/178PBEXP	Size 2 Horizon Support Frame - Pink	E1399	\$2,670
<input type="checkbox"/>	LHZ/3/178BBEXP	Size 3 Horizon Support Frame - Blue	E1399	\$3,302
<input type="checkbox"/>	LHZ/3/178PBEXP	Size 3 Horizon Support Frame - Pink	E1399	\$3,302

Base

Selection Required. Must match Frame size. Includes Worldwide Adapters.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/177-EXP	Size 2 Horizon Base/Chassis	E0641	\$3,283
<input type="checkbox"/>	LHZ/3/177-EXP	Size 3 Horizon Base/Chassis	E0641	\$4,097

Headrest

Head Support Assembly or Extended Head Support Assembly selection required if Head Lateral Supports is selected. Head Support Assembly Size must match Head Lateral Supports size.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/169-BB	Size 2 and 3 Head Support Assembly - Blue	E1399	\$267
<input type="checkbox"/>	LHZ/2/169-PB	Size 2 and 3 Head Support Assembly - Pink	E1399	\$267
<input type="checkbox"/>	LHZ/2/117-BB	Size 2 and 3 Head Lateral Supports - Blue (Pair)	E1399	\$230
<input type="checkbox"/>	LHZ/2/117-PB	Size 2 and 3 Head Lateral Supports - Pink (Pair)	E1399	\$230
<input type="checkbox"/>	LHZ/2/187-BB	Size 2 and 3 Extended 4" Additional Depth Head Support Assembly - Blue	E1399	\$358
<input type="checkbox"/>	LHZ/2/187-PB	Size 2 and 3 Extended 4" Additional Depth Head Support Assembly - Pink	E1399	\$358
<input type="checkbox"/>	LHZ/1/180	Size 1-3 Contoured Head Support	E1399	\$585

Additional Chest Supports

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/168-BB	Size 2 Flexible Chest Laterals - Blue (Pair)		\$480
<input type="checkbox"/>	LHZ/2/168-PB	Size 2 Flexible Chest Laterals - Pink (Pair)		\$480
<input type="checkbox"/>	LHZ/3/168-BB	Size 3 Flexible Chest Laterals - Blue (Pair)		\$536
<input type="checkbox"/>	LHZ/3/168-PB	Size 3 Flexible Chest Laterals - Pink (Pair)		\$536
<input type="checkbox"/>	LHZ/3/154-BB	Size 3 Flip Away Chest Laterals - Blue (Pair)		\$589
<input type="checkbox"/>	LHZ/3/154-PB	Size 3 Flip Away Chest Laterals - Pink (Pair)		\$589

Hip Laterals

Selection required if Hip Belt is not selected

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/164-BB	Size 2 Horizon Hip Laterals - Blue (Pair)		\$343
<input type="checkbox"/>	LHZ/2/164-PB	Size 2 Horizon Hip Laterals - Pink (Pair)		\$343
<input type="checkbox"/>	LHZ/3/164-BB	Size 3 Horizon Hip Laterals - Blue (Pair)		\$356
<input type="checkbox"/>	LHZ/3/164-PB	Size 3 Horizon Hip Laterals - Pink (Pair)		\$356
<input type="checkbox"/>	LHZ/2/162-BB	Size 3 Horizon Hip Laterals for Size 2 Frame - Blue (Pair)		\$356
<input type="checkbox"/>	LHZ/2/162-PB	Size 3 Horizon Hip Laterals for Size 2 Frame - Pink (Pair)		\$356
<input type="checkbox"/>	LHZ/3/153-BB	Size 3 Flip Away Hip Laterals - Blue (Pair)		\$477
<input type="checkbox"/>	LHZ/3/153-PB	Size 3 Flip Away Hip Laterals - Pink (Pair)		\$477

Hip Belts

Selection required if Hip Laterals are not selected . Multiple selections can be made.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/181-BB	Size 2 Supine Belt - Blue		\$205
<input type="checkbox"/>	LHZ/2/181-PB	Size 2 Supine Belt - Pink		\$205
<input type="checkbox"/>	LHZ/3/181-BB	Size 3 Supine Belt - Blue		\$226
<input type="checkbox"/>	LHZ/3/181-PB	Size 3 Supine Belt - Pink		\$226
<input type="checkbox"/>	LHZ/2/182-BB	Size 2 Prone Belt - Blue		\$205
<input type="checkbox"/>	LHZ/2/182-PB	Size 2 Prone Belt - Pink		\$205
<input type="checkbox"/>	LHZ/3/182-BB	Size 3 Prone Belt - Blue		\$226
<input type="checkbox"/>	LHZ/3/182-PB	Size 3 Prone Belt - Pink		\$226

Additional Knee Supports

Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Extended Knee Assembly has 5" - 7" adjustment and Long Extended Knee Assembly has 7" - 9" adjustment.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/167-BB	Size 2 and 3 Standard Knee Assembly - Blue (Pair)		\$480
<input type="checkbox"/>	LHZ/2/167-PB	Size 2 and 3 Standard Knee Assembly - Pink (Pair)		\$480
<input type="checkbox"/>	LHZ/2/183-BB	Size 2 and 3 Medium Extended Knee Assembly - Blue (Pair)		\$632
<input type="checkbox"/>	LHZ/2/183-PB	Size 2 and 3 Medium Extended Knee Assembly - Pink (Pair)		\$632
<input type="checkbox"/>	LHZ/2/184-BB	Size 2 and 3 Long Extended Knee Assembly - Blue (Pair)		\$669
<input type="checkbox"/>	LHZ/2/184-PB	Size 2 and 3 Long Extended Knee Assembly - Pink (Pair)		\$669
<input type="checkbox"/>	LHZ/2/185	Size 1-3 Medium Extended Knee Brackets Only (Pair)		\$532

Lower Leg Positioning

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)		\$254
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)		\$277
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)		\$299
<input type="checkbox"/>	LHZ/2/161	Size 2 and 3 Flip Up Footplate Assembly (Pair)		\$758
<input type="checkbox"/>	LHZ/3/160	Size 1-3 3D Footplate Assembly (Pair) <i>One size fits all</i>		\$1,239

Accessories

Must match Horizon Stander Size

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/143	Size 2 and 3 Supine Tray		\$402
<input type="checkbox"/>	LHZ/2/158	Size 2 and 3 Prone Tray		\$532
<input type="checkbox"/>	LHZ/2/192	Size 2 and 3 Horizon Grab Rail		\$264
<input type="checkbox"/>	LHZ/2/193	Size 2 and 3 Horizon Grab Posts (Pair)		\$264

Technical Specifications - Imperial

Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	50.4" x 25.2" x 31.7"	50.4" x 25.2" x 31.7"
Weight of Product	110.6 lbs	130 lbs
Maximum User Weight	176 lbs	220 lbs
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	7"	7"

Technical Specifications - Metric

Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm
Weight of Product	50 kgs	95 kgs
Maximum User Weight	80 kgs	100 kgs
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm



Sunrise Medical (US) LLC2842 Business Park Ave. · Fresno, CA 93727 · USA