



Horizon Stander

July 2022

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax Email**Submitting for:****Quote****Order**

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____

ST: _____

ZIP: _____

Ship To Phone: _____

Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Horizon Stander**HCPCS: E0641****Frame**

Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals (select Knees and Footplate seperately)

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/198-GRB-EUR	Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E1399	\$1,677
<input type="checkbox"/>	LHZ/3/198-GRB-EUR	Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E1399	\$1,732

Base

Selection Required. Must match Frame size. Includes Worldwide Adapters.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/177-EXP	Size 2 Horizon Base/Chassis	E0641	\$3,283
<input type="checkbox"/>	LHZ/3/177-EXP	Size 3 Horizon Base/Chassis	E0641	\$4,097

Headrest

Head Support Assembly selection required if Head Lateral Support/Extension Bracket is selected.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/169	Size 2 and 3 Head Support Assembly - Black	E1399	\$267
<input type="checkbox"/>	LHZ/2/117	Size 2 and 3 Head Lateral Supports - Black (Pair)	E1399	\$230
<input type="checkbox"/>	101-1959	Headrest 4" Extention Bracket	E1399	\$129
<input type="checkbox"/>	LHZ/1/180	Size 1-3 Contoured Head Support	E1399	\$585

Additional Chest Supports

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/168	Size 2 Flexible Chest Laterals - Black (Pair)		\$480
<input type="checkbox"/>	LHZ/3/168	Size 3 Flexible Chest Laterals - Black (Pair)		\$536

Hip Laterals*Selection required if Hip Belt is not selected*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/164	Size 2 Horizon Hip Laterals - Black (Pair)		\$343
<input type="checkbox"/>	LHZ/3/164	Size 3 Horizon Hip Laterals - Black (Pair)		\$356

Hip Belts*Selection required if Hip Laterals are not selected . Multiple selections can be made.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/181	Size 2 Supine Belt - Black		\$205
<input type="checkbox"/>	LHZ/3/181	Size 3 Supine Belt - Black		\$226

Additional Knee Supports*Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Knee Bracket has 5" - 7" adjustment and Long Knee Bracket has 7" - 9" adjustment.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/167	Size 2 and 3 Standard Knee Assembly - Black (Pair)		\$480
<input type="checkbox"/>	101-2224	Medium Knee Adjustment Bracket (Pair)		\$94
<input type="checkbox"/>	101-2225	Long Knee Adjustment Bracket (Pair)		\$97

Lower Leg Positioning

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)		\$254
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)		\$277
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)		\$299
<input type="checkbox"/>	LHZ/2/161	Size 2 and 3 Flip Up Footplate Assembly (Pair)		\$758
<input type="checkbox"/>	LHZ/3/160	Size 1-3 3D Footplate Assembly (Pair) <i>One size fits all</i>		\$1,239

Accessories*Must match Horizon Stander Size*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/143	Size 2 and 3 Supine Tray		\$402
<input type="checkbox"/>	LHZ/2/158	Size 2 and 3 Prone Tray		\$532

Technical Specifications - Imperial

Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	50.4" x 25.2" x 31.7"	50.4" x 25.2" x 31.7"
Weight of Product	110.6 lbs	130 lbs
Maximum User Weight	176 lbs	220 lbs
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	7"	7"

Technical Specifications - Metric

Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm
Weight of Product	50 kgs	95 kgs
Maximum User Weight	80 kgs	100 kgs
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm



Sunrise Medical (US) LLC 2842 Business Park Ave. • Fresno, CA 93727 • USA