

MWD Crossover Compact Electric Wheelchair

EIPW60




Customer				Ship To			
Account				Name			
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		State	
Buyer				ZIP	Phone		
Mark For				Email / Fax			

Please submit orders to: magicmobility@sunmed.com

Seating Specifications			
HCP	K0848 K0856 K0861 K0850 K0858 K0862 K0849 K0857 K0851 K0859	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23" Depth: 13-24"	STF Height Range	17", 18", 19"
Drive Base Specifications			
Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	350 lbs.	Drive Motors	4-pole 6mph
Base Size: Off-Road	Width 26" Length 39" Turning Radius 22.5"	Drive Wheels	14" Off-Road Drive Wheels 14" Urban Drive Wheels
Base Size: Urban	Width 24" Length 39" Turning Radius 22.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24
		Ground Clearance	3"
		Color Accents	10 colors

Client Measurements

a) Top of head to seat	_____	in
b) Top of shoulder to seat	_____	in
c) Underarm to seat	_____	in
d) Elbow to seat	_____	in
e) Top of leg to seat	_____	in
f) Shoulder width	_____	in
g) Chest width	_____	in
i) Hip width	_____	in
j) Actual sitting depth	_____	in
m) Foot length	_____	in
n) Knee to foot	_____	in
o) Seat width	_____	in



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Seating Options			\$
☐ 600SF4		REHAB Seat and Back - Back Canes and adjustable flat seat pan	NC
☐ 600SF5		MPS Seat and Back - Captain style seating	NC
☐ 600SF6		MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan	NC
Base			
Magic 360 Base - Includes Transit and lights		End Item Number: EIPW60	\$
☐ 600MW61	K0848	Magic Mobility - 360 No Power Seating Options	\$ 18,587.00
☐ 600MW62	K0856	Magic Mobility - 360 Single Power Seating Options	\$ 18,587.00
☐ 600MW63	K0861	Magic Mobility - 360 Multi Power Seating Options	\$ 18,587.00
☐ 600MW64	K0849	Magic Mobility - 360 No Power Seating Options - MPS (Captain)	\$ 18,587.00
☐ 600MW65	K0857	Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	\$ 18,587.00
☐ 600MW66	K0862	Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	\$ 18,587.00
Drive Wheels			\$
☐ 600DW80		Off-Road Drive Wheels and Tires	NC
☐ 600DW81		Urban Drive Wheels and Tires - Pneumatic	NC
☐ 600DW82		Urban Drive Wheels and Tires - Solid	\$ 163.00
Optional Drive Tire Options			\$
☐ 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	\$ 40.00
Base Accessories			\$
☐ 6008A91		Jack for Drive Tire	\$ 213.00
Frame Color			\$
☐ 600C185		Midnight Black	NC
☐ 600C186		Lucent Lime	NC
☐ 600C187		Sunset Red	NC
☐ 600C188		Horizon Blue	NC
☐ 600C189		Dusk Purple	NC
☐ 600C190		Vivid Pink	NC
☐ 600C191		Desert Orange	NC
☐ 600C192		Camouflage - Woods	NC
☐ 600C193		Camouflage - Desert	NC
☐ 600C171		Carbon Fiber	NC
☐ 600C194		American Flag	NC
Seat To Floor Height			\$
Seat to Floor Heights start at 17" +/- 0.25"			
☐ 600HT44		Lowest Seat To Floor Height	NC
☐ 600HT45		Seat To Floor Height +1"	NC
☐ 600HT46		Seat To Floor Height +2"	NC
Batteries			\$
☐ 600AG25	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 1,138.00
Battery Charger			\$
☐ 600BJ28		10A Battery Charger	NC
☐ 600BJ38		12A Battery Charger	\$ 495.00

Seating: Seat Functions			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
No Power Seat Options			\$
☐ 600ST87		No Power Seating System with Power Functionality	\$ 5,586.00
Power Tilt			\$
Utilized without Power Recline.			
☐ 600PA64	E1002	50° Power Tilt in Space	\$ 8,458.00
Power Seat Elevate			\$
Only available when Power Tilt or Power Tilt and Recline is selected.			
☐ 600PI16	E2298	12" Adjustable Seat Elevator	\$ 4,750.00
Power Recline			\$
Utilized without Power Tilt. Requires No Power Seating with Power Functionality.			
☐ 600PA65	E1003	150° Power Recline Only	\$ 7,225.00
Power Tilt and Recline			\$
☐ 600PA66	E1006	50° Tilt & 150° Power Recline	\$ 13,709.00

Seating: Seat Sizes			
MPS (Captain) Seat Cushion (width x depth)			\$
MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.			
☐ 600SI24		16" x 18" Vinyl Seat	NC
☐ 600SI25		18" x 18" Vinyl Seat	NC
☐ 600SI26		18" x 20" Vinyl Seat	NC
☐ 600SI27		20" x 18" Vinyl Seat	NC
☐ 600SI28		20" x 20" Vinyl Seat	NC
☐ 600SI29		22" x 20" Vinyl Seat	\$ 486.00
REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

☐ 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
☐ 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 486.00
☐ 600PN15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 858.00

REHAB Seat Width

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

☐ 13"	☐ 14"	☐ 15"	☐ 16"	☐ 17"	☐ 18"	☐ 19"	☐ 20"	☐ 21"	☐ 22"	☐ 23"	☐ 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options											\$
-------------------------------	--	--	--	--	--	--	--	--	--	--	----

Note cushion size will be determined by Rehab Seat size selected above

☐ 600CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 460.00
☐ 600CU60	E2608	JAY® Ion® (Size 22" Wide)	\$ 578.00
☐ 600CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$ 520.00
☐ 600CU62	E2608	JAY® Union® (Size 22" Wide)	\$ 635.00
☐ 600CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 514.00
☐ 600CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$ 613.00

Seating: Seat Sizes			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)			\$
Lateral Thigh Support Pads			
☐ 600KF22	E0953	Lateral Thigh Support Pad Short - Right (5x5)	\$ 147.00
☐ 600KF23	E0953	Lateral Thigh Support Pad Short - Left (5x5)	\$ 147.00
Lateral Thigh Support Mounts			
☐ 600KF26	E1028	Lateral Thigh Removable Mount - Right	\$ 292.00
☐ 600KF27	E1028	Lateral Thigh Removable Mount - Left	\$ 292.00
Seating: Seatback Options			
MPS Backrest Options - Comes standard with Manual Recline			\$
Backrest Width			
☐ 600BN46		Standard Backrest Width - 16-20" Wide	NC
☐ 600BN47		Wide Backrest Width - 22-23" Wide	\$ 353.00
Backrest Height			
☐ 600B80		Standard Backrest Height - 22" Tall	NC
☐ 600B81		Tall Backrest Height - 26" Tall	\$ 586.00
Backrest Material			
☐ 600BU106		Cloth Backrest	NC
☐ 600BU107		Vinyl Backrest	\$ 319.00
MPS Back Accessories			\$
☐ 600PH33		Push Handle - Needed for O2 holder	\$ 167.00
☐ 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 127.00
REHAB Seating Backrest Options			\$
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
JAY J3 Back Contour			
☐ 600JP01	E2313	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 625.00
☐ 600JP02	E2614	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 928.00
☐ 600JP11	E2615	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 729.00
☐ 600JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 875.00
REHAB Seating Backrest Options			\$
JAY® J3® Back Shell Options and Height Ranges			
☐ 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
☐ 600JR15		18" tall shell (Upper Thoracic Med)	NC
☐ 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
☐ 600JR17		21" tall shell (Shoulder Height)	NC
☐ 600JR27		22.5" tall shell (Shoulder Height)	NC
☐ 600JR28		24" tall shell (Shoulder Height)	NC
REHAB Flip Forward Backrest Option			\$
☐ 600OD23		REHAB Flip Forward Backrest Option	\$ 311.00

Seating: Seatback Options			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
☐ 600JW35	E0956	Lateral Thoracic (L x H) 4" x 5" Pad - Right	\$ 147.00
☐ 600JW38	E0956	Lateral Thoracic (L x H) 4" x 5" Pad - Left	\$ 147.00
☐ 600JW37	E0956	Lateral Thoracic (L x H) 6" x 4" Pad - Right	\$ 147.00
☐ 600JW40	E0956	Lateral Thoracic (L x H) 6" x 4" Pad - Left	\$ 147.00
JAY® Lateral Thoracic Swing-Away Mount			
☐ 600JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$ 292.00
☐ 600JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$ 292.00
JAY® Anterior Thoracic Support			\$
☐ 600ZR10	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 158.00
☐ 600ZR11	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 158.00
☐ 600K10	E0960	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 177.00
☐ 600K11	E0960	JAY® Anterior Thoracic Support - Large ATS8303	\$ 177.00
Pelvic Positioning Options - Select One			\$
☐ 600A848		Positioning Belt, Auto Style Buckle, 2" x 60"	NC
☐ 600A849	E0978	Positioning Belt, Aircraft Buckle, Padded 2" x 50"	\$ 88.00
☐ 600A850	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 203.00
☐ 600A851	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 224.00
Headrest Options - MPS Style			\$
☐ 600NU12	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 432.00
☐ 600NU13	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 432.00
Headrest Options - REHAB Style			\$
Headrest Pad			
☐ 600HE39	E0955	Whitmyer® Plush Headrest Pad, 10"	\$ 239.00
☐ 600HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 367.00
Headrest Mount			
☐ 600MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$ 192.00
☐ 600MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 234.00
Armrest Style Options			\$
☐ 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 462.00
☐ 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 928.00
☐ 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,232.00
Armrest Pad Options			\$
☐ 600AB81		10" Waterfall Armrest Pad - Right	NC
☐ 600AE81		10" Waterfall Armrest Pad - Left	NC
☐ 600AB82		14" Waterfall Armrest Pad - Right	NC
☐ 600AE82		14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
☐ 600H210	K0040	Center Mount (14.5" - 17" Leg Length) with Flip-Up Footplate (11" x 10")	\$ 477.00
Center Mount Fixed Legrest - Accessories			
☐ 600DI59	E0995	Fixed Center Mount Calf Pads	\$ 604.00
☐ 600LS16	E0952	Foot Strap - 1.5" wide strap	\$ 93.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
☐ 600H213	E1012	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 4,384.00
Center Mount Power Elevating Legrest Extension			
☐ 600E42		Center Mount Power Elevating Extension - Standard (15" - 17" Leg Length)	NC
☐ 600E43		Center Mount Power Elevating Extension - Short (10" - 15" Leg Length)	NC
Center Mount Power Elevating Legrest Footplates			
☐ 600F117		Center Mount Power Elevating Footplate - Standard (5" x 12" each)	NC
☐ 600F118		Center Mount Power Elevating Footplate - Short (4.5" x 7" each)	NC
Swing-Away Legrest			\$
Swing-Away Legrest			
☐ 600H211	K0040	Swing-Away Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 835.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
☐ 600H214	E1010	Swing-Away Power Elevating Legrest (12" - 17" Leg Length) with flip-up footplates (6" x 8.5")	\$ 4,783.00
S/A Power Legrest Positioning Options			
☐ 600HL15	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 132.00
Accessories			\$
☐ 6008A93		Waterproof Wheelchair Cover	\$ 173.00
☐ 6008A94		Waterproof Joystick Cover	\$ 54.00
☐ 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 40.00
☐ 600A221		Deluxe Backpack	\$ 113.00
☐ 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 187.00
☐ 600A719	E2208	Articulating Oxygen Holder for MPS Seating System	\$ 1,230.00
☐ 600A612	E2207	Cane Holder for MPS Seating System - Includes push handles	\$ 791.00
☐ 6008A96		Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	\$ 393.00
☐ 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 630.00
☐ 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 345.00
☐ 600A614		Fishing Rod Holder	\$ 345.00
Be Adaptive Products			
☐ 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 1,760.00
☐ 6008A99		Be Adaptive Trigger Mechanism	\$ 975.00
☐ 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 1,760.00

Electronics									
✓ OPTION		HCPCS		DESCRIPTION				RETAIL	
Joystick									\$
☐ 600J121				Standard LED Joystick (w/o Screen Display				NC	
☐ 600J116				Color Joystick (2.3" Color LCD Screen Display, CJSM1)				NC	
☐ 600J118				High Resolution Color Joystick (CJSM2+BT) 3.5" Color LCD Display				\$ 1,230.00	
Joystick Mount									\$
☐ 600JM34				Fixed on Armrest - Depth adjustable 8.25"				NC	
☐ 600JM45		E1028		Swing-Away Joystick Mount - Includes latch to lock in the forward position				\$ 665.00	
☐ 600JM46		E1028		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)				\$ 1,608.00	
Joystick Mount - Location									\$
☐ 600HM16				Right Hand Armrest				NC	
☐ 600HM17				Left Hand Armrest				NC	
Display Language Selection R-net									\$
☐ 600L1	English	NC	☐ 600L3	French	NC	☐ 600L7	Spanish	NC	
Joystick Handle Type									\$
Comes standard with PG Cone-Shaped Joystick Handle									
☐ 600JH13		E2323		Ball Handle				\$ 147.00	
☐ 600JH26		E2323		Bodypoint Grooved Mushroom Joystick Handle				\$ 191.00	
☐ 600JH24		E2323		Bodypoint – Rubber Dome Joystick Handle				\$ 191.00	
☐ 600JH25		E2323		Bodypoint – U-Shaped Joystick Handle (4 inch)				\$ 191.00	
Expandable Power Module									\$
Must select both items									
☐ 600EX05		E2377		Expandable 120A Controller - R-net				\$ 786.00	
☐ 600EX06		E2313		Harness for 120A Expandable Controller				\$ 710.00	
Power Seat Function									\$
☐ 600EO37		E2310		Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function				\$ 2,306.00	
☐ 600EO38		E2311		Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions				\$ 4,174.00	
Power Seat Control Option									\$
☐ 600QG11				Ctrl+5 Switch Box with Voice Activation				\$ 1,459.00	
☐ 600QG12				Ctrl+5 Switch Box with Voice Activation (Multi-power only)				NC	
Ctrl+5 Switch Type									\$
☐ 600EY22				Push Button				NC	
☐ 600EY23		K0108		Toggle Switches				\$ 350.00	
Power Seat Control Option - Ctrl+5 Mounting location									\$
☐ 600EY19	Right Armrest	NC	☐ 600EY20	Left Armrest	NC				
Attendant Controls									\$
☐ 6000A05		E2331		Dual Control Module for Attendant - Includes accessory mount and push handle				\$ 1,561.00	
Attendant Control Mounting Location									\$
☐ 600DN5	Right Back Cane	NC	☐ 600DN6	Left Back Cane	NC				
R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth									\$
☐ 600ON14				R-net Omni2 Only Input Control Module				\$ 2,072.00	
☐ 600ON15				R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices				NC	
R-net Omni2 - Mounting Location									\$
☐ 600DS23	Right Hand Armrest	NC	☐ 600DS24	Left Hand Armrest	NC				
Head Control Packages									\$
☐ 600EU25		E2328		Switch-It DualPro® Proportional Head Array - Facial Pads				\$ 6,668.00	
☐ 600EU26		E2328		Switch-It DualPro® Proportional Head Array - Spot Pads				\$ 6,668.00	