



ORDER FORM

Squiggles Stander

The adjustable 3-in-1 stander designed for growing children.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Squiggles Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Packages			\$
Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps. Substitutions are not available for Packages.			
<input type="checkbox"/> PSSTB	Squiggles Stander Package - Blue	E0641	\$5,928
<input type="checkbox"/> PSSTBV	Squiggles Stander Package - Blue Vinyl	E0641	\$5,928
<input type="checkbox"/> PSSTG	Squiggles Stander Package - Green	E0641	\$5,928
<input type="checkbox"/> PSSTP	Squiggles Stander Package - Pink	E0641	\$5,928
Standard Model			\$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.

<input type="checkbox"/> 118-600	Standard Model with Support Shell	E0641	\$2,218
<input type="checkbox"/> 118-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	E0641	\$1,620

Positioning Pads			\$
-------------------------	--	--	----

Selection Required. Includes chest and hip pad covers and knee support covers.

<input type="checkbox"/> 118-762	Squiggles Stander Positioning Pads - Blue	\$626
<input type="checkbox"/> 118-762-03VIN	Squiggles Stander Positioning Pads - Blue Vinyl	\$626
<input type="checkbox"/> 118-760	Squiggles Stander Positioning Pads - Green	\$626
<input type="checkbox"/> 118-763	Squiggles Stander Positioning Pads - Pink	\$626

Accessories			\$
--------------------	--	--	----

<input type="checkbox"/> 102-1640	QTY	Size 1 Knee Cup	\$128
<input type="checkbox"/> 102-3640	QTY	Size 3 Knee Cup	\$128
<input type="checkbox"/> 118-619		Tray - Gray	\$420
<input type="checkbox"/> 188-1854-05		Padded Tray Insert - Black	\$340
<input type="checkbox"/> 117-769		Grab Rail / Toy Bar for Tray	\$305
<input type="checkbox"/> 152-1600		Small Sandals w/ Straps (Pair)	\$273
<input type="checkbox"/> 152-2600		Medium Sandals w/ Straps (Pair)	\$298
<input type="checkbox"/> 120-816	QTY	Footplate Raiser 1" (Each)	\$92



Sandals



Grab Rail

Squiggles Stander

✓ OPTION	DESCRIPTION	RETAIL
Headrest		\$

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	\$385
<input type="checkbox"/>	137-800-03	Contoured Headrest Cushion - Blue	\$176
<input type="checkbox"/>	137-800-03VIN	Contoured Headrest Cushion - Blue Vinyl	\$176
<input type="checkbox"/>	137-800-01	Contoured Headrest Cushion - Green	\$176
<input type="checkbox"/>	137-800-04	Contoured Headrest Cushion - Pink	\$176
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	\$176
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	\$236
<input type="checkbox"/>	117-850-03	Flat Headrest Cushion - Blue	\$121
<input type="checkbox"/>	117-850-03VIN	Flat Headrest Cushion - Blue Vinyl	\$121
<input type="checkbox"/>	117-850-01	Flat Headrest Cushion - Green	\$121
<input type="checkbox"/>	117-850-04	Flat Headrest Cushion - Pink	\$121
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	\$121
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) Requires Flat Headrest Support Hardware and Cushion	\$217
<input type="checkbox"/>	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	\$197
<input type="checkbox"/>	137-807-03VIN	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	\$197
<input type="checkbox"/>	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	\$197
<input type="checkbox"/>	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	\$197
<input type="checkbox"/>	137-807-05	Flat Headrest Lateral Cushions - Black (Pair)	\$197



Headrest Support



Flat Headrest



Lateral Cushions



Contoured Headrest

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5" Size 3: 3.5"	Size 1: 63 mm Size 3: 89 mm

Sunrise Medical (US) LLC

2842 Business Park Avenue
 Fresno, CA 93727 USA
 Tel.: 1-800-333-4000
 Fax: 1-800-300-7502
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129857 Rev. F ©01/2024



800.333.4000

www.SunriseMedical.com