



# Squiggles Stander

June 2022

**Mark For:**

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via:  Fax  Email

Submitting for:  Quote  Order

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**Squiggles Stander****HCPCS: E0641***Weight Capacity: 48.4 lbs***Packages**

*Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps and Product Assembly. Substitutions are not available for Packages.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	PSSTB	Squiggles Stander Package - Blue	E0641	\$5,642
<input type="checkbox"/>	PSSTBV	Squiggles Stander Package - Blue Vinyl	E0641	\$5,642
<input type="checkbox"/>	PSSTG	Squiggles Stander Package - Green	E0641	\$5,642
<input type="checkbox"/>	PSSTO	Squiggles Stander Package - Orange	E0641	\$5,642
<input type="checkbox"/>	PSSTP	Squiggles Stander Package - Pink	E0641	\$5,642

**Standard Model**

*Selection Required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	118-600	Standard Model with Support Shell		\$2,064
	118-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning		\$1,508

**Positioning Pads**

*Selection Required. Includes chest and hip pad covers and knee support covers.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	118-762	Squiggles Stander Positioning Pads - Blue		\$583
<input type="checkbox"/>	118-762-03VIN	Squiggles Stander Positioning Pads - Blue Vinyl		\$583
<input type="checkbox"/>	118-760	Squiggles Stander Positioning Pads - Green		\$583
<input type="checkbox"/>	118-761	Squiggles Stander Positioning Pads - Orange		\$583
<input type="checkbox"/>	118-763	Squiggles Stander Positioning Pads - Pink		\$583

## Accessories

Item Code	Description	Quantity:	HCPCS	Retail
<input type="checkbox"/> 102-1640	Size 1 Knee Cup	(1 or 2)		\$119
<input type="checkbox"/> 102-3640	Size 3 Knee Cup	(1 or 2)		\$119
<input type="checkbox"/> 118-619	Tray - Gray		E0950	\$390
<input type="checkbox"/> 120-270	Padded Tray Insert - Gray		K0108	\$316
<input type="checkbox"/> 117-769	Grab Rail / Toy Bar for Tray		K0108	\$284
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)		E0951, E0952	\$254
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)		E0951, E0952	\$277
<input type="checkbox"/> 120-816	Footplate Raiser 1" (Each)	(1 to 3)	E0951, E0952	\$86

## Headrest

*Cushion selection required with hardware*

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> <b>120-859</b>	<b>Contoured Headrest Support Hardware</b>	E1028	\$358
<input type="checkbox"/> 137-800-03	Contoured Headrest Cushion - Blue	E0955	\$164
<input type="checkbox"/> 137-800-03VIN	Contoured Headrest Cushion - Blue Vinyl	E0955	\$164
<input type="checkbox"/> 137-800-01	Contoured Headrest Cushion - Green	E0955	\$164
<input type="checkbox"/> 137-800-02	Contoured Headrest Cushion - Orange	E0955	\$164
<input type="checkbox"/> 137-800-04	Contoured Headrest Cushion - Pink	E0955	\$164
<input type="checkbox"/> <b>117-697</b>	<b>Flat Headrest Support Hardware</b>	E1028	\$220
<input type="checkbox"/> 117-850-03	Flat Headrest Cushion - Blue	E0955	\$112
<input type="checkbox"/> 117-850-03VIN	Flat Headrest Cushion - Blue Vinyl	E0955	\$112
<input type="checkbox"/> 117-850-01	Flat Headrest Cushion - Green	E0955	\$112
<input type="checkbox"/> 117-850-02	Flat Headrest Cushion - Orange	E0955	\$112
<input type="checkbox"/> 117-850-04	Flat Headrest Cushion - Pink	E0955	\$112
<input type="checkbox"/> <b>118-667</b>	<b>Lateral Bracket Assembly Hardware (Pair)</b> <i>Requires Flat Headrest Support Hardware and Cushion</i>	E1028	\$202
<input type="checkbox"/> 137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$183
<input type="checkbox"/> 137-807-03VIN	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	E0955	\$183
<input type="checkbox"/> 137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$183
<input type="checkbox"/> 137-807-02	Flat Headrest Lateral Cushions - Orange (Pair)	E0956	\$183
<input type="checkbox"/> 137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$183

## Assembly

Item Code	Description	Retail
<input type="checkbox"/> ASSY-LKY	Product Assembly Fee	\$126

## Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5"    Size 3: 3.5"	Size 1: 63 mm    Size 3: 89 mm



Sunrise Medical (US) LLC2842 Business Park Ave. · Fresno, CA 93727 · USA