



# Squiggles Stander

February 2023

**Mark For:**

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via:  Fax  EmailSubmitting for:  Quote  Order

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**Squiggles Stander****HCPCS: E0641***Weight Capacity: 48.4 lbs***Packages**

*Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps. Substitutions are not available for Packages.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	PSSTB	Squiggles Stander Package - Blue	E0641	\$5,792
<input type="checkbox"/>	PSSTBV	Squiggles Stander Package - Blue Vinyl	E0641	\$5,792
<input type="checkbox"/>	PSSTG	Squiggles Stander Package - Green	E0641	\$5,792
<input type="checkbox"/>	PSSTP	Squiggles Stander Package - Pink	E0641	\$5,792

**Standard Model**

*Selection Required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	118-600	Standard Model with Support Shell		\$2,167
<input type="checkbox"/>	118-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning		\$1,583

**Positioning Pads**

*Selection Required. Includes chest and hip pad covers and knee support covers.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	118-762	Squiggles Stander Positioning Pads - Blue		\$612
<input type="checkbox"/>	118-762-03VIN	Squiggles Stander Positioning Pads - Blue Vinyl		\$612
<input type="checkbox"/>	118-760	Squiggles Stander Positioning Pads - Green		\$612
<input type="checkbox"/>	118-763	Squiggles Stander Positioning Pads - Pink		\$612

**Accessories**

Item Code	Description	Quantity:		HCPCS	Retail
<input type="checkbox"/> 102-1640	Size 1 Knee Cup		(1 or 2)		\$125
<input type="checkbox"/> 102-3640	Size 3 Knee Cup		(1 or 2)		\$125
<input type="checkbox"/> 118-619	Tray - Gray			E0950	\$410
<input type="checkbox"/> 188-1854-05	Padded Tray Insert - Black			K0108	\$332
<input type="checkbox"/> 117-769	Grab Rail / Toy Bar for Tray			K0108	\$298
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)			E0951, E0952	\$267
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)			E0951, E0952	\$291
<input type="checkbox"/> 120-816	Footplate Raiser 1" (Each)		(1 to 3)	E0951, E0952	\$90

**Headrest**

*Cushion selection required with hardware*

Item Code	Description	HCPCS	Retail
<input type="checkbox"/> <b>120-859</b>	<b>Contoured Headrest Support Hardware</b>	E1028	\$376
<input type="checkbox"/> 137-800-03	Contoured Headrest Cushion - Blue	E0955	\$172
<input type="checkbox"/> 137-800-03VIN	Contoured Headrest Cushion - Blue Vinyl	E0955	\$172
<input type="checkbox"/> 137-800-01	Contoured Headrest Cushion - Green	E0955	\$172
<input type="checkbox"/> 137-800-04	Contoured Headrest Cushion - Pink	E0955	\$172
<input type="checkbox"/> 137-800-05	Contoured Headrest Cushion - Black	E0955	\$172
<input type="checkbox"/> <b>117-697</b>	<b>Flat Headrest Support Hardware</b>	E1028	\$231
<input type="checkbox"/> 117-850-03	Flat Headrest Cushion - Blue	E0955	\$118
<input type="checkbox"/> 117-850-03VIN	Flat Headrest Cushion - Blue Vinyl	E0955	\$118
<input type="checkbox"/> 117-850-01	Flat Headrest Cushion - Green	E0955	\$118
<input type="checkbox"/> 117-850-04	Flat Headrest Cushion - Pink	E0955	\$118
<input type="checkbox"/> 117-850-05	Flat Headrest Cushion - Black	E0955	\$118
<input type="checkbox"/> <b>118-667</b>	<b>Lateral Bracket Assembly Hardware (Pair)</b> <i>Requires Flat Headrest Support Hardware and Cushion</i>	E1028	\$212
<input type="checkbox"/> 137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$192
<input type="checkbox"/> 137-807-03VIN	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	E0955	\$192
<input type="checkbox"/> 137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$192
<input type="checkbox"/> 137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$192
<input type="checkbox"/> 137-807-05	Flat Headrest Lateral Cushions - Black (Pair)	E0956	\$192

## Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5"    Size 3: 3.5"	Size 1: 63 mm    Size 3: 89 mm



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA