



Quickie 2 August 2026

EIQ2X

Mark For:

Date: _____
Dealer Acct #: _____
Dealer: _____
Dealer Contact: _____
Dealer Address: _____
Dealer City: _____ ST: _____ ZIP: _____
Dealer Phone: () _____ Fax: () _____
Confirmation Email: _____
Confirm Via: ☐ Fax ☐ Email

Submitting for:☐ Quote☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
Attention: _____
Address: _____
Address: _____
Ship To City: _____ ST: _____ ZIP: _____
Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME**Base model**

E1236 coding still applies for EIQ2X with sizes 14" or less.

☐	002WM32	K0005	Quickie 2	\$2,467.00
☐	002WM109	K0005	Quickie 2 - Transit included	\$2,842.00

Transit Option

N/A on Stroller Handles (Removable Push Handles) or Xtender

☐	002TO46	E1022	Transit System WC19	\$374.00
--------------------------	---------	-------	---------------------	----------

Frame Style

☐	002FS31		Standard Swing Away	NC
☐	002FS30		Hemi Swing Away	NC

Rear Frame Type (300 lbs. Weight Capacity)

☐	002RS1		Rear Side Frame Straight	NC
☐	002RS2		Rear Side Frame Rounded	NC
☐	002RS3		Depth Adjustable Back (DAB) 3.0" of rearward adjustment. N/A w/ swing-away or flip-back armrest.	\$227.00

Heavy Duty Package

☐	002HV1		Heavy Duty (18"-24" Wide only) (350 lbs. weight capacity). Includes Reinforced Back Upholstery & Stabilizer Bar.	\$583.00
--------------------------	--------	--	--	----------

2. SIZE / SEATING**Frame Width**

☐	002W3		12"	NC
☐	002W4		13"	NC
☐	002W5		14"	NC
☐	002W6		15"	NC
☐	002W7		16"	NC
☐	002W8		17"	NC
☐	002W9		18"	NC
☐	002W10		19"	NC
☐	002W11	E2201	20"	\$583.00
☐	002W12	E2201	21" (Only Avail. w/ Heavy Duty)	\$583.00
☐	002W13	E2201	22" (Only Avail. w/ Heavy Duty)	\$583.00
☐	002W57	E2201	B4ME - 23" Width (Heavy Duty)	\$816.00
☐	002W58	E2202	B4ME - 24" Width (Heavy Duty)	\$816.00

Frame Depth

Seat Depth can be ordered up to 3" smaller for growth. Std frame available with 13"-22" frame Depth & HD frame available with 16"-22" frame Depth. Hemi available with 16"-18" & 20"-22" frame depth. Please refer to matrix in page 8 of Appendix for frame Depth Vs seat Depth matrix. Growth is achieved to the front of the depth.

☐	002FD4		13"	NC
☐	002FD5		14"	NC
☐	002FD6		15"	NC
☐	002FD7		16"	NC
☐	002FD8		17"	NC
☐	002FD9		18"	NC
☐	002FD10		19" Steel reinforced	NC
☐	002FD36		20" Steel reinforced	NC
☐	002FD51		B4ME - 21" Frame Depth	\$700.00
☐	002FD52		B4ME - 22" Frame Depth	\$700.00

Seat Depth

12" seat depth not available with Hemi frame

☐	002SD5		12"	NC
☐	002SD6		13"	NC
☐	002SD7		14"	NC
☐	002SD8		15"	NC
☐	002SD9		16"	NC
☐	002SD10		17"	NC
☐	002SD11		18"	NC
☐	002SD72		19"	\$88.00
☐	002SD13	E2203	20"	\$641.00
☐	002SD105	E2203	B4ME - 21" Depth	\$700.00
☐	002SD106	E2204	B4ME - 22" Depth	\$700.00
☐	002SD107	E2203	B4ME - 21" Depth (Heavy Duty)	\$816.00
☐	002SD108	E2204	B4ME - 22" Depth (Heavy Duty)	\$816.00

Power Assist

To order any Empulse Power Assist productd please see the Order Form for Part #'s

Option	HCPCS	DESCRIPTION	RETAIL
☐ 209MW91	E0986	Empulse R90 Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 2.5Ahr Removable Battery	\$7,300.00
☐ 209MW92	E0986	Empulse M90 M90 Wheels, Dual Controller, Two M90 Batteries, Mounting Hardware, Quick Release Axles, and Device Charger	\$7,950.00
☐ 147MW95	E0986	Empulse F35 Includes: F35 Stick, Subframe mount, Frame clamps kit, and Charger	\$6,500.00



✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

2. SIZE / SEATING (CONT)

Seat Sling

☐	002SL1		Standard Seat Sling	NC
☐	002SL6		Omit Seat Sling	-\$15.00

JAY Cushions (Optional)

To order any other size JAY cushions on-chair, please see the JAY price list.

☐	002JC20	Cushion Name or Part #	Price List
		Width: Depth:	

NEW!

JAY seat cushions (On-chair)

By opting in the below options the cushions will match the chair width/depth.

☐	002JC15	E2622	JAY X2	\$501.00
☐	002JC17	E2622	JAY Fusion Reduced Profile / Fluid	\$515.00
☐	002JC18	E2607	JAY Union Reduced Profile	\$520.00
☐	002JC10	E2607	JAY ION	\$462.00
☐	002JC06	E2622	JAY LitePro	\$525.00

3. CASTERS / SEAT HEIGHT

Front Seat Height

Please see matrix on Page 9 for additional information. 13.5"-14" Hemi only

21.5"-22" not available in Hemi Frame

☐	002SH40	13.5"	NC	☐	002SH18	18.0"	NC
☐	002SH41	14.0"	NC	☐	002SH20	18.5"	NC
☐	002SH42	14.5"	NC	☐	002SH21	19.0"	NC
☐	002SH12	15.0"	NC	☐	002SH23	19.5"	NC
☐	002SH13	15.5"	NC	☐	002SH24	20.0"	NC
☐	002SH14	16.0"	NC	☐	002SH26	20.5"	NC
☐	002SH15	16.5"	NC	☐	002SH27	21.0"	NC
☐	002SH16	17.0"	NC	☐	002SH29	21.5"	NC
☐	002SH17	17.5"	NC	☐	002SH36	22.0"	NC

NEW!

Caster Wheels

See matrix on page 9.

☐	002CW1		3" Micro	NC
☐	002CW3		4" Low Profile Polyurethane N/A w/ HD package	NC
☐	002CW31		4" Micro Lighted	\$82.00
☐	002CW72		4" x 1" Aluminum Polyurethane	\$129.00
☐	002CW77		4" x 1.25" Semi Pneumatic	NC
☐	002CW64		4" x 1.5" Aluminum Soft Roll	\$216.00
☐	002CW5		5" Low Profile Polyurethane	NC
☐	002CW51		5" Micro Lighted	\$100.00
☐	002CW6	E2219	5" Soft Roll (Grey) N/A w/ Frog Legs	\$105.00
☐	002CW65		5" x 1" Aluminum Polyurethane	\$129.00
☐	002CW78		5" x 1.5" Semi-Pneumatic	NC
☐	002CW66		5" x 1.5" Aluminum Soft Roll	\$227.00
☐	002CW34		6" Polyurethane	NC
☐	002CW33	E2214 E2215	6" Pneumatic	\$88.00
☐	002CW79		6" x 1.5" Semi Pneumatic	NC
☐	002CW68		6" x 1.5" Aluminum Soft Roll	\$239.00
☐	002CW16		7" Polyurethane	NC
☐	002CW32		8" Polyurethane	NC
☐	002CW21	E2214 E2215	8" Pneumatic	\$100.00
☐	002CW81		8" x 1.5" Semi Pneumatic	NC
☐	002CW36	E2214 E2215	8" x 2" Pneumatic	\$169.00
☐	002CW37	E2214 K0108	8" x 2" Pneumatic w/Airless Insert	\$263.00

Caster Housing

Integral bearing forks only available with Integral caster housing

☐	002CP2		Standard	NC
☐	002CP14		Integral Caster Housing N/A with caster Pin locks	\$88.00

Caster Fork


Please see matrix on Page 9 for additional information. Frog Legs not

available with Integral caster housing

☐	002CF46		Multi Position Fork-Short	NC
☐	002CF47		Multi Position Fork- Tall	NC
☐	002CF62		3" Fork	NC
☐	002CF12		5 1/4" Fork	NC
☐	002CF17		6 " Fork	NC
☐	002CF23		7" Fork	NC
☐	002CF29		8" x 2" Fork	NC
☐	002CF65		Multiple Position Integral caster housing only	NC
☐	002CF69		Multiple Position Single Sided Fork Integral caster housing only	\$197.00
☐	002CF53	E1015	Frog Legs Suspension (0-100 lbs)	\$507.00
☐	002CF54	E1015	Frog Legs Suspension (101-190 lbs)	\$507.00
☐	002CF55	E1015	Frog Legs Suspension (191-265 lbs)	\$507.00

Caster Options

☐	002OC1	K0073	Caster Pin Locks (N/A w/8"x2") NA w/ Xtender, Frog Legs or Transit	\$227.00
--------------------------	--------	-------	---	----------



✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

4. REAR WHEELS / SEAT HEIGHT

Rear Seat Height							
☐	002RH11	13.0"	NC	☐	002RH26	17.5"	NC
☐	002RH12	13.5"	NC	☐	002RH20	18.0"	NC
☐	002RH13	14.0"	NC	☐	002RH27	18.5"	NC
☐	002RH14	14.5"	NC	☐	002RH21	19.0"	NC
☐	002RH15	15.0"	NC	☐	002RH37	19.5"	NC
☐	002RH16	15.5"	NC	☐	002RH22	20.0"	NC
☐	002RH17	16.0"	NC	☐	002RH38	20.5"	NC
☐	002RH18	16.5"	NC	☐	002RH23	21.0"	NC
☐	002RH19	17.0"	NC	☐	002RH39	21.5"	NC
				☐	002RH40	22.0"	NC

Axle Plate			
☐	002AL102	Standard Adjustable (All in-one)	NC
☐	002AL6	E0959 Amputee Adjustable N/A w/ Padded Armrest	\$205.00
Axle			
☐	002AX1	Quick Release Steel	NC
☐	002AX11	Quad Quick Release 5 Spoke Mag & Mountain Wheel only	\$147.00

One Arm Drive									
N/A w/ Xtender, 12", 23" & 24" wide, amputee axle plates, and small tray table. N/A w/ combination of Padded SA arms and Anti tips. Only available w/ 20", 22", 24" mag wheels. Max weight capacity of 300 lbs. Only available in 6 spoke Mags									
☐	002AL16		E0958	One Arm Drive Adj. Axle Plate					\$991.00
	☐	002OA1		Small Handrim (N/A w/ 20" wheels) ☐ RH position ☐ LH position					NC
	☐	002OA2		Standard Handrim ☐ RH position ☐ LH position					NC
Select Wheel Size:	☐	20"	☐	22"	☐	24"			

Xtender Power Assist	
RSH of 14.5"-20" for 22" rear wheels and RSH 15.5"-21" for 24" wheels. Front seat height must be greater than rear. N/A w/ Transit Option Max weight capacity of 265 lbs. Available 12"-22" widths.	

Xtender Power Assist				
☐	002XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	\$10,404.00
☐	002XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	\$11,390.00
☐	002XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	\$10,404.00
☐	002XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	\$11,390.00
Please select one of the following handrim options				
☐	002XH5		Plastic Coated	NC
☐	002XH6		Stainless Steel	NC

Xtender - Other Options				
☐	002XA10		Spare NiMH Power Pack - NICKEL	\$2,128.00
☐	002XA11		Spare Li-ion Power Pack - LITHIUM	\$2,594.00

Xtender Wheel Options				
☐	002XA8		SPARE 22" Lite Spoke Pneu Alum	NC
☐	002XA9	E2211 E2213	SPARE 22" Lite Spoke Pneu Plastic	\$234.00
☐	002XA2		SPARE 24" Lite Spoke Pneu Alum	NC
☐	002XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic	\$234.00

Rear Wheel Type			
NEW!			
Not applicable w/ Xtender Power Assist Option. Ultra Lightweight Wheels only available w/ Alum Anodized Handrim. Handrims cannot be omitted on Lite Spoke Wheels. * Omit Rear wheel will default chair configuration set up to a 24" wheel			
☐	002RW103	20" Mag Wheel 5 Spoke	NC
☐	002RW104	22" Mag Wheel 5 Spoke	NC
☐	002RW102	24" Mag Wheel 5 Spoke	NC
☐	002RW158	20" Lite Spoke Wheel	NC
☐	002RW159	22" Lite Spoke Wheel	NC
☐	002RW92	24" Lite Spoke Wheel	NC
☐	002RW93	25" Lite Spoke Wheel	NC
☐	002RW160	26" Lite Spoke Wheel	NC
☐	002RW77	22" Ultra Lightweight Spoke Wheel	\$292.00
☐	002RW78	24" Ultra Lightweight Spoke Wheel	\$292.00
☐	002RW28	22" Spinergy SPOX Black	\$985.00
☐	002RW95	24" Spinergy SPOX Black	\$985.00
☐	002RW96	25" Spinergy SPOX Black	\$985.00
☐	002RW97	26" Spinergy SPOX Black	\$985.00
☐	002RW32	22" Spinergy LX Black	\$1,237.00
☐	002RW98	24" Spinergy LX Black	\$1,237.00
☐	002RW99	25" Spinergy LX Black	\$1,237.00
☐	002RW40	26" Spinergy LX Black	\$1,237.00
☐	002RW50	24" Mountain	\$292.00
☐	002RW25	Omit Rear Wheel and Axle *	-\$55.00
☐	002RW49	Omit Rear Wheel (Use w/24" Xtender)	NC
☐	002RW111	Omit Rear Wheel (Use w/22" Xtender)	NC

Rear Wheel Tire				
☐	002RT1		Pneumatic	NC
☐	002RT2	E2211 E2213	Pneumatic Airless Insert <i>N/A with Xtender.</i>	\$100.00
☐	002RT24	E2211 E2213	Pneumatic Airless Insert - Infinergy® Lightweight - <i>E-TPU material</i>	\$114.00
☐	002RT3		Full Poly	NC
☐	002RT38	E2218	SHOX Performance Poly	\$110.00
☐	002RT25	E2218	Primo Orion Poly	\$126.00
☐	002RT39	E2211 E2212	High-Pressure Clincher Black	\$110.00
☐	002RT27	E2211 E2212	Knobby V-TRAK	\$110.00
☐	002RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$164.00
☐	002RT12	E2211 E2212	Mountain Tire <i>Only available on 24"</i> <i>Mountain Rear Wheel</i>	\$110.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
4. REAR WHEELS / SEAT HEIGHT (CONT)									


Rear Seat Height Availability Matrices

Rear tire, side frame, and axle plate selections dictate seat height availability. *Not all combinations are possible within acceptable seat height ranges

Rear Tire	Min RSH	Max RSH
20" Pneumatic, Poly	13"	19.0"
22" Pneumatic, Poly	14.0"	20.0"
22" HP Clincher Blk	13.5"	19.5"
22" Schwalbe	13.5"	19.5"
24" Pneumatic, Poly	15"	21.0"
24" HP Clincher Blk	14.5"	20.5"
24" Knobby	15"	21.0"
24" Schwalbe	14.5"	20.5"
24" SHOX Poly	14.5"	20.5"
24" Mountain	15.0"	21.0"
25" HP Clincher Blk	15.0"	21.0"
25" Schwalbe	15.0"	21.0"
25" SHOX Poly	14.5"	20.5"
26" Pneumatic	16.0"	22.0"
26" HP Clincher Blk	15.5"	21.5"
26" Schwalbe	15.5"	21.5"

Axle Plate	Min RSH	Max RSH
Std Adj (All in one)	13"	22.0"
Amputee	14"	22.0"

Rear Side Frame	Min RSH	Max RSH
Curved	13"	21.0"
Straight	13	22.0"
DAB	13	22.0"



Wheel	5 Spoke Mag	Mag	Lite Spoke	Ultra Lt wt	Spinergy
Tire					
Pneumatic	20,22,24	20,22,24	20,22,24,26		24,26
Pneumatic w/Airless Insert	20,22,24	20,22,24	20,22,24,26		24,26
Full Poly Tire	20,22,24	20,22,24	20,22,24		
High Pressure Clinch Blk			22,24,26	22,24	22,24,26
Knobby V-TRAK			24		24
Schwalbe Marathon			24,25,26	22,24	24,25,26
SHOX Perform. Poly	20,22,24		24,25		24,25
20" = 451 mm 22" = 501 mm 24" = 540 mm 25" = 559 26" = 590					

Handrim

See availability matrix in Appendix pages.

☐	002HR1	Aluminum Anodized	NC
☐	002HR4	Plastic Coated - N/A on every rim	\$129.00
☐	002HR22	Natural Fit - Standard Grip	\$549.00
☐	002HR23	Natural Fit - Super Grip	\$611.00
☐	002HR28	Q-Grip Neoprene Coated	\$432.00
☐	002HR39	Ellipse 3R	\$500.00
☐	002HR40	Carbolife - Gekko	\$535.00
☐	002HR41	Carbolife - Curve L Tetra grip	\$575.00
☐	002HR8	8 Oblique N/A w/ Xtender. Only available with 24" or 26" 6 spoke Mag or Lite Spoke rims.	\$379.00
☐	002HR7	Omit Handrim Option may not be available on all wheels	NC

Camber

☐	002CI1	0°	NC
☐	00217CI	2° (default)	NC
☐	002CI3	4°	NC

Center of Gravity Preset

Amputee Axle Plate will come preset at 1.0" Center of Gravity (CG), range of adjustment is -3.0" to -1.75" and 1.0" to 3.5" CG.

☐	002CC2	1.0" (default)	NC
☐	002CC3	1.5"	NC
☐	002CC4	2.0"	NC
☐	002CC5	2.5"	NC
☐	002CC6	3.0"	NC
☐	002CC7	3.5"	NC

One Arm Drive CG

☐	002CC30	1AD Forward position	NC
☐	002CC31	1AD Reward position	NC
☐	002CC32	1AD Center position	NC

5. WHEEL LOCKS

Wheel Lock

Ergo & compact N/A w/ 13"-14" frame depth or amputee axle plate. Additional No Gos may exist depending on wheel size and position

☐	002WL46	High Mount Push	NC
☐	002WL50	High Mount Pull	NC
☐	002WL75	Ergo Scissor	\$129.00
☐	002WL64	Compact	\$164.00
☐	002WL24	Unilateral Wheel Lock Right	\$350.00
☐	002WL23	Unilateral Wheel Lock Left	\$350.00
☐	002WL10	Omit	NC

Wheel Lock Option

☐	002WL12	Do Not Mount	NC
--------------------------	---------	--------------	----

Extension Handle

N/A w/ ergo scissor or compact wheel locks

☐	002EH7	E0961	Extension Handle (Pair)	\$88.00
--------------------------	--------	-------	-------------------------	---------

Grade Aid

Only available w/ high mount push or pull to lock wheel locks on pneumatic or pneumatic w/ airless tire options

☐	002GA1	E0974	Grade Aid	\$175.00
--------------------------	--------	-------	-----------	----------

COMING SOON

Anti Roll-Back

new	E2228	Automatic Bracking system - On chair	\$2,075.00
-----	-------	--------------------------------------	------------

PRATIKO Automatic Braking System. Installed by Sunrise Medical. Wheel Lock selection required - If no selection, Push to Lock will be provided. Weight capacity of 250 lbs (113 kg). Not available with transit, drum brake, seat pan, OneArmDrive, Reclining Back, Heavy Duty, or seat widths wider than 21". Available with the following rear wheels: 24" Lite Spoke and XCore Mags. Folding width maybe be affected.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price			
6. HANGER & FOOTREST												
Hangers					Lower Leg Length Setting (Extension Mount Only)							
Front mount hangers only available w/ aluminum angle adj footplates. ALR/ELR Hangers N/A w/ Push To Lock wheel locks.					Recommend 2" less than FSTFH. (*)denotes frame, 9" & 9.5" Hemi Only							
☐	002H145		60° Swing In/Out Hanger	NC	☐	002LO01	9.0"(*)	NC	☐	002LO13	15.0"	NC
☐	002H139		70° Swing In/Out Hanger	NC	☐	002LO02	9.5"(*)	NC	☐	002LO14	15.5"	NC
☐	002H141		80° Swing In/Out Hanger	NC	☐	002LO03	10.0"	NC	☐	002LO15	16.0"	NC
☐	002H154		60° Swing In/Out Hanger Hemi	NC	☐	002LO04	10.5"	NC	☐	002LO16	16.5"	NC
☐	002H155		70° Swing In/Out Hanger Hemi	NC	☐	002LO05	11.0"	NC	☐	002LO17	17.0"	NC
☐	002H140		70° Swing In/Out Front Mount	NC	☐	002LO06	11.5"	NC	☐	002LO18	17.5"	NC
☐	002H142		80° Swing In/Out Front Mount	NC	☐	002LO07	12.0"	NC	☐	002LO19	18.0"	NC
☐	002H158		90° Swing In/Out Front Mount	NC	☐	002LO08	12.5"	NC	☐	002LO20	18.5"	NC
☐	002H159	K0053	Articulating Hanger - Ext Mount N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates, 14"-22" width	\$431.00	☐	002LO09	13.0"	NC	☐	002LO21	19.0"	NC
☐	002H160	E0990	Elevating Hanger - Ext Mount N/A w/ Platform Flip Up and Alum Angle Adj Front Mount footplates, 14"-22" width	\$367.00	☐	002LO10	13.5"	NC	☐	002LO22	19.5"	NC
☐	002H161	E0990	Elevating Hanger - Front Mount Only available w/ front mount alum angle adj footplates and 14"-22" width	\$367.00	☐	002LO11	14.0"	NC	☐	002LO23	20.0"	NC
☐	002H32		Omit Hanger & Footplate	-\$15.00	☐	002LO12	14.5"	NC	☐	002LO24	20.5"Std. Only	NC
Articulating / Elevating Hanger Cover					Z-finity Hangers/Footrest System							
					(Maximum weight capacity for Z-Finity footrests is 165 lbs)							
					Some interference may exist in certain configurations. N/A with Heavy Duty							
☐	002DI47		Neoprene Cover	\$397.00	☐	002H162		Z-finity Hanger	\$408.00			
☐	002DI48		Gel Pad Cover	\$467.00								
Footplates					Z-finity Footrest Mount							
Heel loop or leg strap not included, See Footplate Options Section					☐	002DI22		Direct Hanger Mt/Fixed Frame - High	NC			
☐	002F1		Composite (13"-22"w)	NC	☐	002DI23		Direct Hanger Mt/Fixed Frame - Low	NC			
☐	002F271	K0040	Composite Angle Adj-(16"-22"w)	\$211.00	☐	002DI24		Extension Tube Mount	NC			
☐	002F3	K0040	Alum Angle Adj - Ext Mt (12"-22"w)	\$221.00								
☐	002F8	K0040	Alum Angle Adj - Front Mt(12"-22"w)	\$221.00								
☐	002F5		Alum Footplate - Fixed (14"-22"w) Available w/ Ext Mount ELR or ALR	NC								
☐	002F48	K0040	Platform Flip Up Angle Adj (12"-22"w) N/A w/ Hemi Hangers,ELR/ALR,Front Mnt	\$221.00								
☐	002F7		Omit Footplate & Ext Tube	-\$15.00								
Footplate Options (Ext Tube & Front Mount)					Z-finity Footrest Pivot Arm							
☐	002HL1	E0951	Heel Loop	\$47.00	☐	002DI20		Short 2"	NC			
☐	002HL5	K0038	Leg Strap	\$47.00	☐	002DI21		Long 3 1/4"	NC			
Footrest Extension Tubes					Z-finity Footplate							
Required w/ extension tube mounted hangers. See leg length matrix					Some interference may exist in certain configurations.							
☐	002E61		Short Hemi	NC	☐	002F320	K0040	Multi Angle Adjustable (Pair) (not avail. w/8"-9" widths)	\$216.00			
☐	002E1		Short	NC	☐	002F321	K0040	Angle Adjustable Platform Footboard	\$250.00			
☐	002E2		Medium	NC	☐	002F6		Omit Footplate	NC			
☐	002E5		Long	NC								
☐	002E16		High Mount N/A w/ one-arm drive option or ELR/ALR or 90° front mount. Only available w/angle adj footplates	\$88.00								
☐	002E23		2" - Front Mount Only N/A w/ 6", 7" or 8" casters with 90°hangers, or ELR's	\$63.00								
☐	002E24		4" - Front Mount Only N/A w/6", 7" or 8" casters w/ 90°hangers or ELR/ALR	\$63.00								
					Z-finity ExtensionTube Length (When Choosing Ext. Tube Mount)							
					☐	002E61		Hemi Short	NC			
					☐	002E1		Short	NC			
					☐	002E2		Medium	NC			
					☐	002E5		Long	NC			

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST / SEATING									

Back Type			
☐	002BT10	Standard	NC
☐	002BT20	Angle Adjustable Back (3-5" Pivot)	\$216.00
☐	002BT22	Angle Adjustable Back (5"-7" Pivot)	\$216.00
☐	002BT6	Depth Adjustable (Required with DAB)	NC

 Back Type - Build 4 Me Options			
Must choose DAB Back type with B4ME options			
☐	002BT127	E2398	B4ME - Dynamic Back
		☐	Fixed Back post
		☐	Std (19") ☐ Med (22") ☐ Tall (25")
		☐	Height Adjustable Stroller handles
☐	002BT128	E1226	B4ME - Reclining Backrest
			☐ Std (18"-23") ☐ Tall (19"-26")
			☐ Med (19") ☐ Tall (22")

 Back Post Options			
Omit push handles N/A w/ stabilizer bar. Folding HA push handle N/A w/8 deg or Low back height.			
Standard Back type options:			
☐	002BP1	Push Handle Straight	NC
☐	002BP24	Push Handle 8° Med or Tall only. N/A w/ Angle Adjustable Backrest	NC
☐	002BP2	Omit Push Handle Straight (N/A w/HD)	NC
☐	002BP76	Omit Push Handle 8° (N/A w/HD) Med or Tall only	NC
☐	002PH17	Folding Push Handle with Ht. Adj. Back N/A w/Angle Adj Back & Heavy D.	\$263.00

Angle Adj Back type options:			
☐	002BP5	Angle Adj. Straight Push Handles	NC
☐	002BP6	Angle Adj. Omit Push Handles, N/A w/HD	NC

Back Height (Std or Angle Adjustable*)			
*For Angle Adj Back see the angle Adj Backrest chart in the right column of the page.			
☐	002B1	Low (9"-12" Std Back height range)	NC
☐	002B5	Medium (12"-15" Std Back height range)	NC
☐	002B9	Tall (16"-19" Std Back height range)	NC

Back Height Setting (Std or Angle Adjustable)					
EXO upholstery not available with 21" back height					
☐	002UX4	9"	NC	☐	002UX18
☐	002UX6	10"	NC	☐	002UX20
☐	002UX8	11"	NC	☐	002UX22
☐	002UX10	12"	NC	☐	002UX24
☐	002UX12	13"	NC	☐	002UX26
☐	002UX14	14"	NC	☐	002UX28
☐	002UX16	15"	NC	(*)20" & 21" only w/ Ang Adj Back	

Stroller Handles (Removable Push Handles)			
Only available w/ omit push handles. N/A w/ DAB, amputee axle plate, push handles, or when angle adjustable back is selected, flip-back or flip-back height adjustable armrests, and JAY backs.			
☐	002SR2	Regular Stroller Handles	\$158.00
☐	002SR3	Tall Stroller Handles	\$158.00

Stabilizer Bar			
N/A with Folding Push Handle, or with Ht. Adj. Back or DAB Stroller handles			
☐	002SZ1	Stabilizer Bar Standard with HeavyDuty Package and Recline	\$123.00


Backrest Upholstery (Standard or Angle Adjustable)			
Std EXO upholstery and Tension adj binding is always black.EXO Tension Adj Uph N/A with Heavy Duty. Max back height 20"for EXO.			
☐	002BU1	Standard Upholstery	NC
☐	002BU88	Standard Upholstery EXO	\$58.00
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue			
☐	002BU70	E2611 Backrest Tension Adjust Ballistic *Reference 186BU70 for billing	\$367.00
☐	002BU89	E2611 Backrest Tension Adj EXO *Reference XXXBU89 for billing	\$387.00
Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue			
☐	002BU71	E2611 Backrest Tension Adjust 3DX Vented *Reference 186BU71 for billing	\$387.00
☐	002BU82	E2611 Jay GO Backrest N/A on Heavy Duty	\$475.00
☐	002BU14	Omit Back Upholstery	-\$15.00

Standard Backrest Option			
☐	002OD4	Reinforced Upholstery Only avail w/ Std Back Upholstery. Std w/HD package	\$58.00

Angle Adjustable Backrest Pivot Height			
Standard push handles only. Some interference w/ Padded SA arms. N/A w/ flip-back or flip-back height adjustable armrests (classic).			
☐	002AH1	Angle Adjustable Pivot 3"	NC
☐	002AH2	Angle Adjustable Pivot 4"	NC
☐	002AH3	Angle Adjustable Pivot 5"	NC
☐	002AH4	Angle Adjustable Pivot 6"	NC
☐	002AH5	Angle Adjustable Pivot 7"	NC

Preset Angle Adjustable Backrest (-15° to +3°) (Optional)			
☐	002BD5	Back Angle - 15°	NC
☐	002BD4	Back Angle - 12°	NC
☐	002BD3	Back Angle - 6°	NC
☐	002BD2	Back Angle - 0°	NC
☐	002BD1	Back Angle +3°	NC

Angle Adjustable Backrest Chart			
Back Tube Height	Low	Medium	Tall
Pivot Regular set at 3"	12" 13" 14"	15" 16" 17"	17" 18" 19"
Pivot Regular set at 4"	13" 14" 15"	16" 17" 18"	18" 19" 20"
Pivot Regular set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 5"	14" 15" 16"	17" 18" 19"	19" 20" 21"
Pivot Tall set at 6"	15" 16" 17"	18" 19" 20"	N/A
Pivot Tall set at 7"	16" 17" 18"	19" 20" 21"	N/A
Adjustable Pivot Point Measurement from Seat Sling			
If back height option is not chosen, middle option will be present			



✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
7. BACKREST / SEATING (cont)									
DAB Back Height					Military Service Patches (Optional)				
☐	002B66		Depth Adjustable 12"	NC	☐	002VU2		Army Patch	\$116.00
☐	002B67		Depth Adjustable 14"	NC	☐	002VU4		Air Force Patch	\$116.00
☐	002B68		Depth Adjustable 16"	NC	☐	002VU12		Coast Guard Patch	\$116.00
☐	002B69		Depth Adjustable 18"	NC	☐	002VU1		Marines Patch	\$116.00
☐	002B70		Depth Adjustable 20"	NC	☐	002VU3		Navy Patch	\$116.00
DAB Back Post Options					☐	002VU5		Purple Heart Patch	\$116.00
<i>Diagram of DAB in Appendix on page 9. SA Stroller Handle height range is 21.25"-28.5"</i>					☐	002VU9		Purple Heart with Army Patch	\$181.00
☐	002BP33		DAB - Swing Away Stroller Handles	\$116.00	☐	002VU11		Purple Heart with Air Force Patch	\$181.00
☐	002BP27		DAB - Post Standard (19")	NC	☐	002VU13		Purple Heart with Coast Guard Patch	\$181.00
☐	002BP28		DAB - Post Medium (22")	NC	☐	002VU8		Purple Heart with Marines Patch	\$181.00
☐	002BP29		DAB - Post Tall (25")	NC	☐	002VU10		Purple Heart with Navy Patch	\$181.00
DAB Backrest options					JAY Back (Optional)				
<i>DAB Solid Back NA w/23" & 24" frame widths</i>					☐	002JB01		J3 Back (See J3 Order Form)	Price List
☐	002BU28		Omit DAB Solid Back and Hardware	NC					
☐	002BU29		Omit DAB Solid Back Include Hardware	\$94.00					
☐	002BU23		Depth Adjustable Solid Back	\$227.00					

8. ARMREST / SIDE GUARDS

NEW!

Padded Swing Away

N/A w/ DAB, Amputee axle. Some limitations may exist with integral push handles, some RSTFHs and wheel sizes. **Now available w/Transit.

☐	002AA18	E0973	Padded Swing-Away STD (9" to 13")	\$263.00
☐	002AA176	E0973	Padded Swing-Away with Side guard Standard (9" to 13")	\$432.00
	Setting	☐ 9"	☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"	

Flip Back

N/A w/ DAB, Xtender, 12" - 14" frame depth, Classic Style armpad in highest/lowest axle RSH on rounded rear frame, or angle adjustable back. 4 Point (FS) J3 mounting bracket will cause interference with Flip Back Height Adj armrest. Flip back N/A with amputee axle plate for seat heights requiring the highest axle plate position

☐	002AA6	E0973	Flip Back Height Adj(Classic 8.5"-13.5")	\$367.00
☐	002AA5		Flip Back (incl. Desk Length Pad)	NC
☐	002AA65	E0973	Flip Back Height Adj (LITE 8.5"-13.5")	\$367.00

Flip Up (For use with Depth Adjustable Back)

Requires an Armpad or Foam from the Arm Pad section below

☐	002AA7	E0973	Length Adjustable Locking Flip-up	\$321.00
☐	002AA1	E0973	Adj. Locking Flip-up (Std. w/ foam)	\$309.00

Height Adjustable

☐	002AA11	E0973	Single Post Hgt. Adj. Std (8.5"-13.5")	\$333.00
☐	002AA29	E0973	Single Post Hgt. Adj. Low (7"-11")	\$333.00

Omit Armrest

☐	002AA14		Omit Armrest	NC
--------------------------	---------	--	--------------	----

Armpad for Height Adjustable

Ergo armpad only available with Single post Height adj & Flip back armrest

Ergo not available with Classic Flip Back armrest

☐	002AM23		Desk Length - Classic (10")	NC
☐	002AM24		Full Length - Classic (14")	NC
☐	002AM9		Armpad Foam - Short (Flip Up armrests)	NC
☐	002AM10		Armpad Foam - Long (Flip Up armrests)	NC
☐	002AM27		Ergo Pad (9.5"-13.5")	NC

Side Guard

Plastic side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist.

☐	002SG33	K0108	Plastic Standard (8.5" High)	\$169.00
☐	002SG29	K0108	Plastic Low (6.5" High)	\$169.00

9. KEY OPTIONS & ACCESSORIES									
Anti-Tip					Additional Cover & Pack Accessories				
☐	002AT1	E0971	Rear Anti-Tip (NC w/ Xtender)	\$175.00	☐	0028A02		Gel Cover Hanger/Frame	\$129.00
Backpack					☐	0028A01		Gel Cover Swing Away Arm	\$175.00
☐	002A1		Quickie Backpack Black	\$88.00	Positioning Belts				
☐	002A2		Quickie Kids Backpack Black	\$70.00	☐	002A141	E0978	Belt 1.5" Auto Buckle w/Clamps	\$82.00
Seat Pouch					☐	002A117	E0978	Belt 2" Aircraft Buckle w/Clamps	\$129.00
☐	0021A01		Seat Pouch Black /Charcoal	\$58.00	☐	002A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$139.00
JAY Positioning					☐	002A120		Positioning Belt Slide Clamps	\$58.00
☐	002K1	E1020	JAY Lower Extremity S/A support	\$460.00	Spoke Guards				
			☐ Left ☐ Right Position		<i>Only available w/Mountain, and Lite Spoke wheels,</i>				
☐	002K2		Omit Hanger	-\$10.00	<i>Not available with 5 spoke MAG wheels or Projection handrims</i>				
			☐ Left ☐ Right		☐	0024A02	K0065	Quickie Clear Spoke Guards	\$134.00
☐	002K7		Pad Size	NC	☐	0024A01	K0065	Quickie Black Spoke Guards	\$134.00
			☐ 6" x 8" ☐ 8" x 10"						

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
9. KEY OPTIONS & ACCESSORIES (cont)									
Tray Tables					Miscellaneous				
☐	002A30	E0950	Small (13"-14"W)	\$298.00	<i>O2 holder N/A with Xtender or 12" - 13" frame width. Crutch holder N/A with EXO Tension Adj upholstery and Omit push handle</i>				
☐	002A31	E0950	Medium (15"-17"W)	\$298.00					
☐	002A32	E0950	Large (18"-20"W)	\$311.00	☐	0028A03		Caddy (N/A w/ front mount hangers)	\$58.00
☐	002A33	E0950	X-Large (21"-22"W)	\$311.00	☐	002A2100	E2207	Crutch Holder	\$227.00
Hanger Options					☐	002LS14	K0038	Extra Leg Strap Black	\$42.00
<i>Aluminum material. This options comes uninstalled.</i>					☐	002IG2		Impact Guard Neoprene	\$105.00
<i>Note: During swing In/Out process may wear out the paint on the frame over time.</i>					☐	002A73	K0105	IV Holder	\$321.00
☐	002IE05	K0108	Enhanced hanger latch	\$79.00	☐	002A145	E2208	O2 Holder	\$362.00
					☐	0028A21		Tool Kit	\$63.00
					☐	002A999		Quickie Water Bottle & Holder	\$63.00

10. COLOR

Standard Frame Colors				
☐	002C214		Black (Gloss)	NC
☐	002C215		Black (Matte)	NC
☐	002C211		Celestial Sapphire	NC
☐	002C212		Desert Sage	NC
☐	002C213		Desert Tan	NC
☐	002C216		Metallic Shadow	NC
☐	002C217		Midnight Blue	NC
☐	002C218		Mystic Purple	NC
☐	002C219		Petal Pink	NC
☐	002C220		Pink Power	NC
☐	002C221		Red Pepper	NC
☐	002C222		Royal Blue	NC
☐	002C223		Salute Green	NC
☐	002C224		Scarlet Noir	NC
☐	002C225		Silver Ore	NC
☐	002C226		Solar Orange	NC
☐	002C227		Starlight White	NC
☐	002C228		Storm Gray	NC
☐	002C229		Turbo Green	NC
☐	002C230		Velocity Yellow	NC

NEW!

Legacy - Frame Colors				
☐	002C141		Black Cherry	NC
☐	002C10		Black Opal	NC
☐	002C14		Candy Purple	NC
☐	002C142		Electric Blue	NC


Personalized Frame Color				
☐	002C112		Desert Camouflage	\$460.00
☐	002C43		Mossy Oak Camouflage	\$460.00
☐	002C42		Stars & Stripes	\$460.00
☐	002C108		Zebra <i>(Choose above base color) Std with White unless above base paint is chosen N/A w/ Black</i>	\$460.00

Anodized Color Package				
Standard Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog Legs.				
☐	002EB07		Black - Caster only	\$58.00
☐	002EB08		Blue	\$187.00
☐	002EB13		Gold	\$187.00
☐	002EB09		Lime Green	\$187.00
☐	002EB10		Purple	\$187.00
☐	002EB11		Red	\$187.00
☐	002EB12		Ti Grey	\$187.00

Additional Options Through Built for Me				
To inquire about any unique personalized Quickie products, contact: Phone: (800) 333 - 4000, Fax: (800) 333 9011				

</

11. APPENDIX

Frame Depth matrix					Depth Adjustable Back Frame				
Seat Depth									
Frame Depth	12"	13"	14"	15"					
13"	S	S							
14"	S	S	S						
15"	S	S	S	S					
16"		S	S	S	S				
17"			Hemi*	S	S	S	S		
18"				S	S	S	S		
19"					Hemi*	S	S	S	
20"						S	S	S	S
21"							Hemi*	S	S
22"								S	S

"S" = Standard Frame Style
"H" = Hemi Frame Style
*Only (6) Frame Depths avail. on Hemi Style Frame

Seat Depth and Frame Depth selections based on Standard Backrest position

11. APPENDIX

Lower Leg Length Matrix

2" & 4" extensions for front mount hangers shipped separate, not mounted. Front seat height must be at least 2" taller than lower leg length on all swing in/out hangers. Front seat height must be at least 2" taller than maximum lower leg length on all front mount hangers. All values shown apply to Standard frame style, subtract 1.5" to get Hemi equivalent. More than 2" will require on field adjustment.

Hanger	Extension	Seat Sling to Top of Footrest in inches		Min Allowable FSH
		Min	Max	
60° Hemi Swing In/Out	Short Hemi	10	13	14
	Short	-	-	-
	Med	12.5	15.5	15.5
	Long	14.5	17.5	17.5
70° Hemi Swing In/Out	Short Hemi	11	14	14
	Short	-	-	-
	Med	13.5	16.5	16
	Long	15.5	18.5	18
60° Swing In/Out	Short	11.5	15	15
	Med	11.5	17	15
	Long	14	18.5	17
70° Swing In/Out	Short	12	16.5	15
	Med	12.5	18.5	15.5
	Long	14.5	20.5	17.5
80° Swing In/Out	Short	13.5	17	15.5
	Med	13.5	19	15.5
	Long	15.5	20.5	17.5
60° Swing In/Out w/Platform Flip Up	Short	11.5	14.5	14
	Long	14	18.5	16
70° Swing In/Out w/Platform Flip Up	Short	13	15.5	15
	Long	15.5	20	17.5
80° Swing In/Out w/Platform Flip Up	Short	13.5	14.5	15.5
	Long	15	19	17
70° Front Mount	No Ext	6	11.5	14
	2"	6	14	16
	4"	6	16	18
80° Front Mount	No Ext	6	12	14
	2"	6	14.5	16.5
	4"	6	16.5	18.5
90° Front Mount	No Ext	5.5	11.5	14
ALR	Short	15	16.5	17
	Med	15	18.5	17
	Long	16.5	20.5	18.5
ELR Ext. Mount	Short	12	16.5	15.5
	Med	12	18.5	15.5
	Long	13.5	20.5	17
ELR Front Mount	No Ext	5.5	9.5	13

Bckrest Height Matrix

Back Type	Ang-Adj Pivot Ht Setting	Push Handles	Backrest Height	Avail. Ht. Range
Standard - (Fixed)	N/A	Integral	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"
		Omit	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"
Angle Adjustable (3"-5" Pivot)	3"	Integral	Low	12"-14"
			Medium	15"-17"
			Tall	17"-19"
	4"	Integral	Low	13"-15"
			Medium	16"-18"
			Tall	18"-20"
	5"	Integral	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
	3"	Omit	Low	12"-14"
			Medium	15"-17"
			Tall	17"-19"
	4"	Omit	Low	13"-15"
			Medium	16"-18"
			Tall	18"-20"
	5"	Omit	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
Angle Adjustable (5"-7" Pivot)	5"	Integral	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
	6"	Integral	Low	15"-17"
			Medium	18"-20"
			Tall	N/A
	7"	Integral	Low	16"-18"
			Medium	18"-20"
			Tall	N/A
	5"	Omit	Low	14"-16"
			Medium	17"-19"
			Tall	19"-20"
	6"	Omit	Low	15"-17"
			Medium	18"-20"
			Tall	N/A
	7"	Omit	Low	16"-18"
			Medium	18"-20"
			Tall	N/A
Folding Push Handle - (Stroller Handle w/ Knuckle)	N/A	N/A	Low	9"-12"
			Medium	12"-15"
			Tall	16"-19"

11. APPENDIX

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"

* Values Given are Estimates. Actual seat to floor height may vary based on seat depth and front seat to floor height selection.

All 3" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	17.5	18.0
Multi Position Tall	1 1/2"	-	-	-	18.0	18.5
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
3"	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	-	15.0
5 1/4"	1/2"	-	-	-	-	15.5
	1"	-	-	-	-	16.0
	1 1/2"	-	-	-	-	16.5
6"	Std	-	-	-	-	17.5
	1/2"	-	-	-	-	18.0
	1"	-	-	-	-	18.5
7"	1 1/2"	-	-	-	-	19.0
	Std	-	-	-	-	18.0
	1/2"	-	-	-	-	18.5
8" x 2"	1"	-	-	-	-	19.0
	1 1/2"	-	-	-	-	19.5
	Std	-	-	-	-	19.0
Frog Legs	1/2"	-	-	-	-	16.5
	1"	-	-	-	-	17.0
	1 1/2"	-	-	-	-	17.5
Integral Multi Pos Forks (Low)	Std	-	-	-	-	15.5
	1/2"	-	-	-	-	16.0
	1"	-	-	-	-	16.5
Integral Multi Pos Forks (Std)	1 1/2"	-	-	-	-	17.0
	Std	-	-	-	-	16.5
	1/2"	-	-	-	-	17.0
	1"	-	-	-	-	17.5

All 4" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	17.5	18.0
Multi Position Tall	1 1/2"	-	-	-	18.0	18.5
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
3"	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	20.5	21.0
5 1/4"	Std	-	-	-	17.0	17.5
	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5
6"	1 1/2"	-	-	-	18.5	19.0
	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	19.5	20.0
7"	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
8" x 2"	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	20.5	21.0
	1/2"	-	-	-	-	-
Frog Legs	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
Integral Multi Pos Forks (Low)	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

All 5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	16.5	17.0
	1/2"	-	-	-	17.0	17.5
	1"	-	-	-	17.5	18.0
Multi Position Tall	1 1/2"	-	-	-	18.0	18.5
	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
3"	1"	-	-	-	19.5	20.0
	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	20.5	21.0
5 1/4"	1/2"	-	-	-	17.5	18.0
	1"	-	-	-	18.0	18.5
	1 1/2"	-	-	-	18.5	19.0
6"	Std	-	-	-	19.0	19.5
	1/2"	-	-	-	19.5	20.0
	1"	-	-	-	20.0	20.5
7"	1 1/2"	-	-	-	20.5	21.0
	Std	-	-	-	21.0	21.5
	1/2"	-	-	-	-	-
8" x 2"	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-

All 6" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	17.5	18.0
	1/2"	-	-	-	18.0	18.5
	1"	-	-	-	18.5	19.0
Multi Position Tall	1 1/2"	-	-	-	19.0	19.5
	Std	-	-	-	19.5	20.0
	1/2"	-	-	-	20.0	20.5
3"	1"	-	-	-	20.5	21.0
	1 1/2"	-	-	-	21.0	21.5
	Std	-	-	-	21.5	22.0
5 1/4"	Std	-	-	-	18.5	19.0
	1/2"	-	-	-	19.0	19.5
	1"	-	-	-	19.5	20.0
6"	1 1/2"	-	-	-	20.0	20.5
	Std	-	-	-	20.5	21.0
	1/2"	-	-	-	21.0	21.5
7"	Std	-	-	-	21.5	22.0
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
8" x 2"	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
Frog Legs	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
Integral Multi Pos Forks (Low)	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

All 7" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	19.0
	1/2"	-	-	-	-	19.5
	1"	-	-	-	-	20.0
Multi Position Tall	1 1/2"	-	-	-	-	20.5
	Std	-	-	-	-	21.0
	1/2"	-	-	-	-	21.5
3"	1"	-	-	-	-	22.0
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
5 1/4"	1/2"	-	-	-	-	19.0
	1"	-	-	-	-	19.5
	1 1/2"	-	-	-	-	20.0
6"	Std	-	-	-	-	20.5
	1/2"	-	-	-	-	21.0
	1"	-	-	-	-	21.5
7"	1 1/2"	-	-	-	-	22.0
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
8" x 2"	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-

All 8" x 1.5" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Multi Position Tall	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
3"	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
5 1/4"	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
6"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
7"	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
8" x 2"	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
Frog Legs	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	1 1/2"	-	-	-	-	-
	Std	-	-	-	-	-
	1/2"	-	-	-	-	-

11. APPENDIX

Quickie 2 Seat To Floor Height Matrix - Standard Frame *Hemi options will result in a reduction of approx 2 1/2"
** Values Given are Estimates. Actual seat to floor height may vary based on seat depth and front seat to floor height selection.*

All 8" x 2" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Multi Position Tall	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
6"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
7"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
8" x 2"	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
	1 1/2"	-	-	-	21.5	22.0
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

All Other 8" Casters						
Fork Size	Stem Size	Bolt Position				
		5	4	3	2	1
Multi Position Short	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Multi Position Tall	Std	19.5	20.0	20.5	21.0	21.5
	1/2"	-	-	-	21.5	22.0
	1"	-	-	-	22.0	-
	1 1/2"	-	-	-	-	-
3"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
5 1/4"	Std	-	-	-	-	20.0
	1/2"	-	-	-	-	20.5
	1"	-	-	-	-	21.0
	1 1/2"	-	-	-	-	21.5
6"	Std	-	-	-	20.0	20.5
	1/2"	-	-	-	20.5	21.0
	1"	-	-	-	21.0	21.5
	1 1/2"	-	-	-	21.5	22.0
7"	Std	-	20.0	20.5	21.0	21.5
	1/2"	-	-	-	21.5	22.0
	1"	-	-	-	22.0	-
	1 1/2"	-	-	-	-	-
8" x 2"	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Frog Legs	Std	-	-	-	-	-
	1 1/2"	-	-	-	-	-
Integral Multi Pos Forks (Low)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-
Integral Multi Pos Forks (Std)	Std	-	-	-	-	-
	1/2"	-	-	-	-	-
	1"	-	-	-	-	-

Handrims	Wheels				
	5 Spoke Mag	Lite Spoke	Ultra Lightweight	Spinergy	Mountain
Aluminum Anodized	20,22,24	20,22,24,26	22,24,25	22,24,25,26	24
Plastic Coat	20,22,24	22,24	-	24,26	-
Natural Fit - Standard Grip	24	22,24,25,26	-	22,24,25,26	-
Natural Fit - Super Grip	24	22,24,25,26	-	22,24,25,26	-
Q Grip Neoprene Coated	-	24,25	-	24,25	-
Ellipse 3R	24	24,25	-	24,25	-
Gekko	-	24,25	-	24,25	-
Curve L - Tetra Grip	-	24,25	-	24,25	-
Omit Handrim	24	24	24	24	24



Sunrise Medical (US) LLC
12002 Volunteer Blvd. Mount Juliet, TN 37122