




## ORDER FORM

To see QUICKIE Power Options  
or Accessories, scan here or visit:  
[https://www.sunrisemedical.com/  
power-wheelchair-options-accessories](https://www.sunrisemedical.com/power-wheelchair-options-accessories)




# Q300M Mini with CAPTAIN'S SEAT EIPW40

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		State	Zip
Mark For		Phone	

## Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



## Notes

## Chair Specifications

HCPC	K0848 / K0849 / K0856 / K0857
Weights starting at	273 lbs.
Power Seating	30° CG Tilt / 12" seat elevator
Seat Type	Captain's Seat w/25° Adjustable Recline
Seat Size Range	Width 16" to 20" Depth 16" to 22"
Seat-to-Floor Height Range	17.5" to 19"

## Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 20.5", Length 33"
Minimum Turning Radius	17" (excludes footplates)
Drive Control	VR2 90A / R-net PM 90A
Driving Range	Up to 17 miles
Drive Motors	4-pole 6 mph
Drive Wheels	12" w/flat-free inserts (black) 14" w/flat-free inserts (black)
Casters	7" front / 7" rear (black)
Charger	8A Dual Mode (Compact)
Battery Size	50 Ah Group 22
Color	9 shroud colors and 7 color accents
Maximum Ground Clearance	3"

BASE			
OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW40</i>			\$
400PC296	K0849	<b>Quickie Q300 M NCC No Power Captain's Seat and Back</b> Included in Base Pricing : Q300 M Base, VR2 90 amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8 amp Charger. Max Speed = 6 mph/10 kph, Weight Capacity = 300 lbs., STFH = 17.5".	\$8,374.00
400PC297	K0857	<b>Quickie Q300 M SCC Single Power Option Captain's Seat and Back</b> Included in Base Pricing : Q300 M Base, VR2 90 amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8 amp Charger. Max Speed = 6 mph/10 kph, Weight Capacity = 300 lbs., STFH = 17.5".	\$8,374.00
400PC298	K0848	<b>Quickie Q300 M NCC No Power Captain's Seat Back and Solid Seat</b> Included in Base Pricing : Q300 M Base, VR2 90 amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8 amp Charger. Max Speed = 6 mph/10 kph, Weight Capacity = 300 lbs., STFH = 17.5".	\$8,374.00
400PC299	K0856	<b>Quickie Q300 M SCC Single Power Option Captain's Seat Back and Solid Seat</b> Included in Base Pricing : Q300 M Base, VR2 90 amp Electronics for No Power Seat Functions, 6-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8 amp Charger. Max Speed = 6 mph/10 kph, Weight Capacity = 300 lbs., STFH = 17.5".	\$8,374.00

**User Weight** \$

400WG23	–	Up to 300 lbs.	NC
---------	---	----------------	----

**Motor Package** \$

400PF99	–	4-pole 10 kph / 6 mph	NC
---------	---	-----------------------	----

**Drive Wheels**  P. 10 \$

400DW13	–	Drive Wheel, 12", Solid, Knobby, Black	NC
400DW58	–	Drive Wheel, 14", Dual Density, Knobby, Black	NC

**Caster Wheels**  P. 10 \$

400CW59	–	Caster Wheels 7", Solid, Black	NC
---------	---	--------------------------------	----

**Shroud Color** \$

400C119	–	Matte Black	NC
400C165	–	Candy Apple Red	NC
400C166	–	Candy Blue	NC
400C167	–	Candy Purple	NC
400C176	–	Green Apple	NC
400C178	–	Pearl Pink	NC
400C182	–	Silver	NC
400C125	–	White	NC
400C170	–	American Flag	\$590.00

**BASE**

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

**Color Accents** \$

*Note: Color Accents and Shroud Colors may not match.*

400EZ20	–	Apple Green	NC
400EZ22	–	Black	NC
400EZ21	–	Blue	NC
400EZ23	–	Orange	NC
400EZ24	–	Red	NC
400EZ25	–	Silver	NC
400EZ19	–	White	NC

**ISO Transit Tie Downs** \$

400T050	–	ISO Transit Package	NC
---------	---	---------------------	----

**LED Lights**  P. 10 \$

400LI04	–	LED Lights and Indicators	\$1,373.00
400LI08	–	LED Fender Lights	\$555.00
400LI02	–	Omit LED Lights	NC

**Seat Heights** \$

*Note: Available seat-to-floor heights are dependent on 2 options: 1- Battery Size and 2- Powered Seating. Use the matrix below to determine available ranges for your desired options. Measured from front of seat pan to floor with 0° frame angle.*

		Seat Height	
400HT104	–	17.5"	NC
400HT56	–	18"	NC
400HT105	–	18.5"	NC
400HT57	–	19"	NC

**Batteries and Chargers**  P. 10 \$

**Group 22 Battery**


400AG35	E2363	GP22 / 50Ah Flat Top Pair	\$667.00
---------	-------	---------------------------	----------

**Omit Batteries**

400AG24	–	Battery Harnesses Only ( <i>Omit Batteries</i> )	NC
---------	---	--	----

**Chargers**

400BJ30	–	10A Dual Mode (Compact)	\$108.00
400BJ19	–	8A Dual Mode (Compact)	NC

SEATING				
OPTION	HCPC	DESCRIPTION		RETAIL
Power Seating Options  P. 10				\$
<b>Seat Elevator</b>				
400PI07	E2300	Seat Elevator, Powered, 300 mm / 12"		\$3,717.00
400PI08	E2300	C-Me Seat Elevator, Powered, 300 mm / 12" (C-Me speed up to 3.0 mph) (N/A with VR-2)		\$3,717.00
<b>Seat Tilt</b>				
400PA46	E1002	Posterior 30° Tilt		\$7,720.00
Pre Set Seat Frame Angle (No Power Seating Options, Increases Seat to Floor Height)				\$
400FA10	–	0° (Adds 0" to STFH)		NC
400FA3	–	3° (Adds 1" to STFH)		NC
Backrest Frame Type				\$
400BF08	–	Captain Seat Frame (Includes Back Cushion)		NC
Seat Cushion Type				\$
400ST77	–	Captain's		NC
400ST80	–	JAY®		NC
400ST75	–	Omit		NC
Seat Size				\$
<b>Seat Pan Width</b>			<b>Available Depths</b>	
400R53	–	16"	16" to 18"	NC
400R55	–	18"	16" to 20"	NC
400R57	–	20"	18" to 22"	NC
<b>Seat Pan Depth (Front of Seat Pan to Back Structure. Back Cushions Can Cause 2" Loss of Effective Depth)</b>				
400RF51	–	16"		NC
400RF45	–	17"		NC
400RF46	–	18"		NC
400RF47	–	19"		NC
400RF48	–	20"		NC
400RF49	–	21"		NC
400RF50	–	22"		NC
JAY® Seat Cushion  P. 10				\$
400CU59	E2607	JAY® Ion® (Sizes 16" - 20" Wide)		\$451.00
400CU61	E2607	JAY® Union® (Sizes 16" - 20" Wide)		\$510.00
400CU64	E2622	JAY® Fusion® (Sizes 16" - 20" Wide)		\$504.00
400CU104	E2622	JAY® Fusion® with Cryo Technology (Sizes 16" - 20" Wide)		\$618.00

## POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
<b>Armrest Structure</b>			\$
400EG37	–	Flip Back <i>(Includes Waterfall Arm Pads) (Height Adjustable 11" to 14")</i>	NC
<b>Fixed Centermount Legrest</b>			\$
400H184	–	Fixed Centermount Legrest, Short <i>(Includes Calf Pads and Single Footplate) (11.5" - 14.5" Leg Length)</i>	NC
400H185	–	Fixed Centermount Legrest, Long <i>(Includes Calf Pads and Single Footplate) (13.5" - 18.5" Leg Length)</i>	NC
400H136	–	Angle Adjustable Centermount Legrest	\$515.00
<b>Power Centermount Legrest  P. 10</b>			\$
400H109	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i>	\$3,717.00
<b>Centermount Legrest Accessories</b>			\$
<b>Power Centermount Footplates</b>			
400F127	K0040	Dual, Height, Width and Angle Adjustable, Standard (14"W x 10"D)	\$193.00
400F128	K0040	Dual, Height, Width and Angle Adjustable, Large (16"W x 12"D)	\$193.00
<b>Angle Adjustable / Fixed Centermount Footplates</b>			
400F45	–	Single / Dual Footplate	NC
<b>Power Centermount Calf Pads</b>			
400DI34	–	Calf Pads, Pair, Standard Integral Foam	NC
400DI35	E0995	Calf Pads, Pair, Padded <i>(N/A with standard footplate)</i>	\$400.00
<b>Angle Adjustable Centermount Calf Pads</b>			
400DI49	–	Calf Pads, Pair, Standard	NC
400DI50	–	Extra Large Single Calf Pad	\$172.00
<b>SEDEO PRO Swing-Away Legrest (Pair)</b>			\$
400H177	–	Legrests, Swing-Away, 70°, SEDEO PRO <i>(14" to 21" Leg Length)</i>	NC
400H178	E0990	Legrests, Manual Elevating, Swing-Away, SEDEO PRO <i>(Includes Padded Calf Pads) (14.5" to 21.5" Leg Length)</i>	\$400.00
400H179	E1010	Legrests, Power Elevating, Swing-Away, SEDEO PRO <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$2,459.00
400H176	–	Omit Legrests	NC
<b>Swing-Away Legrest Footplates</b>			
400F369	K0040	Angle Adjustable Flip-Up Footplates	\$193.00
<b>Swing-Away Legrest Calf Supports</b>			
400DI42	–	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending on chair width, x 9.5" Tall)</i>	\$451.00
<b>Swing-Away Legrest Options</b>			
400DI44	E0951	Heel Loops	\$57.00
400KG4	–	Knee Pads <i>(2.5" Wide x 4.5" Tall)</i>	\$229.00
<b>SEDEO Lower Extremity Swing-Away Support</b>			
400VD01	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO, Right	\$451.00
400VD02	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO, Left	\$451.00

## POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
<b>Pelvic Positioning Belts</b>  P. 10			
400A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
400A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
400A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$86.00
400A845	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$199.00
400A846	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$220.00

### Accessories

\$

#### **O2 Tank Option** (Note: Only one can be mounted on chair)

400A718	E2208	O2 Tank Holder Mounted to Backrest (D and E Styles)	\$372.00
---------	-------	---	----------

#### **Captain Accessories** (Note: Only one can be mounted on chair)

400A608	E2207	Cane / Crutch Holder	\$297.00
400A610	K0108	Walker Holder	\$297.00

#### **Misc. Accessories** (Mounts in front of armpad and opposite of joystick)

4008A52	–	Ram Mount Phone	\$98.00
4008A53	–	Ram Mount Tablet	\$400.00
4008A54	–	Ram Mount Cup Holder	\$98.00
4008A58	–	Rear View Mirror	\$98.00
4009A10	–	Utility Hooks	\$59.00

#### **Medial Thigh Support** (Available with S/A Hangers)

4008A80	E0957	Medial Thigh Support	\$193.00
4008A81	E1028	Pommel Swing-Away Hardware	\$286.00

**JOYSTICK / ELECTRONICS**

OPTION	HCPC	DESCRIPTION	RETAIL
<b>VR-2 Joystick (Remote)</b>			<b>\$</b>
400J127	–	VR-2 Joystick	NC
<b>R-net Joystick (Remote)  P. 11</b>			<b>\$</b>
400J121	–	Standard LED Joystick ( <i>w/o Screen Display</i> )	NC
400J116	–	Color Joystick ( <i>2.3" Color LCD Screen Display, CJSM1</i> )	NC
400J118	–	High Resolution Color Joystick, (CJSM2+BT) <i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>	\$1,206.00
400J119	–	Omit Remote Control	NC
<b>Joystick Mounting</b>			<b>\$</b>
400JM34	–	Fixed on Armrest	NC
400JM35	E1028	Parallel Swing-Away	\$338.00
400JM36	E1028	Link-It Mid Line Mount on Armrest ( <i>Swing-Away</i> )	\$972.00
<b>Joystick Mount (Location)</b>			<b>\$</b>
400HM16	–	Right Hand Armrest	NC
400HM17	–	Left Hand Armrest	NC
<b>Display Language Selection R-net</b>			<b>\$</b>
400L1	–	English	NC
400L3	–	French	NC
400L7	–	Spanish	NC
<b>Joystick Handle (Optional)  P. 11</b>			<b>\$</b>
400JH13	E2323	Ball Handle	\$144.00
400JH15	E2323	T-Handle	\$144.00
400JH24	E2323	Bodypoint - Rubber Dome Joystick Handle	\$187.00
400JH25	E2323	Bodypoint - U-Shaped Joystick Handle ( <i>4 in</i> )	\$187.00
400JH26	E2323	Bodypoint - Grooved Mushroom Joystick	\$187.00
<b>Electronic Programmers  P. 11</b>			<b>\$</b>
400PR10	–	R-net PC Programmer w/ Dongle	\$601.00
<b>R-net Power Modules  P. 11</b>			<b>\$</b>
400ER12	–	Power Module, VR-2, 90A, Non-Expandable	NC
400ER13	–	Power Module, R-net, 90A, Expandable <i>Required for specialty controls or 3+ actuators</i>	\$1,309.00
400EX11	E2377	Expandable 90A Controller - R-net	\$692.00
400EX12	E2313	Harness for 90A Expandable Controller	\$617.00



**JOYSTICK / ELECTRONICS**

OPTION	HCPC	DESCRIPTION	RETAIL
--------	------	-------------	--------

**Seating Control (Button Switches will not be mounted)**  P. 11

\$

**Powered Seating Control**

400E040	–	Omit Thru-Drive <i>(Requires Switch Control)</i>	NC
400E037	E2310	Thru-Drive Control for 1 Actuator	\$2,196.00
400E038	E2311	Thru-Drive Control for 2+ Actuators	\$3,975.00

**Switch Control Only**

400E033	–	Dual Toggle for 1 Actuator	NC
400E034	–	2.5" Single Button Switch for 1 Actuator	NC
400E041	–	Ctrl+5 Switch Box for 2 Actuators or More	NC

**Additional Power Seating Control (Available with Thru-Drive only)**

400QG11	–	Ctrl+5 Switch Box with Free Assignable Functions	\$1,430.00
---------	---	--	------------

**Dual Toggle Mounting Location**

400TN1	–	Toggle Mount, Right	NC
400TN2	–	Toggle Mount, Left	NC

**Ctrl+5 Mounting Location**

400EY19	–	Right Armrest	NC
400EY20	–	Left Armrest	NC

**Ctrl+5 Switch Type**

400EY22	–	Push Buttons	NC
400EY23	K0108	Toggle Switches	\$343.00

**R-net Programming**

\$

**Assignable Buttons**

*Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.*

400PM12	–	Include Assignable Buttons	NC
---------	---	----------------------------	----

**Power Centermount Footplate Extension Option**

400IE03	K0108	Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only available with Powered Articulation.</i>	\$1,430.00
---------	-------	---	------------



Q300 M with Captain's Seat

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year

