



Quickie QS5 X

July 2026

EIF7



Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

Submitting for:

☐ Quote

☐ Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME

E1236 coding still applies for EIF7 with sizes 14" or less. 300 lbs weight capacity

☐	170WM112	K0005	Quickie QS5 X	\$2,550.00
☐	170WM113	K0005	Quickie QS5 X - Transit included	\$2,924.00

Frame Style

See figure 1. Measure from front frame to floor

☐	170FS31	Swing Away Frame	NC
☐	170FS43	Fixed Front 75°	NC
☐	170FS44	Fixed Front 85°	NC

Frame Inset (Fixed Front Only)

See **FIGURE 2** Frame inset is the distance from the outside of the seat tube to the outside of the footrest. Footrest width = frame width minus 2. N/A for 14" and 15" Wide using composite angle adj footrest. **75° only available with 0" inset.**

☐	170FF1	0"	NC
☐	170FF2	1"	NC

Transit Option

ISO 7176-19 Transit Approved

☐	170TO49	E1022	Transit System	\$374.00
--------------------------	---------	-------	----------------	----------

Frame Options

Ergo Touch Points is a rubber cover on the knee bend area of the frame tubes

☐	1709A11	Ergo Touch Points Swing Away hangers or Fixed front	\$150.00
☐	1709A12	Frame Fold-N-Lock mechanism	NC
☐	1709A06	Twin LED lights	\$248.00



QS5 X - Swing Away



QS5 X - Fixed Front

Power Assist

To order any Empulse Power Assist product please see the Order Form for Part #'s

Option	HCPCS	DESCRIPTION	RETAIL
☐ 209MW91	E0986	Empulse R90 Includes: Device Charger, Wireless Control Box, V-Lock Latch, Travel Bag, 2.5Ahr Removable Battery	\$7,300.00
☐ 209MW92	E0986	Empulse M90 M90 Wheels, Dual Controller, Two M90 Batteries, Mounting Hardware, Quick Release Axles, and Device Charger	\$7,950.00
☐ 147MW95	E0986	Empulse F35 Includes: F35 Stick, Subframe mount, Frame clamps kit, and Charger	\$6,500.00



2. SIZE / SEATING

Frame Width

See **FIGURE 2** Measured from outside of seat tube to the outside of opposite seat tube.

☐	170W4	13"	NC
☐	170W5	14"	NC
☐	170W6	15"	NC
☐	170W7	16"	NC
☐	170W8	17"	NC
☐	170W9	18"	NC
☐	170W10	19"	NC
☐	170W11	E2201 20"	\$583.00

Seat Sling

Standard seat sling includes utility bag sewn at the bottom of upholstery.

☐	170SL1	Standard Seat Sling	NC
☐	170SL6	Omit Seat Sling	-\$15.00

Figure 1




Figure 2



Seat Depth

See **Seat Squeeze Matrix** in Appendix section. Measured from back canes to front edge of upholstery. Frame Fold-N-Lock mechanism NA with 13"-14".

☐	170SD6	13"	NC
☐	170SD7	14"	NC
☐	170SD8	15"	NC
☐	170SD9	16"	NC
☐	170SD10	17"	NC
☐	170SD11	18"	NC
☐	170SD12	19"	NC
☐	170SD13	E2203 20"	\$641.00

NEW!

JAY seat cushions (On-chair)

By opting in the below options the cushions will match the chair width/depth.

☐	170JC15	E2622	JAY X2	\$501.00
☐	170JC17	E2622	JAY Fusion Reduced Profile / Fluid	\$515.00
☐	170JC18	E2607	JAY Union Reduced Profile	\$520.00
☐	170JC10	E2607	JAY ION	\$462.00
☐	170JC06	E2622	JAY LitePro	\$525.00

JAY Cushions (Optional)

To order any other size JAY cushions, please see the JAY price list.

☐	170JC20	Cushion Name or Part #	Price List
		Width: _____ Depth: _____	

3. CASTERS / SEAT HEIGHT

Front Seat Height

☐	170SH41	14.0"	NC	☐	170SH18	18.0"	NC
☐	170SH42	14.5"	NC	☐	170SH20	18.5"	NC
☐	170SH12	15.0"	NC	☐	170SH21	19.0"	NC
☐	170SH13	15.5"	NC	☐	170SH23	19.5"	NC
☐	170SH14	16.0"	NC	☐	170SH24	20.0"	NC
☐	170SH15	16.5"	NC	☐	170SH26	20.5"	NC
☐	170SH16	17.0"	NC	☐	170SH27	21.0"	NC
☐	170SH17	17.5"	NC	☐	170SH29	21.5"	NC
				☐	170SH36	22.0"	NC


Caster Fork

☐	170CF5	Multiple Position	NC
☐	170CF69	Multiple Position Single Sided Fork	\$197.00
☐	170CF53	E1015 Frog Legs 0 weight to 100 lbs.	\$507.00
☐	170CF54	E1015 Frog Legs 101 lbs. weight to 190 lbs.	\$507.00
☐	170CF55	E1015 Frog Legs 191 lbs. weight to 265 lbs.	\$507.00

Multi Position

Single side

Frog Leg



Caster Wheels (option A)

☐	170CW3	4" Low Profile Polyurethane	NC
☐	170CW77	4" x 1.25" Semi Pneumatic	NC
☐	170CW5	5" Low Profile Polyurethane	NC
☐	170CW78	5" x 1.5" Semi-Pneumatic	NC
☐	170CW34	6" Polyurethane	NC
☐	170CW33	E2214 E2215 6" Pneumatic	\$88.00
☐	170CW79	6" x 1.5" Semi Pneumatic	NC
☐	170CW16	7" Polyurethane	NC
☐	170CW32	8" Polyurethane	NC

Caster Wheels colors (option B)

Casters are silver color unless you choose color option. Refer to Anodize package section on Page 7 for Anodize color selection. Match Ano color pkg pick is N/A with Black option.

☐	170CW64	4" x 1.5" Aluminum Soft Roll	\$216.00
☐	170CW66	5" x 1.5" Aluminum Soft Roll	\$227.00
☐	170CW68	6" x 1.5" Aluminum Soft Roll	\$239.00
Colors:	☐	Match anodize color package	Package price
	☐	Match frame color	\$57.00
	☐	Black anodize caster only	\$58.00

Seat Squeeze Matrix

Seat Depth	Max Seat Squeeze
13"-15"	2.0"
16"-19"	2.5"
20"	3.0"

4. REAR WHEELS / SEAT HEIGHT

Rear Seat Height							
☐	170RH13	14.0"	NC	☐	170RH19	17.0"	NC
☐	170RH14	14.5"	NC	☐	170RH26	17.5"	NC
☐	170RH15	15.0"	NC	☐	170RH20	18.0"	NC
☐	170RH16	15.5"	NC	☐	170RH27	18.5"	NC
☐	170RH17	16.0"	NC	☐	170RH21	19.0"	NC
☐	170RH18	16.5"	NC	☐	170RH37	19.5"	NC
				☐	170RH22	20.0"	NC

Axle			
☐	170AX1	Quick Release Steel	NC
☐	170AX11	Quad Quick Release	\$144.00

NEW! Rear Wheel Type			
Max weight capacity on Proton wheels is 220 lbs.			
☐	170RW103	20" Mag Wheel 5 Spoke	NC
☐	170RW104	22" Mag Wheel 5 Spoke	NC
☐	170RW102	24" Mag Wheel 5 Spoke	NC
☐	170RW158	20" Lite Spoke Wheel	NC
☐	170RW159	22" Lite Spoke Wheel	NC
☐	170RW92	24" Lite Spoke Wheel	NC
☐	170RW93	25" Lite Spoke Wheel	NC
☐	170RW160	26" Lite Spoke Wheel	NC
☐	170RW95	24" Spinergy SPOX Black <i>(silver hub)</i>	\$985.00
☐	170RW96	25" Spinergy SPOX Black <i>(silver hub)</i>	\$985.00
☐	170RW98	24" Spinergy LX Black <i>(silver hub)</i>	\$1,237.00
☐	170RW99	25" Spinergy LX Black <i>(silver hub)</i>	\$1,237.00
☐	170RW148	24" Proton Wheel	\$234.00
☐	170RW149	25" Proton Wheel	\$234.00
☐	170RW25	Omit Rear Wheel and Axle *	-\$55.00

Spinergy Spoke Colors (Optional)	
Please select one of the following Spinergy wheel options first. Includes black hub.	
Spinergy Spox spokes color:	
☐ Black	☐ Blue
☐ Red	☐ White
☐ Yellow	
Spinergy LX spokes color:	
☐ Black	☐ Blue
☐ Green	☐ Orange
☐ Pink	
☐ Purple	☐ Red
☐ White	☐ Yellow

Rear Wheel Tire			
☐	170RT1	Pneumatic	NC
☐	170RT24	E2211 Pneumatic Airless Insert - Infinergy® E2213 Lightweight - <i>E-TPU material</i>	\$114.00
☐	170RT3	Full Poly	NC
☐	170RT38	E2218 SHOX Performance Poly	\$110.00
☐	170RT25	E2218 Primo Orion Poly	\$126.00
☐	170RT39	E2211 High-Pressure Clincher Black E2212	\$110.00
☐	170RT21	E2211 Schwalbe Marathon Plus Evolution E2212	\$164.00

5. WHEEL LOCKS

Wheel Lock			
No Gos may exist depending on wheel size and position. Ergo & compact N/A w/ 13"-14" frame depth			
☐	170WL1	High Mount Push	NC
☐	170WL2	High Mount Pull	NC
☐	170WL45	Ergo Scissor	NC
☐	170WL88	Fold Down Push (Knee lever)	\$86.00
☐	170WL30	Compact	\$164.00
☐	170WL10	Omit	NC

Wheel Lock Option			
☐	170WL12	Do Not Mount	NC

Wheel Dia.	RSH Range	
	Min Squeeze (0° seat angle)	Max Squeeze (9° seat angle)
20"	14.0"-18.0"	14.0"-17.5"
22"	14.0"-19.0"	14.0"-18.5"
24"	15.0"-19.5"	14.5"-19.0"
25"	15.5"-20.0"	15.0"-19.5"
26"	16.0"-20.0"	15.5"-20.0"

Handrim			
See page 6 for matrices.			
☐	170HR1	Aluminum Anodized	NC
☐	170HR29	Aluminum Powder coated colored	\$82.00
☐	170HR4	Plastic Coated	\$129.00
☐	170HR22	Natural Fit - Standard Grip	\$549.00
☐	170HR23	Natural Fit - Super Grip	\$611.00
☐	170HR28	Q-Grip Neoprene Coated	\$408.00
☐	170HR39	Ellipse 3R	\$500.00
☐	170HR40	Carbolife - Gekko	\$535.00
☐	170HR41	Carbolife - Curve L Tetra grip	\$575.00
☐	170HR7	Omit Handrim <i>Option may not be available on all wheels</i>	NC

Camber			
☐	170CI1	0°	NC
☐	170CI13	2° <i>(default)</i>	NC
☐	170CI3	4°	NC

Center of Gravity Preset			
☐	170CC20	1.0" <i>(default)</i>	NC
☐	170CC15	1.5"	NC
☐	170CC17	2.0"	NC
☐	170CC18	2.5"	NC
☐	170CC19	3.0"	NC

Extension Handle			
N/A w/ ergo scissor or compact wheel locks			
☐	170EH7	E0961 Extension Handle (Pair)	\$88.00

Anti Roll-Back			
	new	E2228 Automatic Bracking system - On chair	\$2,075.00
PRATIKO Automatic Braking System. Installed by Sunrise Medical. Wheel Lock selection required - If no selection, Push to Lock will be provided. Weight capacity of 250 lbs (113 kg). Not available with transit. Available with the following rear wheels: 24" Lite Spoke and XCore Mags. Folding width maybe be affected.			

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price											
6. HANGER & FOOTREST																				
Swing Away Hangers					Footrest Extension Tubes															
Swing Awat frame only					See Extension matrix for leg length range															
☐	170H145		60° Swing In/Out Hanger	NC	☐	170E1		Short	NC											
☐	170H139		70° Swing In/Out Hanger	NC	☐	170E2		Medium	NC											
☐	170H141		80° Swing In/Out Hanger	NC	☐	170E5		Long	NC											
☐	170H160	E0990	Elevating Hanger - Ext Mount.	\$367.00																
☐	170H159	K0053	Articulating Hanger - Ext Mount	\$431.00																
☐	170H32		Omit Hanger and Footplate	-\$15.00																
Footplates																				
Heel loop or leg strap not included, See Footplate Options Section. Platform footplates flip up (sideways).																				
☐	170F1		Composite	STD																
☐	170F271	K0040 E0951	Composite Angle Adj-(16"-22"w)	\$211.00																
☐	170F5		Aluminum (ELR & ALR)	NC																
☐	170F3	K0040	Alum Adult Angle Adjustable - Ext Mour	\$221.00																
☐	170F259	K0040	Multi Angle Adustable ☐ Black NC ☐ Std Frame color\$51	\$221.00																
☐	170F50	K0040	Platform Footrest ☐ Black NC ☐ Std Frame color\$51	\$110.00																
☐	170F335	K0040	Platform Ultra Lite Footrest	\$134.00																
☐	170F336		Platform Carbon Fiber footrest	\$158.00																
☐	170F7		Omit Footplate and Ext. Tube	-\$15.00																
Lower Leg Length (LLL)																				
Lower Leg Length is measured from seat sling to the top rear of the footrest . See Appendix on Page 7. A minimum 2" ground clearance is required below the footrest. If the selected configuration results in less than 2" of ground clearance the footrest will be raised to 2" ground clearance, which may alter the lower leg length measurement. Available lower leg length is dependent on frame selection including adapter spacers. Reference LLL matrix in page 7																				
☐	10"	☐	11"	☐	12"	☐	13"	☐	14"	☐	15"	☐	16"	☐	17"	☐	18"	☐	19"	
7. BACKREST / SEATING																				
Back Type										Backrest Upholstery										
☐	170BT123		Standard Fix (90°)	NC	EXO standard and Tension adjustable upholstery binding is always black .															
☐	170BT47		Angle Adjustable Back (80° - 105°)	NC	☐	170BU95		Standard Upholstery EXO	NC											
Select desired setting				☐ 80° ☐ 85° ☐ 90° ☐ 100°	Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue \$26.00															
☐	170BT124		Fixed Half Fold down (to the front)	\$216.00	☐	170BU70	E2611	Backrest Tension Adjust Ballistic <i>*Reference 186BU70 for billing</i>	\$367.00											
☐	170BT49		Lumbar curve bend	NC	☐	170BU89	E2611	Backrest Tension Adj EXO <i>*Reference XXXBU89 for billing</i>	\$387.00											
Select desired setting				☐ 5° forward ☐ 5° backward	Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue \$26.00															
Back Height Setting										☐	170BU90	E2611	Backrest Tension Adj EXO Pro 3DX <i>*Reference 186BU71 for billing</i>	\$387.00						
Half Fold Down N/A with the 10" - 13" back heights. Lumbar curve backrest with integral push handles available with 12" thru 20" back heights. Lumbar curve backrest with Fold Down or Omit push handles available with 14" thru 20" back heights.										Color Binding: ☐ Silver ☐ Orange ☐ Red ☐ Blue \$26.00										
☐	170B161	Setting	10 "- 12" Back Height ☐ 10" ☐ 11" ☐ 12"			NC	☐ 170BU14 Omit Back Upholstery -\$15.00													
☐	170B162	Setting	12 "- 14" Back Height ☐ 12" ☐ 13" ☐ 14"			NC														
☐	170B163	Setting	14 "- 16" Back Height ☐ 14" ☐ 15" ☐ 16"			NC														
☐	170B164	Setting	16 "- 18" Back Height ☐ 16" ☐ 17" ☐ 18"			NC														
☐	170B165	Setting	18 "- 20" Back Height ☐ 18" ☐ 19" ☐ 20"			NC														
Push Handles																				
Std EXO Uph N/A with Fold Down. Ht Adj Push Handles only available with tension adjustable EXO / EXO Pro upholstery.																				
☐	170PH6		Integral Push Handles	NC																
☐	170PH5		Fold Down Push Handles <i>Available w/12"-20" back height</i>	\$116.00																
☐	170PH21		Height adjustable push Handles <i>Available w/ 12"-20" back height. N/A w/Half fold down back</i>	\$175.00																
☐	170PH25		Height adjustable push Handles for JAY Backs <i>N/A w/Half fold down back</i>	\$175.00																
☐	170PH22		Omit Push Handles - <i>Straight Post</i>	NC																
JAY Back (Optional)																				
To order a J3 Back please select Omit Back Upholstery. See J3 Back Order Form for Part #'s																				
☐	170JB01		J3 Back (See J3 Order Form)	Price List																

8. ARMREST / SIDE GUARDS

NEW!

Padded Swing Away

Some limitations when used w/ angle adj back in the rear most setting. 10"-12" height ranges N/A w/Swing away anti tip. Please select a setting, will default to 11" if no selection is made.

☐	170AA162		Padded Swing-Away Fixed height (13")	NC
☐	170AA18	E0973	Padded Swing-Away STD (9" to 14")	\$263.00
☐	170AA177	E0973	Padded Swing-Away with Side guard Standard (9" to 14")	\$432.00
Setting ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14"				

Armrest - Height Adjustable

☐	170AA11	E0973	Single Post Hgt. Adj. Std (8.5"-13.5")	\$333.00
☐	170AA29	E0973	Single Post Hgt. Adj. Low (7"-11")	\$333.00
☐	170AA65	E0973	Flip Back Height Adj <i>LITE</i> (8.5"-13.5")	\$367.00

Armpad - For Height Adjustable

Ergo armpad only available with Single post Height adj & Flip back armrest

☐	170AM23		Desk Length - Classic (10")	NC
☐	170AM24		Full Length - Classic (14")	NC
☐	170AM27		Ergo Pad (9.5"-13.5")	NC

Side Guard

Side guards N/A w/ height adj. arms. Side guard interference w/ Padded SA arms at some heights may exist. Aluminum & carbon side guards some limitations with low RSH and large rear wheels

☐	170SG33	K0108	Plastic Standard	\$169.00
☐	170SG29	K0108	Plastic Low	\$169.00
☐	170SG35	K0108	Aluminum	\$286.00
☐	170SG36	K0108	Aluminum w/ Fender	\$321.00
☐	170SG12	K0108	Carbon Fiber	\$537.00
☐	170SG41	K0108	Carbon Fiber Long w/LG Fender	\$607.00

Side Guard - Color Option

Available with Aluminum or Aluminum with Fender. Not applicable to carbon fiber

☐	170GE01		Black	NC
☐	170GE02		Matching Frame- Standard Colors	NC
☐	170GE03		Matching Frame- Personalized Color	\$58.00

Omit Armrest

☐	170AA14		Omit Armrest	NC
--------------------------	---------	--	--------------	----

9. KEY OPTIONS & ACCESSORIES

Anti-Tip

☐	170AT1	E0971	Anti Tips (Pair)	\$175.00
☐	170AT	E0971	Anti Tip swing away <i>N/A Transit wheels</i>	
			☐ 170AT88 Right Mount	\$102.00
			☐ 170AT89 Left Mount	\$102.00
			☐ 170AT90 Pair	\$204.00

Additional Accessories

☐	170A1		Backpack - Adult	\$88.00
☐	170A2		Backpack - Kids	\$70.00
☐	1701A01		Seat Pouch Black	\$58.00

Additional Cover Accessories

☐	1708A02		Gel Cover Hanger/Frame	\$129.00
☐	1708A01		Gel Cover Swing Away Arm	\$175.00

Positioning Belts

☐	170ZZ15	E0978	JAY 2PT Belt 1.5" Push Btn 7.5" Pad	\$101.00
☐	170ZZ26	E0978	JAY 2PT Belt 2" Push Btn 9.5" Pad	\$101.00
☐	170ZZ14	E0978	JAY 2PT Belt 1.5" Push Btn No Pad	\$48.00
☐	170A118	E0978	Belt 2" Aircraft Buckle Pad w/Clamp	\$139.00
☐	170A120		Positioning Belt Slide Clamps	\$58.00

Spoke Guards

Only available w/ Lite Spoke wheels.

☐	1704A01	K0065	Quickie Black Spoke Guards	\$134.00
☐	1704A02	K0065	Quickie Clear Spoke Guards	\$134.00

Miscellaneous

Transit wheels N/A with Swing away anti tip.

☐	1708A03		Caddy	\$58.00
☐	170A2100	E2207	Crutch Holder	\$227.00
☐	1708A21		Tool Kit	\$63.00
☐	170A999		Quickie Water Bottle & Holder	\$63.00
☐	1708A47		Mobile phone pocket	\$30.00
☐	170FO4	E0971	Transit wheels (Pair)	\$122.00
☐	170AT	E0971	Tip Assist	
			☐ 170AT83 Right Mount	\$36.00
			☐ 170AT84 Left Mount	\$36.00
			☐ 170AT91 Pair	\$71.00

10. COLOR

Personalized Crossbrace Color

Option has en extended leadtime. Frame colors only, NA with Personalized frame colors.
If no color pick crossbrace will be Gloss Black.

☐	170GS15	Crossbrace Color package	\$250.00
--------------------------	---------	--------------------------	----------

Standard Frame Colors

Handrim color option only available if Aluminum Coated colored option was selected.

			Frame	X-brace	Handrim
170C214		Black (Gloss)	☐	STD	☐
170C215		Black (Matte)	☐	☐	☐
170C211		Celestial Sapphire	☐	☐	☐
170C212		Desert Sage	☐	☐	☐
170C213		Desert Tan	☐	☐	☐
170C216		Metallic Shadow	☐	☐	☐
170C217		Midnight Blue	☐	☐	☐
170C218		Mystic Purple	☐	☐	☐
170C219		Petal Pink	☐	☐	☐
170C220		Pink Power	☐	☐	☐
170C221		Red Pepper	☐	☐	☐
170C222		Royal Blue	☐	☐	☐
170C223		Salute Green	☐	☐	☐
170C224		Scarlet Noir	☐	☐	☐
170C225		Silver Ore	☐	☐	☐
170C226		Solar Orange	☐	☐	☐
170C227		Starlight White	☐	☐	☐
170C228		Storm Gray	☐	☐	☐
170C229		Turbo Green	☐	☐	☐
170C230		Velocity Yellow	☐	☐	☐

Built For Me

To inquire about any unique personalized for Quickie products, contact:
Phone: (800) 333 - 4000, Option #3, Fax: (800) 333 9011



Personalized Frame Color

Only frame will have the personalized color; Crossbrace, handrim, footplate and Aluminum side guard not included (unless selected).

☐	170C112	Desert Camouflage	\$460.00
☐	170C43	Mossy Oak Camouflage	\$460.00
☐	170C42	Stars & Stripes	\$460.00
☐	170C108	Zebra (Choose above base color) Std frame colors with White unless above base paint is chosen N/A w/ Black & Matte Black)	\$460.00

NEW! Legacy - Frame Colors

☐	170C141	Black Cherry	NC
☐	170C10	Black Opal	NC
☐	170C14	Candy Purple	NC
☐	170C142	Electric Blue	NC

Frame Decal Colors

Applies to front caster link decals and aluminum & carbon fiber sideguards.

☐	170D16	Black	NC
☐	170D17	Blue	NC
☐	170D18	Orange	NC
☐	170D19	Red	NC
☐	170D20	Silver	NC

Anodized Color Package

Axle plate, Multi Forks and Aluminum Soft Roll casters only. N/A with Frog suspension Legs.

☐	170EB08	Blue	\$187.00
☐	170EB13	Gold	\$187.00
☐	170EB09	Lime Green	\$187.00
☐	170EB10	Purple	\$187.00
☐	170EB11	Red	\$187.00
☐	170EB12	Ti Grey	\$187.00

11. APPENDIX

Wheel - Handrims Matrix				
Handrims	Wheels			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Aluminum Ano / Painted	20,22,24	20,22,24,25,26	24,25	24,25
Plastic Coat	20,22,24	22,24	24,25	-
Natural Fit - Standard Grip	24	22,24,25,26	24,25	-
Natural Fit - Super Grip	24	22,24,25,26	24,25	-
Q Grip Neoprene Coated	-	24,25	24,25	-
Ellipse 3R	24	24,25	24,25	24,25
Carbolife - Gekko	-	24,25	24,25	24,25
Carbolife - Curve L Tetra grip	-	24,25	24,25	24,25
Omit Handrim	20,22,24	-	24,25	-

Wheel - Tire Matrix				
Tire	Rear Wheel			
	5 Spoke Mag	Lite Spoke	Spinergy	Proton
Pneumatic	20,22,24	20,22,24,26	24	24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	24	-
Full Poly	20,22,24	20,22,24	24	-
SHOX Performance Poly	20,22,24	20,22,24,25	24,25	-
Primo Orion Poly	20,22,24	20,22,24	24	-
High-Pressure Clincher Blk	24	22,24,26	24	24
Schwalbe Marathon	22, 24	24,25,26	24,25	24,25
20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm				

11. APPENDIX

Note 1: These Front Seat Height (FSH) matrices contain rounded estimates and assume flat seat heights (0 deg seat angle). Increasing seat angle can affect FSH setting.

Note 2: Highlighted cells indicate all valid combinations the Configurator will possibly choose based on FSH and Caster selection.

Note 3: Blank spaces (denoted with a dash "-") indicate at least one of two things: A) The caster won't fit in that fork hole, and/or B) The fork hole may not exist.

Note 4: For all forks except Frog Legs, Hole Position 1+ indicates that a 10mm spacer is used above the fork bearings to achieve additional FSH. For all other hole positions, the 10mm spacer sits below the fork bearings, except for the Frog Legs which can only have the spacer on the top side of the fork.

Note 5: The lowest FSH 14.0" is only available with the 4" Low Profile Poly Caster (Multi 90 hole 2 is not obtainable with 4"x1.25" Semi Pneu or 4"x1.5" Alum Soft Roll).

FSH for 4" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	14.0	14.5	15.0
	Multi 115	-	14.5	15.0	15.5	16.0
	Multi 131	-	15.0	15.5	16.0	16.5
	Multi 158	-	16.0	16.5	17.0	17.5
	SS Low	-	-	14.0	14.5	15.0
	SS Std	-	14.5	15.0	15.5	16.0
	Frog Legs	-	15.5	16.0	16.5	-
Short	Multi 90	-	-	15.5	16.0	16.0
	Multi 115	-	16.0	16.5	17.0	17.0
	Multi 131	-	16.5	17.0	17.5	18.0
	Multi 158	-	17.5	18.0	18.5	19.0
	SS Low	-	-	15.5	16.0	16.0
	SS Std	-	16.0	16.5	17.0	17.0
	Frog Legs	-	17.0	17.5	18.0	-
Tail	Multi 90	-	-	16.5	17.0	17.5
	Multi 115	-	17.0	17.5	18.0	18.5
	Multi 131	-	17.5	18.0	18.5	19.0
	Multi 158	-	18.5	19.0	19.5	20.0
	SS Low	-	-	16.5	17.0	17.5
	SS Std	-	17.0	17.5	18.0	18.5
	Frog Legs	-	18.0	18.5	19.0	-

FSH for 5" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	15.5	16.0	16.0
	Multi 131	15.0	15.5	16.0	16.5	17.0
	Multi 158	16.0	16.5	17.0	17.5	18.0
	SS Low	-	-	-	-	-
	SS Std	-	15.0	15.5	16.0	16.0
	Frog Legs	-	-	16.5	17.0	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	16.5	17.0	17.5
	Multi 131	16.5	17.0	17.5	18.0	18.0
	Multi 158	17.5	18.0	18.5	19.0	19.5
	SS Low	-	-	-	-	-
	SS Std	-	16.0	16.5	17.0	17.5
	Frog Legs	-	-	17.5	18.0	-
Tail	Multi 90	-	-	-	-	-
	Multi 115	-	-	18.0	18.5	19.0
	Multi 131	17.5	18.0	18.5	19.0	19.5
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
	SS Std	-	17.5	18.0	18.5	19.0
	Frog Legs	-	-	19.0	19.5	-

FSH for 6" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	16.0	16.5
	Multi 131	-	16.0	16.5	17.0	17.0
	Multi 158	16.5	17.0	17.5	18.0	18.5
	SS Low	-	-	-	-	-
	SS Std	-	-	15.5	16.0	16.5
	Frog Legs	-	-	-	17.0	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	17.5	18.0
	Multi 131	-	17.5	18.0	18.5	18.5
	Multi 158	18.0	18.5	19.0	19.5	19.5
	SS Low	-	-	-	-	-
	SS Std	-	-	17.0	17.5	18.0
	Frog Legs	-	-	-	18.5	-
Tall	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	19.0	19.5
	Multi 131	-	18.5	19.0	19.5	20.0
	Multi 158	19.0	19.5	20.0	20.5	21.0
	SS Low	-	-	-	-	-
	SS Std	-	-	18.5	19.0	19.5
	Frog Legs	-	-	-	20.0	-

FSH for 7" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	17.0	17.5	18.0
	Multi 158	17.0	17.5	18.0	18.5	19.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	18.5	19.0	19.0
	Multi 158	18.5	19.0	19.5	20.0	20.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Tall	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	19.5	20.0	20.5
	Multi 158	19.5	20.0	20.5	21.0	21.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-

FSH for 8" Casters						
Caster Link	Fork Size	Hole Position				
		4	3	2	1	1+
Hemi	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	17.5	18.0	18.5
	Multi 158	-	18.0	18.5	19.0	19.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Short	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	19.0	19.5	19.5
	Multi 158	-	19.5	20.0	20.5	20.5
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-
Tall	Multi 90	-	-	-	-	-
	Multi 115	-	-	-	-	-
	Multi 131	-	-	20.0	20.5	21.0
	Multi 158	-	20.5	21.0	21.5	22.0
	SS Low	-	-	-	-	-
	SS Std	-	-	-	-	-
	Frog Legs	-	-	-	-	-

Lower Leg Length Matrix				
Note for Swing Away Hangers only: subtract 1" from the max achievable LLL for long extension tubes if the Multi Angle Adjustable footplate is selected.				
Hanger	Extension	LLL (inches) (measured from top of seat sling to top of footrest)		Min Allowable FSH (due to 2" ground clearance requirement)
		Min	Max	
60° Swing In/Out	Short	10	14	No lower FSH limit
	Med	11	16	
	Long	12	18	
70° Swing In/Out	Short	10	16	No lower FSH limit
	Med	12	17	
	Long	13	19	
80° Swing In/Out	Short	11	15	No lower FSH limit
	Med	12	17	
	Long	14	19	
ELR	Short	11	16	No lower FSH limit
	Med	11	18	
	Long	13	19	
ALR	Short	11	16	No lower FSH limit
	Med	13	18	
	Long	15	19	
75° Fixed Front	Short	13	18	15.5
	Med	13	20	
	Long	14	20	
85° Fixed Front 0" inset	Short	13	18	15.5
	Med	13	20	
	Long	15	20	
85° Fixed Front 1" inset	Short	13	18	15.5
	Med	15	20	
	Long	17	20	

Back Height Go/No-Go Table																			
Back Type	Push Handle Type	Back Height Range																	
		10-12			12-14			14-16			16-18			18-20					
		10	11	12	12	13	14	14	15	16	16	17	18	18	19	20			
Standard / Angle Adj	Integral	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO			
	Fold Down	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Ht Adj (reg & JAY)	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Omit PH	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO	GO			
Lumbar (FWD or REV)	Integral	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Fold Down	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Ht Adj (reg & JAY)	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Omit PH	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
Half Fold	Integral	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Fold Down	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Ht Adj (reg & JAY)	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			
	Omit PH	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO	NO-GO			

Back Upholstery Go/No-Go Table				
Push Handle Type	Back Upholstery Type			
	STD EXO (BU95)	TA EXO (BU89 & BU90)	TA Ballistic (BU70)	Omit (BU14)
Integral	GO	GO	GO	GO
Fold Down	NO-GO	GO	GO	GO
Ht Adj reg (PH21)	NO-GO	GO	NO-GO	GO
Ht Adj JAY (PH25)	NO-GO	NO-GO	NO-GO	GO
Omit PH	GO	GO	GO	GO