



Squiggles Plus Stander

September 2020

Mark For:

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: _____ Fax: _____

Confirmation Email: _____

Confirm Via: Fax Email**Submitting for:** Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: _____ Fax: _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Squiggles + Stander**HCPCS: E0641***Weight Capacity: 48.4 lbs***Packages**

Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps and Product Assembly.

	Item Code	Description	HCPCSs	Retail
<input type="checkbox"/>	PSPSTB	Squiggles Plus Stander Package - Blue	E0641	\$6,102
<input type="checkbox"/>	PSPSTP	Squiggles Plus Stander Package - Purple	E0641	\$6,102
<input type="checkbox"/>	PSPSTR	Squiggles Plus Stander Package - Red	E0641	\$6,102
<input type="checkbox"/>	PSPSTRV	Squiggles Plus Stander Package - Red Vinyl	E0641	\$6,102

Standard Model

Selection Required. Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplates and tray attachment.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	184-600	Standard Model with Support Shell		\$2,805
<input type="checkbox"/>	184-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning		\$1,372

Positioning Pads

Selection Required. Includes chest and hip pad covers and knee support covers.

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	184-1800-02	Squiggles Plus Stander Positioning Pads - Blue		\$623
<input type="checkbox"/>	184-1800-04	Squiggles Plus Stander Positioning Pads - Purple		\$623
<input type="checkbox"/>	184-1800-09	Squiggles Plus Stander Positioning Pads - Red		\$623
<input type="checkbox"/>	184-1800-09VIN	Squiggles Plus Stander Positioning Pads - Red Vinyl		\$623

Accessories

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	118-619	Tray - Gray	E0950	\$371
<input type="checkbox"/>	120-270	Padded Tray Insert - Gray	K0108	\$300
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray	K0108	\$270
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$241
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$263

Headrest

Cushion selection required with hardware

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	E1028	\$209
<input type="checkbox"/>	184-1812-02	Flat Headrest Cushions - Blue	E0955	\$110
<input type="checkbox"/>	184-1812-04	Flat Headrest Cushions - Purple	E0955	\$110
<input type="checkbox"/>	184-1812-09	Flat Headrest Cushions - Red	E0955	\$110
<input type="checkbox"/>	184-1812-09VIN	Flat Headrest Cushions - Red Vinyl (Pair)	E0955	\$110
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) <i>Requires Flat Headrest Support Hardware and Cushion</i>	E1028	\$192
<input type="checkbox"/>	184-1814-02	Flat Headrest and Lateral Cushions - Blue	E0956	\$193
<input type="checkbox"/>	184-1814-04	Flat Headrest and Lateral Cushions - Purple	E0956	\$193
<input type="checkbox"/>	184-1814-09	Flat Headrest and Lateral Cushions - Red	E0956	\$193
<input type="checkbox"/>	184-1814-09VIN	Flat Headrest and Lateral Cushions - Red Vinyl (Pair)	E0956	\$193
<input type="checkbox"/>	184-1838-02	Spine Cap Cover - Blue		\$34
<input type="checkbox"/>	184-1838-04	Spine Cap Cover - Purple		\$34
<input type="checkbox"/>	184-1838-09	Spine Cap Cover - Red		\$34
<input type="checkbox"/>	184-1838-09VIN	Spine Cap Cover - Red Vinyl		\$34

Assembly

	Item Code	Description	Retail
<input type="checkbox"/>	ASSY-LKY	Product Assembly Fee	\$120

Technical Specifications

Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Distance Between Midline of Knees	5.5" - 8.4"	140 mm - 215 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.5" - 31.5"	480 mm - 800 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 mm L x 390 mm W x 390 mm H
Pivot Base Weight	14.3 lbs	6.5 kgs
Pivot Base Dimensions (Unfolded)	37.4" L x 23.2" W x 20.3" H	950 mm L x 590 mm W x 515 mm H
Pivot Base Dimensions (Folded)	26" L x 23.2" W x 20.3" H	660 mm L x 590 mm W x 515 mm H
Pivot Base Angle Range	Prone 90-70° / Supine 90-160°	Prone 90-70° / Supine 90-160°



Sunrise Medical (US) LLC2842 Business Park Ave. • Fresno, CA 93727 • USA