



# MyWay

## August 2020

**Mark For:**

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via:  Fax  Email

**Submitting for:**
 Quote  Order

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**MyWay**
**HCPCS: E8000**
**Packages**

*Includes Frame, Harness, User Handles, Guide Handles, Headrest, Saddle and Ankle Guides. Packages for frame sizes 1 and 2 include Size 1 Saddle; Packages for frame size 3 include Size 2 Saddle.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	PMWS1P	Size 1 MyWay Package with Size 1 Harness - Purple		\$5,171
<input type="checkbox"/>	PMWS1G	Size 1 MyWay Package with Size 2 Harness - Green		\$5,175
<input type="checkbox"/>	PMWS1B	Size 1 MyWay Package with Size 3 Harness - Blue		\$5,204
<input type="checkbox"/>	PMWS2B	Size 2 MyWay Package with Size 3 Harness - Blue		\$5,765
<input type="checkbox"/>	PMWS2R	Size 2 MyWay Package with Size 4 Harness - Red		\$5,799
<input type="checkbox"/>	PMWS2O	Size 2 MyWay Package with Size 5 Harness - Orange		\$5,836
<input type="checkbox"/>	PMWS3R	Size 3 MyWay Package with Size 4 Harness - Red		\$5,982
<input type="checkbox"/>	PMWS3O	Size 3 MyWay Package with Size 5 Harness - Orange		\$6,019
<input type="checkbox"/>	PMWS3G	Size 3 MyWay Package with Size 6 Harness - Gray		\$6,035

**Frame**

*Selection Required. Includes multi-functional casters.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	173-1600	Size 1 MyWay Frame		\$3,628
<input type="checkbox"/>	173-2600	Size 2 MyWay Frame		\$4,184
<input type="checkbox"/>	173-3600	Size 3 MyWay Frame		\$4,357

**Harness**

*Selection Required.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	173-4100-04	Size 1 Harness - Purple		\$329
<input type="checkbox"/>	173-4200-07	Size 2 Harness - Green		\$333
<input type="checkbox"/>	173-4300-02	Size 3 Harness - Blue		\$362
<input type="checkbox"/>	173-4400-09	Size 4 Harness - Red		\$397
<input type="checkbox"/>	173-4500-06	Size 5 Harness - Orange		\$433
<input type="checkbox"/>	173-4600-08	Size 6 Harness - Gray		\$450

**Handles***Must match Frame Size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	173-1667	Size 1 User Handles (Pair)		\$169
<input type="checkbox"/>	173-2667	Size 2 User Handles (Pair)		\$173
<input type="checkbox"/>	173-3667	Size 3 User Handles (Pair)		\$180
<input type="checkbox"/>	173-1691	Size 1 Guide Handles		\$290
<input type="checkbox"/>	173-2691	Size 2 Guide Handles		\$290
<input type="checkbox"/>	173-3691	Size 3 Guide Handles		\$290

**Positioning Accessories***Ankle Guides must match Frame Size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	173-1680	Headrest (One Size)		\$205
<input type="checkbox"/>	173-1671	Size 1 Saddle		\$321
<input type="checkbox"/>	173-2671	Size 2 Saddle		\$324
<input type="checkbox"/>	173-1700	Size 1 Ankle Guides (Pair)		\$229
<input type="checkbox"/>	173-2700	Size 2 Ankle Guides (Pair)		\$229
<input type="checkbox"/>	173-3700	Size 3 Ankle Guides (Pair)		\$229

**Other Accessories***Connection belt available only with Size 2 or Size 3 Frames. Selection of Guide Handles required with Guide Connection Belt.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	173-1693	Guide Connection Belt		\$104
<input type="checkbox"/>	173-1800	Hoist Straps (Pair)		\$55

## Technical Specifications

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years
User Shoulder Height (to Floor)	25" - 36"	34" - 47"	42" - 59"
User Saddle Height (Inside Leg)	8.5" - 17.5"	16" - 24"	21.5" - 32"
Maximum User Weight	55 lbs	110 lbs	176 lbs
Frame Weight	26.5 lbs	14.5 lbs	21 lbs
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do no separate		23 lbs
Frame Width	22.5"	26.5"	28.5"
Frame Length	28.5"	33.5"	43.5"
Saddle	9"L x 3.25"W	10"L x 3.25"W	

  

Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	18" - 22"	19.5" - 23.5"	21.5" - 25.5"	23.5" - 27.5"	27.5" - 31.5"	31.5" - 35.5"
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	110 lbs	110 lbs	110 lbs	176 lbs	176 lbs	176 lbs

## Technical Specifications (Metric)

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 11 years	9 - 16 years
User Shoulder Height (to Floor)	635 mm - 920 mm	870 mm - 1190 mm	1070 mm - 1490 mm
User Saddle Height (Inside Leg)	220 mm - 440 mm	400 mm - 610 mm	550 mm - 810 mm
Maximum User Weight	25 kgs	50 kgs	80 kgs
Frame Weight	12 kgs	6.5 kgs	9.5 kgs
Upper Support Unit Weight	Size 1 Frame and Upper Support Unit do no separate		10.5 kgs
Frame Width	570 mm	670 mm	720 mm
Frame Length	720 mm	850 mm	1100 mm
Saddle	225 mm L x 80 mm W	255 mm L x 80 mm W	

  

Harness Sizing Specification	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6
Age (approximate)	1 - 4 years	2 - 6 years	5 - 9 years	8 - 12 years	12 - 15 years	16 - 18 years
Waist at Navel Level	460 mm - 560 mm	500 mm - 600 mm	550 mm - 650 mm	600 mm - 700 mm	700 mm - 800 mm	800 mm - 900 mm
Color	Purple	Green	Blue	Red	Orange	Gray
Maximum User Weight	50 kgs	50 kgs	50 kgs	80 kgs	80 kgs	80 kgs



Sunrise Medical (US) LLC 2842 Business Park Ave. • Fresno, CA 93727 • USA